

Two-Bin Mobile Phlebotomy Stand

MODEL 6902

- ClintonClean™ solid plastic top
- Chrome plated steel pole with knob-lock adjustment
- 5-leg, cast aluminum base
- 2" dual wheel, soft roll casters
- Two 10-3/4" L x 8-1/4" W x 3-1/2" D plastic bins (27 cm L x 21 cm W x 9 cm D)

NOTE: *Supplies pictured not included.*


All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



MODEL 6902 SPECIFICATIONS

Top	13.5"L x 21"W
Base	23" D
Height Range	30"— 38"

Prop 65 Warning

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1