

Half Round, Stationary, Padded Phlebotomy Stand

MODEL 6940

- Vinyl upholstered, padded surface
- Easy-clean, 10-3/4" x 8-1/4" W x 3-1/2" D plastic bin (27cm L x 21cm W x 9cm D)
- Chrome plated steel pole with knob-lock adjustment
- 5-leg, cast aluminum base with black, non-marking, nylon glides

NOTE: *Supplies pictured not included.*

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



MODEL 6940 SPECIFICATIONS

Top	22" L x 11" Radius
Height Range	30" — 38"
Base	23" D

Prop 65 Warning

! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Upholstery Colors

