

Half Round, Stationary, ClintonClean™ Phlebotomy Stand

MODEL 6950

- ClintonClean™ solid plastic top
- Easy-clean, 10-3/4" x 8-1/4" W x 3-1/2" D plastic bin (27cm L x 21cm W x 9cm D)
- Chrome plated steel pole with knob-lock adjustment
- 5-leg, cast aluminum base with black, non-marking, nylon glides

NOTE: *Supplies pictured not included.*

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



MODEL 6950 SPECIFICATIONS

Top	22" L x 11.5" Radius
Height Range	30" — 38"
Base	23" D

Prop 65 Warning



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1