


## Discover Series, Sweet Dreams Candy Factory, Treatment Table



### MODEL 7370-71 (Adjustable Backrest)

### MODEL 7370-71-088 (Flat Top)

- No upcharge for adjustable backrest
- Clinton original themes
- Easy to clean laminate surfaces
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- 2" of firm foam padding
- 4 door storage with 2 adjustable shelves
- Scratch and mar resistant surfaces
- Chip resistant PVC edges
- Premium, stain-resistant, woven, knit-backed vinyl upholstery
- Soft-close, Euro style door hinges
- Graphics are on one long side with gray ends and back
- Top features seamless rounded outside corners
- Paper dispenser included
- Available with or without backrest

#### SPECIFICATIONS (BOTH MODELS)

72"L x 27"W x 33"H


California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1



#### Options



066 Door & Inside Latch Combo



040 Paper Cutter

#### Standard Upholstery



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black