

Imagination Series, Cool Pals, Pediatric Scale Table

MODEL 7831

TABLE FEATURES

- Clinton original illustrations
- Easy to clean laminate surfaces
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Graphics are on one long side and two short ends. Back is slate gray. Only manufactured as pictured.
- Comfortable work height of 36"
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- Choose to access concealed paper dispenser from either end
- Recessed front toe kick
- Easy-clean all laminate base
- Ample 2 door and 1 drawer storage
- Soft-close, door hinges



MODEL 7831 SPECIFICATIONS

43"L x 24"W x 38.5"H

SCALE FEATURES

- Dependable, accurate Doran® scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Minimum registered weight is 5 lbs.—Maximum weight for scale is 40 lbs.


California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Options



066 Door & Inside Latch Combo



055 Drawer Locks



078-S Seal Scale Pal



078-B Bear Scale Pal- Cool Pals

Upholstery



Gunmetal

Features



Infantometer

Built-in infantometer with folding ends for measuring height



Paper Dispenser

Choose to access concealed paper dispenser from either end



Rimmed Top with Beveled Edges

Top is rimmed on three sides for safety and the front and back rim edges are beveled and the corners rounded for safety



Padded Top

Easy-clean, padded, upholstered top for comfort



Scale Head

Easy-read, digital display