

Imagination Series, Rainforest Follies, Scale Table OS

MODEL 78322

TABLE FEATURES

- Clinton original illustrations
- Easy-to-lean laminate surfaces
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Graphics on front only
- Two drawer storage with full open shelf
- Slate gray legs and shelf
- Comfortable work height of 36"
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- Paper dispenser included
- Optional drawer locks available
- Ships knocked-down to save on freight



MODEL 78322 SPECIFICATIONS

43"L x 24"W x 38.5"H

SCALE FEATURES

- Dependable, accurate Doran® scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Minimum registered weight is 5 lbs.—Maximum weight for scale is 40 lbs.

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Options



055 Drawer Locks




78-P Parrot Scale Pal

Upholstery



Gunmetal

Features



Infantometer

Built-in infantometer with folding ends for measuring height



Padded Top

Easy-clean, padded, upholstered top for comfort



Rimmed Top with Beveled Edges

Top is rimmed on three sides for safety and the front and back rim edges are beveled and the corners rounded for safety



Scale Head

Easy-read, digital display