

## Fun Series, Zoo Bus, Pediatric Scale Table



**MODEL 7822**

### TABLE FEATURES

- Clinton original illustrations
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- Built-in concealed paper dispenser and cutter
- Adjustable built-in levelers provide a stable surface
- Easy-clean all laminate base
- Optional storage area locks available
- Ample 2-door storage
- Soft-close, door hinges

### SCALE FEATURES

- Dependable, accurate Doran® scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Minimum registered weight is 5 lbs.—Maximum weight for scale is 40 lbs.



### MODEL 7822 SPECIFICATIONS

Top:	43"W x 24"D
Base:	47"W x 21"D x 36"H

California residents please be advised of the following, pursuant to Proposition 65:

**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1




### Upholstery



Gunmetal

## Features



### Infantometer

Built-in infantometer with folding ends for measuring height



### Padded Top

Easy-clean, padded, upholstered top for comfort



### Rimmed Top with Beveled Edges

Top is rimmed on three sides for safety and the front and back rim edges are beveled and the corners rounded for safety



### Scale Head

Easy-read, digital display