

## Select Series, Pediatric Scale Table

### MODEL 7840

#### TABLE FEATURES

- Comfortable work height of 36"
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- Two drawers and two-door storage
- Built-in concealed paper dispenser and cutter
- Recessed front toe kick
- Easy-clean all laminate base

#### SCALE FEATURES

- Dependable, accurate Doran® scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Minimum registered weight is 5 lbs.—Maximum weight for scale is 40 lbs.



MODEL 7840 SPECIFICATIONS

43"L x 24"W x 38.5"H

California residents please be advised of the following, pursuant to Proposition 65:

**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).



#### Options



066 Door & Inside Latch Combo



055 Drawer Locks



078-S Seal Scale Pal



078-B Bear Scale Pal



078-K Kitty Scale Pal



078-D Doggie Scale Pal




Select Series Graphics

#### Laminate Colors



Maple



Dark Cherry



Gray



Slate Gray

#### Standard Upholstery



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

## Features




### Infantometer

Built-in infantometer with folding ends for measuring height



### Paper Dispenser

Choose to access concealed paper dispenser from either end



### Rimmed Top with Beveled Edges

Top is rimmed on three sides for safety and the front and back rim edges are beveled and the corners rounded for safety



### Padded Top

Easy-clean, padded, upholstered top for comfort



### Scale Head

Easy-read, digital display