



Completely Coordinated.
Completely Cost Effective.
Completely Easy.

SOUTHWESTFEST SCALE COMPLETE



Imagination Series, Southwestfest, Pediatric Scale Table

MODEL 7852

TABLE FEATURES

- Clinton original illustrations
- Easy to clean laminate surfaces
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Graphics are on one long side and two short ends. Back is slate gray. Only manufactured as pictured.
- Comfortable work height of 36"
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- Choose to access concealed paper dispenser from either end
- Recessed front toe kick
- Easy-clean all laminate base
- Ample 2 door and 1 drawer storage
- Soft-close, door hinges




MODEL 7852 SPECIFICATIONS

43"L x 24"W x 38.5"H

SCALE FEATURES

- Dependable, accurate Doran® scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Minimum registered weight is 5 lbs.—Maximum weight for scale is 40 lbs.


California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Options



066 Door & Inside Latch Combo



055 Drawer Locks

Upholstery



Gunmetal

Features



Infantometer

Built-in infantometer with folding ends for measuring height



Paper Dispenser

Choose to access concealed paper dispenser from either end



Rimmed Top with Beveled Edges

Top is rimmed on three sides for safety and the front and back rim edges are beveled and the corners rounded for safety



Padded Top

Easy-clean, padded, upholstered top for comfort



Scale Head

Easy-read, digital display

Southwestfest Cabinets

MODEL 6152-BW

"Southwestfest" Base Cabinets

4 doors with adjustable shelves and 3 drawers.
Shipped in one unit.

"Southwestfest" Wall Cabinets

4 doors with adjustable shelves. Shipped in one unit.

Base & wall cabinets are only sold in sets

- Clinton original illustrations
- Pre-assembled units save time
- Easy-clean all laminate construction
- All graphic are permanently embedded in the high pressure laminate and will not peel off
- Tough, chip-resistant PVC edges
- Assembled with steel fasteners
- 4 adjustable leveling feet on each unit
- 3-way adjustable, soft-close 110° hinges
- Flat panel drawer fronts with metal sides that hold up to 75 lbs.
- Recessed toe kick
- Easy hook and rail installation for upper cabinets
- Slate Gray edgebanding on all edges
- Customizable with 4 options
- Easily moves for relocation
- Drawers can be added to any base unit
- Doors and drawer can be moved to accommodate sink location
- Shown with optional sink and faucet




SPECIFICATIONS

BASE: 60"L x 18"D x 35"H

WALL: 60"L x 12"D x 24"H

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.


1



Options



066 Door & Inside Latch Combo



055 Door Locks




044 Add a Cabinet Drawer



22 Sink & Faucet

10' X 10' FLOOR PLAN

Optional Items in Red



1 Square = 1 foot

NOTE: These are suggested floor plans. Please check local building and fire codes before beginning any construction or renovations.




Walls
314 Imperial Yellow
2009-20 Salsa



Flooring
Rejuvenations™ Timberline®
Maple 37350 you better

OPTIONAL ACCESSORIES



- 15-CC Cactus Wall Sticker (26" x 54")
- 16-CC Joshua Tree Wall Sticker (50" x 64")
- 2135 5-Leg. Pneumatic Stool with Gunmetal upholstery (-3GM)
- C-40 Side Chair in Gunmetal upholstery (-3GM)
- T-70 Clear Labeled Sundry Jars (5 in set)



800-441-9131 • 717-848-3519
Toll-Free Fax 866-522-1233
Fax 717-843-5871
www.clinton-ind.com