

IMAGINATION SERIES TREATMENT TABLES

Dino Days Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7934	length	width	height
	68"	27"	33"
	173 cm	69 cm	84 cm

7934-1 Same as 7934 without adjustable backrest

DINO DAYS COMPLETE—discount package pricing when combining table with coordinating cabinets (see page 45 for cabinets and features)

7934-X Must order using this model # for package discount

7934-1-X Must order using this model # for package discount



Ocean Commotion Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7936	length	width	height
	68"	27"	33"
	173 cm	69 cm	84 cm

7936-1 Same as 7936 without adjustable backrest

OCEAN COMMOTION COMPLETE—discount package pricing when combining table with coordinating cabinets (see page 42 for cabinets and features)

7936-X Must order using this model # for package discount

7936-1-X Must order using this model # for package discount



Alley Cats & Dogs Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7937	length	width	height
	68"	27"	33"
	173 cm	69 cm	84 cm

7937-1 Same as 7937 without adjustable backrest

ALLEY CATS & DOGS COMPLETE—discount package pricing when combining table with coordinating cabinets (see page 46 for cabinets and features)

7937-X Must order using this model # for package discount

7937-1-X Must order using this model # for package discount



Serengeti Sunrise Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7951	length	width	height
	68"	27"	33"
	173 cm	69 cm	84 cm

7951-1 Same as 7951 without adjustable backrest

SERENGETI SUNRISE COMPLETE—discount package pricing when combining table with coordinating cabinets (see page 46 for cabinets and features)

7951-X Must order using this model # for package discount

7951-1-X Must order using this model # for package discount



Southwestfest Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7952	length	width	height
	68"	27"	33"
	173 cm	69 cm	84 cm

7952-1 Same as 7952 without adjustable backrest

SOUTHWESTFEST COMPLETE—discount package pricing when combining table with coordinating cabinets (see page 46 for cabinets and features)

7952-X Must order using this model # for package discount

7952-1-X Must order using this model # for package discount



Kangaroo Country Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7953	length	width	height
	68"	27"	33"
	173 cm	69 cm	84 cm

7953-1 Same as 7953 without adjustable backrest

KANGAROO COUNTRY COMPLETE—discount package pricing when combining table with coordinating cabinets (see page 46 for cabinets and features)

7953-X Must order using this model # for package discount

7953-1-X Must order using this model # for package discount



Alphabet Animals Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7955	length	width	height
	68"	27"	33"
	173 cm	69 cm	84 cm

7955-1 Same as 7953 without adjustable backrest



DOWNLOAD A FUN FACT SHEET ABOUT THE ANIMALS PICTURED ON THE TABLE AT:
<https://www.clinton-ind.com/catalogs#alphabet-animals-fun-facts>

IMAGINATION SERIES
TREATMENT TABLES

Cool Pals Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7931	length	width	height
	68"	27"	33"
	172.72 cm	68.58 cm	83.82 cm

7931-1 Same as 7931 without adjustable backrest



COOL PALS COMPLETE includes: 7931 Cool Pals Treatment Table & 6131-BW Coordinating Cabinets (see page 44 for cabinet features)

ADD:

Coordinating Graphics (page 60) shown here for a total "Cool Pals" exam room designed environment

For special package discount pricing, ORDER USING MODEL # 7931-X

For special package discount pricing with table 7831-1, USE # 7931-1-X

Cool Pals Complete Exam Room uses four colors of floor covering to simulate a 3-dimensional iceberg effect.

IMAGINATION SERIES SCALE TABLES



The Clinton Kids™ Imagination Series Scale Tables take office decor to a new dimension with bright illustrations and the ultimate in comfort, performance and accuracy.

Imagination Series Scale Tables feature:

- Comfortable work height of 36" (91.44 cm)
- Top rimmed on 3 sides for safety
- Beveled edges & rounded corners on top
- Easy-clean, replaceable center pad
- Built-in infantometer
- Easy access to weight measurement and storage
- Built-in concealed paper dispenser
- Easy-clean laminate surfaces
- Ample 2-door and 1 drawer storage
- Optional storage locks available (-055)
- All models feature Clinton Complete™ pricing when combined with coordinating cabinets

See page 10 for Doran® scale features.

Dino Days Scale Table

- 2-door and 1 drawer storage


7834	length	width	height
	43"	24"	36"
	109.22 cm	60.96 cm	91.44 cm



right side detail



left side detail



DINO DAYS SCALE COMPLETE includes: 7834 Dino Days Scale Table & 6134-BW Coordinating Cabinets (see page 45 for cabinet features)

ADD:

Coordinating Graphics (page 61) shown here for a total "Dino Days" exam room designed environment

For special package discount pricing, ORDER USING MODEL # **7834-X**

IMAGINATION SERIES SCALE TABLES

Cool Pals Scale Table

- 2-door and 1 drawer storage

7831	length	width	height
	43"	24"	36"
	109.22 cm	60.96 cm	91.44 cm

COOL PALS SCALE COMPLETE—discount package pricing when combining 7831 table with coordinating cabinets (see page 44 for cabinets and features)

7831-X Must order using this model # for package discount



right side detail

left side detail



Rainforest Follies Scale Table

- 2-door and 1 drawer storage

7832	length	width	height
	43"	24"	36"
	109.22 cm	60.96 cm	91.44 cm

RAINFOREST FOLLIES SCALE COMPLETE—discount package pricing when combining 7832 table with coordinating cabinets (see page 44 for cabinets and features)

7832-X Must order using this model # for package discount



right side detail

left side detail



Space Place Scale Table

- 2-door and 1 drawer storage

7835	length	width	height
	43"	24"	36"
	109.22 cm	60.96 cm	91.44 cm

SPACE PLACE SCALE COMPLETE—discount package pricing when combining 7835 table with coordinating cabinets (see page 45 for cabinets and features)

7835-X Must order using this model # for package discount



right side detail

left side detail



Ocean Commotion Scale Table

- 2-door and 1 drawer storage

7836	length	width	height
	43"	24"	36"
	109.22 cm	60.96 cm	91.44 cm

OCEAN COMMOTION SCALE COMPLETE—discount package pricing when combining 7836 table with coordinating cabinets (see page 42 for cabinets and features)

7836-X Must order using this model # for package discount



right side detail

left side detail



IMAGINATION SERIES SCALE TABLES

Alley Cats & Dogs Scale Table

- 2-door and 1 drawer storage

7837	length	width	height
	43"	24"	36"
	109.22 cm	60.96 cm	91.44 cm



right side detail



left side detail



ALLEY CATS & DOGS SCALE COMPLETE includes 7837 Alley Cats & Dogs Scale Table & 6137-BW Coordinating Cabinets (see page 46 for cabinet features)

ADD:
Coordinating Graphics (page 59) shown here for a total "Alley Cats & Dogs" exam room designed environment



For special package discount pricing, ORDER USING MODEL # 7837-X

IMAGINATION SERIES SCALE TABLES

Serengeti Sunrise Scale Table

- 2-door and 1 drawer storage

7851	length	width	height
	43"	24"	36"
	109.22 cm	60.96 cm	91.44 cm

SERENGETI SUNRISE SCALE COMPLETE—discount package pricing when combining 7851 table with coordinating cabinets (see page 46 for cabinets and features)

7851-X Must order using this model # for package discount



right side detail

left side detail



Southwestfest Scale Table

- 2-door and 1 drawer storage

7852	length	width	height
	43"	24"	36"
	109.22 cm	60.96 cm	91.44 cm

SOUTHWESTFEST SCALE COMPLETE— discount package pricing when combining 7852 table with coordinating cabinets (see page 47 for cabinets and features)

7852-X Must order using this model # for package discount



right side detail

left side detail



Kangaroo Country Scale Table

- 2-door and 1 drawer storage

7853	length	width	height
	43"	24"	36"
	109.22 cm	60.96 cm	91.44 cm

KANGAROO COUNTRY SCALE COMPLETE—discount package pricing when combining 7853 table with coordinating cabinets (see page 47 for cabinets and features)

7853-X Must order using this model # for package discount



right side detail

left side detail



Merry Mermaids Scale Table

- 2-door and 1 drawer storage
- Paper dispenser from either side

7854	length	width	height
	43"	24"	36"
	109.22 cm	60.96 cm	91.44 cm



right side detail

left side detail



ADVENTURE SERIES TREATMENT TABLES



Take young patients on safari, down under or on a visit to America's icons. Clinton's Adventure Series Treatment Tables offer vivid photo collages to engage school-aged children.

Adventure Series Treatment Tables feature:

- Clinton original adventures
- Easy to clean laminate surfaces
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- 2" (5.08 cm) of firm foam padding
- 4 door storage with 2 adjustable shelves
- Slate Gray ends and backs
- Scratch and mar resistant surfaces
- Chip resistant PVC edges
- Premium, stain-resistant, woven, knit-backed vinyl upholstery in Gunmetal

Safari Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7940	length	width	height
	68"	27"	33"
	172.72 cm	68.58 cm	83.82 cm

7940-1 Same as 7940 without adjustable backrest



Down Under Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7942	length	width	height
	68"	27"	33"
	172.72 cm	68.58 cm	83.82 cm

7942-1 Same as 7942 without adjustable backrest



Americana Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7945	length	width	height
	68"	27"	33"
	172.72 cm	68.58 cm	83.82 cm

7945-1 Same as 7945 without adjustable backrest



Wild West Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7943	length	width	height
	68"	27"	33"
	172.72 cm	68.58 cm	83.82 cm

7943-1 Same as 7943 without adjustable backrest



Alaskan Wilderness Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7946	length	width	height
	68"	27"	33"
	172.72 cm	68.58 cm	83.82 cm

7946-1 Same as 7946 without adjustable backrest



Cosmos Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7944	length	width	height
	68"	27"	33"
	172.72 cm	68.58 cm	83.82 cm

7944-1 Same as 7944 without adjustable backrest



Equestrian Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest

7947	length	width	height
	68"	27"	33"
	172.72 cm	68.58 cm	83.82 cm

7947-1 Same as 7947 without adjustable backrest





Discover four brown eggs! Discover four cats, four birds or four little mice hiding in the pictures of our Discover Series Tables. Keeping children engaged in discovery while waiting for the doctor makes every visit easy on staff and parents.

Discover Series Treatment Tables feature:

- Clinton original discovery illustrations
- Easy to clean laminate surfaces
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- 2" (5.08 cm) of firm foam padding
- 4-door storage with 2 adjustable shelves
- Gray ends and backs
- Scratch and mar resistant surfaces
- Chip resistant PVC edges
- Premium, stain-resistant, woven, knit-backed vinyl upholstery in your choice of color

Here's how it works.

Hidden in each scene there are **4 Brown Eggs, 4 Cats, 4 Birds, and 4 Mice**. Keep the number of objects a secret. Ask the young patient how many of any one of the objects they can find. Whether it's all or a few, every child's answer is correct.

Next visit ask them to find one of the other three remaining objects and keep them amused for several visits to come.




DISCOVER SERIES SCALE TABLES

Sweet Dreams Candy Factory Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest
- Picture object map included

7370-71	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.74 cm

7370-71-088 Same as 7370-71 without adjustable backrest



Wings 'N Things Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest
- Picture object map included

7370-72	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.74 cm

7370-72-088 Same as 7370-72 without adjustable backrest



Secret Garden Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest
- Picture object map included

7370-73	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.74 cm

7370-73-088 Same as 7370-73 without adjustable backrest



Fairy Tale Dale Treatment Table

- 4-door storage with 2 adjustable shelves and adjustable backrest
- Picture object map included

7370-74	length	width	height
	72"	27"	31"
	182.88 cm	68.58 cm	78.74 cm

7370-74-088 Same as 7370-74 without adjustable backrest





OCEAN COMMOTION COMPLETE includes 6136-BW Ocean Commotion Cabinets & Coordinating Treatment or Scale Table (see pages 22 & 29 for tables and features)

ADD:

Coordinating Graphics (pages 59 & 60) shown here for a total "Ocean Commotion" exam room designed environment

For package price with Treatment Table, ORDER USING MODEL #s 7936-X or 7936-1-X

For package price with Scale Table, ORDER USING MODEL # 7836-X

Clinton Imagination Series Cabinets are the perfect complement for Imagination Series treatment and scale tables.

Imagination Series Cabinets feature:

- Clinton original illustrations
- Pre-assembled
- Easy hook and rail wall cabinet installation
- Recessed toe kick
- Easy-clean interior & exterior laminate surfaces
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Chip resistant PVC edges
- Steel fastener assembly
- Slate Gray ends and tops
- Scratch and mar resistant surfaces
- 4 adjustable, leveling feet
- Euro-style drawer fronts with metal sides and 75 lbs. (34.01 kg) capacity
- Black door and drawer pulls
- 3-way adjustable, 110°, soft close hinges
- All units shown with optional sink
- *Customize by changing sink location or adding other options (See page 49 for all options.)*
- **All models feature Clinton Complete™ pricing when combined with coordinating tables.**

Ocean Commotion Base & Wall Cabinets

- Wall unit has 4 doors/2 adjustable shelves
- Base unit has 4 doors/3 adjustable shelves and 3 drawers
- Base & wall units not sold separately

6136-BW (base & wall cabinets are only sold in sets)

base	length	depth	height
	61"	18"	35"
	154.94 cm	45.72 cm	88.9 cm

wall	length	depth	height
	61"	12"	24"
	154.94 cm	30.48 cm	60.96 cm

IMAGINATION SERIES CABINETS

Cool Pals Base & Wall Cabinets

- Wall unit has 4 doors/2 adjustable shelves
- Base unit has 4 doors/3 adjustable shelves and 3 drawers
- Base & wall units not sold separately
- Base unit shown with optional sink

6131-BW (base & wall cabinets are only sold in sets)

base	length	depth	height
	61"	18"	35"
	154.94 cm	45.72 cm	88.9 cm

wall	length	depth	height
	61"	12"	24"
	154.94 cm	30.48 cm	60.96 cm

COOL PALS COMPLETE—discount package pricing when combining 6131-BW cabinets with coordinating tables (see pages 24 & 28 for tables and features)

For package price with Treatment Table, ORDER USING MODEL #s 7931-X or 7931-1-X

For package price with Scale Table, ORDER USING MODEL # 7831-X



Rainforest Follies Base & Wall Cabinets

- Wall unit has 4 doors/2 adjustable shelves
- Base unit has 4 doors/3 adjustable shelves and 3 drawers
- Base & wall units not sold separately
- Base unit shown with optional sink

6132-BW (base & wall cabinets are only sold in sets)

base	length	depth	height
	61"	18"	35"
	154.94 cm	45.72 cm	88.9 cm

wall	length	depth	height
	61"	12"	24"
	154.94 cm	30.48 cm	60.96 cm

RAINFOREST FOLLIES COMPLETE—discount package pricing when combining 6132-BW cabinets with coordinating tables (see pages 20 & 28 for tables and features)

For package price with Treatment Table, ORDER USING MODEL #s 7932-X or 7932-1-X

For package price with Scale Table, ORDER USING MODEL # 7832-X



Dino Days Base & Wall Cabinets

- Wall unit has 4 doors/2 adjustable shelves
- Base unit has 4 doors/3 adjustable shelves and 3 drawers
- Base & wall units not sold separately
- Base unit shown with optional sink

6134-BW (base & wall cabinets are only sold in sets)

base	length	depth	height
	61"	18"	35"
	154.94 cm	45.72 cm	88.9 cm

wall	length	depth	height
	61"	12"	24"
	154.94 cm	30.48 cm	60.96 cm

DINO DAYS COMPLETE—discount package pricing when combining 6134-BW cabinets with coordinating tables (see pages 26 & 22 for tables and features)

For package price with Treatment Table, ORDER USING MODEL #s 7934-X or 7934-1-X

For package price with Scale Table, ORDER USING MODEL # 7834-X



Space Place Base & Wall Cabinets

- Wall unit has 4 doors/2 adjustable shelves
- Base unit has 4 doors/3 adjustable shelves and 3 drawers
- Base & wall units not sold separately
- Base unit shown with optional sink

6135-BW (base & wall cabinets are only sold in sets)

base	length	depth	height
	61"	18"	35"
	154.94 cm	45.72 cm	88.9 cm

wall	length	depth	height
	61"	12"	24"
	154.94 cm	30.48 cm	60.96 cm

SPACE PLACE COMPLETE—discount package pricing when combining 6135-BW cabinets with coordinating tables (see pages 18 & 29 for tables and features)

For package price with Treatment Table, ORDER USING MODEL #s 7935-X or 7935-1-X

For package price with Scale Table, ORDER USING MODEL # 7835-X

IMAGINATION SERIES CABINETS

Alley Cats & Dogs Base & Wall Cabinets

- Wall unit has 4 doors/2 adjustable shelves
- Base unit has 4 doors/3 adjustable shelves and 3 drawers
- Base & wall units not sold separately
- Base unit shown with optional sink

6137-BW (base & wall cabinets are only sold in sets)

base	length	depth	height
	61"	18"	35"
	154.94 cm	45.72 cm	88.9 cm

wall	length	depth	height
	61"	12"	24"
	154.94 cm	30.48 cm	60.96 cm

ALLEY CATS & DOGS COMPLETE—discount package pricing when combining 6137-BW cabinets with coordinating tables (see pages 22 & 30 for tables and features)

For package price with Treatment Table, ORDER USING MODEL #s 7937-X or 7937-1-X

For package price with Scale Table, ORDER USING MODEL # 7837-X



Serengeti Sunrise Base & Wall Cabinets

- Wall unit has 4 doors/2 adjustable shelves
- Base unit has 4 doors /3 adjustable shelves and 3 drawers
- Base & wall units not sold separately
- Base unit shown with optional sink

6151-BW (base & wall cabinets are only sold in sets)

base	length	depth	height
	61"	18"	35"
	154.94 cm	45.72 cm	88.9 cm

wall	length	depth	height
	61"	12"	24"
	154.94 cm	30.48 cm	60.96 cm

SERENGETI SUNRISE—discount package pricing when combining 6151-BW cabinets with coordinating tables (see pages 23 & 32 for tables and features)

For package price with Treatment Table, ORDER USING MODEL #s 7951-X or 7951-1-X

For package price with Scale Table, ORDER USING MODEL # 7851-X



Southwestfest Base & Wall Cabinets

- Wall unit has 4 doors/2 adjustable shelves
- Base unit has 4 doors/3 adjustable shelves and 3 drawers
- Base & wall units not sold separately
- Base unit shown with optional sink

6152-BW (base & wall cabinets are only sold in sets)

base	length	depth	height
	61"	18"	35"
	154.94 cm	45.72 cm	88.9 cm

wall	length	depth	height
	61"	12"	24"
	154.94 cm	30.48 cm	60.96 cm

SOUTHWESTFEST COMPLETE—discount package pricing when combining 6152-BW cabinets with coordinating tables (see pages 23 & 33 for tables and features)

For package price with Treatment Table, ORDER USING MODEL #s 7952-X or 7952-1-X

For package price with Scale Table, ORDER USING MODEL # 7852-X

Kangaroo Country Base & Wall Cabinets

- Wall unit has 4 doors/2 adjustable shelves
- Base unit has 4 doors/3 adjustable shelves and 3 drawers
- Base & wall units not sold separately
- Base unit shown with optional sink

6153-BW (base & wall cabinets are only sold in sets)

base	length	depth	height
	61"	18"	35"
	154.94 cm	45.72 cm	88.9 cm

wall	length	depth	height
	61"	12"	24"
	154.94 cm	30.48 cm	60.96 cm

KANGAROO COUNTRY COMPLETE—discount package pricing when combining 6153-BW cabinets with coordinating tables (see pages 23 & 32 for tables and features)

For package price with Treatment Table, ORDER USING MODEL #s 7953-X or 7953-1-X

For package price with Scale Table, ORDER USING MODEL # 7853-X



Kangaroo Country Cabinets also coordinate with Fun Series Outback Buggy models 7844 & 7040 (pages 8 & 12)

FUNSERIES CABINETS



Fun Series Cabinets have the same features as the Imagination Series Cabinets.

Schoolhouse Base & Wall Cabinets

- Wall unit has 4 doors/2 adjustable shelves
- Base unit has 4 doors/3 adjustable shelves and 3 drawers
- Base & wall units not sold separately
- Base unit shown with optional sink

6120-BW (base & wall cabinets are only sold in sets)

base	length	depth	height
	61"	18"	35"
	154.94 cm	45.72 cm	88.9 cm

wall	length	depth	height
	61"	12"	24"
	154.94 cm	30.48 cm	60.96 cm

ZOO BUS COMPLETE—discount package pricing when combining 6120-BW cabinets with coordinating tables (see pages 6 & 10 for tables and features)

For package price with Treatment Table, ORDER USING MODEL # 7020-X

For package price with Scale Table, ORDER USING MODEL # 7822-X



Firehouse Base & Wall Cabinets

- Wall unit has 4 doors/2 adjustable shelves
- Base unit has 4 doors/3 adjustable shelves and 3 drawers
- Base & wall units not sold separately
- Base unit shown with optional sink

6130-BW (base & wall cabinets are only sold in sets)

base	length	depth	height
	61"	18"	35"
	154.94 cm	45.72 cm	88.9 cm

wall	length	depth	height
	61"	12"	24"
	154.94 cm	30.48 cm	60.96 cm

ENGINE K-9 COMPLETE—discount package pricing when combining 6130-BW cabinets with coordinating tables (see pages 4 & 12 for tables and features)

For package price with Treatment Table, ORDER USING MODEL # 7030-X

For package price with Scale Table, ORDER USING MODEL # 7833-X



Cool Park Campgrounds Base & Wall Cabinets

- Wall unit has 4 doors with 2 adjustable shelves
- Base unit has 4 doors/3 adjustable shelves and 3 drawers
- Base & wall units not sold separately
- Base unit shown with optional sink

6150-BW (base & wall cabinets are only sold in sets)

base	length	depth	height
	61"	18"	35"
	154.94 cm	45.72 cm	88.9 cm

wall	length	depth	height
	61"	12"	24"
	154.94 cm	30.48 cm	60.96 cm

COOL CAMPER COMPLETE—discount package pricing when combining 6150-BW coordinating cabinets with coordinating table (see page 9 for table features)

For package price with Treatment Table, ORDER USING MODEL # 7050-X

CABINET OPTIONS

Stainless Steel Sink and Gooseneck Faucet with Wing Levers

(shown on all cabinet photos)

- Choose your location in middle, on right or on left end of the base

022

Add a Drawer

- Available for any base model
- Maximum 5 drawer limit in each vertical cabinet section

044

Door and Drawer Locks

- Individual cam lock for doors and drawers
- Independently operated
- Orders with multiple locks will be all keyed alike unless specified otherwise

055

Door Lock and Inside Latch Combo

- Used on double door units
- One door will have a lock
- The other door will have a slide latch on the inside back

066

OFFICE ACCESSORIES

No office is complete without accessories. Clinton offers a variety of pediatric and non-pediatric specific accessories to fit your needs.

Koala Laminate Foot Stool

- Heavy laminate wood construction
- Non-skid tread
- Your choice of laminate
- Koala graphic stickers/2 per stool

10-K	length	depth	height
	20"	16"	7"
	50.8 cm	40.64 cm	17.78 cm

Dalmatian Laminate Foot Stool

- Heavy laminate wood construction
- Non-skid tread
- Your choice of laminate
- Dalmatian graphic stickers/2 per stool

10-D	length	depth	height
	20"	16"	7"
	50.8 cm	40.64 cm	17.78 cm

Hyena Laminate Foot Stool

- Heavy laminate wood construction
- Non-skid tread
- Your choice of laminate
- Hyena graphic stickers/2 per stool

10-H	length	depth	height
	20"	16"	7"
	50.8 cm	40.64 cm	17.78 cm

More stool designs available.
Visit www.clinton-ind.com for complete selection



Black Frame Chair/No Arms

- Nylon glides
- Contoured seat
- 2" (5.08 cm) thick firm foam padded seat
- Now with wall bumper
- Your choice of 15 colors of upholstery (See page 65 for upholstery colors)
- Box Dimensions 22" x 21" x 27" (55.88 cm x 53.34 cm x 68.58 cm)

C-40	seat	seat height
	17 1/2" x 18"	17 3/4"
	44.45 cm x 45.72 cm	45.08 cm



Black Frame Chair with Arms

- Nylon glides
- Contoured seat
- 2" (5.08 cm) thick firm foam padded seat
- Now with wall bumper
- Your choice of 15 colors of upholstery (See page 65 for upholstery colors)
- Box Dimensions 22" x 21" x 27" (55.88 cm x 53.34 cm x 68.58 cm)

C-50	seat	seat height
	17 1/2" x 18"	18"
	44.45 cm x 45.72 cm	45.72 cm



5-Leg Pneumatic Stool/Multi-Color Top

- Single lever pneumatic height adjustment
- 24" (60.96 cm) diameter nylon base
- Black hooded dual wheel casters
- 4" (10.16 cm) thick padded seat
- Easy-clean, Perma-Block 3 upholstery in 15 colors, mix or match (See page 65 for upholstery colors)

2135-M	seat diameter	height range
	16"	18" - 23"
	40.64 cm	45.72 cm - 58.42 cm

Treatment Cabinet with 2 Doors & 1 Drawer

- 1 adjustable shelf
- Drawer with metal sides that hold up to 75 lbs. (34.01 kg)
- All laminate construction

8821	length	depth	height
	24"	18"	35"
	60.96 cm	45.72 cm	88.9 cm



Treatment Cabinet with 2 Doors & 2 Drawers

- 1 adjustable shelf
- Drawers with metal sides that hold up to 75 lbs. (34.01 kg)
- All laminate construction

8822	length	depth	height
	24"	18"	35"
	60.96 cm	45.72 cm	88.9 cm



Treatment Cabinet with 5 Drawers

- Drawers with metal sides that hold up to 75 lbs. (34.01 kg)
- All laminate construction

8805	length	depth	height
	24"	18"	35"
	60.96 cm	45.72 cm	88.9 cm



Mobile Treatment Cabinet with 2 Doors & 2 Drawers

- 1 adjustable shelf
- Drawers with metal sides that hold up to 75 lbs. (34.01 kg)
- All laminate construction
- Mounted on 3" dual wheel, swivel casters (2 locking)
- Plastic bumpers on each side of the base

8922	length	depth	height
	24"	18"	35"
	60.96 cm	45.72 cm	88.9 cm



Mobile Equipment Cart with 2 Drawers

- Mounted on a stable 5-leg black nylon base
- 2 Drawers
- Top rimmed on 3 sides
- All laminate cabinet

8920	length	depth	height
	22"	17"	33"
	55.88 cm	43.18 cm	83.82 cm



Mobile Equipment Cart with 1 Drawer & 1 Shelf

- Mounted on a stable 5-leg black nylon base
- 1 Drawer
- Top and back have wire-feed cutouts
- Top rimmed on 3 sides
- All laminate cabinet

8940	length	depth	height
	22"	17"	33"
	55.88 cm	43.18 cm	83.82 cm



Optional Basket

- Wire basket with wire epoxy finish

089	length	depth	height
	17 1/2"	6 3/4"	5"
	44.45 cm	17.14 cm	12.7 cm

