


## Imagination Series, Alphabet Animals, Treatment Table



### MODEL 7955 (Adjustable Backrest)

### MODEL 7955-1 (Flat Top)

- No upcharge for adjustable backrest
- Clinton original themes
- Easy to clean laminate surfaces
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- 2" of firm foam padding
- 4 door storage with 2 adjustable shelves
- Scratch and mar resistant surfaces
- Chip resistant PVC edges
- Premium, Gunmetal, stain-resistant, woven, knit-backed vinyl upholstery
- Soft-close, door hinges
- Graphics are on one long side with slate gray ends and back
- Top features seamless rounded outside corners
- Paper dispenser included
- Available with or without backrest
- Fun Facts about the animals pictured included with the table

### SPECIFICATIONS (BOTH MODELS)

68"L x 27"W x 33"H

California residents please be advised of the following, pursuant to Proposition 65:

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).


1



### Options




066 Door & Inside Latch Combo




055 Door Locks



040 Paper Cutter



20 Small Pillow



60 Large Pillow

### Upholstery



Gunmetal