

General, Flat Top, Ultrasound Table



MODEL 80061

- One piece flat top with firm foam padding
- Trendelenburg/Reverse Trendelenburg plus or minus 15 °
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser, cutter and positioning straps included
- Optional side rails & casters available
- 34" wide model available
- Model 80061: 500 lbs load capacity under normal use
- Model 80061-X: 600 lbs load capacity under normal use



Trendelenburg




Reverse Trendelenburg



Optional Side Rails & Casters

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov. 1

80061 SPECIFICATIONS

	length	width	height range
overall	72"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "

80061-X Extra-Wide 34" width model



Upholstery Colors



OPTIONS

Side Rails (pair)

- 3 height positions
- Plunger height adjustment

098



Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087

