



Clinton takes power tables to new heights with a complete line of power exam, ultrasound and imaging tables designed and engineered to fit the needs of physicians, technicians and therapists.

Clinically proven in the field, Clinton power tables, are among the best in their class. With quality construction, dependable smooth, operation and enough power to satisfy nearly every need, Clinton Power Tables are practical and, best of all, affordable.

Welded steel construction makes them tough enough for years of worry-free service, and two inches of firm foam padding topped with premium, knit-backed vinyl makes patient comfort a priority in all our models. Choose from Shrouded Power Exam Tables, Open Base Power Exam Tables, Open Base Power Imaging Tables or Power Ultrasound and Imaging Trendelenburg Tables.

Make Clinton your one source for affordable, convenient power tables — and plug into the power of Clinton.

#### INDEX

Shrouded, Power Tables .....	2-5
Open Base, Power Tables.....	6-9
Open Base, Power, Imaging Tables.....	10-11
Power, Ultrasound, Trendelenburg Tables....	12-15
Power, Imaging, Trendelenburg Tables.....	16-19
Power Table Options .....	20
Color Selector.....	21



Clinton power tables are among the most practical and affordable power exam tables on the market. They offer smooth dependable power, stylish design and years of dependable service.

Clockwise from left to right: 81360—Shrouded Power Table with Adjustable Backrest, Drop Section & Stirrups; 81330—Shrouded Power Table with Adjustable backrest & Drop Section; 80079—Multi-Use, Power Imaging, Trendelenburg Table with Adjustable Backrest, Drop Section & Stirrups; 80330 Open Base Power Imaging Table with Adjustable Backrest & Drop Section.

## SHROUDED, POWER TABLES

### Shrouded, Power Tables feature:

- Ultra-strong, uni-frame constructed base and top frame provides a rigid patient platform
- All welded steel base and top frame
- 3/4" plywood top for extra support (1.9 cm)
- Textured ABS plastic encloses lift mechanism and makes cleaning easier
- Smooth electric height adjustment moves straight up and down
- Pneumatic, backrest cylinder is internally mounted on the frame for a cleaner appearance
- Hardened bolts and self lubricating pivot points for years of service
- Dependable 110 volt motor
- Hands-free foot control
- Optional casters, safety rails and more
- Made in the USA by Clinton Industries
- Low heights may help achieve ADA accessibility tax credit
- Dependable scissors lift mechanism
- Many different top configurations and options available
- 5-year limited warranty on frame and electrical components
- Ships fully assembled with no upcharge



### Shrouded, Power Table with One Piece Top

- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (page 20)

81100	length	width	height range
overall	72"	27"	21"-34"
	183 cm	69 cm	53 cm - 86 cm

### Shrouded, Power Table with Adjustable Backrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (page 20)

81200	length	width	height range
overall	72"	27"	21"-34"
	183 cm	69 cm	53 cm - 86 cm
backrest	38"	27"	
	97 cm	69 cm	



### Shrouded, Power Table with Adjustable Backrest & Drop Section

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section is adjustable from either side and raises up to 30° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (specify when ordering)
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (page 20)

81330	length	width	height range
overall	76"	27"	21"-34"
	193 cm	69 cm	53 cm - 86 cm
backrest	38"	27"	
	97 cm	69 cm	
drop section	13"	27"	
	33 cm	69 cm	

## SHROUDED, POWER TABLES

### Shrouded, Space Saver Power Table with Drop Section

- Drop section is adjustable from either side and raises up to 30° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (*specify when ordering*)
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (*page 20*)

81220	length	width	height range
overall	76"	27"	21"-34"
	193 cm	69 cm	53 cm - 86 cm
drop section	13"	27"	
	33 cm	69 cm	



81220

### Shrouded, Power Adjustable, Casting Table with ClintonClean™ Leg Rest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- ClintonClean™ Solid Plastic Surface leg rest is adjustable from either side and raises up to 30° from flat position
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (*page 20*)

81350	length	width	height range
overall	76"	27"	21"-34"
	193 cm	69 cm	53 cm - 86 cm
backrest	38"	27"	
	97 cm	69 cm	
leg rest	13"	26"	
	33 cm	66 cm	



81350

### Shrouded, Power Table with Stirrups, Adjustable Backrest & Footrest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and raises up to 30° from flat position
- Retractable, 2-position, adjustable stirrups
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (*page 20*)

81360	length	width	height range
overall	76"	27"	21"-34"
	193 cm	69 cm	53 cm - 86 cm
backrest	38"	27"	
	97 cm	69 cm	
footrest	13"	19"	
	33 cm	48 cm	



81360

## SHROUDED, EXTRA-WIDE, BARIATRIC, POWER TABLES

### Shrouded Extra-Wide Bariatric Straight Power Table

- All welded, uni-frame steel base and top frame
- 600 lbs. load capacity under normal use (272 kg)
- Optional side rails, locking casters & more (*page 20*)

84108-34	length	width	height range
overall	72"	34"	21"-34"
	183 cm	86 cm	53 cm - 86 cm

84108-40 40" (101.60 cm) width model



84108-40

### Shrouded, Extra-Wide, Bariatric, Power Table with Adjustable Backrest

- Dual cylinder, pneumatic backrest, adjustable from either side and raises up to 80° from flat
- 600 lbs. load capacity under normal use (272 kg)
- Optional side rails, locking casters & more (*page 20*)

84208-34	length	width	height range
overall	72"	34"	21"-34"
	183 cm	86 cm	53 cm - 86 cm
backrest	38"	34"	
	97 cm	86 cm	

84208-40 40" (101.60 cm) width model



84208-34

### Shrouded, Extra-Wide, Bariatric, Power Table with Adjustable Backrest & Drop Section

- Dual cylinder pneumatic backrest, adjustable from either side and raises up to 80° from flat
- Drop section is adjustable from either side and raises up to 30° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (*specify when ordering*)
- 600 lbs. load capacity under normal use (272 kg)
- Optional side rails, locking casters & more (*page 20*)

84438-34	length	width	height range
overall	76"	34"	21"-34"
	193 cm	86 cm	53 cm - 86 cm
backrest	38"	34"	
	97 cm	86 cm	
drop section	13"	34"	
	33 cm	86 cm	

84438-40 40" (101.60 cm) width model



84438-34

**Open Base, Power Tables feature:**

- Ultra-strong uni-frame constructed base and top frame provides a rigid patient platform
- 3/4" plywood top for extra support (1.9 cm)
- All welded steel base, top frame and scissors lift
- Smooth electric height adjustment moves straight up and down
- Pneumatic, backrest cylinder are internally mounted on the frame for a cleaner appearance
- Hardened bolts and self lubricating pivot points for years of service
- Dependable 110 volt motor
- Hands-free foot control
- Optional casters, safety rails and more
- Made in the USA by Clinton Industries
- Low heights may help achieve ADA accessibility tax credit
- Dependable scissors lift mechanism
- Many different top configurations and options available
- 5-year limited warranty on frame and electrical components
- Ships fully assembled with no upcharge

**Open Base, Power Table with Stirrups, Adjustable Backrest & Footrest**

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Retractable, 2-position, adjustable stirrups
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (page 20)

80360	length	width	height range
overall	76" 193 cm	27" 69 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	27" 69 cm	
footrest	13" 33 cm	19" 48 cm	



**Open Base, Power Table with One Piece Top**

- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (page 20)

80100	length	width	height range
	72" 183 cm	27" 69 cm	18"-35" 46 cm - 89 cm



**Open Base, Power Table with Adjustable Backrest**

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (page 20)

80200	length	width	height range
overall	76" 193 cm	27" 69 cm	18"-35" 46 cm - 89 cm
Backrest	38" 97 cm	27" 69 cm	



**Open Base, Power Table with Adjustable Backrest & Drop Section**

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Drop section, adjustable from either side and lowers up to 80° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (specify when ordering)
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (page 20)

80330	length	width	height range
overall	76" 193 cm	27" 69 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	27" 69 cm	
drop section	13" 33 cm	27" 69 cm	



## OPEN BASE, POWER TABLES

### Space Saver Power Table with Drop Section

- Drop section is adjustable from either side and lowers up to 80° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (*specify when ordering*)
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (*page 20*)

80220	length	width	height range
overall	76"	27"	18"-35"
	193 cm	69 cm	46 cm - 89 cm
drop section	13"	27"	
	33 cm	69 cm	



### Power Adjustable Casting Table with ClintonClean™ Surface Leg Rest

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and lowers up to 80° from flat position
- ClintonClean™ Solid Plastic Surface leg rest is adjustable from either side and raises up to 30° from flat position
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (*page 20*)

80350	length	width	height range
overall	76"	27"	18"-35"
	193 cm	69 cm	46 cm - 89 cm
backrest	38"	27"	
	97 cm	69 cm	
leg rest	13"	26"	
	33 cm	66 cm	



## OPEN BASE, EXTRA-WIDE, BARIATRIC, POWER TABLES

### Open Base, Extra-Wide, Bariatric, Straight Power Table

- All welded, uni-frame steel base and top frame
- 600 lbs. load capacity under normal use (272 kg)
- Optional side rails, locking casters & more (*page 20*)

84100-34	length	width	height range
overall	72"	34"	18"-35"
	183 cm	86 cm	46 cm - 89 cm

84100-40 Same as above with 40" width (102 cm)



### Open Base, Power Table with Adjustable Backrest

- Dual cylinder, pneumatic backrest, adjustable from either side and raises up to 80° from flat
- All welded, uni-frame steel base and top frame
- 600 lbs. load capacity under normal use (272 kg)
- Optional side rails, locking casters & more (*page 20*)

84200-34	length	width	height range
overall	72"	34"	18"-35"
	183 cm	86 cm	46 cm - 89 cm
backrest	38"	34"	
	97 cm	86 cm	

84200-40 Same as above with 40" widths (102 cm)



### Open Base, Power Table with Adjustable Backrest & Drop Section

- Dual cylinder pneumatic backrest, adjustable from either side and raises up to 80° from flat
- Drop section is adjustable from either side and lowers up to 80° from flat position
- Standard drop section has nose hole, but it may be ordered without nose hole (*specify when ordering*)
- 600 lbs. load capacity under normal use (272 kg)
- Optional side rails, locking casters & more (*page 20*)

84430-34	length	width	height range
overall	76"	34"	18"-35"
	193 cm	86 cm	46 cm - 89 cm
backrest	38"	34"	
	97 cm	86 cm	
drop section	13"	34"	
	33 cm	86 cm	

84430-40 Same as above with 40" widths (102 cm)



## OPEN BASE, POWER IMAGING TABLES

### Open Base, Power Imaging Tables feature:

- Ultra-strong uni-frame constructed base and top frame provides a rigid patient platform
- 3/4" plywood top for extra support (1.9 cm)
- All welded steel base, top frame and scissors lift
- Smooth electric height adjustment moves straight up and down
- Pneumatic, backrest cylinder are internally mounted on the frame for a cleaner appearance
- Hardened bolts and self lubricating pivot points for years of service
- Built-in hinged window drop. 10" x 10"
- Metal slide bar supports window drop
- Dependable 110 volt motor
- Hands-free foot control
- Optional casters, safety rails and more
- Made in the USA by Clinton Industries
- Low heights may help achieve ADA accessibility tax credit
- Dependable scissors lift mechanism that does not need to roll on the floor
- Many different top configurations and options available
- 5-year limited warranty on frame and electrical components
- Ships fully assembled with no upcharge

### Open Base, Multi-Use, Power Imaging Table with Stirrups

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Retractable, 2-position, adjustable stirrups
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (page 20)

85309	length	width	height range
overall	76" 193 cm	27" 69 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	27" 69 cm	
footrest	13" 33 cm	19" 48 cm	



Flat position with window drop engaged

### Open Base, Power Imaging Table with Window Drop

- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (page 20)

85100	length	width	height range
overall	72" 183 cm	27" 69 cm	18"-35" 46 cm - 89 cm



85100

### Open Base, Power Imaging Table with Adjustable Backrest and Window Drop

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- All welded steel base, top frame and scissors lift
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (page 20)

85200	length	width	height range
overall	72" 183 cm	27" 69 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	27" 69 cm	



85200

### Power Imaging Table with 3 Section Top

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and lowers up to 80° from flat position
- 500 lbs. load capacity under normal use (227 kg)
- Optional side rails, locking casters & more (page 20)

85300	length	width	height range
overall	76" 193 cm	27" 69 cm	18"-35" 46 cm - 89 cm
backrest	38" 97 cm	27" 69 cm	
footrest	13" 33 cm	27" 69 cm	



85300

**E-Series Power Tables feature:**

- Electric power height adjustment with foot control
- Ultra-strong, uni-frame base and top construction
- Heavy duty, all welded, lift mechanism
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Hands-free low voltage foot control
- Easy service, self contained, electrical components
- Seamless, rounded corner, upholstered top
- 1.5" (4 cm) ultra-firm padding
- Low height may help achieve ADA tax credit
- Customize with options
- Made in the USA by Clinton Industries
- Backed by a 5-year limited warranty
- 500 lbs. (227 kg) load capacity under normal use



**Open Base, Power Table with One Piece Top**

86100	length	width	height range
	72"	27"	18"-35"
	183 cm	69 cm	46 cm - 89 cm



86100

**Open Base, Power Table with Adjustable Backrest**

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position

86200	length	width	height range
overall	72"	27"	18"-35"
	183 cm	69 cm	46 cm - 89 cm
Backrest	38"	27"	
	97 cm	69 cm	



86200

**Open Base, Power Table with Adjustable Backrest & Drop Section**

- Drop section, adjustable from either side and lowers up to 80° from flat position

86220	length	width	height range
overall	76"	27"	18"-35"
	193 cm	69 cm	46 cm - 89 cm
drop section	13"	27"	
	33 cm	69 cm	



86220

**Open Base, Power Table with Adjustable Backrest & Footrest**

- Pneumatic backrest, adjustable from either side and raises up to 80° from flat position
- Footrest is adjustable from either side and lowers up to 80° from flat position

86300	length	width	height range
overall	76"	27"	18"-35"
	193 cm	69 cm	46 cm - 89 cm
backrest	38"	27"	
	97 cm	69 cm	
footrest	13"	19"	
	33 cm	48 cm	




86300



**Power, Ultrasound, Trendelenburg Tables feature:**

- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Power supplies and motors meet UL 60601-1 and CAN/CSA C22.2 601.1
- Industry leading 5-year, limited warranty
- Standard models: 500 lbs load capacity under normal use (227 kg)
- Extra-wide models: 600 lbs load capacity under normal use (272 kg)



80061



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

**General, Flat Top, Ultrasound Table**

- One piece flat top with firm foam padding
- Optional side rails, locking casters & more (page 20)

80061	length	width	height range
	72"	27"	21 1/2" – 37 1/2"
	183 cm	69 cm	55 cm – 95 cm

**80061-X** Extra-wide: 34" (86 cm)

**General Ultrasound Table with Three-Section Top**

- Three-piece top with pneumatic adjustable backrest and manual footrest
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Optional side rails, locking casters & more (page 20)

80063	length	width	height range
overall	76"	27"	21 1/2" – 37 1/2"
	193 cm	68.58 cm	55 cm – 95 cm
backrest	38"	27"	
	97 cm	68.58 cm	
footrest	13"	27"	
	33 cm	69 cm	

**80063-X** Extra-wide: 34" (86 cm)



80063



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

**Multi-Use, Ultrasound Table with Stirrups**

- Three-piece top with pneumatic adjustable backrest, manual footrest, and stirrups
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Optional side rails, locking casters & more (page 20)

80069	length	width	height range
overall	76"	27"	21 1/2" – 37 1/2"
	193 cm	68.58 cm	55 cm – 95 cm
backrest	38"	27"	
	97 cm	69 cm	
footrest	12"	18.5"	
	30 cm	47 cm	

**80069-X** Extra-wide: 34" (86 cm)



80069



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

POWER, ULTRASOUND,  
TRENDELENBURG TABLES

**General, Ultrasound Table with  
Adjustable Backrest**

- Two-piece top with easy-to-adjust, pneumatic backrest
- Optional side rails, locking casters & more (page 20)

**80062**

overall	72"	27"	21 <sup>1</sup> / <sub>2</sub> " – 37 <sup>1</sup> / <sub>2</sub> "
	183 cm	69 cm	55 cm – 95 cm

backrest	38"	27"	<b>POWER</b> Backrest Available
	97 cm	69 cm	

**80062-X** Extra-wide: 34" (86 cm)



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters



**Power Imaging, Trendelenburg Tables feature:**

- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Power supplies and motors meet UL 60601-1 and CAN/CSA C22.2 601.1
- Industry leading 5-year, limited warranty
- Standard models: 500 lbs load capacity under normal use (227 kg)
- Extra-wide models: 600 lbs load capacity under normal use (272 kg)



**Flat Top, Power Imaging Table with Drop Window**

- One piece flat top with drop window
- Optional side rails, locking casters & more (page 20)

80071	length	width	height range
overall	72" 183 cm	27" 69 cm	21 1/2" – 37 1/2" 55 cm – 95 cm
window	10" 25 cm	10" 25 cm	

**80071-X** Extra-wide: 34" (86 cm)



80061



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

**Power Imaging Table with Fowler Back and Drop Window**

- Two-piece top with pneumatic adjustable backrest and window drop
- Optional side rails, locking casters & more (page 20)

80072	length	width	height range
overall	72" 183 cm	27" 69 cm	21 1/2" – 37 1/2" 55 cm – 95 cm
backrest	38" 97 cm	27" 69 cm	
window	10" 25 cm	10" 25 cm	

**80072-X** Extra-wide: 34" (86 cm)



80072



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

**Power Imaging Table with Three-Section Top and Drop Window**

- Three-piece top with pneumatic adjustable backrest, manual footrest and window drop
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Optional side rails, locking casters & more (page 20)

80073	length	width	height range
overall	76" 193 cm	27" 68.58 cm	21 1/2" – 37 1/2" 54.61 cm – 95.25 cm
backrest	38" 96.52 cm	27" 68.58 cm	
footrest	12" 30.48 cm	27" 68.58 cm	
window	10" 25.4 cm	10" 25.4 cm	

**80073-X** Extra-wide: 34" (86 cm)



80073



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

POWER, IMAGING,  
TRENDELENBURG TABLES

**Multi-Use, Power Imaging Table with Stirrups and Drop Window**

- Three-section table with pneumatic adjustable backrest, manual footrest, stirrups and window drop
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Optional side rails, locking casters & more (page 20)

80079	length	width	height range
overall	76"	27"	21 1/2" – 37 1/2"
	193 cm	69 cm	55 cm – 95 cm
backrest	38"	27"	
	97 cm	69 cm	
footrest	12"	18 1/2"	
	30 cm	47 cm	
window	10"	10"	
	25 cm	25 cm	

**POWER**  
Backrest & Footrest  
Available

**80079-X** Extra-wide: 34" (86 cm)



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters



**5 Year**  
Limited  
Warranty

## POWER TABLE OPTIONS

### Safety Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 safety straps per set

080	length	width
	92"	2"
	233.67 cm	5 cm

### Safety Rails (pair)

- 3 height position
- Plunger height adjustment

082

### Power Backrest

*(Trendelenburg Tables only)*

- Smooth power control raises backrest to 80° and lowers to flat

011

### Power Footrest

*(Trendelenburg Tables only)*

- Smooth power control lowers footrest to 30° and raises to flat

090

### Hand-Held Height Control

*(Not available on E-Series Tables)*

- Alternative to foot control
- Stores easily under top when not in use

084

### Locking Casters

*(Not available on E-Series Tables)*

- Heavy duty, rubber wheel, easy-roll casters

087

### 30" Width Top

*(Trendelenburg Tables only)*

037



## E-SERIES OPTIONS

### Nose Hole

- Nose hole in drop section of 86220 or 86300

009

### Paper Dispenser

- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030

### Paper Cutter

- Clear plastic with elastic straps and buckle

009

### Paper Dispenser/Cutter Combo

- Paper dispenser and paper cutter ordered as one option

009





**One Complete Exam Room, One Order Number, One Price**

A complete office, well equipped and ready for patients. Clinton's Power Ready Room or E-Series Ready Room package offers big advantages in function, one source convenience and a pricepoint that fits your budget.

**81360-RR**

*(no model substitutions, your choice of Fashion Finish cabinetry and upholstery colors)*

- 81360 in your choice of upholstery color
- 8348 Wall Cabinet in your choice of Fashion Finish
- 8648 Base Cabinets with sink in your choice of Fashion Finish
- 48P your choice of Postform Countertop
- 2135 in your choice of upholstery color
- C-40B Premium, Padded, Side Chair in your choice of upholstery color
- GS-1000 Stainless Steel Glove Box holder
- TR-20W 20 Quart, Medium Round, White Enamel Steel, Waste Receptacle




**86300-RR**

*(no model substitutions, your choice of Fashion Finish cabinetry and upholstery colors)*

- 86300 in your choice of upholstery color
- 8348 Wall Cabinet in your choice of Fashion Finish
- 8648 Base Cabinets with sink in your choice of Fashion Finish
- 48P your choice of Postform Countertop
- 2135 in your choice of upholstery color
- C-40B Premium, Padded, Side Chair in your choice of upholstery color
- GS-1000 Stainless Steel Glove Box holder
- TR-20W 20 Quart, Medium Round, White Enamel Steel, Waste Receptacle

**COUNTERTOP COLORS**



*fashionfinish*  
CABINETRY



800-441-9131 • 717-848-3519

Toll-Free Fax 866-522-1233

Fax 717-843-5871

[www.clinton-ind.com](http://www.clinton-ind.com)