

Multi-Use, Ultrasound Table with Stirrups



MODEL 80069

- Three-piece top with pneumatic adjustable backrest, manual footrest, and stirrups
- Footrest is adjustable from either side and lowers up to 80° from flat position
- Trendelenburg/Reverse Trendelenburg plus or minus 15 °
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser, cutter and positioning straps included
- Optional side rails & casters available
- 34" wide model available
- Model 80069: 500 lbs load capacity under normal use
- Model 80069-X: 600 lbs load capacity under normal use



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

80069 SPECIFICATIONS

	length	width	height range
overall	76"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
backrest	38"	27"	
footrest	12"	18 ¹ / ₂ "	

80069-X Extra-Wide 34" width model

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black

OPTIONS

Side Rails (pair)

- 3 height positions
- Plunger height adjustment

098



Power Backrest

(Trendelenburg Tables only)

- Smooth, power control raises to 80° and lowers to flat

011



Power Footrest

- Smooth, power control lowers to 30° and raises to flat

090



Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087

