


Power, Imaging Table with Three-Section Top and Drop Window



MODEL 80073

- Three-piece top with pneumatic adjustable backrest, manual footrest and window drop
- Footrest is adjustable from either side and raises up to 30° from flat position
- Optional side rails & casters (see options on page 2)
- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser, cutter and positioning straps included
- Optional side rails & casters available
- 34" wide model available
- Model 80073: 500 lbs load capacity under normal use
- Model 80073-X: 600 lbs load capacity under normal use

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

80073 SPECIFICATIONS

	length	width	height range
overall	76"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
backrest	38"	27"	
footrest	12"	27"	
window	10"	10"	

80073-X Extra-Wide 34" width model



Upholstery Colors



OPTIONS

Side Rails (pair)

- 3 height positions
- Plunger height adjustment

098




Power Backrest

(Trendelenburg Tables only)

- Smooth, power control raises to 80° and lowers to flat

011



Power Footrest

- Smooth, power control lowers to 30° and raises to flat

090



Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087

