


Multi-Use, Imaging Table with Stirrups and Drop Window



MODEL 80079

- Three-section table with pneumatic adjustable backrest, manual footrest, stirrups and window drop
- Footrest is adjustable from either side and raises up to 30° from flat position
- Optional side rails & casters
- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser, cutter and positioning straps included
- Optional side rails & casters available
- 34" wide model available
- Model 80079: 500 lbs load capacity under normal use
- Model 80079-X: 600 lbs load capacity under normal use

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

80079 SPECIFICATIONS

	length	width	height range
overall	72"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
backrest	38"	27"	
footrest	12"	18 ¹ / ₂ "	
window	10"	10"	

80079-X Extra-Wide 34" width model (86.36 cm)



Upholstery Colors



OPTIONS

Side Rails (pair)

- 3 height positions
- Plunger height adjustment

098




Power Backrest

(Trendelenburg Tables only)

- Smooth, power control raises to 80° and lowers to flat

011



Power Footrest

- Smooth, power control lowers to 30° and raises to flat

090



Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087

