

## Mobile, Molded Top Bedside Cabinet

CLASSIC CLINTON  
• LAMINATE CASEWORK •



### MODEL 8720-A

- 1 self closing drawer and door
- 1 adjustable shelf
- 1-piece, vacuum formed gray haircell textured top rimmed on all sides
- Laminate drawer glides on nylon rollers
- Easy-clean laminate base
- Concealed, adjustable, soft-close, door hinges
- Cabinet is mounted on four, 2" dual wheel casters
- Satin chrome pulls
- Choice of solid color or woodgrain laminates

### 8720-A SPECIFICATIONS

18 W x 16 D x 29.25" H



California residents please be advised of the following, pursuant to Proposition 65:

**!** **WARNING:** This product can expose you to chemicals including 1, 3-Butadiene which is known to the State of California to cause cancer and reproductive harm. For more information go to: [www.P65Warning.ca.gov](http://www.P65Warning.ca.gov).

### Options




055 Door & Drawer Locks

### Laminate



Maple



Dark Cherry



Gray



Slate Gray