

Power-Glide, Hi-Lo, Treatment Table with Adjustable Backrest, Footrest & Stirrups

MODEL 87360

- One of the most affordable power tables in its class
- Features an adjustable backrest, friction-lock, footrest and stirrups
- Pneumatic backrest is adjustable from either side and raises to a maximum of approximately 80°
- Footrest is adjustable from either side and raises to approximately 75° from flat
- Retractable, stirrups swing out to a maximum of 35°
- Smooth-glide, electric power height adjustment
- Ultra heavy-duty base with solid, one-piece, steel ends
- Easy-clean, ABS plastic center shroud hides lift components and provides a clean look
- Rigid, one-piece, laser-cut steel top frame
- Durable gray powder-coated steel frame
- Self-lubricating pivot points
- Smooth 110 volt electric height adjustment
- Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1
- Fast easy motor access via ABS side panels
- Hands-free low voltage foot control
- Seamless, rounded corner, upholstered top
- 2" ultra-firm padding
- Customize with options
- Made in the USA by Clinton Industries
- 500 lbs. (227 kg) load capacity under normal use
- 3-year limited factory warranty from Clinton Industries

power-glide™
Hi-Lo Treatment Tables by Clinton



California residents please be advised of the following, pursuant to Proposition 65:

1
WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

87330 SPECIFICATIONS

Overall Size	76" L x 27" W x 18.5-34" H
Backrest	38" L x 27" W
Footrest	13" L x 18" W
Travel Space	97"

Upholstery Colors



OPTIONS

Positioning Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 straps per set (92" x 2")

080



Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

082



Locking Casters

- Four, heavy duty, rubber wheel, easy-roll, locking casters

087



030 Paper Dispenser


- Accommodates paper rolls from 18" (45.72cm) to 21" (53.34cm) wide

030

Paper Cutter

- Clear plastic with elastic strap and buckle

040



Paper Dispenser/Cutter Combo

- Paper dispenser and paper holder ordered as one option

034