

Molded Top, Mobile Treatment Cabinet with 2 Doors and 2 Drawers



CLASSIC CLINTON
• LAMINATE CASEWORK •

MODEL 8922-A


- Mobile cabinet with 2 doors, 2 drawers and 1 adjustable shelf
- Easy-clean, all-laminate construction
- One-piece, vacuum-formed, gray-haircell, textured, top rimmed on all sides
- Fully adjustable, soft-close, door hinges
- Smooth glide, drawers with one-piece metal sides and a load capacity of 75 lbs.
- Tough chip resistant edgebanding
- Tough, plastic bumpers on each side of the base
- Cabinet is mounted on four, 3" swivel casters (two locking)
- Satin chrome pulls
- Choice of solid color or woodgrain laminates
- Optional locks are available

8922-A SPECIFICATIONS

24.875 W x 18.375 D x 35.25" H




California residents please be advised of the following, pursuant to Proposition 65:

7
 **WARNING:** This product can expose you to chemicals including 1, 3-Butadiene which is known to the State of California to cause cancer and reproductive harm. For more information go to: www.P65Warning.ca.gov.

Options



055 Door & Drawer Locks



066 Door and Inside Latch Combo

Laminate



Maple



Dark Cherry



Gray



Slate Gray