

## Cabinet Style, Space Saver Table with Adjustable Backrest and Leg Rest



### MODEL 91002-27

- No additional charge for adjustable backrest
- All laminated construction
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Leg rest has 4 settings of approximately: 0°, 12°, 24° & 36°
- Two, closed-door, storage compartments with access from both sides
- 1 adjustable shelf behind each door
- Seamless rounded corner top
- 2" firm foam padding
- Heavyweight, knit-backed upholstery
- Ships assembled
- 350 lbs. load capacity under normal use (159 kg)

### SPECIFICATIONS

91002-27	54"/72" Extended L x 27" W x 31" H
Top	54" L x 27" W
Leg Rest	18" L x 22" W


### Prop 65 Warning



**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

### Options



### Laminate Colors



### Upholstery Colors

