

## Cabinet Table with Center Drawers & Closed-Door Storage



### MODEL 91009

- No upcharge for adjustable backrest
- All laminated construction
- Choose standard adjustable backrest or flat top option
- No additional charge for adjustable backrest or flat top
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Center drawer section has three drawers that hold up to 75 lbs.
- Drawer section is easily reversible for access from either side
- Drawers have smooth glides for easy open and close
- Two, closed-door, storage compartments with access from both sides
- 1 adjustable shelf behind each door
- Seamless, rounded corner top
- 2" firm foam padding
- Heavyweight, knit-backed upholstery
- Ships assembled
- 425 lbs. load capacity under normal use (193 kg)



### SPECIFICATIONS

91009-27	72" L x 27" W x 31" H
91009-30	72" L x 30" W x 31" H
91009-28	78" L x 27" W x 31" H
91009-38	78" L x 30" W x 31" H

### Prop 65 Warning

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

### Options



20-Small Pillow

60-Large Pillow

040-Paper Cutter

030-Paper Dispenser

034-Paper Dispenser & Cutter Combo

009-Nose Hole

055-Door & Drawer Locks

088-Flat Top

### Laminate Colors



Dark Cherry



Maple



Gray

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black