

Panel Leg Series, Treatment Table with Sliding Doors



MODELS: 91040-27 & 91040-30

- No upcharge for adjustable backrest
- All laminate construction
- Sleek panel leg design
- Adjustable backrest with 6 settings of approximately: 27°, 36°, 43°, 50°, 55° & 60°
- Sliding door storage with access from either side
- Scratch resistant topcoat finish
- Ships flat for easy handling and extra savings
- Fast and easy assembly
- Seamless, rounded corner top
- 2" firm foam padding (5 cm)
- Heavyweight, knit-backed upholstery
- 400 lbs. load capacity under normal use (181 kg)

SPECIFICATIONS

91040-27	72" L x 27" W x 31" H
91040-30	72" L x 30" W x 31" H



Prop 65 Warning

! **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

2

Options



034– Paper Dispenser & Holder Combo

040– Paper Holder

030– Paper Dispenser

088–Flat Top


Laminate Colors



Gray



Maple



Dark Cherry

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black