

Cabinet Style, Treatment Table with Drawers & Doors



MODEL 9105-30

- No upcharge for adjustable backrest
- Features 3 drawers, adjustable backrest and 1 storage compartment with 1 adjustable shelf
- Easy-clean laminate surfaces
- Flat panel drawers with all steel sides
- Adjustable doors with concealed soft-close hinges
- Ships assembled
- Seamless, rounded corner top
- Heavyweight, knit-backed upholstery
- 2" firm foam padding (5 cm)
- 425 lbs. load capacity under normal use

SPECIFICATIONS

9105-30

72" L x 30" W x 31" H

California residents please be advised of the following, pursuant to Proposition 65:

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Options



20- Small Pillow

60- Large Pillow

034- Paper Dispenser & Holder Combo

040- Paper Holder

030- Paper Dispenser

009- Nose Hole

055- Drawer Locks

066- Door & Inside Latch

088- Flat Top

Laminate Colors



Dark Cherry



Maple



Gray



Slate Gray

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black