

## Extra Long, Step-Up Table with 2 Sliding Doors



### MODEL 9777-38

- No upcharge for adjustable backrest
- Sturdy, all laminate and steel construction
- Extra-long 78" length
- Handy, steel stool with slip-resistant top pulls out for use, and neatly stores out of the way when not in use
- Stool can be easily changed on site for use on either side of the table
- Adjustable backrest with 2" for firm foam padding
- 2 sliding door storage, accessible from either side of the table
- Non-marking glides for easy positioning
- 400 lbs. weight capacity under normal use

### 9777-38 SPECIFICATIONS

78 L x 30" W x 31" H




California residents please be advised of the following, pursuant to Proposition 65:

**!** **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Options



055- Door Locks  
(Must select two)



034- Paper Dispenser & Holder Combo



040- Paper Holder



030- Paper Dispenser



088-Flat Top

### Laminate Colors



Dark Cherry



Maple



Gray

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black