

Cabinet Table with Center Drawers & Open Storage



MODEL 98013

- Choose standard adjustable backrest or flat top option at no additional charge
- All laminated construction
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Center drawer section has three drawers that hold up to 75 lbs.
- Drawer section is easily reversible for access from either side
- Drawers have smooth glides for easy open and close
- Open storage with shelf on each side of center section
- Seamless rounded corner top
- 2" firm foam padding
- Heavyweight, knit-backed upholstery
- Ships assembled
- 425 lbs. load capacity under normal use (193 kg)



SPECIFICATIONS

98013-27	72" L x 27" W x 31" H
98013-30	72" L x 30" W x 31" H
98013-28	78" L x 27" W x 31" H
98013-38	78" L x 30" W x 31" H

Prop 65 Warning

! **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Options



20- Small Pillow

60- Large Pillow

040- Paper Cutter

030- Paper Dispenser

034- Paper Dispenser & Cutter Combo

009- Nose Hole

055- Drawer Locks

088- Flat Top

Laminate Colors



Maple



Gray

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black