

## Cabinet Table with Open Access Storage



### MODEL 98022

- Choose standard adjustable backrest or flat top option at no additional charge
- All laminated construction
- Backrest has 6 settings of approximately: 60°, 55°, 50°, 43°, 36° & 27°
- Center section has an adjustable shelf on both sides with pass-through access
- Open storage on each side of center section
- Seamless rounded corner top
- 2" firm foam padding
- Heavyweight, knit-backed upholstery
- Ships assembled
- 425 lbs. load capacity under normal use (193 kg)

### SPECIFICATIONS

98022-27	72" L x 27" W x 31" H
98022-30	72" L x 30" W x 31" H
98022-28	78" L x 27" W x 31" H
98022-38	78" L x 30" W x 31" H

### Prop 65 Warning

**WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

### Options



20-Small Pillow

60-Large Pillow

040-Paper Cutter

030-Paper Dispenser

034-Paper Dispenser & Cutter Combo

009-Nose Hole

088-Flat Top

### Laminate Colors



Maple



Gray

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black