


Heavy Base, 4-Hook IV Pole



MODEL: IV-314

- 4-hook top
- Chrome plated steel pole/knob-lock adjustment
- 4-leg epoxy coated steel base with bumper guards
- Detachable ram's horn hooks with plastic tips
- 2" rubber wheel ball bearing casters
- See available options on next page
- 45 lbs. load capacity under normal use


Shipping Weight: 18 lbs.

SPECIFICATIONS

Base Diameter—22"

Height Range: 52"—91"

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

OPTIONS:



IV-55 Extra Hooks & Organizer Hooks



IV-41 Universal Clamp



IV-50 Pole Outlet Strip



Heavy Duty Black Epoxy Wire Basket
IV-51B - 14" x 8" x 8.5"
IV-52B - 12" x 6" x 6.5"



Heavy Duty Stainless Steel Wire Basket
IV-51S - 14" x 8" x 8.5"
IV-52S - 12" x 6" x 6.5"



IV-57 Single Ram's Horn



IV-43 Steering Handle



IV-44 Oxygen Cylinder Holder



IV-56 Drainage Hook

