


## Economy 4-Leg, 2-Hook IV Pole



### MODEL: IV-40

- 2-hook top
- Chrome plated steel pole/twist-lock adjustment
- Epoxy coated, 4-leg base
- Detachable ram's horn hooks
- 4, easy-roll, ball bearing casters
- 30 lbs. load capacity under normal use

*Shipping Weight: 17 lbs.*

#### SPECIFICATIONS

Base Diameter—22"

Height Range: 46"—83"

#### OPTIONS




*IV-41 Universal Clamp*



*IV-43 Steering Handle*

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1