


Six-Leg, 4-Hook, Heavy Base Infusion Pump Stand



MODEL: IV-604

- 4-hook top
- Chrome plated steel pole
- Powder-coated, ultra heavy, tip resistant steel 6-leg base
- Detachable ram's horn hooks with plastic tips
- 4" Heavy duty, rubber wheel, ball bearing casters
- Bent, low center of gravity style base
- Extra rigid, pole-in-tube mounting system

See available options

Maximum load capacity 65 lbs. under normal use


Shipping Weight: 37 lbs.

SPECIFICATIONS

Base Diameter—23"

Height Range: 51"—92"

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

OPTIONS:



IV-46 Pump Tray Support



IV-42 Infusion Pump Frame



IV-55 Extra Hooks and Organizer Hooks




IV-41 Universal Clamp



Heavy Duty Black Epoxy Wire Basket
IV-51B – 14" x 8" x 8.5"
IV-52B – 12" x 6" x 6.5"



Heavy Duty Stainless Steel Wire Basket
IV-51S – 14" x 8" x 8.5"
IV-52S – 12" x 6" x 6.5"



IV-57 Single Ram's Horn



IV-44 Oxygen Cylinder Holder



IV-43 Steering Handle



IV-56 Drainage Hook