


## Stainless Steel\*, Five-Leg, 2-Hook IV Pole



### MODEL: IVS-352

- 2-hook top
- 5-leg cast aluminum base
- Detachable ram's horn hooks
- Soft roll, dual wheel casters
- Stainless steel\* pole/knob-lock height adjustment
- See available accessories on next page
- 45 lbs. load capacity under normal use

\*Clinton stainless steel items are not M.R.I. compatible.


*Shipping Weight: 15 lbs.*

#### SPECIFICATIONS

Base Diameter—23"

Height Range: 60"—100"

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

**OPTIONS:**



*IV-46 Pump Tray Support*



*IV-41 Universal Clamp*



*IV-50 Pole Outlet Strip*



*Heavy Duty Black Epoxy Wire Basket*  
*IV-51B – 14" x 8" x 8.5"*  
*IV-52B – 12" x 6" x 6.5"*



*Heavy Duty Stainless Steel Wire Basket*  
*IV-51S – 14" x 8" x 8.5"*  
*IV-52S – 12" x 6" x 6.5"*



*IV-57 Single Ram's Horn*



*IV-44 Oxygen Cylinder Holder*



*IV-56 Drainage Hook*



*IV-43 Steering Handle*



*IV-55 Extra Hooks and Organizer Hooks*



*IV-48 2 lb. Base Stability Weight / V-485 2 lb. Base Stability Weight*