

Premium Mobile Instrument Stand



MODEL: M-20

- Stable, 5-leg cast aluminum base
- Easy positioning with no need to tip style base
- Removable 19" x 12-3/4" x 3/4" stainless steel tray
- Dual wheel soft roll casters
- 1" chrome-plated pole and tray frame
- Knob height adjustment
- Maximum load capacity 15 lbs.


Shipping Weight: 17 lbs.

SPECIFICATIONS

Base: 23" D

Height Range: 31"—50"

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1