

## Double Post Mayo Stand



### MODEL: M-22

- Removable stainless steel tray
- 2 ball bearing, rubber casters
- 1" chrome-plated pole & base
- Tip style position base
- Knob height adjustment
- Removable 19" x 12-3/4" x 3/4" stainless steel tray

15 lbs. load capacity under normal use


*Shipping Weight: 19 lbs.*

#### SPECIFICATIONS

Base: 21" W x 20" D

Height Range: 34"—54"

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1