

Value Instrument Stand



MODEL: M-29

- Stable, 5-leg nylon base
- Easy positioning with no need to tip style base
- Removable 19" x 12-3/4" x 3/4" stainless steel tray
- Dual wheel casters
- 1" chrome plated pole and frame
- Knob height adjustment
- 15 lbs. load capacity under normal use


Shipping Weight: 17 lbs.

SPECIFICATIONS

Base: 21" D

Height Range: 31"—50"

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1