

Mobile *Stainless Steel Instrument Stand



MODEL: MS-29

- Stainless steel frame and base
- Easy positioning with no need to tip style base
- Removable, 19" x 12-3/4" x 3/4" stainless steel tray
- Moves easily on 4 ball bearing casters
- Stable heavy base design
- Knob lock height adjustment

*Clinton stainless steel items are not M.R.I. compatible.


Shipping Weight: 18 lbs.

SPECIFICATIONS

Base: 23" D

Height Range: 33"—49"

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1