

***Stainless Steel Stool with Casters**

MODEL SS-2162

- 15" diameter stainless steel seat
- Extra wide base for stability
- Polished stainless steel seat with recessed black vinyl insert
- 2" Rubber wheel ball bearing casters
- Round foot ring
- All welded stainless steel frame


**Clinton stainless steel items are not M.R.I. compatible*



SPECIFICATIONS

Seat Diameter	15"
Height Range	21.5"-28.5"

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1