

*Stainless Steel Stool with Casters

MODEL SS-2172

- 4" thick padded seat for comfort
- Smooth machine screw height adjustment
- 2" Rubber wheel ball bearing casters
- Round foot ring
- All welded stainless steel frame
- Easy clean, premium, knit-backed upholstery
- Extra wide base for stability


**Clinton stainless steel items are not M.R.I. compatible*



SPECIFICATIONS

| | |
|---------------|-----------|
| Seat Diameter | 14" |
| Height Range | 25"-29.5" |

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Upholstery Colors

