

Sharps Disposal

Phlebotomy and
Multi-Purpose Containers

Multi-Purpose Containers

In-Room™ System

Large Volume Containers

ChemoSafety™ Containers

PGII Large Volume and
ChemoSafety™ and
Transportable Sharps
Container

Transportable Sharps
Containers

Brackets and Holders

Wall Enclosures and Carts

Foot Pedal Carts and Floor
Brackets

For general enquiries call +353 14381600
or fax +353 14073174






[Click here for a full list
of International contacts](#)

Phlebotomy and Multi-Purpose Containers


Key

● Red	● Yellow	● Transparent Red	○ Transparent Yellow	○ Transparent	○ Moulded with recycled medical waste plastic bin
--	--	---	--	--	---

Phlebotomy Sharps Containers

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit	
8901SA	Phlebotomy Sharps Container adapts to present blood drawing trays and techniques	0.5	10.8 x 10.8 x 11.4cm	1	●	100	
8900SA 8906	Phlebotomy Sharps Container adapts to present blood drawing trays and technique	1	15.9 x 10.8 x 11.4cm	1 11	● ●	100 100	
1316676202 1316676228	Small Monoject™ Sharps Container, clear top for easier assessment of fill level	1.5	13.3 x 20 x 13.3cm	-	● ●	40	
1522SA	Phlebotomy Sharps Container with rotor opening	2	18.4 x 12 x 16.5cm	2	●	60	
1525SA	Sharps Container for all sharps waste with rotor opening	7.4	25.4 x 18.4 x 26.7cm	3	●	20	

Brackets & Holders






Code 8910

Code 1523SA

Code 8974

Code 8912SA

Multi-Purpose Sharps Containers

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit	
1100676376 1100676517	Monoject™ Sharps Container with chimney drop top lid. Clear top on translucent container for easier assessment of fill level. Can be decontaminated in an Autoclave.	4	13.9 x 17.3 x 28.6cm	12	● ●	40	
1100676384 1100676533	Monoject™ Sharps Container with chimney drop top lid. Clear top on translucent container for easier assessment of fill level. Can be decontaminated in an Autoclave.	8	23.7 x 17.3 x 28.6cm	12	● ●	20	
1100676392 1100676558	Monoject™ Sharps Container with chimney drop top lid. Clear top on translucent container for easier assessment of fill level. Can be decontaminated in an Autoclave.	11	40.8 x 17.3 x 28.6	12	● ●	10	

Brackets & Holders



Code 8881676012

Multi-Purpose Containers

Multi-Purpose Sharps Containers

	Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
	8950SA	Multi-purpose round Sharps Container with rotor opening	4.6	17 (H) x 22.2cm (D)	6	●	40
	8970 8970Y	Multi-purpose container for all sharps waste with rotor opening	7.4	25.4 x 18.4 x 26.6cm	3 8 9	● ●	20
	8990SA	Multi-purpose container with clear hinged lid	7.4	25.4 x 18.4 x 26.6cm	3 8 9	●	20
	8990Y 8996	Multi-purpose container with special gasket hinged lid	7.4	25.4 x 18.4 x 26.6cm	3 8 9	● ●	20
	89671 8967Y	Multi-purpose container with horizontal drop opening lid	7.4	32.4 x 18.4 x 26.6cm	3 8 9	○	20
	8525Y 8527R	Multi-purpose container with rotor opening and hinged lid to accommodate a variety of sharps sizes and provide temporary and permanent closure lid	11	35 x 15.2 x 35cm	5 10	● ●	10

Brackets & Holders



Code 8974



Code 8529



Code 8975



Code 8528C

In-Room™ System

SharpSafety™ Lid Containers

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
8506Y	Sharps Container with SharpSafety™ lid and counter-balanced door	4.6	31.7 x 14 x 27.3cm	4 7 15	○	20
8533SA 8534SA	Sharps Container with SharpSafety™ lid and counter-balanced door	7.4	30 x 16.5 x 35cm	10 16	○ ●	10
8536SA 8537SA 8537Y	Sharps Container with SharpSafety™ lid and counter-balanced door	11	42 x 15.2 x 35cm	5 10 16	○ ● ○	10



Mailbox Style Lid Containers

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
85021	Mailbox style lid Sharps Container can be kept closed when container is not in use	1.8	15.9 x 12 x 27.3cm	4 7 15	○	20
85121	Mailbox style lid Sharps Container can be kept closed when container is not in use	4.6	28 x 12 x 27.3cm	4 7 15	○	20




Always-Open Lid Containers

Code	Description	Litres Capacity	Dimensions HxDxW	Brackets & Holders	Colour	Sales Unit
850201 850301	Sharps Container with always-open lid provides convenient, quick and easy point-of-use sharps disposal	1.8	15.9 x 12 x 27.3cm	4 7 15	○ ●	20
851201 85120Y	Sharps Container with always-open lid provides convenient, quick and easy point-of-use sharps disposal	4.6	28 x 12 x 27.3cm	4 7 15	○ ○	20
85321R	Sharps Container with always-open lid provides convenient, quick and easy point-of-use sharps disposal	7.4	29.2 x 15.2 x 35cm	10 16	○ ●	10
85221 85221R 8522Y	Sharps Container with always-open lid provides convenient, quick and easy point-of-use sharps disposal	11	41.3 x 15.2 x 35cm	5 10 16	○ ● ○	10



Brackets & Holders



Code 8515

Code 8529

Code 8518X

Code 8528C

Code 85301H