

ENDO STITCH™ SUTURING DEVICES

SILS™ Stitch

75° articulation, 360° tip rotation



ORDER CODE	DESCRIPTION	BOX QTY
SILSSTITCH	10 mm Single Use Articulating Endoscopic Suturing Device Compatible with the current EndoStitch™ Loading Units 75° articulation, 360° tip rotation, shaft length 43 cm	3

Endo Stitch™ 10 mm

Suturing Device



ORDER CODE	DESCRIPTION	BOX QTY
173016	Endo Stitch™ 10 mm Suturing Device, 34 cm Shaft Length	3

Endo Stitch™ Triple-Stitch Loading Unit

Single Use Loading Units for use with SILS™ Stitch and Endo Stitch™ 10 mm Suturing Devices. Each Loading Unit contains three strands of suture, each attached to a 9 mm needle.



ORDER CODE	SUTURE MATERIAL	USP SIZE	METRIC SIZE	LENGTH	BOX QTY
170043	Surgidac™ Ctd, Brd Polyester non absorbable	0	3.5	120 cm / 48"	3
170040	Surgidac™ Ctd, Brd Polyester non absorbable	0	3.5	18 cm / 7"	3
170041	Surgidac™ Ctd, Brd Polyester non absorbable	2-0	3	18 cm / 7"	3
170044	Surgidac™ Ctd, Brd Polyester non absorbable	2-0	3	120 cm / 48"	3