



8507SA



8540SA



85401H



8537MW



85301H



8506SA



85161H



8555SA



8556H

Item	Description	Qty/Case	Item	Description	Qty/Case
3 AND 4 GALLON SHARPSTAR™ IV CONTAINERS			5 QUART SHARPSTAR™ CONTAINERS		
8540SA	4 Gallon, Clear	10	8506SA	5 Quart, Clear	20
8541SA	4 Gallon, Trans Red	10	8507SA	5 Quart, Trans Red	20
8544SA	3 Gallon, Clear	10	8507MW	5 Quart, Recycled Plastics	20
8545SA	3 Gallon, Trans Red	10	8519C	Non-Locking Wire Bracket	5
8543C	Locking Wire Bracket	5	8518X	Locking Wire Bracket	5
85401H	3 & 4 Gallon Wall Enclosure	1	85161H	5 Quart Wall Enclosure	1
8549SA	3 & 4 Gallon Tabletop Holder	5	85165H	5 Quart Wall Enclosure	5
2 AND 3 GALLON SHARPSTAR™ CONTAINERS			8556H	5 Quart Wall Enclosure with Container and Glove Box	1
8536SA	3 Gallon, Clear	10	8515	5 Quart Tabletop Holder	5
8537SA	3 Gallon, Trans Red	10	GLOVE BOXES		
8537MW	3 Gallon, Recycled Plastics	10	8555SA	Stackable Glove Box, Side Load	10
8533SA	2 Gallon, Clear	10	8550LG	Top Loading Glove Box	10
8534SA	2 Gallon, Trans Red	10	8554B	Glove Box Mounting Bracket	5
8524C	Non-Locking Wire Bracket	5			
8528C	Locking Wire Bracket	5			
85301H	2 & 3 Gallon Wall Enclosure	1			
8529	3 Gallon Tabletop Holder	5			
8539SA	2 Gallon Tabletop Holder	5			



SharpStar™ In Room™ Disposal System



COVIDIEN, COVIDIEN with logo and all ™ marked brands are trademarks of Covidien AG or an affiliate.
© 2009 Covidien. All rights reserved.
H6364 XX 0209

For Ordering Information Contact CHARTECH INC.
1-877-658-1849 Fax 770-442-2843
sales@chartech.net

Anti-kickback lid is designed to prevent sharps from accidentally exiting the container.

No other manufacturer offers this important engineering control.



SharpStar™ In-Room™ Disposal System

The Covidien SharpSafety Division is committed to providing real solutions for all of your sharps disposal needs. The SharpStar™ family of containers is engineered with controls to prevent overfilling and limit accidental or intentional access to the container contents. These containers are offered in multiple sizes from 5 quarts to 4 gallons to accommodate the increasing volume and size of safety engineered sharps devices used today.



- Counterbalanced lid is designed to allow for easy one-handed disposal of sharps.

- Safety lid design restricts hand access to contents.

- Automatic shut off prevents overfilling and a bold, red "FULL" indicator alerts clinician to change the container.

- If the full indicator door is forced open, the anti-kickback feature is designed to prevent sharps from being accidentally ejected from the container.

- Locking tabs keep hands away from contents during final closure.

How will the SharpStar™ In-Room™ System help me to comply with the Needlestick Safety & Prevention Legislation?
OSHA (29 CFR 1910.1030)

