



CLEANROOM

C O N N E C T I O N

For Pricing and Ordering Please Call 800.616.5319



Needles and Syringes
Prefilled I.V. Flush Syringes, Blood Collection
and Sharps Disposal Products

1-800-962-9888



Table of Contents

Needles and Syringes, Prefilled I.V. Flush Syringes, Blood Collection

MAGELLAN™ SAFETY PRODUCTS

Magellan™ Safety Needles	5
Magellan™ Safety Needle & Syringe Combinations	6
Magellan™ Insulin & Tuberculin Safety Syringes.....	7
Safety Blood Collection Device	7
3mL, 6mL, and 12mL Monoject™ Safety Syringes.....	8
Insulin & Tuberculin Monoject™ Safety Syringes.....	9

STANDARD NEEDLES & SYRINGES

Hypodermic Needles	10-12
Blunt Cannula	13
Low Dead Space Needles.....	13
Insulin Syringes	14
Tuberculin Syringes.....	15
Allergy Trays	16
1mL, 3mL, 6mL, 12mL, 20mL, 35mL and 60mL Syringes.....	16-20
Pharmacy Trays	21
Non-Sterile Syringes.....	21
Oral Medication Syringes.....	22
Piston Syringes.....	22
Control Syringes	22

PREFILLED FLUSH SYRINGES

Saline & Heparin I.V. Flush Syringes	23
--	----

SAFETY I.V. ACCESS PRODUCTS

Monoject SmarTip™ Needleless Med-Prep Cannula.....	24
Lifeshield™* Blunt Cannula & Syringes	24
I.V. Access Cannula.....	25
I.V. Security Clips	25

PHARMACY PRODUCTS

Pharmacy Trays	26
Filter Connector	26
Medication Transfer Needle.....	26
Filter Needles	27
Filter Aspirator	27

INFUSION PRODUCTS

Argyle™ EZ-FLO Stopcocks.....	28
Argyle™ Intermittent Infusion Plug.....	28

MONOJECT™ ACCESSORIES

Syringe Tip Caps.....	29
Safety Syringe Tip Caps	29
Finger Grip Extenders.....	29

SAFETY PHLEBOTOMY PRODUCTS

Magellan™ Safety Blood Collection Device.....	30
Defender™ Safety Needle Holder	30
Safety Winged Sets with Female Luer.....	31
Safety Winged Sets with Multi-Sample Luer.....	31
Safety Winged Sets with Holder Attached.....	31

BLOOD COLLECTION AND TRANSFER SETS

Luer Lock Collection Set	33
Safety Transfer Device	32
Safety Collection Device	32
Tube Holder with Adapter	33
Blood Collection Assembly with Monoject BlunTip™ Safety I.V. Access Cannula.....	33

STANDARD PHLEBOTOMY PRODUCTS

Blood Collection Needles.....	34
Luer Adapter	34
Blood Collection Needle & Tube Holders.....	34

Table of Contents CONTINUED

Sharps Disposal

SHARPSTAR™ IN-ROOM™ SYSTEMS & CONTAINERS	36	PG II-RATED LARGE VOLUME CONTAINERS	44
IN-ROOM™ CONTAINERS WITH ALWAYS-OPEN LIDS	37	TRACE CHEMOTHERAPY SHARPS CONTAINERS	45
IN-ROOM™ CONTAINERS WITH MAILBOX STYLE LIDS	37	NON-BIOHAZARDOUS WASTE CONTAINERS	46
SHARPS-A-GATOR™ TORTUOUS PATH SHARPS CONTAINERS	38	DIALYSAFETY CONTAINERS	46
GATORGUARD™ IN-ROOM™ SYSTEMS & CONTAINERS	38	TRANSPORTABLE SHARPS CONTAINERS	46
RENEWABLES™ SHARPS CONTAINERS	39	MULTI-PURPOSE SHARPS CONTAINERS	47-48
IN-ROOM™ SYSTEM WALL ENCLOSURES & GLOVE BOXES	40	PHLEBOTOMY CONTAINERS	49
GATORGUARD™ WALL CABINETS, GLOVE BOXES & HOLDERS	40	RCRA HAZARDOUS WASTE CONTAINERS	50
LARGE VOLUME SHARPS CONTAINERS	41	PHARMACEUTICAL WASTE CONTAINERS	52
LARGE VOLUME SHARPS CART™ FOOT PEDAL CARTS & FLOOR BRACKETS	42	BRACKETS, HOLDERS & ACCESSORIES	52-54
SHARPS-A-GATOR™ LARGE VOLUME CONTAINERS, FLOOR STANDS, BRACKETS & CARTS	43	HARDWARE & CONTAINER CROSS REFERENCE GUIDE	55-62

Magellan™ Safety Products LATEX FREE



All Magellan™ Safety Products include many innovative features:

- A choice of easy, one-handed activation methods that minimizes technique change and in-service time
- A fixed safety shield that eliminates the ability to remove the safety device which increases OSHA compliance
- Bevel orientation so that if the safety shield is up, the needle bevel is up
- An audible and tactile confirmation click that confirms the safety device has been activated for ultimate protection
- Robust safety shield that completely encases the needle and remains in effect after disposal
- Color coded hub to correspond with ISO standards
- A full range of gauges and lengths
- Manufactured in the United States



Magellan™ Hypodermic Safety Needles - Sterile

- Needle-based safety technology
- Intuitive, one-handed activation
- Bevel orientation, so that if safety shield is up, the needle bevel is up
- Self-leveling needle sheath
- Compatible with any standard luer lock syringe

Code	Description	Box/Case
8881850810	Safety Needle 18 x 1"	50/500
8881850815	Safety Needle 18 x 1 1/2"	50/500
8881850910	Safety Needle 19 x 1"	50/500
8881850915	Safety Needle 19 x 1 1/2"	50/500
8881850010	Safety Needle 20 x 1"	50/500
8881850015	Safety Needle 20 x 1 1/2"	50/500
8881850158	Safety Needle 21 x 5/8"	50/500
8881850110	Safety Needle 21 x 1"	50/500
8881850115	Safety Needle 21 x 1 1/2"	50/500
8881850210	Safety Needle 22 x 1"	50/500
8881850215	Safety Needle 22 x 1 1/2"	50/500
8881850358	Safety Needle 23 x 5/8"	50/500
8881850310	Safety Needle 23 x 1"	50/500
8881850558	Safety Needle 25 x 5/8"	50/500
8881850510	Safety Needle 25 x 1"	50/500
8881850515	Safety Needle 25 x 1 1/2"	50/500

Safety Needle Gauge Color Chart

18 Gauge Pink	19 Gauge Cream	20 Gauge Yellow	21 Gauge Green	22 Gauge Black	23 Gauge Blue	25 Gauge Orange	27 Gauge Gray	28 Gauge Teal
------------------	-------------------	--------------------	-------------------	-------------------	------------------	--------------------	------------------	------------------

All Magellan Safety Products are ISO Color Coded

Magellan™ Safety Products LATEX FREE



Magellan™ Hypodermic Safety Needle & Syringe Combinations - Sterile

- Needles “pre-attached” to syringes
- Needle-based safety technology
- Intuitive, one-handed activation
- Bevel orientation, so that if safety shield is up, the needle bevel is up
- Self-leveling needle sheath
- Low dead space

Code	Description	Box/Case
8881811310	1mL Syringe 23 x 1" Safety Needle	50/200
8881811558	1mL Syringe 25 x 5/8" Safety Needle	50/200
8881811510	1mL Syringe 25 x 1" Safety Needle	50/200
8881833810	3mL Syringe 18 x 1" Safety Needle	50/400
8881833010	3mL Syringe 20 x 1" Safety Needle	50/400
8881833015	3mL Syringe 20 x 1 1/2" Safety Needle	50/400
8881833110	3mL Syringe 21 x 1" Safety Needle	50/400
8881833115	3mL Syringe 21 x 1 1/2" Safety Needle	50/400
8881833210	3mL Syringe 22 x 1" Safety Needle	50/400
8881833215	3mL Syringe 22 x 1 1/2" Safety Needle	50/400
8881833310	3mL Syringe 23 x 1" Safety Needle	50/400
8881833558	3mL Syringe 25 x 5/8" Safety Needle	50/400
8881833510	3mL Syringe 25 x 1" Safety Needle	50/400
8881833515	3mL Syringe 25 x 1 1/2" Safety Needle	50/400
8881866810	6mL Syringe 18 x 1" Safety Needle	50/200
8881866015	6mL Syringe 20 x 1 1/2" Safety Needle	50/200
8881866110	6mL Syringe 21 x 1" Safety Needle	50/200
8881866115	6mL Syringe 21 x 1 1/2" Safety Needle	50/200
8881866215	6mL Syringe 22 x 1 1/2" Safety Needle	50/200
8881822810	12mL Syringe 18 x 1" Safety Needle	50/200
8881822015	12mL Syringe 20 x 1 1/2" Safety Needle	50/200
8881822110	12mL Syringe 21 x 1" Safety Needle	50/200
8881822115	12mL Syringe 21 x 1 1/2" Safety Needle	50/200
8881822215	12mL Syringe 22 x 1 1/2" Safety Needle	50/200

Magellan™ Safety Products LATEX FREE



Magellan™ Insulin and Tuberculin Safety Syringes - Sterile

- Permanently attached needle
- Needle-based safety technology
- Intuitive, one-handed activation
- Clear, bold unobstructed graduations
- Fixed safety shield that is not an accessory
- Bevel orientation, so that if safety shield is up, the needle bevel is up
- Low dead space

Code	Description	Box/Case
INSULIN		
8881892930	3/10mL 29 x 1/2"	50/500
8881892950	1/2mL 29 x 1/2"	50/500
8881892910	1mL 29 x 1/2"	50/500
8881893030	3/10mL 30 x 5/16"	50/500
8881893050	1/2mL 30 x 5/16"	50/500
8881893010	1mL 30 x 5/16"	50/500
TUBERCULIN		
8881882712	1mL 27 x 1/2"	50/500
8881882812	1mL 28 x 1/2"	50/500

Note: 3/10mL = 30 units

1/2mL = 50 units

1mL = 100 units



Magellan™ Safety Blood Collection Device - Sterile

- Ready to use 3-in-1 safety design
- One piece safety blood collection needle and blood tube holder combination that is permanently attached
- Intuitive, one-handed activation
- Fixed safety shield that is not an accessory
- Bevel orientation, so that if safety shield is up, the needle bevel is up
- Single use to help comply with OSHA regulations
- Contains latex

Code	Description	Case
8881227121	21 x 1 1/4" Safety Blood Collection Device	300
8881227122	22 x 1 1/4" Safety Blood Collection Device	300

Monoject™ Safety Syringes LATEX FREE

All Monoject™ Safety Syringes offer the following features and benefits:

- Barrel-based safety technology
- Robust safety shield that covers the entire needle during transport and after disposal
- Transport position that allows clinicians to safely transport medication and eliminates the need to recap the needle
- Rigid pack packaging that may be autoclaved or gas sterilized (Insulin and TB are available only in soft pack packaging)
- Luer lock tip is compatible with any standard needle (up to 1 1/2" length) making the combination OSHA compliant



3mL Safety Syringes - Sterile

- Rigid pack packaging
- Luer lock syringe tip

Code	Description	Box/Case
8881533056	3mL Safety Syringe 20 x 1 1/2"	100/500
8881533130	3mL Safety Syringe 21 x 1"	100/500
8881533155	3mL Safety Syringe 21 x 1 1/2"	100/500
8881533239	3mL Safety Syringe 22 x 1"	100/500
8881533254	3mL Safety Syringe 22 x 1 1/2"	100/500
8881533338	3mL Safety Syringe 23 x 1"	100/500
8881533510	3mL Safety Syringe 25 x 5/8"	100/500
8881533908	3mL Safety Syringe Luer Lock Tip	100/500



6mL Safety Syringes - Sterile

- Rigid pack packaging
- Luer lock syringe tip

Code	Description	Box/Case
8881566015	6mL Safety Syringe 20 x 1 1/2"	50/300
8881566023	6mL Safety Syringe 21 x 1 1/2"	50/300
8881566007	6mL Safety Syringe Luer Lock Tip	50/300



12mL Safety Syringes - Sterile

- Rigid pack packaging
- Luer lock syringe tip

Code	Description	Box/Case
8881522018	12mL Safety Syringe 20 x 1 1/2"	50/300
8881522026	12mL Safety Syringe 21 x 1 1/2"	50/300
8881522000	12mL Safety Syringe Luer Lock Tip	50/300

Monoject™ Safety Syringes LATEX FREE



Insulin Safety Syringes - Sterile

- Soft pack packaging
- Low dead space
- Permanently attached needle
- For U-100 Insulin

Code	Description	Box/Case
8881511144	3/10mL Insulin Safety Syringe 29 x 1/2"	100/500
8881511136	1/2mL Insulin Safety Syringe 29 x 1/2"	100/500
8881511110	1mL Insulin Safety Syringe 29 x 1/2"	100/500
8881511344	3/10mL Insulin Safety Syringe 30 x 5/16"	100/500
8881511336	1/2mL Insulin Safety Syringe 30 x 5/16"	100/500
8881511310	1mL Insulin Safety Syringe 30 x 5/16"	100/500



Tuberculin Safety Syringes - Sterile

- Soft pack packaging
- Low dead space
- Permanently attached needle

Code	Description	Box/Case
8881511235	1mL Tuberculin Safety Syringe 25 x 5/8"	100/500
8881511201	1mL Tuberculin Safety Syringe 28 x 1/2"	100/500



Safety Syringe Tip Caps - Sterile

- Individually packaged
- May be used with Monoject Safety Syringe
- Integral luer plug fits regular and luer lock syringes

Code	Description	Box/Case
8881682010	Safety Syringe Tip Cap	50/250

Monoject™ Standard Needles & Syringes LATEX FREE



SoftPack Hypodermic Needle with Polypropylene Hub - Sterile

- Soft pack packaging
- Polypropylene luer lock hub

Code	Description	Box/Case
1188818100	Hypo Needle 18 x 1" A	100/1000
1188818112	Hypo Needle 18 x 1 1/2" A	100/1000
1188818113	Hypo Needle 18 x 1 1/2" B	100/1000
1188819100	Hypo Needle 19TW* x 1" A	100/1000
1188819112	Hypo Needle 19TW* x 1 1/2" A	100/1000
1188820100	Hypo Needle 20 x 1" A	100/1000
1188820112	Hypo Needle 20 x 1 1/2" A	100/1000
1188821100	Hypo Needle 21 x 1" A	100/1000
1188821112	Hypo Needle 21 x 1 1/2" A	100/1000
1188822100	Hypo Needle 22 x 1" A	100/1000
1188822112	Hypo Needle 22 x 1 1/2" A	100/1000
1188823012	Hypo Needle 23 x 1/2" IN	100/1000
1188823100	Hypo Needle 23 x 1" A	100/1000
1188825058	Hypo Needle 25 x 5/8" A	100/1000
1188825100	Hypo Needle 25 x 1" A	100/1000
1188825112	Hypo Needle 25 x 1 1/2" A	100/1000
1188826112	Hypo Needle 26 x 1 1/2" A	100/1000
1188827012	Hypo Needle 27 x 1/2" A	100/1000
1188827112	Hypo Needle 27 x 1 1/2" A	100/1000
1188830340	Hypo Needle 30 x 3/4" A	100/1000


A Bevel- Regular bevel for intramuscular and subcutaneous injections

B Bevel- Short bevel for intravenous or intra-arterial applications and nerve blocks

IN Bevel- Intermediate bevel for intradermal injections/skin testing

*Thin wall needle

Standard Needle Gauge Color Chart



Monoject™ Standard Needles & Syringes LATEX FREE



Rigid Pack Hypodermic Needles with Polypropylene Hub - Sterile

- Rigid pack packaging
- Polypropylene luer lock hub
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Code	Description	Box/Case
8881250040	Hypo Needle 18 x 1" A	100/1000
8881250057	Hypo Needle 18 x 1" B	100/1000
8881250016	Hypo Needle 18 x 1 1/2" A	100/1000
8881250024	Hypo Needle 18 x 1 1/2" B	100/1000
8881250081	Hypo Needle 19 TW* x 1" A	100/1000
8881250099	Hypo Needle 19 TW* x 1" B	100/1000
8881250065	Hypo Needle 19 TW* x 1 1/2" A	100/1000
8881250073	Hypo Needle 19 TW* x 1 1/2" B	100/1000
8881250123	Hypo Needle 20 x 1" A	100/1000
8881250107	Hypo Needle 20 x 1 1/2" A	100/1000
8881250115	Hypo Needle 20 x 1 1/2" B	100/1000
8881250172	Hypo Needle 21 x 1" A	100/1000
8881250149	Hypo Needle 21 x 1 1/2" A	100/1000
8881250222	Hypo Needle 22 x 1" A	100/1000
8881250230	Hypo Needle 22 x 1" B	100/1000
8881250206	Hypo Needle 22 x 1 1/2" A	100/1000
8881250214	Hypo Needle 22 x 1 1/2" B	100/1000
8881250263	Hypo Needle 23 x 1/2" IN	100/1000
8881250289	Hypo Needle 23 x 3/4" A	100/1000
8881250255	Hypo Needle 23 x 1" A	100/1000
8881250305	Hypo Needle 25 x 1" A	100/1000
8881250313	Hypo Needle 25 x 5/8" A	100/1000
8881250545	Hypo Needle 25 x 1 1/2" A	100/1000
8881250321	Hypo Needle 26 x 1/2" IN	100/1000
8881250362	Hypo Needle 27 x 1/2" IN	100/1000
8881250032	Hypo Needle 30 x 3/4" A	100/1000

A Bevel- Regular bevel for intramuscular and subcutaneous injections

B Bevel- Short bevel for intravenous or intra-arterial applications and nerve blocks

IN Bevel- Intermediate bevel for intradermal injections/skin testing

*Thin wall needle

Monoject™ Standard Needles & Syringes LATEX FREE



Rigid Pack Hypodermic Needles with Aluminum Hub - Sterile

- Rigid pack packaging
- Aluminum luer lock hub
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Code	Description	Box/Case
8881200805	Hypo Needle 14 x 1" IN	100/1000
8881200011	Hypo Needle 14 x 1 1/2" IN	100/1000
8881200573	Hypo Needle 14 x 2" IN	100/1000
8881200029	Hypo Needle 15 x 1 1/2" B	100/1000
8881200755	Hypo Needle 16 x 5/8" B	100/1000
8881200060	Hypo Needle 16 x 3/4" IN	100/1000
8881200052	Hypo Needle 16 x 1" A	100/1000
8881200037	Hypo Needle 16 x 1 1/2" A	100/1000
8881200045	Hypo Needle 16 x 1 1/2" B	100/1000
8881200110	Hypo Needle 18 x 1" B	100/1000
8881200714	Hypo Needle 18 x 1" A	100/1000
8881200078	Hypo Needle 18 x 1 1/2" A	100/1000
8881200086	Hypo Needle 18 x 1 1/2" B	100/1000
8881200169	Hypo Needle 19 TW*x 1" A	100/1000
8881200151	Hypo Needle 19 TW*x 1 1/4" A	100/1000
8881200136	Hypo Needle 19 TW*x 1 1/2" A	100/1000
8881200144	Hypo Needle 19 TW* x 1 1/2" B	100/1000
8881200219	Hypo Needle 20 x 1" A	100/1000
8881200227	Hypo Needle 20 x 1" B	100/1000
8881200185	Hypo Needle 20 x 1 1/2" A	100/1000
8881200193	Hypo Needle 20 x 1 1/2" B	100/1000
8881200235	Hypo Needle 20 x 3" A	10/1000
8881200318	Hypo Needle 21 x 2" A	100/1000
8881200342	Hypo Needle 22 x 1" A	100/1000
8881200326	Hypo Needle 22 x 1 1/2" A	100/1000
8881200383	Hypo Needle 23 x 1" A	100/1000
8881200466	Hypo Needle 25 x 5/8" A	100/1000
8881200433	Hypo Needle 25 x 1 1/4" A	100/1000
8881200441	Hypo Needle 25 x 2" A	100/1000
8881200516	Hypo Needle 27 x 1/2" IN	100/1000
8881200508	Hypo Needle 27 x 1 1/4" A	100/1000

A Bevel- Regular bevel for intramuscular and subcutaneous injections

B Bevel- Short bevel for intravenous or intra-arterial applications and nerve blocks

IN Bevel- Intermediate bevel for intradermal injections/skin testing

*Thin wall needle

Monoject™ Standard Needles & Syringes LATEX FREE



Rigid Pack Blunt Cannula - Sterile

- Rigid pack packaging
- Aluminum luer lock hub
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Code	Description	Box/Case
8881202314	Blunt Cannula 15 x 1 1/2"	25/100
8881202322	Blunt Cannula 16 x 1 1/2"	25/100
8881202330	Blunt Cannula 17 x 1 1/2"	25/100
8881202348	Blunt Cannula 18 x 1"	25/100
8881202355	Blunt Cannula 19 x 1 1/2"	25/100
8881202363	Blunt Cannula 20 x 1 1/2"	25/100
8881202371	Blunt Cannula 21 x 1"	25/100
8881202389	Blunt Cannula 22 x 1 1/2"	25/100
8881202397	Blunt Cannula 23 x 1"	25/100



Low Dead Space (LDS) Needles - Sterile

- Ultra-low dead space

Code	Description	Box/Case
1188001112	21 x 1 1/2" LDS Needle	10/1000
1188005058	25 x 5/8 " LDS Needle	10/1000

Monoject™ Standard Insulin Syringes LATEX FREE



SoftPack Insulin Syringes - Sterile

- Soft pack packaging
- Permanently attached needle
- Low dead space
- For U-100 Insulin

Code	Description	Box/Case
1188125058	1mL Insulin Syringe 25 x 5/8" Det*	100/500
8881701406	1mL Insulin Syringe 27 x 1/2"	100/300
1188127012	1mL Insulin Syringe 27 x 1/2" Det*	100/500
8881600004	1/2mL Insulin Syringe 28 x 1/2"	100/300
1188528012	1/2mL Insulin Syringe 28 x 1/2"	100/500
8881601101	1mL Insulin Syringe 28 x 1/2"	100/300
1188128012	1mL Insulin Syringe 28 x 1/2"	100/500
8881600145	3/10mL Insulin Syringe 29 x 1/2"	100/300
8881600350	1/2mL Insulin Syringe 29 x 1/2"	100/300
8881601358	1mL Insulin Syringe 29 x 1/2"	100/300
8881600800	3/10mL Insulin Syringe 30 x 5/16"	100/300
8881600700	1/2mL Insulin Syringe 30 x 5/16"	100/300
8881601600	1mL Insulin Syringe 30 x 5/16"	100/300
1188100555	1mL Insulin Syringe Regular Tip	100/500

* Detachable needle



Rigid Pack Insulin Syringes - Sterile

- Rigid pack packaging
- Permanently attached needle
- Low dead space
- For U-100 Insulin
- May be autoclaved or gas sterilized for use in O.R./procedure paks

Code	Description	Box/Case
8881501822	1mL Insulin Syringe 25 x 5/8" Det*	100/500
8881501970	1mL Insulin Syringe 27 x 1/2" Det*	100/500
8881500014	1/2mL Insulin Syringe 28 x 1/2"	100/500
8881501210	1mL Insulin Syringe 28 x 1/2"	100/500
8881501384	1mL Insulin Syringe Regular Tip	100/500

* Detachable needle

Monoject™ Standard Tuberculin Syringes LATEX FREE



SoftPack Tuberculin Syringes - Sterile

- Soft pack packaging
- Permanently attached needle
- Low dead space

Code	Description	Box/Case
1180125158	1mL TB Syringe 25 x 5/8"	100/500
1180125058	1mL TB Syringe 25 x 5/8" Det*	100/500
1180126038	1mL TB Syringe 26 x 3/8" Det*	100/500
1180127012	1mL TB Syringe 27 x 1/2" Det*	100/500
1180128012	1mL TB Syringe 28 x 1/2"	100/500
1180100555	1mL TB Syringe Regular Tip	100/500

*Detachable needle



Rigid Pack Tuberculin Syringes - Sterile

- Rigid pack packaging
- Permanently attached needle
- Low dead space
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Code	Description	Box/Case
8881501640	1mL TB Syringe 25 x 5/8"	100/500
8881501160	1mL TB Syringe 25 x 5/8" Det*	100/500
8881501178	1mL TB Syringe 26 x 3/8" Det*	100/500
8881501368	1mL TB Syringe 27 x 1/2" Det*	100/500
8881500105	1/2mL TB Syringe 28 x 1/2"	100/500
8881501608	1mL TB Syringe 28 x 1/2"	100/500
8881501400	1mL TB Syringe Regular Luer Tip	100/500

*Detachable needle

Monoject™ Standard Insulin & Tuberculin Syringes LATEX FREE



1mL Luer Lock Syringe - Sterile

- True 1mL size syringe ensures clinicians have the right size syringe to minimize confusion and promote accuracy
- Luer lock tip provides a secure connection that reduces risk of leakage and splatter
- Easy to read graduations with big bold scale markings
- Insulin syringe comes with a pre-attached pharmacy grade tip cap
- Soft pack packaging

Code	Description	Box/Case
1188100777	1mL Insulin Syringe Luer Lock	60/240
1180100777	1mL Tuberculin Syringe Luer Lock	60/240




Allergy Trays - Sterile

- Each tray contains 25 allergy syringes
- Ultra-Comfort bevel for injection
- Graduations: 0.01mL
- Accu-Tip flat plunger tip

Code	Description	Box/Case
8881501962	1mL Tuberculin 27 x 1/2" Det*	25/1000
8881500501	1/2mL Tuberculin Syringe 28 x 1/2"	25/1000
8881501236	1mL Tuberculin Syringe 28 x 1/2"	25/1000

*Detachable needle

Syringe Tip Styles



Regular

Luer Lock

Eccentric

Catheter

Toomey

Monoject™ Standard Needles & Syringes LATEX FREE



SoftPack 3mL Syringes - Sterile

- Soft pack packaging
- Bold graduations: 0.1mL
- Ultra-sharp, tri-beveled, anti-coring needle

Code	Description	Box/Case
1180320100	3mL Syringe 20 x 1"	100/800
1180320112	3mL Syringe 20 x 1 1/2"	100/800
1180321100	3mL Syringe 21 x 1"	100/800
1180321112	3mL Syringe 21 x 1 1/2"	100/800
1180322100	3mL Syringe 22 x 1"	100/800
1180322112	3mL Syringe 22 x 1 1/2"	100/800
1180323100	3mL Syringe 23 x 1"	100/800
1180325058	3mL Syringe 25 x 5/8"	100/800
1180325100	3mL Syringe 25 x 1"	100/800
1180325114	3mL Syringe 25 x 1 1/4"	100/800
1180327114	3mL Syringe 27 x 1 1/4"	100/800
1180300555	3mL Syringe Regular Tip	100/800
1180300777	3mL Syringe Luer Lock Tip	100/800



Rigid Pack 3mL Syringes - Sterile

- Rigid pack packaging
- Bold graduations: 0.1mL
- Ultra-sharp, tri-beveled, anti-coring needle
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Code	Description	Box/Case
8881513025	3mL Syringe 20 x 3/4"	100/1000
8881513033	3mL Syringe 20 x 1"	100/1000
8881513058	3mL Syringe 20 x 1 1/2"	100/1000
8881513132	3mL Syringe 21 x 1"	100/1000
8881513157	3mL Syringe 21 x 1 1/2"	100/1000
8881513231	3mL Syringe 22 x 1"	100/1000
8881513256	3mL Syringe 22 x 1 1/2"	100/1000
8881513330	3mL Syringe 23 x 1"	100/1000
8881513512	3mL Syringe 25 x 5/8"	100/1000
8881513538	3mL Syringe 25 x 1"	100/1000
8881513744	3mL Syringe 27 x 1 1/4"	100/1000
8881513918	3mL Syringe Regular Tip	100/1000
8881513934	3mL Syringe Luer Lock Tip	100/1000

Monoject™ Standard Needles & Syringes LATEX FREE



SoftPack 6mL Syringes - Sterile

- Soft pack packaging
- Bold graduations: 0.2mL
- Ultra-sharp, tri-beveled, anti-coring needle

Code	Description	Box/Case
1181620112	6mL Syringe 20 x 1 1/2"	100/400
1181621100	6mL Syringe 21 x 1"	100/400
1181621112	6mL Syringe 21 x 1 1/2"	100/400
1181622112	6mL Syringe 22 x 1 1/2"	100/400
1180600555	6mL Syringe Regular Tip	100/400
1180600777	6mL Syringe Luer Lock Tip	100/400



Rigid Pack 6mL Syringes - Sterile

- Rigid pack packaging
- Bold graduations: 0.2mL
- Ultra-sharp, tri-beveled, anti-coring needle
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Code	Description	Box/Case
8881516051	6mL Syringe 20 x 1 1/2"	50/500
8881516135	6mL Syringe 21 x 1"	50/500
8881516150	6mL Syringe 21 x 1 1/2"	50/500
8881516911	6mL Syringe Regular Tip	50/500
8881516937	6mL Syringe Lock Tip	50/500

Monoject™ Standard Needles & Syringes LATEX FREE



SoftPack 12mL Syringes - Sterile

- Soft pack packaging
- Bold graduations: 0.2mL
- Ultra-sharp, tri-beveled, anti-coring needle

Code	Description	Box/Case
1181218100	12mL Syringe 18 x 1"	80/480
1181220112	12mL Syringe 20 x 1 1/2"	80/480
1181221100	12mL Syringe 21 x 1"	80/480
1181221112	12mL Syringe 21 x 1 1/2"	80/480
1181200555	12mL Syringe Regular Tip	80/480
1181200777	12mL Syringe Luer Lock Tip	80/480



Rigid Pack 12mL Syringes - Sterile

- Rigid pack packaging
- Bold graduations: 0.2mL
- Ultra-sharp, tri-beveled, anti-coring needle
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Code	Description	Box/Case
8881512597	12mL Syringe 18 x 1"	80/480
8881512662	12mL Syringe 20 x 1 1/2"	80/480
8881512738	12mL Syringe 21 x 1"	80/480
8881512746	12mL Syringe 21 x 1 1/2"	80/480
8881512811	12mL Syringe 22 x 1 1/2"	80/480
8881512852	12mL Syringe Regular Tip	80/480
8881512878	12mL Syringe Luer Lock Tip	80/480
8881512860	12mL Syringe Eccentric Tip	80/480



SoftPack 20mL Syringes - Sterile

- Soft pack packaging
- Bold graduations: 1mL

Code	Description	Box/Case
1182000555	20mL Syringe Regular Tip	40/160
1182000777	20mL Syringe Luer Lock Tip	40/160



Rigid Pack 20mL Syringes - Sterile

- Rigid pack packaging
- Bold graduations: 1mL
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Code	Description	Box/Case
8881520673	20mL Syringe Regular Tip	50/300
8881520657	20mL Syringe Luer Lock Tip	50/300
8881520665	20mL Syringe Eccentric Tip	50/300

Monoject™ Standard Needles & Syringes LATEX FREE



SoftPack 35mL Syringes - Sterile

- Soft pack packaging
- Bold dual graduations: 1mL increments and 1 oz. in 1/4 oz. increments
- Pharmacy grade tip cap included

Code	Description	Box/Case
1183500555	35mL Syringe Regular Tip	40/160
1183500777	35mL Syringe Luer Lock Tip	40/160
1183500888	35mL Syringe Catheter Tip	40/160



Rigid Pack 35mL Syringes - Sterile

- Rigid pack packaging
- Bold dual graduations: 1mL increments and 1/4 oz. increments
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Code	Description	Box/Case
8881535796	35mL Syringe Regular Tip	30/180
8881535762	35mL Syringe Luer Lock Tip	30/180
8881535788	35mL Syringe Eccentric Tip	30/180
8881535770	35mL Syringe Catheter Tip	30/180



SoftPack 60mL Syringes - Sterile

- Soft pack packaging
- Bold dual graduations: 60mL in 5mL increments and 2 oz. in 1/4 oz. increments

Code	Description	Box/Case
1186000555	60mL Syringe Regular Tip	30/120
1186000777	60mL Syringe Luer Lock Tip	30/120
1186000444	60mL Syringe Catheter Tip	30/120



Rigid Pack 60mL Syringes - Sterile

- Rigid pack packaging
- Bold dual graduations: 5mL increments and 1/4 oz. increments
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Code	Description	Box/Case
8881560224	60mL Syringe Regular Tip	20/100
8881560125	60mL Syringe Luer Lock Tip	20/100
8881560182	60mL Syringe Eccentric Tip	20/100
8881560141	60mL Syringe Catheter Tip	20/100
8881560265	60mL Syringe Toomey Tip	20/100

Monoject™ Syringes LATEX FREE



Pharmacy Trays - Sterile

- Bulk packaged syringes for pre-filling and batch filling procedures
- Tamper-evident trays

Code	Description	Box/Case
8881501459	1mL Regular Tip	25/1000
8881513207	3mL Luer Lock Tip	25/200
8881516200	6mL Luer Lock Tip	25/200
8881512258	12mL Luer Lock Tip	25/200
8881520251	20mL Luer Lock Tip	25/200
8881535101	35mL Luer Lock Tip	10/120
8881560232	60mL Luer Lock Tip	10/120



Non-Sterile Syringes

- Bulk packed
- May be autoclaved or gas sterilized
- Standard graduations

Code	Description	Box/Case
8881103025	3mL Syringe Regular Tip	1000
8881103066	3mL Syringe Luer Lock Tip	2000
8881106028	6mL Syringe Regular Tip	500
8881106010	6mL Syringe Luer Lock Tip	500
8881112059	12mL Syringe Regular Tip	500
8881112083	12mL Syringe Luer Lock Tip	800
8881112075	12mL Syringe Eccentric Tip	500
8881120037	20mL Syringe Regular Tip	250
8881120193	20mL Syringe Luer Lock Tip	400
8881120029	20mL Syringe Eccentric Tip	250
8881135084	35mL Syringe Regular Tip	100
8881135019	35mL Syringe Luer Lock Tip	300
8881135068	35mL Syringe Eccentric Tip	100
8881160405	60mL Syringe Regular Tip	100
8881160629	60mL Syringe Luer Lock Tip	155
8881160157	60mL Syringe Catheter Tip	100

Monoject™ Syringes LATEX FREE



Oral Medication Syringes - Non-Sterile

- Bulk packaged
- Permanent, heat-etched dual graduations
- Scale markings in milliliters and teaspoons
- Tip will not accept hypodermic needle
- Ribbed tip caps packaged separately
- Amber color available for light-sensitive drugs

Code	Description	Graduations	Box/Case
8881901014	Clear 1mL	0.01ml/0.5 minims	100/500
8881903002	Clear 3mL	0.1ml/0.25 tsp	100/500
8881906104	Clear 6mL	0.2ml/0.25 tsp	100/500
8881907102	Clear 10mL	0.2ml/0.25 tsp	100/500
8881901006	Amber 1mL	0.01ml/0.5 minims	100/500
8881903010	Amber 3mL	0.1ml/0.25 tsp	100/500
8881906005	Amber 6mL	0.2ml/0.25 tsp	100/500
8881907003	Amber 10mL	0.2ml/0.25 tsp	100/500



140mL Piston Syringes

- Bold dual graduations: 5mL increments and 1/4 oz. increments
- May be autoclaved or gas sterilized

Code	Description	Box/Case
STERILE		
8881114063	140mL Syringe Luer Lock Tip	20
NON-STERILE		
8881114014	140mL Syringe Regular Tip	20
8881114030	140mL Syringe Luer Lock Tip	20
8881114055	140mL Syringe Catheter Tip	20



Control Syringes - Sterile

- Ring controls permanently bonded to syringe
- Attached flexible formulation tip packing

Code	Description	Box/Case
8881512977	12mL Control Syringe Luer Lock Tip (0.2mL increments)	40/160
8881520178	20mL Control Syringe Luer Lock Tip (1mL increments)	40/160

Monoject™ I.V. Flush Syringes LATEX FREE



I.V. Flush Syringes

- Saline and heparin flush syringes
- Terminally sterile
- Pharmaceutical grade manufacturing
- Drug stopper quality plunger tip
- Virtual elimination of syringe reflux
- Low breakforce
- Preservative free, DEHP free



Code	HRI/NDC Number	Description	Box/Case
Saline Flush Syringes (0.9% Sodium Chloride)			
8881570300	08881-5703-00	3mL Syringe 2.5mL Fill	30/180
8881570123	08881-5701-23	12mL Syringe 3mL Fill	30/180
8881570125	08881-5701-25	12mL Syringe 5mL Fill	30/180
8881570121	08881-5701-21	12mL Syringe 10mL Fill	30/180
8881571300	08881-5713-00	3mL Syringe 2.5mL Fill with Cannula	30/180
8881571123	08881-5711-23	12mL Syringe 3mL Fill with Cannula	30/180
8881571125	08881-5711-25	12mL Syringe 5mL Fill with Cannula	30/180
8881571121	08881-5711-21	12mL Syringe 10mL Fill with Cannula	30/180
8881570128	08881-5701-28	12mL Syringe 10mL Fill for Sterile Delivery (1 syringe per pouch)	20/60
8881570129	08881-5701-29	2mL Syringe 10mL Fill for Sterile Delivery (2 syringes per pouch)	20/60
Heparin Flush Syringes (10 units/mL)			
8881580123	08881-5801-23	12mL Syringe 3mL Fill	30/180
8881580125	08881-5801-25	12mL Syringe 5mL Fill	30/180
8881580121	08881-5801-21	12mL Syringe 10mL Fill	30/180
Heparin Flush Syringes (100 units/mL)			
8881590123	08881-5901-23	12mL Syringe 3mL Fill	30/180
8881590125	08881-5901-25	12mL Syringe 5mL Fill	30/180
8881590121	08881-5901-21	12mL Syringe 10mL Fill	30/180

Monoject SmarTip™ Needleless Vial Access Cannula LATEX FREE



Needleless Vial Access Cannula - Sterile

- Plastic cannula specifically designed to access rubber stopper medication vials
- 16g flow rate
- Distinctive color and appearance
- Anti-coring/anti-particulate design

Code	Description	Box/Case
8881540111	Needleless Vial Access Cannula	100/1000
8881540133	3mL Needleless Vial Access Combination	100/800
8881540166	6mL Needleless Vial Access Combination	100/400
8881540122	12mL Needleless Vial Access Combination	80/480

Lifeshield™ Safety I.V. Access Products LATEX FREE



Lifeshield Blunt Cannula & Syringes - Sterile

- For use with LifeShield system
- Stainless steel cannula
- Polypropylene luer lock hub

Code	Description	Box/Case
8881202017	Blunt Cannula 18 x 1" Cannula	100/1000
8881202025	1mL Syringe, Regular Tip 18 x 1" Cannula	100/500
8881202033	3mL Syringe, Luer Lock Tip 18 x 1" Cannula	100/1000
8881202116	6mL Syringe, Luer Lock Tip 18 x 1" Cannula	50/500
8881202132*	12mL Syringe, Luer Lock Tip 18 x 1" Cannula	80/480

Monoject™ Blunt Fill Cannula LATEX FREE



Monoject Blunt Fill Cannula- Sterile

- Used to draw up medication from a vial
- Compatible with a standard luer lock syringe
- Polypropylene luer lock hub

Code	Description	Box/Case
8881202018	Blunt Fill Cannula	100/1000

Monoject™ Safety I.V. Access Products LATEX FREE



Safety I.V. Access Cannula - Sterile

- Compatible with a variety of needleless I.V. systems
- Distinctive color
- Rigid, autoclavable hardpack
- Product identification clearly printed on rigid pack
- Easy-grip design
- 15 gauge x 1/2" cannula

Code	Description	Box/Case
8881540101	I.V. Access Cannula	100/1000
8881541034	3mL Syringe with Monoject BlunTip™ Cannula	100/1000
8881541067	6mL Syringe with Monoject BlunTip™ Cannula	50/500
8881541125	12mL Syringe with Monoject BlunTip™ Cannula	80/480



Safety I.V. Access Cannula with Vial Access Pin - Sterile

- Designed to provide needleless access to drug vials with rubber closures
- SoftPack packaging
- Includes a vial-piercing pin and I.V. access cannula

Code	Description	Box/Case
8881540120	Monoject BlunTip™ with Vial Access Pin	100/1000
8881541307	3mL Syringe with Monoject BlunTip™ & Vial Access Pin	100/800
8881541604	6mL Syringe with Monoject BlunTip™ & Vial Access Pin	100/400
8881541208	12mL Syringe with Monoject BlunTip™ & Vial Access Pin	80/480

Pharmacy Products



Monoject™ Pharmacy Trays - Sterile

- Bulk packaged syringes for pre-filling procedures
- Tamper-evident trays

Code	Description	Box/Case
8881501459	1mL Regular Tip	25/1000
8881513207	3mL Luer Lock Tip	25/200
8881516200	6mL Luer Lock Tip	25/200
8881512258	12mL Luer Lock Tip	25/200
8881520251	20mL Luer Lock Tip	25/200
8881535101	35mL Luer Lock Tip	10/120
8881560232	60mL Luer Lock Tip	10/120



Argyle™ Filter Connector - Sterile

- Microporous 5µm depth filter
- Individually tested for proper filter seal, efficiency and flow rate
- Sufficient flow for gravity infusions
- Male-female luer connector
- Latex free

Code	Description	Box/Case
8888118059	Filter Connector	100/1000



Monoject™ Medication Transfer Needle - Sterile

- Double-ended stainless steel needle with special anti-coring bevels
- Smooth, rounded aluminum hub
- Short butt end for evacuation
- Latex free

Code	Description	Box/Case
8881204005	Medication Transfer Needle 20 x 1"	500/1000

Pharmacy Products LATEX FREE



Monoject™ Filter Needles - Sterile

- Polypropylene hub
- Microporous 5µm depth filter
- Individually tested for proper filter seal, efficiency and flow rate
- Anti-coring stainless steel needle

Code	Description	Box/Case
8881305117	Filter Needle 18 x 1 1/2"	100/1000
8881305018	Filter Needle 20 x 1 1/2"	100/1000



Monoject™ Filter Aspirator - Sterile

- Polypropylene hub
- Microporous 5µm depth filter
- Individually tested for proper filter seal, efficiency and flow rate
- Specially designed stainless steel blunt cannula

Code	Description	Box/Case
8881305109	Filter Aspirator 18 x 3"	25/100

Infusion Products LATEX FREE



1 -Way Stopcock



3 -Way Stopcock



4 -Way Stopcock



3 -Way Stopcock with Tube



Argyle™ EZ-FLO Stopcocks - Sterile

- Blue tint polycarbonate construction
- Low torque, smooth turning control
- Unique threaded female luer connection

Code	Description	Box/Case
8888173005	1-way with male luer lock adapter port covers	50
8888173104	1-way with port cover	50
8888173203	3-way with injection site and 20" extension tube	50
8888173336*	3-way with injection site and 33" extension tube	50
8888173500	3-way only	50
8888173518	3-way with port covers	50
8888173526	3-way with 20" extension tube and port covers	50
8888173534	3-way with 33" extension tube and port covers	50
8888173542	Double 3-way with male luer lock adapter port covers	50
8888173690	3-way with male luer lock adapter and port covers	50
8888174201	4-way with 20" extension tube and port covers	50
8888174698	4-way with male luer lock adapter and port covers	50

DISPOSABLE STOPCOCKS WITH & WITHOUT EXTENSION SETS, CLEAR POLYSTYRENE

170030	3-way stopcock (no port covers)	50
170031	3-way stopcock with male luer lock adapter	50
170042	4-way stopcock (male luer slip) w/20" extension set	50

OPAQUE NYLON

170060	3-way stopcock (no port covers)	50
170061	3-way stopcock with male luer lock adapter	50

EXTENSION TUBES (SETS) WITH MALE/FEMALE LUER SLIP ADAPTERS, STERILE

155922	20" Length	50
155923	30" Length	50
155924	40" Length	50

*Product contains latex

Argyle™ Intermittent Infusion Plug - Sterile

- Male luer lock design also accommodates luer slip fitting
- Special synthetic diaphragm retained by secure plastic shrink sleeve
- Protective sleeve over connector end

Code	Description	Box/Case
8888115006	Heparin lock for I.V. device	50/250

Monoject™ Accessories LATEX FREE



Syringe Tip Caps - Sterile

- Flexible
- Integral luer plug fits luer tip and luer lock syringes
- Caps positioned upright in tray wells
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Code	Description	Box/Case
8881682101	Tray pack of 10	1000
8881682119	Tray pack of 25	1000
8881682085	Bag pack of 100	1000



Finger Grip Extenders

- Bulk packed
- Slip over barrel of syringe for administration of viscous fluids where extra pressure is necessary
- May be autoclaved or gas sterilized for use in O.R./procedure packs

Code	Description	Box/Case
8881820008	Finger grip extender for 20mL	100/500
8881835006	Finger grip extender for 35mL	100/500
8881860004	Finger grip extender for 60mL	100/500

Regional Anesthesia LATEX FREE



Introducer Needle with Aluminum Hub - Sterile

- Rigid pack packaging
- Stainless steel needle
- Ultra-smooth cannula-to-needle point transition

CODE	NEEDLE GAUGE & LENGTH	DESCRIPTION	BOX/CASE
8881205010	18 x 1 1/4"	For use with 22, 25 or 26 x 3 1/2" Spinal Needles	100/1000

Safety Blood Collection Tube Holders



Magellan™ Safety Blood Collection Device - Sterile

- Ready to use 3-in-1 safety design
- One piece safety blood collection needle and blood tube holder combination that is permanently attached
- Intuitive, one-handed activation
- Fixed safety shield that is not an accessory
- Bevel orientation, so that if safety shield is up, the needle bevel is up
- Single use to help comply with OSHA regulations
- Contains latex

Code	Description	Case
8881227121	21 x 1 1/4" Safety Blood Collection Device	300
8881227122	22 x 1 1/4" Safety Blood Collection Device	300



Defender™ Safety Needle Holder

- Bulk packed
- Simultaneous front and back end needle protection
- Reduced sharps disposal cost

Code	Description	Case
8881225110	1" Safety Needle Holder	500
8881225115	1 1/2" Safety Needle Holder	500

Safety Needle Gauge Color Chart

18 Gauge Pink	19 Gauge Cream	20 Gauge Yellow	21 Gauge Green	22 Gauge Black	23 Gauge Blue	25 Gauge Orange	27 Gauge Gray	28 Gauge Teal
------------------	-------------------	--------------------	-------------------	-------------------	------------------	--------------------	------------------	------------------

All Magellan Safety Products are ISO Color Coded

Safety Winged Sets



Blood Collection Sets with Female Luer - Sterile

- One-handed activation of safety mechanism
- Easy-to-handle textured wings
- Female luer connector compatible with all blood drawing equipment
- Latex free, DEHP free
- Color coded wings to coincide with ISO standards

Code	Needle Gauge & Length	Tubing Length	Case
8881225174	19 x 3/4"	12"	50
8881225182	21 x 3/4"	12"	50
8881225390	23 x 3/4"	3"	50
8881225190	23 x 3/4"	12"	50
8881225208	25 x 3/4"	12"	50



Blood Collection Sets with Multi-Sample Luer Adapter - Sterile

- One-handed activation of safety mechanism
- Easy-to-handle textured wings
- Luer adapter for multi-sample collection
- Latex free, DEHP free
- Color coded wings to coincide with ISO standards

Code	Needle Gauge & Length	Tubing Length	Case
8881225281	19 x 3/4"	12"	50
8881225299	21 x 3/4"	12"	50
8881225707	23 x 3/4"	7"	50
8881225307	23 x 3/4"	12"	50
8881225315	25 x 3/4"	12"	50



Blood Collection Sets with Holder Attached - Sterile

- Unique self-sheathing needle system
- One-handed activation of safety mechanism
- Easy-to-handle textured wings
- Superior needle sharpness
- Blood tube holder attached
- Color coded wings to coincide with ISO standards

Code	Needle Gauge & Length	Tubing Length	Case
8881225221	21 x 3/4"	12"	50
8881225265	23 x 3/4"	12"	50
8881225273	25 x 3/4"	12"	50

Blood Collection and Transfer Sets



BacT/Alert™* Safety Transfer Device - Female - Sterile

- Single use for safe transfer of blood specimens from a syringe into a blood culture bottle
- Designed specifically to fit BacT/Alert blood culture bottles
- Adapter insert for evacuated tubes

Code	Description	Box/Case
8881225241	Female Safety Transfer Device	50



Multi-Sample Transfer Set (Dome) - Female - Sterile

- Needle-shielded holder for transferring blood specimens from syringes into blood culture bottles
- Protective cap

Code	Description	Box/Case
8881225240	Holder and Safety Cap (Female)	50



BacT/Alert™* Safety Collection Device - Male - Sterile

- Single use for safe collection of blood specimens
- Male adapter attaches to a female luer for accurate collection sample
- Adapter insert for evacuated tubes
- Designed specifically to fit BacT/Alert blood culture bottles

Code	Description	Box/Case
8881225245	Male Safety Collection Device	50



Multi-Sample Collection Set (Dome) - Male - Sterile

- Needle-shielded holder for use with blood culture and blood collection tubes
- Designed for use with Covidien Blood Collection Sets
- Protective cap

Code	Description	Box/Case
8881225224	Holder and Safety Cap (Male)	50

Blood Collection and Transfer Sets LATEX FREE



Tube Holder with Male Adapter - Sterile

- Needle-shielded holders for use with blood culture and blood collection tubes
- Designed for use with Covidien Blood Collection Sets
- Protective cap

Code	Description	Box/Case
8881225216	Tube Holder (Male) collection tubes	50



Tube Holder with Female Adapter - Sterile

- Needle-shielded holders for transferring blood specimens from syringes into blood culture bottles
- Protective cap

Code	Description	Box/Case
8881225232	Tube Holder (Female)	50



Luer Lock Collection Set - Sterile

- Eliminates the need for a needle when collecting blood from valve-based IV ports
- Pre-assembled combination of blood collection tube holder and luer lock adapter
- Provides a secure connection when performing direct draw from a catheter line into blood collection tube
- Fast, effective and safe

Code	Description	Box/Case
8881225119	Luer Lock Collection Set	50



Blood Collection Assembly with Monoject BlunTip™ Safety I.V. Access Cannula - Sterile

- Pre-assembled combination of blood collection tube holder with a Monoject BlunTip™ Safety I.V. access cannula
- For use with pre-pierced ports
- Safe and effective for patient care
- Provides safe and accurate blood draw
- Reduces needlestick by eliminating a sharp

Code	Description	Box/Case
8881225220	Blood Collection Assembly with BlunTip Cannula	50

Standard Phlebotomy Products



Monoject™ Supra/Hone Blood Collection Needles - Sterile

- Multiple sample needles
- Unique polymer coating on stainless steel needle
- Single piece flow-through cannula design
- Polypropylene hub
- Anti-deflecting butt-end

Code	Needle Gauge & Length	Color Code	Case
8881216033	21 x 1"	Lavender	1000
8881216066	22 x 1"	Blue	1000



Luer Adapter - Sterile

- Multiple sample

Code	Description	Box/Case
8881225257	Luer Adapter	100/1000



Blood Collection Tube Holders

- Translucent polypropylene

Code	Size	Tube Diameter	Case
8881610102	Regular	13mm	100
1500SA	Regular	13mm	1000

SharpStar™ Safety In-Room™ Sharps Containers

The Covidien SharpSafety Division is committed to providing real solutions for all of your sharps disposal needs. The SharpStar™ family of containers is engineered with controls to prevent overfilling and limit accidental or intentional access to the container contents. These containers are offered in multiple sizes from 5 quarts to 4 gallons to accommodate the increasing volume and size of safety engineered sharps devices used today.




Disposal and Processing Recommendations:

- All non-gasket sharps disposal containers manufactured by Covidien are autoclavable.
- All trace chemotherapy, hazardous and non-hazardous pharmaceutical waste must be incinerated.
- When incinerated properly, containers only give off carbon dioxide and water.
- Material Safety Data Sheets (MSDS) are not needed for Covidien sharps disposal containers. They contain no heavy metals or DEHP and they are Latex Free.

Covidien Sharps Disposal containers meet or exceed the OSHA Bloodborne Pathogens Standard and conform to NIOSH/CDC recommendations. Contact your Covidien representative to begin an evaluation.

SharpStar™ Safety In-Room™ Sharps Containers

The SharpStar™ family of containers is engineered with controls to prevent overfilling and limit accidental or intentional access to the container contents. These containers are offered in multiple sizes from 5 quarts to 4 gallons to accommodate the increasing volume and size of safety engineered sharps devices used today.



SharpStar™ In-Room™ Sharps Containers with Counter Balanced Lids

Minimizes the potential for overfilling by stopping in the “FULL” position when it reaches maximum capacity.

- Unique lid design limits accidental or intentional access to container contents for pediatric and high-risk areas.
- Allows for visual monitoring.
- Economical forced horizontal drop maximizes container volume.
- The SharpStar IV™ container has an expanded disposal opening to safely accept the largest sharp used in a patient’s room.
- *OSHA Regulation 29 CFR 1910.1030(d)(4)(iii)*, “During use, containers for contaminated sharps shall be maintained upright throughout use...”

Code	Description	Dimensions	Ship Case
In-Room™ Sharps Containers with Counter Balanced Lid			
8506SA	5 Quart, Clear	12.5”H x 5.5”D x 10.75”W	20
8507SA	5 Quart, Transparent Red	12.5”H x 5.5”D x 10.75”W	20
8533SA	2 Gallon, Clear	11.75”H x 6”D x 13.75”W	10
8534SA	2 Gallon, Transparent Red	11.75”H x 6”D x 13.75”W	10
8536SA	3 Gallon, Clear	16.5”H x 6”D x 13.75”W	10
8537SA	3 Gallon, Transparent Red	16.5”H x 6”D x 13.75”W	10
SharpStar™ IV Sharps Containers with Counter Balanced Lid			
8540SA	4 Gallon, Clear	18.5”H x 6”D x 16.5”W	10
8541SA	4 Gallon, Transparent Red	18.5”H x 6”D x 16.5”W	10

In-Room™ Sharps Containers with Always-Open Lid



851201

In-Room™ System with Always-Open Lid

- Always-Open lids provide convenient, quick and easy point-of-use sharps disposal.
- Accommodates larger syringes.
- Economical horizontal drop maximizes container volume.
- Containers easily lock for final disposal.
- 8.5" wide opening will accept up to 60mL syringe.
- Wall enclosures, brackets and adhesive-backed holders are available.



85221

Code	Description	Dimensions	Ship Case
851201	5 Quart, Clear	11"H x 4.75"D x 10.75"W	20
851301	5 Quart, Transparent Red	11"H x 4.75"D x 10.75"W	20
85321	2 Gallon, Clear	11.5"H x 6"D x 13.75"W	10
85321R	2 Gallon, Transparent Red	11.5"H x 6"D x 13.75"W	10
85221	3 Gallon, Clear	16.25"H x 6"D x 13.75"W	10
85221R	3 Gallon, Transparent Red	16.25"H x 6"D x 13.75"W	10

In-Room™ Sharps Containers with Mailbox-Style Lid



85021

In-Room™ System with Mailbox-Style Lid

- Mailbox-style lid can be kept closed when container is not in use.
- Economical horizontal drop maximizes container volume.
- Containers lock for final disposal.
- Wall enclosures, brackets and adhesive-backed holders are available.



85121

Code	Description	Dimensions	Ship Case
85021	2 Quart, Clear	6.25"H x 4.75"D x 10.75"W	20
85031	2 Quart, Transparent Red	6.25"H x 4.75"D x 10.75"W	20
85121	5 Quart, Clear	11"H x 4.75"D x 10.75"W	20
85131	5 Quart, Transparent Red	11"H x 4.75"D x 10.75"W	20

In-Room™ Sharps Containers with Tortuous Path Lid



In-Room™ System with Tortuous Path Lid

- One-hand, one-step, no-touch sharps depositing.
- Horizontal loading, non-mechanical sharps opening.
- Tortuous-path design impedes hand entry.
- Various mounting brackets and glove box holders are available.

Code	Description	Dimensions	Ship Case
31144010	5 Quart, Transparent Red	12.25"H x 4.25"D x 11"W	30
31143897	5 Quart, Clear	12.25"H x 4.25"D x 11"W	30

GatorGuard™ Safety In-Room™ Sharps Containers



GatorGuard™ In-Room™ Sharps Containers

- Counter-balanced door design allows for disposal of various sized items, from a lightweight butterfly to a 60mL syringe.
- Easy visibility of recommended fill line.
- Signals when container's capacity is reached with a FULL indicator.
- Horizontal drop system maximizes full use of container.
- Thumb tabs lock the container closed for final disposal when full.
- Complete line of brackets and accessories available.

Code	Description	Dimensions	Ship Case
GatorGuard™ Jr. Sharps Containers			
31353603	5 Quart, Transparent Red	14"H x 6"D x 13"W	14
GatorGuard™ Sharps Containers			
31323333	2 Gallon, Transparent Red	17.5"H x 6"D x 14"W	12
31314886	3 Gallon, Transparent Red	20.5"H x 6"D x 14"W	12

Renewables™ Sharps Disposal Containers



8900MW



8980MW



8507MW



8979MW



8985MW

Renewables™ Sharps Disposal Containers

- Manufactured using recycled plastics.
- Reduces waste into landfills and incinerators.
- Promotes a positive image of the Healthcare industry.
- Saves fossil fuels, water, and other precious resources.

Code	Description	Dimensions	Ship Case
Phlebotomy Sharps Container			
8900MW	1 Quart, Clear Lid	6.25"H x 4.25"D x 4.5"W	100
SharpStar™ In-Room™ Sharps Containers with Counter-Balanced Lid			
8507MW	5 Quart	12.5"H x 5.5"D x 10.75"W	20
8537MW	3 Gallon	16.5"H x 6"D x 13.75"W	10
In-Room™ Sharps Containers with Always-Open Lid			
8522MW	3 Gallon	16.25"H x 6"D x 13.75"W	10
Ancillary / Multi-Purpose Containers with Horizontal-Drop Opening Lid			
8967MW	2 Gallon, Clear Lid	12.75"H x 7.25"D x 10.5"W	20
Ancillary / Multi-Purpose Containers with Rotor Opening Lid			
8979MW	2 Gallon, Clear Rotor Lid	10"H x 7.25"D x 10.5"W	20
Large Volume Sharps Containers with Hinged Lid			
8980MW	8 Gallon	17.75"H x 11"D x 15.5"W	10
ChemoSafety™ Container with Hinged Lid			
8985MW	8 Gallon	17.75"H x 11"D x 15.5"W	10

SharpSafety™ In-Room™ Wall Enclosures, Brackets, Glove Boxes & Holders



85301H

In-Room™ System Wall Enclosures and Glove Boxes

- Sturdy, locked sharps container mounting allows for secure point of use disposal of sharps.
- For convenience, all cabinets are universally keyed.
- Sturdy, durable, long lasting ABS plastic provides a secondary barrier of protection.
- Easily monitored fill level through viewing windows on front and sides.
- Keep gloves readily available with optional glove dispensers.
- Beige, unobtrusive style fits any room decor.



85161H



8555SA

Code	Description	Dimensions	Ship Case
85161H	2 & 5 Quart Wall Enclosure	11.75"H x 5.5"D x 13"W	1
85165H	2 & 5 Quart Wall Enclosure	11.75"H x 5.5"D x 13"W	5
8556H	2 & 5 Quart Wall Enclosure, Glove Dispenser (8550LG) and SharpStar™ Container (8507SA)	24.25"H x 5.5"D x 13"W	1
85301H	2 & 3 Gallon Wall Enclosure	17.5"H x 6.5"D x 14.5"W	1
85401H	3 & 4 Gallon SharpStar IV™ Wall Enclosure	26.5"H x 7.25"D x 17.75"W	1
8550B	Glove Dispenser (8550LG) with attached Wall Mounting Bracket (8554B)	7.25"H x 4.5"D x 12.25"W	5
8550LG	Glove Dispenser, Top Loading	6.75"H x 4.5"D x 12"W	10
8554B	Wall Mounting Bracket for 8550LG & 8555SA	8"H x 4.5"D x 12.5"W	5
8555SA	Stackable Multi-Glove Dispenser, Side Loading	8"H x 4.5"D x 12.5"W	10
8516HDL	Replacement Wall Enclosure Door with Window and Lock		1

Sharps-A-Gator™ Wall Enclosures, Glove Boxes & Holders



31353553

In-Room™ System Wall Enclosures and Glove Boxes

Code	Description	Dimensions	Ship Case
31353553	5 Quart GatorGuard™ JR. Cabinet	9.75"H x 6.75"D x 13.5"W	1
31307005	2 & 3 Gallon GatorGuard™ Cabinet	16"H x 6.5"D x 14.25"W	1
31320032	Glove Box Dispenser	6.5"H x 4.5"D x 12"W	10
31369955	5 Quart GatorGuard™ JR. Table Top Holder	2.25"H x 7.25"D x 14.5"W	1
31139762	Single Glove Box Holder	5"H x 4"D x 3.5"W	20



31369955

SharpSafety™ Large Volume Sharps Containers



Large Volume Containers with Sliding, Hinged, and Sealing Gasket Lid

- Adjustable sliding and hinged lids accommodate large sharps.
- Lids can provide temporary closure when container is not in use to reduce exposure to contents.
- Carrying handles and carts with wheels are available for easy transport and lifting of heavy, full containers.
- Optional gasket lids provide leak-resistant containment of contents.
- Easy, simple assembly.
- Containers lock for final disposal.
- Nestable containers save valuable storage space.
- Sliding gasket lids have a round 6" opening to accommodate up to 1000mL bottles.

Code	Description	Dimensions	Ship Case
Large Volume Containers with Sliding Lid			
8980S	8 Gallon, Red	17.75"H x 11"D x 15.5"W	10
8997S	8 Gallon, Red, with Sealing Gasket Lid	17.75"H x 11"D x 15.5"W	10
8935	12 Gallon, Red	18.75"H x 12.75"D x 18.25"W	10
8936SA	12 Gallon, Red, with Sealing Gasket Lid	18.75"H x 12.75"D x 18.25"W	10
8938	18 Gallon, Red	26"H x 12.75"D x 18.25"W	5
8998S	18 Gallon, Red, with Sealing Gasket Lid	26"H x 12.75"D x 18.25"W	5
Large Volume Containers with Hinged Lid			
8980	8 Gallon, Red	17.75"H x 11"D x 15.5"W	10
8997	8 Gallon, Red, with Sealing Gasket Lid	17.75"H x 11"D x 15.5"W	10
8933	12 Gallon, Red	18.75"H x 12.75"D x 18.25"W	10
8932	12 Gallon, Red, with Sealing Gasket Lid	18.75"H x 12.75"D x 18.25"W	10
8991	18 Gallon, Red	26"H x 12.75"D x 18.25"W	5
8998	18 Gallon, Red, with Sealing Gasket Lid	26"H x 12.75"D x 18.25"W	5
8930SA	30 Gallon, Red Perfusion Waste Container with Sealing Gasket Lid	27.5"H x 15.25"D x 21.25"W	3

SharpSafety™ Large Volume SharpsCart™ Foot Pedal Carts and Floor Brackets



Foot Pedal Carts and Floor Brackets

- One-handed disposal.
- Lids lift or slide open easily with foot-operated pedal.
- Lids may remain closed when not in use to reduce exposure to contents.
- Ergonomic handle is telescopic when transporting and retractable when stationary.
- Heavy containers can be removed from the side with minimal lifting.
- *OSHA Regulation 29 CFR 1910.1030(d)(4)(iii), "During use, containers for contaminated sharps shall be maintained upright throughout use..."*
- Covidien offers a complete line of brackets and carts for large volume containers to ensure your safety and compliance.

Code	Description	Dimensions	Ship Case
8984	Locking Wire Wall Bracket for 8 Gallon Containers	18"H x 9.75"D x 15.5"W	5
8992H	Floor cart with wheels for 8, 12 & 18 Gallon Containers	22"H x 12.75"D x 22.5"W	1
8980FP	SharpsCart™ Foot Pedal Cart for 8 Gallon Hinged Lid Containers	22"H x 21.5"D x 16"W	1
8981FP	SharpsCart™ Foot Pedal Cart for 8 Gallon Sliding Lid Containers	22"H x 21.5"D x 16"W	1
8991FP	SharpsCart™ Foot Pedal Cart for 12 & 18 Gallon Hinged Lid Containers	31"H x 22.75"D x 15.25"W	1
8938FP	SharpsCart™ Foot Pedal Cart for 12 & 18 Gallon Sliding Containers	31"H x 22.75"D x 15.25"W	1
8940FP	Molded SharpsCart™ Foot Pedal Cart for 12 & 18 Gallon Sliding Lid Containers	31"H x 22.75"D x 15.25"W	1
8941FP	Molded SharpsCart™ Foot Pedal Cart for 12 & 18 Gallon Hinged Lid Containers	31"H x 22.75"D x 15.25"W	1
8993H	Floor bracket with wheels for 30 Gallon Container	35.5"H x 22.25"D x 15.25"W	1

Sharps-A-Gator™ Large Volume Containers, Floor Stands, Brackets & Carts



Sharps-A-Gator™ Large Volume Sharps Containers

- Designed for the convenient disposal of large instruments and high volumes.
- The versatile lid on the 7 and 10 gallon containers features both a wide opening for large instruments and an integrated sharps port for smaller items.
- Gasketed containers for spill resistance.
- Clear lid keeps contents and fill level visible for increased safety.
- Convenient molded carrying handles to aid in lifting heavy, filled containers for final disposal.



Code	Description	Dimensions	Ship Case
Sharps-A-Gator™ Large Volume Sharps Containers			
31156550	7 Gallon, Red, Split Lid, Sharps Port	14"H x 12"D x 15.5"W	10
31143533	7.5 Gallon, Red, Round	13.5"H x 14.5"D x 14.5"W	6
31143665	10 Gallon, Red, Split Lid, Sharps Port	21.5"H x 12"D x 15.5"W	6
31378089	19 Gallon, Red	25.5"H x 16.5"D x 16.5"W	5

Sharps-A-Gator™ Large Volume Containers, Floor Stands, Brackets & Carts



Code	Description	Dimensions	Ship Case
31140349	Cart for 19 Gallon Containers	30"H x 15"D x 17.5"W	1
31140109	Cart for 7 & 10 Gallon Containers Sharps-A-Gator™ (4-Castered)	30"H x 20"D x 16"W	1
31217949	Cart for 7 & 10 Gallon Containers	21.75"H x 12"D x 15"W	1



PG II-Rated Large Volume Containers

The PG II-Rated containers allow for transport of sharps, cultures, and stocks. This rating eliminates the need for costly secondary packaging and provides a solution for safe disposal of sharps, chemotherapy waste, cultures, and stocks.



PG II Sharps Containers and PG II Chemotherapy D.O.T. Compliant

- Multiple sizes (8, 12, 18 gallon) for containment of disposable contaminated sharps and trace chemotherapy waste.
- Properly labeled ID# “UN3291”, “PG II” and “Regulated Medical Waste”. Eliminates the need for costly D.O.T. compliant secondary packaging during transport of sharps, cultures, and stocks when transported by a private or contract carrier.
- Absorbent pad helps contain excess fluid.
- Gasketed lid provides leak resistance.
- Label differentiates containers from those not D.O.T. compliant.
- Hinged lids or sliding lids allow for temporary and final opening and closing.
- Molded-in carrying handles allow for easier transport.
- SharpsCart™ foot pedal carts and floor brackets provide extra stability and easier transport during use.
- Sliding gasket lids have a round 6” opening to accommodate up to 1000mL bottles.

Code	Description	Dimensions	Ship Case
PG II Sharps Disposal Containers			
8997PG2	8 Gallon, Red, Hinged Lid	17.75”H x 11”D x 15.5”W	10
8932PG2	12 Gallon, Red, Hinged Lid	18.75”H x 12.75”D x 18.25”W	10
8998PG2	18 Gallon, Red, Hinged Lid	26”H x 12.75”D x 18.25”W	5
8997SPG2	8 Gallon, Red, Sliding Lid	17.75”H x 11”D x 15.5”W	10
8936PG2	12 Gallon, Red, Sliding Lid	18.75”H x 12.75”D x 18.25”W	10
8998SPG2	18 Gallon, Red, Sliding Lid	26”H x 12.75”D x 18.25”W	5
PG II Chemotherapy Containers			
8985PG2	8 Gallon, Yellow, Hinged Lid	17.75”H x 11”D x 15.5”W	10
8931PG2	12 Gallon, Yellow, Hinged Lid	18.75”H x 12.75”D x 18.25”W	10
8989PG2	18 Gallon, Yellow, Hinged Lid	26”H x 12.75”D x 18.25”W	5
8985SPG2	8 Gallon, Yellow, Sliding Lid	17.75”H x 11”D x 15.5”W	10
8934PG2	12 Gallon, Yellow, Sliding Lid	18.75”H x 12.75”D x 18.25”W	10
8939PG2	18 Gallon, Yellow, Sliding Lid	26”H x 12.75”D x 18.25”W	5

Trace Chemotherapy Sharps Containers



ChemoSafety™ Containers with Sliding Lid or Hinged Lid

- For the disposal of trace chemotherapy waste.
- Sealing gasket lids provide leak-resistant containment of fluid.
- Absorbent pad in bottom of container absorbs excess liquid.
- Adjustable sliding lid reduces the possibility of aerosolization and accommodates chemotherapy waste of various types and sizes.
- Available carts with wheels and built-in carrying handles for easy transport and lifting of heavy, full containers.
- Containers provide temporary and final locking closure for disposal.
- Nestable containers save valuable storage space.
- Sliding gasket lids have a round 6" opening to accommodate up to 1000mL bottles.

Code	Description	Dimensions	Ship Case
8985S	8 Gallon, Sliding Lid, Yellow	17.75"H x 11"D x 15.5"W	10
8934	2 Gallon, Sliding Lid, Yellow	18.75"H x 12.75"D x 18.25"W	10
8939	18 Gallon, Sliding Lid, Yellow	26"H x 12.75"D x 18.25"W	5
8982	2 Gallon, Hinged Lid, Yellow	10"H x 7.25"D x 10.5"W	20
8985	8 Gallon, Hinged Lid, Yellow	17.75"H x 11"D x 15.5"W	10
8931	12 Gallon, Hinged Lid, Yellow	18.75"H x 12.75"D x 18.25"W	10
8989	18 Gallon, Hinged Lid, Yellow	26"H x 12.75"D x 18.25"W	5

DialySafety Sharps Containers



DialySafety Dialysis Sharps Disposal Containers

- 3-gallon capacity designed for dialysis use to accommodate fistula needles.
- Rotor lid adjusts to accommodate small and large sharps.
- Temporary closure and final lock prevents accidental contact with container contents when not in use; secures for final disposal.
- Locking and non-locking wall brackets are available.

Code	Description	Dimensions	Ship Case
DialySafety Containers with Rotor & Hinged Opening Lid			
8964	3 Gallon DialySafety Container, Transparent Red	13.75"H x 6"D x 13.75"W	20

Transportable Sharps Containers



Transportable Sharps Containers

- Simple, one-handed operation opens and closes the one quart disposal opening.
- Convenient size fits into carry bags for Nursing, EMS and Fire/Rescue.
- Leak-resistant screw cap offers secure closure between uses when properly closed.
- Large disposal opening easily accommodates small syringes and other sharps waste.
- Sharps Shuttle snap closure for one-handed opening and closing.



Code	Description	Dimensions	Ship Case
Transportable Sharps Disposal Containers			
8301	Sharps Shuttle with Locking Mechanism	1.5"D x 6.5" Inside Length	24
8303SA	1 Quart, Red	8.75"H x 2.5"D x 4.5"W	20
8909	1.5 Quart, Red	10.5"H x 3.5"D x 3.5"W	20

Non-Biohazardous Waste Sharps Disposal Containers



EnviroStar™ Sharps Disposal Containers

- Environmental green in color differentiate from other containers.
- Designed to help segregate waste to reduce waste costs.

Code	Description	Dimensions	Ship Case
8790	2 Gallon	10"H x 7.25"D x 10.5"W	20
8781	8 Gallon	17.75"H x 15.5"D x 11"W	10
8791	18 Gallon	26"H x 12.75"D x 18.25"W	5

Multi-Purpose Containers

Rotor Opening Lids



8970

Hinged Opening Lids



8990SA

Rotor & Hinged Opening Lids



8527R

Multi-Purpose Sharps Containers

- Adjustable rotor opening or hinged lids accommodate a variety of sharps sizes and provide temporary and permanent closure.
- Easy, simple assembly.
- Containers lock for final disposal.
- Nestable containers save valuable storage space.
- Available holders and brackets secure containers wherever needed.
- Plastics and pigments are safe to autoclave or incinerate.
- For use on countertops or mounted on the wall.
- Economical horizontal drop maximizes container volume.
- Large horizontal openings accommodate larger syringes.

Code	Description	Dimensions	Ship Case
Multi-Purpose Containers with Rotor Opening Lid			
8920SA	2 Quart, Red	4.58"H x 6.285"D x 4.75"W	60
8950SA	5 Quart, Round, Red	6.75"H x 8.75"D	40
8970	2 Gallon, Red	10"H x 7.25"D x 10.5"W	20
8979C	2 Gallon, Transparent Red	10"H x 7.25"D x 10.5"W	20
Multi-Purpose Containers with Hinged Lid			
8990SA	2 Gallon, Red	10"H x 7.25"D x 10.5"W	20
Multi-Purpose Containers with Rotor & Hinged Opening Lid			
8527R	3 Gallon, Transparent Red	13.75"H x 6"D x 13.75"W	10

Multi-Purpose Containers



Horizontal-Drop Opening Lid

Code	Description	Dimensions	Ship Case
89651	2 Gallon, Red	12.75"H x 7.25"D x 10.5"W	20
89671	2 Gallon, Transparent Red	12.75"H x 7.25"D x 10.5"W	20



Muti-Purpose Sharps Containers with Sliding Lid

Code	Description	Dimensions	Ship Case
31143699	1 Gallon, Red	6.25"H x 7"D x 10.5"W	32
31142222	2 Gallon, Red	10.25"H x 7"D x 10.5"W	20



Monoject™ Sharps Containers

Code	Description	Dimensions	Ship Case
Monoject™ Vertical Drop Containers			
8881676236	4 Quart Red, Small	6.75"W x 10.56"L x 7.08"H	40
8881676285	8 Quart Red, Medium	6.75"W x 10.56"L x 10.89"H	20
8881676434	14 Quart Red, Large	6.75"W x 10.56"L x 17.77"H	10



Phlebotomy Containers



8901SA



8900SA



1522SA

Phlebotomy Sharps Containers

- Adapts to present blood drawing trays and techniques.
- Available holders and brackets.
- A variety of container sizes and openings allow for disposal of other sharps such as lancets, butterfly tubing and small syringes.
- Containers lock for final disposal.

Code	Description	Dimensions	Ship Case
8900SA	1 Quart, Red	6.25"H x 6.3"D x 6.3"W	100
8901SA	1 Pint, Red	4.25"H x 4.25"D x 4.5"W	100
1522SA	2.2 Quart, Red	7.25"H x 4.75"D x 6.5"W	60

RCRA Hazardous Waste Containers



RCRA Hazardous Waste Containers

- For the disposal of hazardous pharmaceutical waste.
- RCRA hazardous waste label meets EPA satellite accumulation labeling requirements.
- White & black color easily distinguished from Sharps, Trace Chemo and other disposal containers.
- Tamper evident label helps deter tampering during final disposal.
- Most are PGII rated, meets D.O.T. requirements for transport without secondary packaging.
- Leak-Resistant gasket and absorbent pad help contain liquid contents should a primary container break.
- Locking floor carts.
- All containers come with an absorbent pad to help contain excess fluid.



8601RC



8602RC



8605RC

Code	Description	Dimensions	Ship Case
8601RC	1.5 Quart RCRA Hazardous Waste Container (Not PGII rated)	10"H x 3.50"D x 3.5"W	20
8602RC	2 Gallon RCRA Hazardous Waste Container	10"H x 7.25"D x 10.5"W	20
8605RC	5 Quart In-Room™ Hazardous Waste Container (Not PGII rated)	14"H x 6"D x 13"W	14
8607RC	8 Gallon RCRA Hazardous Waste Container, Hinged Lid	17.75"H x 11"D x 15.5"W	10
8608RC	8 Gallon RCRA Hazardous Waste Container, Sliding Lid	17.75"H x 11"D x 15.5"W	10
8611RC	12 Gallon RCRA Hazardous Waste Container, Hinged Lid	18.75"H x 12.75"D x 18.25"W	10
8612RC	12 Gallon RCRA Hazardous Waste Container, Sliding Lid	18.75"H x 12.75"D x 18.25"W	10
8617RC	18 Gallon RCRA Hazardous Waste Container, Hinged Lid	26"H x 12.75"D x 18.25"W	5
8618RC	18 Gallon RCRA Hazardous Waste Container, Sliding Lid	26"H x 12.75"D x 18.25"W	5

Pharmaceutical Waste Disposal Containers



PharmaSafety™ Sharps Disposal Containers

- For the disposal of non-hazardous pharmaceutical waste.
- White and blue color easily distinguished from Sharps, Chemotherapy and other disposal containers.
- Leak-resistant gasket and absorbent pad help contain liquid contents.
- Tamper evident label helps deter tampering during final disposal.
- Sturdy, rigid polypropylene construction provides excellent impact and leak resistance (unlike cardboard or fiber) to contain heavy vials, IV bags, and more.
- Two-gallon, three-gallon and eight-gallon capacity sizes are suitable for multiple locations - nursing units, med carts, anesthesia, pharmacies, etc.
- Secure, hinged lid with temporary and final closure accepts larger, bulkier waste; closes temporarily between uses; locks for final disposal.
- Nestable containers save valuable storage space.
- Locking wall-mounting brackets and floor carts available allow for point-of-use placement and one-handed disposal; increases mobility, stability, and security.
- All containers come with an absorbent pad to help contain excess fluid.

PharmaSafety™ Sharps Containers

Code	Description	Dimensions	Ship Case
8820	2 Gallon Gasket Container, Hinged Lid	10"H x 7.25"D x 10.5"W	20
8836SA	3 Gallon PharmaStar™ Container, Tamper-Resistant, Counter-Balanced Lid	16.5"H x 6"D x 13.75"W	10
88361H	3 Gallon PharmaStar™ Wall Enclosure	17.5"H x 6.5"D x 14.5"W	1
8850	8 Gallon Gasket Container, Hinged Lid	17.75"H x 11"D x 15.5"W	10
8851	8 Gallon Gasket Container, Sliding Lid	17.75"H x 11"D x 15.5"W	10
8860	12 Gallon Gasket Container, Hinged Lid	18.75"H x 12.75"D x 18.25"W	10
8870	18 Gallon Gasket Container, Hinged Lid	26"H x 12.75"D x 18.25"W	5
8871	18 Gallon PGII Rated Gasket Container, Sliding Lid	26"H x 12.75"D x 18.25"W	5

SharpSafety™ Brackets, Holders and Accessories



8518X



8528C



8912SA



8975



8984

Brackets, Holders and Accessories

- Locking wall brackets are unobtrusive and securely hold container in place.
- All brackets are keyed alike for your convenience.
- Adhesive-backed holders can be mounted in trays or countertops to help secure container and standardize sharps disposal locations.
- **OSHA Regulation 29 CFR 1910.1030(d)(4)(iii), “During use, containers for contaminated sharps shall be maintained upright throughout use...”**
- Covidien offers a complete line of brackets and holders to ensure your safety and compliance.

Code	Description	Dimensions	Ship Case
8518X	Locking Bracket for 2 & 5 Quart In-Room™ Containers	8.25”H x 6.5”D x 11.25”W	5
8519C	Non-Locking Bracket for 2 & 5 Quart In-Room™ Containers	8.25”H x 6.5”D x 11.25”W	5
8524C	Non-Locking Bracket for 3 Gallon Multi-Purpose and 2 & 3 Gallon In-Room™ Containers	11.75”H x 6.5”D x 13.75”W	5
8528C	Locking Bracket for 3 Gallon Multi-Purpose and 2 & 3 Gallon In-Room Containers	11.75”H x 6.5”D x 13.75”W	5
8543C	Locking Bracket for SharpStar™ IV 3 & 4 Gallon In-Room™ Containers	11.75”H x 6.5”D x 13.75”W	5
8912SA	Wall Bracket for 1 Quart Phlebotomy Containers	5”H x 4.25”D x 5”W	5
8963	Locking Bracket for 2 Gallon Multi-Purpose, RCRA, Pharmaceutical and Chemosafety™ Containers	8.25”H x 7.5”D x 11.5”W	5
8975	Wall Bracket for 2 Gallon Multi-Purpose, RCRA, Pharmaceutical Chemosafety™ Containers	7.5”H x 6.5”D x 10”W	5
8984	Locking Wire Bracket for 8 Gallon Large Volume Containers	18”H x 9.75”D x 15.5”W	5

SharpSafety™ Brackets, Holders and Accessories



8910



8515



8974

Brackets, Holders and Accessories

Code	Description	Dimensions	Ship Case
8910	Holder for 1 Quart Phlebotomy Container	2.75"H x 3.75"D x 3.75"W	10
1523SA	Holder for 2 Quart Phlebotomy Container	2.25"H x 4.25"D x 5.75"W	5
8539SA	Holder for 2 Gallon In-Room™ Container	3.25"H x 7.75"D x 11"W	5
8974	Holder for 2 Gallon Phlebotomy, Multi-Purpose, RCRA, Pharmaceutical and ChemoSafety™ Container	3.25"H x 7.75"D x 11"W	5
8515	Holder for 2 & 5 Quart In-Room™ Containers	3.25"H x 7.75"D x 14.25"W	5
8529	Holder for 3 Gallon, Multi-Purpose SharpStar, and In-Room™ Container	2.5"H x 6.5"D x 13"W	5
8549SA	Holder for SharpStar™ IV 3 & 4 Gallon In-Room Container	3.25"H x 6.75"D x 15.75"W	5

Sharps-A-Gator™ Brackets, Holders and Accessories



Brackets, Holders and Accessories

Code	Description	Dimensions	Ship Case
Monoject™ Sharps Disposal Container Accessories			
8881676012	Aluminum Bracket		10
8881676624	Locking Wall Cabinet	12"H x 13"D x 7"W	2
8881676632	Glove Box		2



Code	Description	Ship Case
Sharps-A-Gator™ Brackets and Accessories		
31307054	Locking Wire Bracket, 5 Quart., 2 and 3 Gallon	10
31307062	Non-Locking 5 Quart, 2 and 3 Gallon GatorGuard	10
31158960	Hidden Bracket and Key	20
31139747	Extra Keys for Hidden Bracket Set	10
31158549	Cabinet, Clear Front, Almond Back, 5 Quart	10
31140109	4-Castered Cart 7 and 10 Gallon	1
31217949	Cart, 7 and 10 Gallon	1
31140349	Frame Cart ,19 Gallon	1
31139762	Glove Box Holder, Single	20



Hardware and Container Cross Reference Guide

CONTAINER NUMBER	CONTAINER VOLUME	CONTAINER WASTE DISPOSAL	HARDWARE NUMBER	HARDWARE DESCRIPTION	LOCKING BRACKET?	HARDWARE CASE QTY
1522SA	2.2 Quart	Biohazard	1523SA	Table Top Holder		5 each
31142222	2 Gallon	Biohazard	8963	Wire Bracket	√	5 each
31142222	2 Gallon	Biohazard	8974	Table Top Holder		5 each
31142222	2 Gallon	Biohazard	8975	Wire Bracket		5 each
31143533	7.5 Gallon	Biohazard	N/A	N/A		N/A
31143665	10 Gallon	Biohazard	31140109	Floor Cart (4-Caster)		1 each
31143665	10 Gallon	Biohazard	31217949	Floor Cart		1 each
31143699	1 Gallon	Biohazard	8963	Wire Bracket	√	5 each
31144010	5 Quart	Biohazard	31158549	Wall Cabinet	√	10 each
31144010	5 Quart	Biohazard	31158960	Hidden Bracket and Key	√	20 each
31156550	7 Gallon	Biohazard	31140109	Floor Cart (4-Caster)		1 each
31156550	7 Gallon	Biohazard	31217949	Floor Cart		1 each
31314886	3 Gallon	Biohazard	31307005	Wall Cabinet	√	1 each
31314886	3 Gallon	Biohazard	31307054	Wire Bracket	√	10 each
31314886	3 Gallon	Biohazard	31307062	Wire Bracket		10 each
31314886	3 Gallon	Biohazard	31320016	Table Top Holder		10 each
31323333	2 Gallon	Biohazard	31158960	Hidden Bracket and Key	√	20 each
31323333	2 Gallon	Biohazard	31307005	Wall Cabinet	√	1 each
31323333	2 Gallon	Biohazard	31307054	Wire Bracket	√	10 each
31323333	2 Gallon	Biohazard	31307062	Wire Bracket		10 each
31323333	2 Gallon	Biohazard	31320016	Table Top Holder		10 each
31353603	5 Quart	Biohazard	31307054	Wire Bracket	√	10 each
31353603	5 Quart	Biohazard	31307062	Wire Bracket		10 each
31353603	5 Quart	Biohazard	31353553	Wall Cabinet	√	1 each
31353603	5 Quart	Biohazard	N/A	Table Top Holder		1 each
31378089	19 Gallon	Biohazard	31140349	Floor Cart		1 each
8301	Sharp Shuttle	Biohazard	N/A	N/A		N/A
8303SA	1 Quart	Biohazard	N/A	N/A		N/A
85021	2 Quart	Biohazard	8515	Table Top Holder		5 each
85021	2 Quart	Biohazard	85161H	Wall Cabinet	√	1 each
85021	2 Quart	Biohazard	85165H	Wall Cabinet	√	5 each
85021	2 Quart	Biohazard	8518X	Wire Bracket	√	5 each
85021	2 Quart	Biohazard	8519C	Wire Bracket		5 each
85021	2 Quart	Biohazard	8556H	Wall Cabinet	√	1 each
85031	2 Quart	Biohazard	8515	Table Top Holder		5 each
85031	2 Quart	Biohazard	8518X	Wire Bracket	√	5 each
85031	2 Quart	Biohazard	8519C	Wire Bracket		5 each
8506SA	5 Quart	Biohazard	8515	Table Top Holder		5 each
8506SA	5 Quart	Biohazard	85161H	Wall Cabinet	√	1 each
8506SA	5 Quart	Biohazard	85165H	Wall Cabinet	√	5 each
8506SA	5 Quart	Biohazard	8518X	Wire Bracket	√	5 each
8506SA	5 Quart	Biohazard	8519C	Wire Bracket		5 each
8506SA	5 Quart	Biohazard	8556H	Wall Cabinet	√	1 each
8507MW	5 Quart	Biohazard	8515	Table Top Holder		5 each

Hardware and Container Cross Reference Guide

CONTAINER CODE	CONTAINER VOLUME	CONTAINER WASTE DISPOSAL	HARDWARE NUMBER	HARDWARE DESCRIPTION	LOCKING BRACKET?	HARDWARE CASE QTY
8507MW	5 Quart	Biohazard	85161H	Wall Cabinet	√	1 each
8507MW	5 Quart	Biohazard	85165H	Wall Cabinet	√	5 each
8507MW	5 Quart	Biohazard	8518X	Wire Bracket	√	5 each
8507MW	5 Quart	Biohazard	8519C	Wire Bracket		5 each
8507MW	5 Quart	Biohazard	8556H	Wall Cabinet	√	1 each
8507SA	5 Quart	Biohazard	8515	Table Top Holder		5 each
8507SA	5 Quart	Biohazard	85161H	Wall Cabinet	√	1 each
8507SA	5 Quart	Biohazard	85165H	Wall Cabinet	√	5 each
8507SA	5 Quart	Biohazard	8518X	Wire Bracket	√	5 each
8507SA	5 Quart	Biohazard	8519C	Wire Bracket		5 each
8507SA	5 Quart	Biohazard	8556H	Wall Cabinet	√	1 each
851201	5 Quart	Biohazard	8515	Table Top Holder		5 each
851201	5 Quart	Biohazard	85161H	Wall Cabinet	√	1 each
851201	5 Quart	Biohazard	85165H	Wall Cabinet	√	5 each
851201	5 Quart	Biohazard	8518X	Wire Bracket	√	5 each
851201	5 Quart	Biohazard	8519C	Wire Bracket		5 each
851201	5 Quart	Biohazard	8556H	Wall Cabinet	√	1 each
85121	5 Quart	Biohazard	8515	Table Top Holder		5 each
85121	5 Quart	Biohazard	85161H	Wall Cabinet	√	1 each
85121	5 Quart	Biohazard	85165H	Wall Cabinet	√	5 each
85121	5 Quart	Biohazard	8518X	Wire Bracket	√	5 each
85121	5 Quart	Biohazard	8519C	Wire Bracket		5 each
85121	5 Quart	Biohazard	8556H	Wall Cabinet	√	1 each
851301	5 Quart	Biohazard	8515	Table Top Holder		5 each
851301	5 Quart	Biohazard	85161H	Wall Cabinet	√	1 each
851301	5 Quart	Biohazard	85165H	Wall Cabinet	√	5 each
851301	5 Quart	Biohazard	8518X	Wire Bracket	√	5 each
851301	5 Quart	Biohazard	8519C	Wire Bracket		5 each
851301	5 Quart	Biohazard	8556H	Wall Cabinet	√	1 each
85131	5 Quart	Biohazard	8515	Table Top Holder		5 each
85131	5 Quart	Biohazard	85161H	Wall Cabinet	√	1 each
85131	5 Quart	Biohazard	85165H	Wall Cabinet	√	5 each
85131	5 Quart	Biohazard	8518X	Wire Bracket	√	5 each
85131	5 Quart	Biohazard	8519C	Wire Bracket		5 each
85131	5 Quart	Biohazard	8556H	Wall Cabinet	√	1 each
85221	3 Gallon	Biohazard		Wire Bracket		5 each
85221	3 Gallon	Biohazard	8528C	Wire Bracket	√	5 each
85221	3 Gallon	Biohazard	8529	Table Top Holder		5 each
85221	3 Gallon	Biohazard	85301H	Wall Cabinet	√	1 each
85221R	3 Gallon	Biohazard		Wire Bracket		5 each
85221R	3 Gallon	Biohazard	8528C	Wire Bracket	√	5 each
85221R	3 Gallon	Biohazard	8529	Table Top Holder		5 each
85221R	3 Gallon	Biohazard	85301H	Wall Cabinet	√	1 each

Hardware and Container Cross Reference Guide

CONTAINER CODE	CONTAINER VOLUME	CONTAINER WASTE DISPOSAL	HARDWARE NUMBER	HARDWARE DESCRIPTION	LOCKING BRACKET?	HARDWARE CASE QTY
8522MW	3 Gallon	Biohazard	8528C	Wire Bracket	√	5 each
8522MW	3 Gallon	Biohazard	8529	Table Top Holder		5 each
8522MW	3 Gallon	Biohazard	85301H	Wall Cabinet	√	1 each
8527R	3 Gallon	Biohazard		Wire Bracket		5 each
8527R	3 Gallon	Biohazard	8528C	Wire Bracket	√	5 each
8527R	3 Gallon	Biohazard	8529	Table Top Holder		5 each
85321	2 Gallon	Biohazard		Wire Bracket		5 each
85321	2 Gallon	Biohazard	8528C	Wire Bracket	√	5 each
85321	2 Gallon	Biohazard	85301H	Wall Cabinet	√	1 each
85321	2 Gallon	Biohazard	8539SA	Table Top Holder		5 each
85321R	2 Gallon	Biohazard		Wire Bracket		5 each
85321R	2 Gallon	Biohazard	8528C	Wire Bracket	√	5 each
85321R	2 Gallon	Biohazard	85301H	Wall Cabinet	√	1 each
85321R	2 Gallon	Biohazard	8539SA	Table Top Holder		5 each
8533SA	2 Gallon	Biohazard		Wire Bracket		5 each
8533SA	2 Gallon	Biohazard	8528C	Wire Bracket	√	5 each
8533SA	2 Gallon	Biohazard	85301H	Wall Cabinet	√	1 each
8533SA	2 Gallon	Biohazard	8539SA	Table Top Holder		5 each
8534SA	2 Gallon	Biohazard		Wire Bracket		5 each
8534SA	2 Gallon	Biohazard	8528C	Wire Bracket	√	5 each
8534SA	2 Gallon	Biohazard	85301H	Wall Cabinet	√	1 each
8534SA	2 Gallon	Biohazard	8539SA	Table Top Holder		5 each
8536SA	3 Gallon	Biohazard		Wire Bracket		5 each
8536SA	3 Gallon	Biohazard	8528C	Wire Bracket	√	5 each
8536SA	3 Gallon	Biohazard	8529	Table Top Holder		5 each
8536SA	3 Gallon	Biohazard	85301H	Wall Cabinet	√	1 each
8537MW	3 Gallon	Biohazard		Wire Bracket		5 each
8537MW	3 Gallon	Biohazard	8528C	Wire Bracket	√	5 each
8537MW	3 Gallon	Biohazard	8529	Table Top Holder		5 each
8537MW	3 Gallon	Biohazard	85301H	Wall Cabinet	√	1 each
8537SA	3 Gallon	Biohazard		Wire Bracket		5 each
8537SA	3 Gallon	Biohazard	8528C	Wire Bracket	√	5 each
8537SA	3 Gallon	Biohazard	8529	Table Top Holder		5 each
8537SA	3 Gallon	Biohazard	85301H	Wall Cabinet	√	1 each
8540SA	4 Gallon	Biohazard	8543C	Wire Bracket	√	5 each
8540SA	4 Gallon	Biohazard	85401H	Wall Cabinet	√	1 each
8540SA	4 Gallon	Biohazard	8549SA	Table Top Holder		5 each
8541SA	4 Gallon	Biohazard	8543C	Wire Bracket	√	5 each
8541SA	4 Gallon	Biohazard	85401H	Wall Cabinet	√	1 each
8541SA	4 Gallon	Biohazard	8549SA	Table Top Holder		5 each
8601RC	1.5 Quart	Hazardous	N/A	N/A		N/A
8602RC	2 Gallon	Hazardous	8963	Wire Bracket	√	5 each
8602RC	2 Gallon	Hazardous	8974	Table Top Holder		5 each

Hardware and Container Cross Reference Guide

CONTAINER CODE	CONTAINER VOLUME	CONTAINER WASTE DISPOSAL	HARDWARE NUMBER	HARDWARE DESCRIPTION	LOCKING BRACKET?	HARDWARE CASE QTY
8602RC	2 Gallon	Hazardous	8975	Wire Bracket		5 each
8605RC	5 Quart	Hazardous	31307054	Wire Bracket	√	10 each
8605RC	5 Quart	Hazardous	31307062	Wire Bracket		10 each
8605RC	5 Quart	Hazardous	31353553	Wall Cabinet	√	1 each
8605RC	5 Quart	Hazardous	31369955	Table Top Holder		1 each
8607RC	8 Gallon	Hazardous	8980FP	Foot Pedal Cart		1 each
8607RC	8 Gallon	Hazardous	8984	Wire Bracket	√	5 each
8607RC	8 Gallon	Hazardous	8992H	Floor Bracket w/Wheels		1 each
8608RC	8 Gallon	Hazardous	8981FP	Foot Pedal Cart		1 each
8608RC	8 Gallon	Hazardous	8984	Wire Bracket	√	5 each
8608RC	8 Gallon	Hazardous	8992H	Floor Bracket w/Wheels		1 each
8611RC	12 Gallon	Hazardous	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8611RC	12 Gallon	Hazardous	8991FP	Foot Pedal Cart		1 each
8611RC	12 Gallon	Hazardous	8992H	Floor Bracket w/Wheels		1 each
8612RC	12 Gallon	Hazardous	8938FP	Foot Pedal Cart		1 each
8612RC	12 Gallon	Hazardous	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8612RC	12 Gallon	Hazardous	8992H	Floor Bracket w/Wheels		1 each
8617RC	18 Gallon	Hazardous	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8617RC	18 Gallon	Hazardous	8991FP	Foot Pedal Cart		1 each
8617RC	18 Gallon	Hazardous	8992H	Floor Bracket w/Wheels		1 each
8618RC	18 Gallon	Hazardous	8938FP	Foot Pedal Cart		1 each
8618RC	18 Gallon	Hazardous	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8618RC	18 Gallon	Hazardous	8992H	Floor Bracket w/Wheels		1 each
8781	8 Gallon	Non-Biohazard	8980FP	Foot Pedal Cart		1 each
8781	8 Gallon	Non-Biohazard	8984	Wire Bracket	√	1 each
8781	8 Gallon	Non-Biohazard	8992H	Floor Bracket w/Wheels		1 each
8790	2 Gallon	Non-Biohazard	8963	Wire Bracket	√	1 each
8790	2 Gallon	Non-Biohazard	8974	Table Top Holder		1 each
8790	2 Gallon	Non-Biohazard	8975	Wire Bracket		1 each
8791	18 Gallon	Non-Biohazard	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8791	18 Gallon	Non-Biohazard	8991FP	Foot Pedal Cart		1 each
8791	18 Gallon	Non-Biohazard	8992H	Floor Bracket w/Wheels		1 each
8800		Biohazard	N/A	N/A		N/A
8820	2 Gallon	Pharmaceutical	8963	Wire Bracket	√	5 each
8820	2 Gallon	Pharmaceutical	8974	Table Top Holder		5 each
8820	2 Gallon	Pharmaceutical	8975	Wire Bracket		5 each
8836SA	3 Gallon	Pharmaceutical		Wire Bracket		5 each
8836SA	3 Gallon	Pharmaceutical	8528C	Wire Bracket	√	5 each
8836SA	3 Gallon	Pharmaceutical	8529	Table Top Holder		5 each
8836SA	3 Gallon	Pharmaceutical	88361H	Wall Cabinet	√	1 each
8850	8 Gallon	Pharmaceutical	8980FP	Foot Pedal Cart		1 each
8850	8 Gallon	Pharmaceutical	8984	Wire Bracket	√	5 each
8850	8 Gallon	Pharmaceutical	8992H	Floor Bracket w/Wheels		1 each

Hardware and Container Cross Reference Guide

CONTAINER CODE	CONTAINER VOLUME	CONTAINER WASTE DISPOSAL	HARDWARE NUMBER	HARDWARE DESCRIPTION	LOCKING BRACKET?	HARDWARE CASE QTY
8851	8 Gallon	Pharmaceutical	8981FP	Foot Pedal Cart		1 each
8851	8 Gallon	Pharmaceutical	8984	Wire Bracket	√	5 each
8851	8 Gallon	Pharmaceutical	8992H	Floor Bracket w/Wheels		1 each
8860	12 Gallon	Pharmaceutical	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8860	12 Gallon	Pharmaceutical	8991FP	Foot Pedal Cart		1 each
8860	12 Gallon	Pharmaceutical	8992H	Floor Bracket w/Wheels		1 each
8870	18 Gallon	Pharmaceutical	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8870	18 Gallon	Pharmaceutical	8991FP	Foot Pedal Cart		1 each
8870	18 Gallon	Pharmaceutical	8992H	Floor Bracket w/Wheels		1 each
8871	18 Gallon	Pharmaceutical	8938FP	Foot Pedal Cart		1 each
8871	18 Gallon	Pharmaceutical	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8871	18 Gallon	Pharmaceutical	8992H	Floor Bracket w/Wheels		1 each
8881676236	4 Quart	Biohazard	8881676012	Aluminum Bracket	√	1 each
8881676236	4 Quart	Biohazard	8881676624	Wall Cabinet	√	2 each
8881676285	8 Quart	Biohazard	8881676012	Aluminum Bracket	√	10 each
8881676285	8 Quart	Biohazard	8881676624	Wall Cabinet	√	2 each
8881676434	14 Quart	Biohazard	8881676012	Aluminum Bracket	√	10 each
8900MW	1 Quart	Biohazard	8910	Table Top Holder		10 each
8900MW	1 Quart	Biohazard	8912SA	Wire Bracket		5 each
8900SA	1 Quart	Biohazard	8910	Table Top Holder		10 each
8900SA	1 Quart	Biohazard	8912SA	Wire Bracket		5 each
8901SA	1 Pint	Biohazard	N/A	N/A		N/A
8909	1.5 Quart	Biohazard	N/A	N/A		N/A
8920SA	2 Quart	Biohazard	8926	Locking Wire Bracket	√	5 each
8930SA	30 Gallon	Biohazard	8993H	Floor Cart		1 each
8931	12 Gallon	Chemotherapy	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8931	12 Gallon	Chemotherapy	8991FP	Foot Pedal Cart		1 each
8931	12 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each
8931PG2	12 Gallon	Chemotherapy	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8931PG2	12 Gallon	Chemotherapy	8991FP	Foot Pedal Cart		1 each
8931PG2	12 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each
8932	12 Gallon	Biohazard	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8932	12 Gallon	Biohazard	8991FP	Foot Pedal Cart		1 each
8932	12 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8932PG2	12 Gallon	Biohazard	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8932PG2	12 Gallon	Biohazard	8991FP	Foot Pedal Cart		1 each
932PG2	12 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8933	12 Gallon	Biohazard	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
88933	12 Gallon	Biohazard	8991FP	Foot Pedal Cart		1 each
8933	12 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8934	12 Gallon	Chemotherapy	8938FP	Foot Pedal Cart		1 each
8934	12 Gallon	Chemotherapy	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8934	12 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each

Hardware and Container Cross Reference Guide

CONTAINER CODE	CONTAINER VOLUME	CONTAINER WASTE DISPOSAL	HARDWARE NUMBER	HARDWARE DESCRIPTION	LOCKING BRACKET?	HARDWARE CASE QTY
8934PG2	12 Gallon	Chemotherapy	8938FP	Foot Pedal Cart		1 each
8934PG2	12 Gallon	Chemotherapy	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8934PG2	12 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each
8935	12 Gallon	Biohazard	8938FP	Foot Pedal Cart		1 each
8935	12 Gallon	Biohazard	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8935	12 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8936PG2	12 Gallon	Biohazard	8938FP	Foot Pedal Cart		1 each
8936PG2	12 Gallon	Biohazard	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8936PG2	12 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8936SA	12 Gallon	Biohazard	8938FP	Foot Pedal Cart		1 each
8936SA	12 Gallon	Biohazard	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8936SA	12 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8938	18 Gallon	Biohazard	8938FP	Foot Pedal Cart		1 each
8938	18 Gallon	Biohazard	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8938	18 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8939	18 Gallon	Chemotherapy	8938FP	Foot Pedal Cart		1 each
8939	18 Gallon	Chemotherapy	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8939	18 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each
8939PG2	18 Gallon	Chemotherapy	8938FP	Foot Pedal Cart		1 each
8939PG2	18 Gallon	Chemotherapy	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8939PG2	18 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each
8950SA	5 Quart	Biohazard	N/A	N/A		N/A
8964	3 Gallon	Biohazard		Wire Bracket		5 each
8964	3 Gallon	Biohazard	8528C	Wire Bracket	√	5 each
8964	3 Gallon	Biohazard	8529	Table Top Holder		5 each
89651	2 Gallon	Biohazard	8963	Wire Bracket	√	5 each
89651	2 Gallon	Biohazard	8974	Table Top Holder		5 each
89651	2 Gallon	Biohazard	8975	Wire Bracket		5 each
89671	2 Gallon	Biohazard	8963	Wire Bracket	√	5 each
89671	2 Gallon	Biohazard	8974	Table Top Holder		5 each
89671	2 Gallon	Biohazard	8975	Wire Bracket		5 each
8967MW	2 Gallon	Biohazard	8963	Wire Bracket	√	5 each
8967MW	2 Gallon	Biohazard	8974	Table Top Holder		5 each
8967MW	2 Gallon	Biohazard	8975	Wire Bracket		5 each
8970	2 Gallon	Biohazard	8963	Wire Bracket	√	5 each
8970	2 Gallon	Biohazard	8974	Table Top Holder		5 each
8970	2 Gallon	Biohazard	8975	Wire Bracket		5 each
8979C	2 Gallon	Biohazard	8963	Wire Bracket	√	5 each
8979C	2 Gallon	Biohazard	8974	Table Top Holder		5 each
8979C	2 Gallon	Biohazard	8975	Wire Bracket		5 each
8979MW	2 Gallon	Biohazard	8963	Wire Bracket	√	5 each
8979MW	2 Gallon	Biohazard	8974	Table Top Holder		5 each
8979MW	2 Gallon	Biohazard	8975	Wire Bracket		5 each

Hardware and Container Cross Reference Guide

CONTAINER CODE	CONTAINER VOLUME	CONTAINER WASTE DISPOSAL	HARDWARE NUMBER	HARDWARE DESCRIPTION	LOCKING BRACKET?	HARDWARE CASE QTY
8980	8 Gallon	Biohazard	8980FP	Foot Pedal Cart		1 each
8980	8 Gallon	Biohazard	8984	Wire Bracket	√	5 each
8980	8 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8980MW	8 Gallon	Biohazard	8980FP	Foot Pedal Cart		1 each
8980MW	8 Gallon	Biohazard	8984	Wire Bracket	√	5 each
8980MW	8 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8980S	8 Gallon	Biohazard	8981FP	Foot Pedal Cart		1 each
8980S	8 Gallon	Biohazard	8984	Wire Bracket	√	5 each
8980S	8 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8982	2 Gallon	Chemotherapy	8963	Wire Bracket	√	5 each
8982	2 Gallon	Chemotherapy	8974	Table Top Holder		5 each
8982	2 Gallon	Chemotherapy	8975	Wire Bracket		5 each
8985	8 Gallon	Chemotherapy	8980FP	Foot Pedal Cart		1 each
8985	8 Gallon	Chemotherapy	8984	Wire Bracket	√	5 each
8985	8 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each
8985MW	8 Gallon	Chemotherapy	8980FP	Foot Pedal Cart		1 each
8985MW	8 Gallon	Chemotherapy	8984	Wire Bracket	√	5 each
8985MW	8 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each
8985PG2	8 Gallon	Chemotherapy	8980FP	Foot Pedal Cart		1 each
8985PG2	8 Gallon	Chemotherapy	8984	Wire Bracket	√	5 each
8985PG2	8 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each
8985S	8 Gallon	Chemotherapy	8981FP	Foot Pedal Cart		1 each
8985S	8 Gallon	Chemotherapy	8984	Wire Bracket	√	5 each
8985S	8 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each
8985SPG2	8 Gallon	Chemotherapy	8981FP	Foot Pedal Cart		1 each
8985SPG2	8 Gallon	Chemotherapy	8984	Wire Bracket	√	5 each
8985SPG2	8 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each
8989	18 Gallon	Chemotherapy	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8989	18 Gallon	Chemotherapy	8991FP	Foot Pedal Cart		1 each
8989	18 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each
8989PG2	18 Gallon	Chemotherapy	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8989PG2	18 Gallon	Chemotherapy	8991FP	Foot Pedal Cart		1 each
8989PG2	18 Gallon	Chemotherapy	8992H	Floor Bracket w/Wheels		1 each
8990SA	2 Gallon	Biohazard	8963	Wire Bracket	√	5 each
8990SA	2 Gallon	Biohazard	8974	Table Top Holder		5 each
8990SA	2 Gallon	Biohazard	8975	Wire Bracket		5 each
8991	18 Gallon	Biohazard	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8991	18 Gallon	Biohazard	8991FP	Foot Pedal Cart		1 each
8991	18 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8997	8 Gallon	Biohazard	8980FP	Foot Pedal Cart		1 each
8997	8 Gallon	Biohazard	8984	Wire Bracket	√	5 each
8997	8 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8997PG2	8 Gallon	Biohazard	8980FP	Foot Pedal Cart		1 each

Hardware and Container Cross Reference Guide

CONTAINER CODE	CONTAINER VOLUME	CONTAINER WASTE DISPOSAL	HARDWARE NUMBER	HARDWARE DESCRIPTION	LOCKING BRACKET?	HARDWARE CASE QTY
8997PG2	8 Gallon	Biohazard	8984	Wire Bracket	√	5 each
8997PG2	8 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8997S	8 Gallon	Biohazard	8981FP	Foot Pedal Cart		1 each
8997S	8 Gallon	Biohazard	8984	Wire Bracket	√	5 each
8997S	8 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8997SPG2	8 Gallon	Biohazard	8981FP	Foot Pedal Cart		1 each
8997SPG2	8 Gallon	Biohazard	8984	Wire Bracket	√	5 each
8997SPG2	8 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8998	18 Gallon	Biohazard	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8998	18 Gallon	Biohazard	8991FP	Foot Pedal Cart		1 each
8998	18 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8998PG2	18 Gallon	Biohazard	8941FP	Molded Plastic Foot Pedal Cart	√	1 each
8998PG2	18 Gallon	Biohazard	8991FP	Foot Pedal Cart		1 each
8998PG2	18 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8998S	18 Gallon	Biohazard	8938FP	Foot Pedal Cart		1 each
8998S	18 Gallon	Biohazard	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8998S	18 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
8998SPG2	18 Gallon	Biohazard	8938FP	Foot Pedal Cart		1 each
8998SPG2	18 Gallon	Biohazard	8940FP	Molded Plastic Foot Pedal Cart	√	1 each
8998SPG2	18 Gallon	Biohazard	8992H	Floor Bracket w/Wheels		1 each
			31139762	S-A-G Single Glove Box		20 each
				Wall Bracket		
			31320032	Glove Box for GatorGuard		10 each
				Wall Enclosures		
			8881676632	Glove Box for 8881676624		2 each
			8555SA	Side Load, Stackable Multi-Glove Dispenser		10 each
			8550B	Top Load Glove Dispenser)		5 each
				(8550LG) w/Wall Mounting Bracket (8554B)		
			8550LG	Top Load Glove Dispenser for SharpStar Wall Enclosures		10 each
			8516HDL	Door for 2 & 5 Qt SharpSafety Wall Enclosure		1 each
			8546	Glove Box Prevention Adapter for 2 & 5 Qt SharpStar Wall Enclosures		5 each
			8554B	Wall Mounting Bracket for Glove Dispensers 8550LG & 8555SA		5 each

INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
1180100555	SoftPack 1mL TB Syringe Regular Tip	15	1188127012	SoftPack 1mL Insulin Syringe 27 x 1/2" Det	14
1180100777	1mL Tuberculin Syringe Luer Lock	16	1188128012	SoftPack 1mL Insulin Syringe 28 x 1/2"	14
1180125058	SoftPack 1mL TB Syringe 25 x 5/8"	15	1188528012	SoftPack 1/2mL Insulin Syringe 28 x 1/2"	14
1180125158	SoftPack 1mL TB Syringe 25 x 5/8"	15	1188818100	SoftPack Hypodermic Needle 18 x 1"	10
1180126038	SoftPack 1mL TB Syringe 26 x 3/8"	15	1188818112	SoftPack Hypodermic Needle 18 x 1 1/2"	10
1180127012	SoftPack 1mL TB Syringe 27 x 1/2"	15	1188818113	SoftPack Hypodermic Needle 18 x 1 1/2" B	10
1180128012	SoftPack 1mL TB Syringe 28 x 1/2"	15	1188819100	SoftPack Hypodermic Needle 19TW x 1"	10
1180300555	SoftPack 3mL Syringe Regular Tip	17	1188819112	SoftPack Hypodermic Needle 19TW x 1 1/2"	10
1180300777	SoftPack 3mL Syringe Luer Lock Tip	17	1188820100	SoftPack Hypodermic Needle 20 x 1"	10
1180320100	SoftPack 3mL Syringe 20 x 1"	17	1188820112	SoftPack Hypodermic Needle 20 x 1 1/2"	10
1180320112	SoftPack 3mL Syringe 20 x 1 1/2"	17	1188821100	SoftPack Hypodermic Needle 21 x 1"	10
1180321100	SoftPack 3mL Syringe 21 x 1"	17	1188821112	SoftPack Hypodermic Needle 21 x 1 1/2"	10
1180321112	SoftPack 3mL Syringe 21 x 1 1/2"	17	1188822100	SoftPack Hypodermic Needle 22 x 1"	10
1180322100	SoftPack 3mL Syringe 22 x 1"	17	1188822112	SoftPack Hypodermic Needle 22 x 1 1/2"	10
1180322112	SoftPack 3mL Syringe 22 x 1 1/2"	17	1188823012	SoftPack Hypodermic Needle 23 x 1/2" IN	10
1180323100	SoftPack 3mL Syringe 23 x 1"	17	1188823100	SoftPack Hypodermic Needle 23 x 1"	10
1180325058	SoftPack 3mL Syringe 25 x 5/8"	17	1188825058	SoftPack Hypodermic Needle 25 x 5/8"	10
1180325100	SoftPack 3mL Syringe 25 x 1"	17	1188825100	SoftPack Hypodermic Needle 25 x 1"	10
1180325114	SoftPack 3mL Syringe 25 x 1 1/4"	17	1188825112	SoftPack Hypodermic Needle 25 x 1 1/2"	10
1180327114	SoftPack 3mL Syringe 27 x 1 1/4"	17	1188826112	SoftPack Hypodermic Needle 26 x 1 1/2"	10
11806000555	SoftPack 6mL Syringe Regular Tip	18	1188827012	SoftPack Hypodermic Needle 27 x 1/2"	10
1180600777	SoftPack 6mL Syringe Luer Lock Tip	18	1188827112	SoftPack Hypodermic Needle 27 x 1 1/2"	10
1181200555	SoftPack 12mL Syringe Regular Tip	19	1188830340	SoftPack Hypodermic Needle 30 x 3/4"	10
1181200777	SoftPack 12mL Syringe Luer Lock Tip	19	1500SA	Blood Collection Needle and Tube Holders Regular 13mm	34
1181218100	SoftPack 12mL Syringe 18 x 1"	19	1522SA	Phlebotomy Sharps Container, 2.2 Quart, Red	49
1181220112	SoftPack 12mL Syringe 20 x 1 1/2"	19	1523SA	Holder for 2 Quart Phlebotomy Container	53
1181221100	SoftPack 12mL Syringe 21 x 1"	19	155922	Extension Tubes (Sets) Luer Slip Adapters, 20" Length	28
1181221112	SoftPack 12mL Syringe 21 x 1 1/2"	19	155923	Extension Tubes (Sets) Luer Slip Adapters, 30" Length	28
1181620112	SoftPack 6mL Syringe 20 x 1 1/2"	18	155924	Extension Tubes (Sets) Luer Slip Adapters, 40" Length	28
1181621100	SoftPack 12mL Syringe 21 x 1"	19	170030	3-way stopcock (no port covers)	28
1181621112	SoftPack 6mL Syringe 21 x 1 1/2"	18	170031	3-way stopcock with male luer lock adapter	28
1181622112	SoftPack 6mL Syringe 22 x 1 1/2"	18	170042	4-way stopcock (male luer slip) with 20" extension set	28
1182000555	SoftPack 20mL Syringe Regular Tip	19	170060	Opaque Nylon 3-way stopcock (no port covers)	28
1182000777	SoftPack 20mL Syringe Luer Lock Tip	19	170061	Opaque Nylon 3-way stopcock with male luer lock adapter	28
1183500555	SoftPack 35mL Syringe Regular Tip	20	31139747	Extra Keys for Hidden Bracket Set 31158960	54
1183500777	SoftPack 35mL Syringe Luer Lock Tip	20	31139762	Single Glove Box Holder	54
1183500888	SoftPack 35mL Syringe Catheter Tip	20	31140109	Cart for 7 & 10 Gallon Sharps-A-Gator™ (4-Castered)	43
1186000444	SoftPack 60mL Syringe Catheter Tip	20	31140349	Frame Cart for 19 Gallon Container	43
1186000555	SoftPack 60mL Syringe Regular Tip	20	31142222	Muti-Purpose Sharps Container with Sliding Lid 2 Gallon	48
1186000777	SoftPack 60mL Syringe Luer Lock Tip	20	31143533	Sharps-A-Gator™ Large Volume Sharps Container, 7.5 Gallon	43
1188001112	21 x 1 1/2" LDS Needle	13	31143665	Sharps-A-Gator™ Large Volume Sharps Container, 10 Gallon	43
1188005058	25 x 5/8" LDS Needle	13	31143699	Muti-Purpose Sharps Container with Sliding Lid, 1 Gallon	48
1188100555	SoftPack 1mL Syringe Regular Tip	20	31143897	5 Quart Tortuous Path Sharps Container, Clear	
1188100777	1mL Insulin Syringe Luer Lock	16	31144010	5 Quart Tortuous Path Sharps Container, Transparent Red	38
1188125058	SoftPack 1mL Insulin Syringe 25 x 5/8"	14	31156550	Sharps-A-Gator™ Large Volume Sharps Container, 7 Gallon	43

INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
31158549	5 Quart Cabinet	54	8539SA	Holder for 2 Gallon & In-Room™ Containers	53
31158960	Hidden Bracket and Key	54	85401H	3 & 4 Gallon SharpStar™ IV Wall Enclosure	40
31217949	Cart for 7 & 10 Gallon Sharps-A-Gator™	43	8540SA	SharpStar™ Container with Counter Balanced Lid, 4 Gallon, Clear	36
31307005	2 & 3 Gallon GatorGuard™ Cabinet	40	8541SA	SharpStar™ Container with Counter Balanced Lid, 4 Gallon, Trans. Red	36
31307054	Locking Wire Bracket, 5 Quart, 2 and 3 Gallon	54	8543C	Non-Locking Bracket for SharpStar™ IV In-Room™ Containers	52
31307062	Non Locking 5 Quart, 2& 3 Gallon	54	8549SA	Holder for SharpStar™ IV 3 & 4 Gallon In-Room Container	53
31314886	Gatorguard™ Sharps Container, 3 Gallon, Transparent Red	38	8550B	Glove Dispenser with attached Wall Mounting Bracket	40
31320032	Glove Box Dispenser	40	8550LG	Glove Dispenser, Top Loading	40
31323333	Gatorguard™ Sharps Container, 2 Gallon, Transparent Red	38	8554B	Wall Mounting Bracket	40
31353553	5 Quart GatorGuard™ Jr. Cabinet	40	8555SA	Stackable Multi-Glove Dispenser, Side Loading	40
31353603	Gatorguard™ Jr. Sharps Container, 5 Quart, Transparent Red	38	8556H	2 & 5 Quart Wall Enclosure, Glove Dispenser & SharpStar™ Container	40
31369955	5 Quart GatorGuard™ Jr. Table Top Holder	40	8601RC	RCRA Hazardous Waste Container, 1.5 Quart	50
31378089	Sharps-A-Gator™ Large Volume Sharps Container, 19 Gallon	43	8602RC	RCRA Hazardous Waste Container, 2 Gallon	50
8301	Sharps Shuttle with Locking Mechanism	46	8605RC	RCRA Hazardous Waste Container, 5 Quart In-Room™	50
85021	In-Room™ System with Mailbox-Style Lid, 2 Quart, Clear	37	8607RC	RCRA Hazardous Waste Container, 8 Gallon	50
85031	In-Room™ System with Mailbox-Style Lid, 2 Quart, Transparent Red	37	8608RC	RCRA Hazardous Waste Container, 8 Gallon	50
8303SA	Transportable Sharps Container, 1 Quart, Red	46	8611RC	RCRA Hazardous Waste Container, 12 Gallon	50
8506SA	In-Room™ Sharps Containers with Counter Balanced Lid, 5 Quart	36	8612RC	RCRA Hazardous Waste Container, 12 Gallon	50
8507MW	In-Room™ Sharps Containers with Counter Balanced Lid, 1 Quart	39	8617RC	RCRA Hazardous Waste Container, 18 Gallon	50
8507SA	In-Room™ Sharps Container with Counter-Balanced Lid, 5 Quart	36	8618RC	RCRA Hazardous Waste Container, 18 Gallon	50
851201	In-Room™ Sharps Containers with Always-Open Lid, 5 Quart, Clear	37	8781	EnviroStar™ Sharps Disposal Container, 8 Gallon	46
85121	In-Room™ System with Mailbox-Style Lid, 5 Quart, Clear	37	8790	EnviroStar™ Sharps Disposal Container, 2 Gallon	46
851301	In-Room™ System with Always-Open Lid, 5 Quart, Transparent Red	37	8791	EnviroStar™ Sharps Disposal Container, 18 Gallon	46
85131	In-Room™ System with Mailbox-Style Lid, 5 Quart, Transparent Red	37	8820	PharmaSafety™ Sharps Containers, 2 Gallon	51
8515	Holder for 2 & 5 Quart In-Room™ Containers	53	88361H	PharmaStar™ Wall Enclosure, 3 Gallon	51
85161H	In-Room™ System Wall Enclosure, 2 & 5 Quart	40	8836SA	PharmaSafety™ Sharps Container, 3 Gallon	51
85165H	In-Room™ System Wall Enclosure, 2 & 5 Quart	40	8850	PharmaSafety™ Sharps Container, 8 Gallon, Hinged Lid	51
8516HDL	Replacement Wall Enclosure Door	40	8851	PharmaSafety™ Sharps Container, 8 Gallon, Sliding Lid	51
8518X	Locking Bracket for 2 & 5 Quart In-Room™ Containers	52	8860	PharmaSafety™ Sharps Container, 12 Gallon	51
8519C	Non-Locking Bracket for 2 & 5 Quart In-Room™ Containers	52	8870	PharmaSafety™ Sharps Container, 18 Gallon, Hinged Lid	51
8522MW	In-Room™ Sharps Containers with Counter Balanced Lid, 3 Quart	39	8871	PharmaSafety™ Sharps Container, 18 Gallon, Sliding Lid	51
85221	In-Room™ System with Always-Open Lid 3 Gallon, Clear	37	8881103025	Non-sterile 3mL Syringe Regular Tip	21
85221R	In-Room™ System with Always-Open Lid 3 Gallon, Transparent Red	37	8881103066	Non-sterile 3mL Syringe Luer Lock Tip	21
8524C	Non-Locking Bracket for 3 Gallon Multi-Purpose Containers	52	8881106010	Non-sterile 6mL Syringe Luer Lock Tip	21
8527R	Multi-Purpose Container with Rotor & Hinged Opening Lid, 3 Gallon	47	8881106028	Non-sterile 6mL Syringe Regular Tip	21
8528C	Locking Bracket for 3 Gallon Multi-Purpose 2 & 3 Gal. In-Room™	52	8881112059	Non-sterile 12mL Syringe Regular Tip	21
8529	Holder for 3 Gallon Multi-Purpose & In-Room™ Containers	53	8881112075	Non-sterile 12mL Syringe Eccentric Tip	21
85301H	2 & 3 Gallon Wall Enclosure	40	8881112083	Non-sterile 12mL Syringe Luer Lock Tip	21
85321	In-Room™ System with Always-Open Lid 2 Gallon, Clear	37	8881114014	Non-sterile 140mL Syringe Regular Tip	22
85321R	In-Room™ System with Always-Open Lid 2 Gallon, Transparent Red	37	8881114030	Non-sterile 140mL Syringe Luer Lock Tip	22
8533SA	In-Room™ Sharps Container with Counter Balanced Lid, 2 Gallon	36	8881114055	Non-sterile 140mL Syringe Catheter Tip	22
8534SA	In-Room™ Sharps Container with Counter Balanced Lid, 2 Gallon	36	8881114063	140mL Syringe Luer Lock Tip Sterile	22
8536SA	In-Room™ Sharps Container with Counter Balanced Lid, 3 Gallon	36	8881120029	Non-sterile 20mL Syringe Eccentric Tip	21
8537MW	In-Room™ Sharps Container with Counter-Balanced Lid, 3 Gallon	36	8881120037	Non-sterile 20mL Syringe Regular Tip	21
8537SA	In-Room™ Sharps Container with Counter Balanced Lid, 3 Gallon	39	8881120193	Non-sterile 20mL Syringe Luer Lock Tip	21

INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
8881135019	Non-sterile 35mL Syringe Luer Lock Tip	21	8881202330	Rigid Pack Blunt Cannula 17 x 1 1/2"	13
8881135068	Non-sterile 35mL Syringe Eccentric Tip	21	8881202348	Rigid Pack Blunt Cannula 18 x 1"	13
8881135084	Non-sterile 35mL Syringe Regular Tip	21	8881202355	Rigid Pack Blunt Cannula 19 x 1 1/2"	13
8881160157	Non-sterile 60mL Syringe Catheter Tip	21	8881202363	Rigid Pack Blunt Cannula 20 x 1 1/2"	13
8881160405	Non-sterile 60mL Syringe Regular Tip	21	8881202371	Rigid Pack Blunt Cannula 21 x 1"	13
8881160629	Non-sterile 60mL Syringe Luer Lock Tip	21	8881202389	Rigid Pack Blunt Cannula 22 x 1 1/2"	13
8881200011	Rigid Pack Hypodermic Needle 14 x 1 1/2" IN	12	8881202397	Rigid Pack Blunt Cannula 23 x 1"	13
8881200029	Rigid Pack Hypodermic Needle 15 x 1 1/2" B	12	8881204005	Medication Transfer Needle 20 x 1"	26
8881200037	Rigid Pack Hypodermic Needle 16 x 1 1/2"	12	8881205010	Introducer Needle with Aluminum Hub 18 x 1 1/4"	29
8881200045	Rigid Pack Hypodermic Needle 16 x 1 1/2" B	12	8881216033	Monoject™ Supra/Hone Blood Collection Needle 21 x 1 1/4"	34
8881200052	Rigid Pack Hypodermic Needle 16 x 1"	12	8881216066	Monoject™ Supra/Hone Blood Collection Needle 22 x 1"	34
8881200060	Rigid Pack Hypodermic Needle 16 x 3/4" IN	12	8881225110	1" Safety Needle Holder	30
8881200078	Rigid Pack Hypodermic Needle 18 x 1 1/2"	12	8881225115	1 1/2" Safety Needle Holder	30
8881200086	Rigid Pack Hypodermic Needle 18 x 1 1/2" B	12	8881225119	Luer Lock Collection Set	30
8881200110	Rigid Pack Hypodermic Needle 18 x 1" B	12	8881225174	Blood Collection Set with Female Luer 19 x 3/4"	31
8881200136	Rigid Pack Hypodermic Needle 19 TW x 1 1/2"	12	8881225182	Blood Collection Set with Female Luer 21 x 3/4"	31
8881200144	Rigid Pack Hypodermic Needle 19 TW x 1 1/2"	12	8881225190	Blood Collection Set with Female Luer 23 x 3/4"	31
8881200151	Rigid Pack Hypodermic Needle 19 TW x 1 1/4" B	12	8881225208	Blood Collection Set with Female Luer 25 x 3/4"	31
8881200169	Rigid Pack Hypodermic Needle 19 TW x 1"	12	8881225216	Tube Holder with Male Adapter	33
8881200185	Rigid Pack Hypodermic Needle 20 x 1 1/2"	12	8881225220	Blood Collection Assembly with BlunTip Cannula	33
8881200193	Rigid Pack Hypodermic Needle 20 x 1 1/2" B	12	8881225221	Blood Collection Set with Holder Attached 21 x 3/4"	31
8881200219	Rigid Pack Hypodermic Needle 20 x 1"	12	8881225224	Multi-Sample Collection Set (Dome) - Male	32
8881200227	Rigid Pack Hypodermic Needle 20 x 1" B	12	8881225232	Tube Holder with Female Adapter	33
8881200235	Rigid Pack Hypodermic Needle 20 x 3"	12	8881225240	Holder and Safety Cap (Female)	32
8881200318	Rigid Pack Hypodermic Needle 21 x 2"	12	8881225241	Female Safety Transfer Device	32
8881200326	Rigid Pack Hypodermic Needle 22 x 1 1/2"	12	8881225245	Male Safety Collection Set	32
8881200342	Rigid Pack Hypodermic Needle 22 x 1"	12	8881225257	Luer Adapter	34
8881200383	Rigid Pack Hypodermic Needle 23 x 1"	12	8881225265	Blood Collection Set with Holder Attached 23 x 3/4"	31
8881200433	Rigid Pack Hypodermic Needle 25 x 1 1/4"	12	8881225273	Blood Collection Set with Holder Attached 25 x 3/4"	31
8881200441	Rigid Pack Hypodermic Needle 25 x 2"	12	8881225281	Blood Collection Set with Multi-Sample Luer Adapter 19 x 3/4"	31
8881200466	Rigid Pack Hypodermic Needle 25 x 5/8"	12	8881225299	Blood Collection Set with Multi-Sample Luer Adapter 21 x 3/4"	31
8881200508	Rigid Pack Hypodermic Needle 27 x 1 1/4"	12	8881225307	Blood Collection Set with Multi-Sample Luer Adapter 23 x 3/4"	31
8881200516	Rigid Pack Hypodermic Needle 27 x 1/2" IN	12	8881225315	Blood Collection Set with Multi-Sample Luer Adapter 25 x 3/4"	31
8881200573	Rigid Pack Hypodermic Needle 14 x 2" IN	12	8881225390	Blood Collection Set with Female Luer 23 x 3/4"	31
8881200714	Rigid Pack Hypodermic Needle 18 x 1"	12	8881225707	Blood Collection Set with Multi-Sample Luer Adapter 23 x 3/4"	31
8881200755	Rigid Pack Hypodermic Needle 16 x 5/8"	12	8881227121	Safety Blood Collection Device 21G x 1 1/4"	30
8881200805	Rigid Pack Hypodermic Needle 14 x 1"	12	8881227122	Safety Blood Collection Device 22G x 1 1/4"	30
8881202017	Blunt Fill Cannula	24	8881250016	Rigid Pack Hypodermic Needle with Polypropylene Hub 18 x 1 1/2"	11
8881202018	1mL Regular 18 x 1"	24	8881250024	Rigid Pack Hypodermic Needle with Polypropylene Hub 18 x 1 1/2" B	11
8881202025	1mL Regular 18 x 1"	24	8881250032	Rigid Pack Hypodermic Needle with Polypropylene Hub 30 x 3/4"	11
8881202033	3mL Luer Lock 18 x 1"	24	8881250040	Rigid Pack Hypodermic Needle with Polypropylene Hub 18 x 1"	11
8881202116	6mL Luer Lock 18 x 1"	24	8881250057	Rigid Pack Hypodermic Needle with Polypropylene Hub 18 x 1" B	11
8881202132	12mL Luer Lock 18 x 1"	24	8881250065	Rigid Pack Hypodermic Needle with Polypropylene Hub 19 x 1 1/2"	11
8881202314	Rigid Pack Blunt Cannula 15 x 1 1/2"	13	8881250073	Rigid Pack Hypodermic Needle with Polypropylene Hub 19 x 1 1/2" B	11
8881202322	Rigid Pack Blunt Cannula 16 x 1 1/2"	13	8881250081	Rigid Pack Hypodermic Needle with Polypropylene Hub 19 TW x 1"	11

INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
8881250099	Rigid Pack Hypodermic Needle with Polypropylene Hub 19 TW x 1"	11	8881512258	Pharmacy Tray 12mL Luer Lock Tip	21, 26
8881250107	Rigid Pack Hypodermic Needle with Polypropylene Hub 20 x 1 1/2"	11	8881512597	Rigid Pack 12mL Syringe 18 x 1"	19
8881250115	Rigid Pack Hypodermic Needle with Polypropylene Hub 20 x 1 1/2"	11	8881512662	Rigid Pack 12mL Syringe 20 x 1 1/2"	19
8881250123	Rigid Pack Hypodermic Needle with Polypropylene Hub 20 x 1"	11	8881512738	Rigid Pack 12mL Syringe 21 x 1"	19
8881250149	Rigid Pack Hypodermic Needle with Polypropylene Hub 21 x 1 1/2"	11	8881512746	Rigid Pack 12mL Syringe 21 x 1 1/2"	19
8881250172	Rigid Pack Hypodermic Needle with Polypropylene Hub 21 x 1"	11	8881512811	12mL Syringe 22 x 1 1/2"	19
8881250206	Rigid Pack Hypodermic Needle with Polypropylene Hub 22 x 1 1/2"	11	8881512852	Rigid Pack 12mL Syringe Regular Tip	19
8881250214	Rigid Pack Hypodermic Needle with Polypropylene Hub 22 x 1 1/2" B	11	8881512860	12mL Syringe Eccentric Tip	19
8881250222	Rigid Pack Hypodermic Needle with Polypropylene Hub 22 x 1"	11	8881512878	12mL Syringe Luer Lock Tip	19
8881250230	Rigid Pack Hypodermic Needle with Polypropylene Hub 22 x 1"	11	8881512977	12mL Control Syringe Luer Lock Tip	22
8881250255	Rigid Pack Hypodermic Needle with Polypropylene Hub 23 x 1"	11	8881513025	3mL Syringe 20 x 3/4"	17
8881250263	Rigid Pack Hypodermic Needle with Polypropylene Hub 23 x 1/2" IN	11	8881513033	3mL Syringe 20 x 1"	17
8881250289	Rigid Pack Hypodermic Needle with Polypropylene Hub 23 x 3/4" IN	11	8881513058	3mL Syringe 20 x 1 1/2"	17
8881250305	Rigid Pack Hypodermic Needle with Polypropylene Hub 25 x 1"	11	8881513132	3mL Syringe 21 x 1"	17
8881250313	Rigid Pack Hypodermic Needle with Polypropylene Hub 25 x 5/8"	11	8881513157	3mL Syringe 21 x 1 1/2"	17
8881250321	Rigid Pack Hypodermic Needle with Polypropylene Hub 26 x 1/2"	11	8881513207	Pharmacy Trays 3mL Luer Lock Tip	21,26
8881250362	Rigid Pack Hypodermic Needle with Polypropylene Hub 27 x 1/2"	11	8881513231	Rigid Pack 12mL Syringe 22 x 1"	17
8881250545	Rigid Pack Hypodermic Needle with Polypropylene Hub 25 x 1 1/2"	11	8881513256	Rigid Pack 12mL Syringe 22 x 1 1/2"	17
8881305018	Filter Needle 20 x 1 1/2"	27	8881513330	Rigid Pack 12mL Syringe 23 x 1"	17
8881305109	Filter Aspirator 18 x 3"	27	8881513512	Rigid Pack 12mL Syringe 25 x 5/8"	17
8881305117	Filter Needle 18 x 1 1/2"	27	8881513538	Rigid Pack 3mL Syringe 25 x 1"	17
8881500014	Rigid Pack 1/2mL Insulin Syringe 28 x 1/2"	14	8881513744	Rigid Pack 3mL Syringe 27 x 1 1/4"	17
8881500105	Rigid Pack 1/2mL TB Syringe 28 x 1/2"	15	8881513918	Rigid Pack 3mL Syringe Regular Tip	17
8881500501	Allergy Trays 1/2mL Tuberculin Syringe 28 x 1/2"	16	8881513934	Rigid Pack 3mL Syringe Luer Lock Tip	17
8881501160	Rigid Pack 1mL TB Syringe 25 x 5/8"	15	8881516051	Rigid Pack 6mL Syringe 20 x 1 1/2"	18
8881501178	Rigid Pack 1mL TB Syringe 26 x 3/8"	15	8881516135	Rigid Pack 6mL Syringe 21 x 1"	18
8881501210	Rigid Pack 1mL Insulin Syringe 28 x 1/2"	14	8881516150	Rigid Pack 6mL Syringe 21 x 1 1/2"	18
8881501236	Allergy Trays 1mL Tuberculin Syringe 28 x 1/2"	16	8881516200	Pharmacy Trays 6mL Luer Lock Tip	21, 26
8881501368	Rigid Pack 1mL TB Syringe 27 x 1/2"	15	8881516911	Rigid Pack 6mL Syringe Regular Tip	18
8881501384	Rigid Pack 1mL Insulin Syringe Regular Tip	14	8881516937	Rigid Pack 6mL Syringe Lock Tip	18
8881501400	Rigid Pack 1mL TB Syringe Regular Luer Tip	15	8881520178	20mL Control Syringe Luer Lock Tip	22
8881501459	Pharmacy Tray 1mL Regular Tip	21, 26	8881520251	Pharmacy Tray 20mL Luer Lock Tip	21, 26
8881501608	Rigid Pack 1mL TB Syringe 28 x 1/2"	15	8881520657	Rigid Pack 20mL Syringe Luer Lock Tip	19
8881501640	Rigid Pack 1mL TB Syringe 25 x 5/8"	15	8881520665	Rigid Pack 20mL Syringe Eccentric Tip	19
8881501822	Rigid Pack 1mL Insulin Syringe 25 x 5/8"	14	8881520673	Rigid Pack 20mL Syringe Regular Tip	19
8881501962	Allergy Trays 1mL TB 27 x 1/2"	16	8881522000	Rigid Pack 12mL Safety Syringe Luer Lock Tip	8
8881501970	Rigid Pack 1mL Insulin Syringe 27 x 1/2"	14	8881522018	Rigid Pack 12mL Safety Syringe 20 x 1 1/2"	8
8881511110	1mL Insulin Safety Syringe 29 x 1/2"	9	8881522026	Rigid Pack 12mL Safety Syringe 21 x 1 1/2"	8
8881511136	1/2mL Insulin Safety Syringe 29 x 1/2"	9	8881533056	Rigid Pack 3mL Safety Syringe 20 x 1 1/2"	8
8881511144	3/10mL Insulin Safety Syringe 29 x 1/2"	9	8881533130	Rigid Pack 3mL Safety Syringe 21 x 1"	8
8881511201	1mL Tuberculin Safety Syringe 28 x 1/2"	9	8881533155	Rigid Pack 3mL Safety Syringe 21 x 1 1/2"	8
8881511235	1mL Tuberculin Safety Syringe 1 1/2"	9	8881533239	Rigid Pack 3mL Safety Syringe 22 x 1"	8
8881511310	1mL Insulin Safety Syringe 30 x 5/16"	9	8881533254	Rigid Pack 3mL Safety Syringe 22 x 1 1/2"	8
8881511336	1/2mL Insulin Safety Syringe 30 x 5/16"	9	8881533338	Rigid Pack 3mL Safety Syringe 23 x 1"	8
8881511344	3/10mL Insulin Safety Syringe 30 x 5/16"	9	8881533510	Rigid Pack 3mL Safety Syringe 25 x 5/8"	8

INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
8881533908	Rigid Pack 3mL Safety Syringe Luer Lock Tip	8	8881600350	SoftPack 1/2mL Insulin Syringe 29 x 1/2"	14
8881535101	Pharmacy Tray 35mL Luer Lock Tip	21, 26	8881600700	SoftPack 1/2mL Insulin Syringe 30 x 5/16"	14
8881535762	Rigid Pack 35mL Syringe Luer Lock Tip	20	8881600800	SoftPack 3/10mL Insulin Syringe 30 x 5/16"	14
8881535770	Rigid Pack 35mL Syringe Catheter Tip	20	8881601101	SoftPack 1mL Insulin Syringe 28 x 1/2"	14
8881535788	Rigid Pack 35mL Syringe Eccentric Tip	20	8881601358	SoftPack 1mL Insulin Syringe 29 x 1/2"	14
8881535796	Rigid Pack 35mL Syringe Regular Tip	20	8881601600	SoftPack 1mL Insulin Syringe 30 x 5/16"	14
8881540101	I.V. Access Cannula	25	8881610102	Blood Collection Needle and Tube Holder Regular 13mm	34
8881540111	Med-Prep Cannula	24	8881676012	Sharps-A-Gator™ Aluminum Bracket	52
8881540120	Monoject BlunTip™ with Vial Access Pin	25	8881676236	Monoject™ Vertical Top Container, 4 Quart	48
8881540122	12mL Med-Prep Combination	24	8881676285	Monoject™ Vertical Top Container, 8 Quart	48
8881540133	3mL Med-Prep Combination	24	8881676434	Monoject™ Vertical Top Container, 14 Quart	48
8881540166	6mL Med-Prep Combination	24	8881676624	Sharps-A-Gator™ Locking Wall Cabinet	52
8881541034	3mL Monoject BlunTip™ Cannula	25	8881676632	Sharps-A-Gator™ Glove Box	52
8881541067	6mL Monoject BlunTip™ Cannula	25	8881676814	Monoject™ Horizontal Drop Container, 4 Quart	48
8881541125	12mL Monoject BlunTip™ Cannula	25	8881676822	Monoject™ Horizontal Drop Container, 8 Quart	48
8881541208	12mL Monoject BlunTip™ with Vial Access Pin	25	8881682010	Safety Syringe Tip Cap	9
8881541307	3mL Monoject BlunTip™ with Vial Access Pin	25	8881682085	Syringe Tip Caps Bag pack of 100	29
8881541604	6mL Monoject BlunTip™ with Vial Access Pin	25	8881682101	Syringe Tip Caps Tray pack of 10	29
8881560125	Rigid Pack 60mL Syringe Luer Lock Tip	20	8881682119	Syringe Tip Caps Tray pack of 25	29
8881560141	Rigid Pack 60mL Syringe Catheter Tip	20	8881701406	SoftPack 1mL Insulin Syringe 27 x 1/2"	14
8881560182	Rigid Pack 60mL Syringe Eccentric Tip	20	8881811310	Magellan™ Safety Needle & Syringe 1mL 23 x 1"	6
8881560224	Rigid Pack 60mL Syringe Regular Tip	20	8881811510	Magellan™ Safety Needle & Syringe 1mL 25 x 1"	6
8881560232	Pharmacy Tray 60mL Luer Lock Tip	21,26	8881811558	Magellan™ Safety Needle & Syringe 1mL 25 x 5/8"	6
8881560265	Rigid Pack 60mL Syringe Toomey Tip	20	8881820008	Finger Grip Extender for 20mL	29
8881566007	Rigid Pack 6mL Safety Syringe Luer Lock Tip	8	8881822015	Magellan™ Safety Needle & Syringe 12mL 20 x 1 1/2"	6
8881566015	Rigid Pack 6mL Safety Syringe 20 x 1 1/2"	8	8881822110	Magellan™ Safety Needle & Syringe 12mL 21 x 1"	6
8881566023	Rigid Pack 6mL Safety Syringe 21 x 1 1/2"	8	8881822115	Magellan™ Safety Needle & Syringe 12mL 21 x 1 1/2"	6
8881570121	Saline I.V. Flush 12mL Syringe 10mL Fill	23	8881822215	Magellan™ Safety Needle & Syringe 12mL 22 x 1 1/2"	6
8881570123	Saline I.V. Flush 12mL Syringe 3mL Fill	23	8881822810	Magellan™ Safety Needle & Syringe 12mL 18 x 1"	6
8881570125	Saline I.V. Flush 12mL Syringe 5mL Fill	23	8881833010	Magellan™ Safety Needle & Syringe 3mL 20 x 1"	6
8881570128	Saline I.V. Flush 12mL Syringe 10mL Fill for Sterile Delivery	23	8881833015	Magellan™ Safety Needle & Syringe 3mL 20 x 1 1/2"	6
8881570129	Saline I.V. Flush 12mLSyringe 10mL Fill for Sterile Delivery	23	8881833110	Magellan™ Safety Needle & Syringe 3mL 21 x 1"	6
8881570300	Saline I.V. Flush 3mL Syringe 2.5mL Fill	23	8881833115	Magellan™ Safety Needle & Syringe 3mL 21 x 1 1/2"	6
8881571121	Saline I.V. Flush 12mL Syringe 2.5mL Fill with Cannula	23	8881833210	Magellan™ Safety Needle & Syringe 3mL 22 x 1"	6
8881571123	Saline I.V. Flush 12mL Syringe 3mL Fill with Cannula	23	8881833215	Magellan™ Safety Needle & Syringe 3mL 22 x 1 1/2"	6
8881571125	Saline I.V. Flush 12mL Syringe 5mL Fill with Cannula	23	8881833310	Magellan™ Safety Needle & Syringe 3mL 23 x 1"	6
8881571300	Saline I.V. Flush 3mL Syringe 2.5mL Fill with Cannula	23	8881833510	Magellan™ Safety Needle & Syringe 3mL 25 x 1"	6
8881580121	Heparin I.V. Flush (10 units/mL) 12mL Syringe 10mL Fill	23	8881833515	Magellan™ Safety Needle & Syringe 3mL 25 x 1 1/2"	6
8881580123	Heparin I.V. Flush (10 units/mL) 12mL Syringe 3mL Fill	23	8881833558	Magellan™ Safety Needle & Syringe 3mL 25 x 5/8"	6
8881580125	Heparin I.V. Flush (10 units/mL) 12mL Syringe 5mL Fill	23	8881833810	Magellan™ Safety Needle & Syringe 3mL 18 x 1"	6
8881590121	Heparin I.V. Flush (100 units/mL) 12mL Syringe 10mL Fill	23	8881835006	Finger Grip Extender for 35mL	29
8881590123	Heparin I.V. Flush (100 units/mL) 12mL Syringe 3mL Fill	23	8881850010	Magellan™ Hypodermic Safety Needle 20 x 1"	5
8881590125	Heparin I.V. Flush (100 units/mL) 12mL Syringe 5mL Fill	23	8881850015	Magellan™ Hypodermic Safety Needle 20 x 1 1/2"	5
8881600004	SoftPack 1/2mL Insulin Syringe 28 x 1/2"	14	8881850110	Magellan™ Hypodermic Safety Needle 21 x 1"	5
8881600145	SoftPack 3/10mL Insulin Syringe 29 x 1/2"	14	8881850115	Magellan™ Hypodermic Safety Needle 21 x 1 1/2"	5

INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
8881850158	Magellan™ Hypodermic Safety Needle 21 x 5/8"	5	8888173542	EZ-FLO Stopcock Double 3-way, male luer lock adapter port covers	28
8881850210	Magellan™ Hypodermic Safety Needle 22 x 1"	5	8888173690	EZ-FLO Stopcock 3-way with male luer lock adapter and port covers	28
8881850215	Magellan™ Hypodermic Safety Needle 22 x 1 1/2"	5	8888174201	EZ-FLO Stopcock 4-way with 20" extension tube and port covers	28
8881850310	Magellan™ Hypodermic Safety Needle 23 x 1"	5	8888174698	EZ-FLO Stopcock 4-way with male luer lock adapter and port covers	28
8881850358	Magellan™ Hypodermic Safety Needle 23 x 5/8"	5	8900MW	Phlebotomy Sharps Container, 1 Quart Container	39
8881850510	Magellan™ Hypodermic Safety Needle 25 x 1"	5	8900SA	Phlebotomy Sharps Container, 1 Quart	49
8881850515	Magellan™ Hypodermic Safety Needle 25 x 1 1/2"	5	8901SA	Phlebotomy Sharps Container, 1 Pint	49
8881850558	Magellan™ Hypodermic Safety Needle 25 x 5/8"	5	8909	Transportable Sharps Container, 1.5 Quart	46
8881850810	Magellan™ Hypodermic Safety Needle 18 x 1"	5	8910	Holder for 1 Quart Phlebotomy Container	53
8881850815	Magellan™ Hypodermic Safety Needle 18 x 1 1/2"	5	8912SA	Wall Bracket for 1 Quart Phlebotomy Container	52
8881850910	Magellan™ Hypodermic Safety Needle 19 x 1"	5	8920SA	Multi-Purpose Container with Rotor Opening Lid, 2 Quart	47
8881850915	Magellan™ Hypodermic Safety Needle 19 x 1 1/2"	5	8930SA	Large Volume Container with Hinged Lid, 30 Gallon	41
8881852010	Magellan™ Hypodermic Safety Needle	6	8931	ChemoSafety™ Container, 12 Gallon, Hinged Lids	45
8881860004	Finger Grip Extender for 60mL	29	8931PG2	PG II Chemotherapy Container, 12 Gallon, Hinged Lid	44
8881866015	Magellan™ Safety Needle & Syringe 6mL 20 x 1 1/2"	6	8932	Large Volume Container with Hinged Lid, 12 Gallon	41
8881866110	Magellan™ Safety Needle & Syringe 6mL 21 x 1"	6	8932PG2	PG II Sharps Disposal Container, 12 Gallon, Hinged Lid	44
8881866115	Magellan™ Safety Needle & Syringe 6mL 21 x 1 1/2"	6	8933	Large Volume Container with Hinged Lid, 12 Gallon	41
8881866215	Magellan™ Safety Needle & Syringe 6mL 22 x 1 1/2"	6	8934	ChemoSafety™ Container, 12 Gallon, Sliding Lid	45
8881866810	Magellan™ Safety Needle & Syringe 6mL 18 x 1"	6	8934PG2	PG II Chemotherapy Container, 12 Gallon, Sliding Lid	44
8881882712	Magellan™ Tuberculin Safety Syringe 1mL 27 x 1/2"	7	8935	Large Volume Container with Sliding Lid, 12 Gallon	41
8881882812	Magellan™ Tuberculin Safety Syringe 1mL 28 x 1/2"	7	8936PG2	PG II Sharps Disposal Container, 12 Gallon, Red, Sliding Lid	44
8881892910	Magellan™ Insulin Safety Syringe 1mL 29 x 1/2"	7	8936SA	Large Volume Container with Sliding Lid, 12 Gallon	41
8881892930	Magellan™ Insulin Safety Syringe 3/10mL 29 x 1/2"	7	8938	Large Volume Container with Sliding Lid, 18 Gallon	41
8881892950	Magellan™ Insulin Safety Syringe 1/2mL 29 x 1/2"	7	8938FP	SharpsCart™ Foot Pedal Cart for 12 & 18 Gallon, Slide	42
8881893010	Magellan™ Insulin Safety Syringe 1mL 30 x 5/16"	7	8938FP	SharpsCart™ Foot Pedal Cart for 12 & 18 Gallon, Slide	42
8881893030	Magellan™ Insulin Safety Syringe 3/10mL 30 x 5/16"	7	8939	ChemoSafety™ Container, 18 Gallon, Sliding Lid	45
8881893050	Magellan™ Insulin Safety Syringe 1/2mL 30 x 5/16"	7	8939PG2	PG II Chemotherapy Container 18 Gallon, Sliding Lid	44
8881901006	Oral Medication Syringe Amber 1mL	22	8940FP	Molded SharpsCart™ Foot Pedal Cart 12 & 18 Gallon, Slide	42
8881901014	Oral Medication Syringe Clear 1mL	22	8941FP	Molded SharpsCart™ Foot Pedal Cart 12 & 18 Gallon, Hinged	42
8881903002	Oral Medication Syringe Clear 3mL	22	8950SA	Multi-Purpose Container with Rotor Opening Lid, 5	47
8881903010	Oral Medication Syringe Amber 3mL	22	8963	Locking Bracket for Multi-Purpose & Chemosafety™ Container	52
8881906005	Oral Medication Syringe Amber 6mL	22	8964	DialySafety Sharps Disposal Container, 3 Gallon	46
8881906104	Oral Medication Syringe Clear 6mL	22	89651	Multi-Purpose Container with Horizontal Drop, 2 Gallon, Red	48
8881907003	Oral Medication Syringe Amber 10mL	22	89671	Multi-Purpose Container with Horizontal Drop, 2 Gal. Transparent Red	48
8881907102	Oral Medication Syringe Clear 10mL	22	8967MW	Ancillary / Multi-Purpose Container, 2 Gallon	39
8888115006	Heparin Lock For I.V. Device	28	8970	Multi-Purpose Container with Rotor Opening Lid, 2 Gallon	47
8888118059	Filter Connector	26	8974	Holder for Phlebotomy, Multi-Purpose & ChemoSafety™ Container	53
8888173005	EZ-FLO Stopcock 1-way with port covers male luer lock adapter	28	8975	Wall / Cart Bracket for Multi-Purpose & Chemosafety™ Container	52
8888173104	EZ-FLO Stopcock 1-way with port cover	28	8979C	Multi-Purpose Containers with Rotor Opening Lid 2 Gallon	47
8888173203	EZ-FLO Stopcock 3-way with injection site and 20" extension tube	28	8979MW	Ancillary / Multi-Purpose Container with Rotor Opening Lid, 2 Gallon	39
8888173336	EZ-FLO Stopcock 3-way with injection site and 33" extension tube	28	8980	Large Volume Container with Hinged Lid, 8 Gallon	41
8888173500	EZ-FLO Stopcock 3-way only	28	8980FP	SharpsCart™ Foot Pedal Cart for Gallon Hinged Lid Containers	42
8888173518	EZ-FLO Stopcock 3-way with port covers	28	8980MW	Large Volume Sharps Container with Hinged Lid, 8 Gallon	39
8888173526	EZ-FLO Stopcock 3-way with 20" extension tube and port covers	28	8980S	Large Volume Container with Sliding Lid, 8 Gallon	41
8888173534	EZ-FLO Stopcock 3-way with 33" extension tube and port covers	28	8981FP	SharpsCart™ Foot Pedal Cart, 8 Gallon, Slide Lid	42

INDEX

CODE	PRODUCT	PAGE
8982	ChemoSafety™ Container, 2 Gallon Hinged Lids	45
8984	Locking Wire Bracket for 8 Gallon Large Volume Container	52
8985	ChemoSafety™ Container, 8 Gallon Hinged Lid	45
8985MW	ChemoSafety™ Container, with Hinged Lid, 8 Gallon	39
8985PG2	PG II Chemotherapy Container, 8 Gallon, Hinged Lid	44
8985S	ChemoSafety™ Container, 8 Gallon, Sliding Lid	45
8985SPG2	PG II Chemotherapy Container, 8 Gallon, Sliding Lid	44
8989	ChemoSafety™ Container, 18 Gallon, Hinged Lid	45
8989PG2	PG II Chemotherapy Container, 18 Gallon, Hinged Lid	44
8990SA	Multi-Purpose Container with Hinged Lid, 2 Gallon	47
8991	Large Volume Container with Hinged Lid, 18 Gallon	41
8991FP	SharpsCart™ Foot Pedal Cart, 12 & 18 Gallon, Hinged Lid	42
8992H	Floor cart with wheels for 8, 12 & 18 Gallon Container	42
8993H	Floor bracket with wheels for 30 Gallon Container	42
8997	Large Volume Container with Hinged Lid, 8 Gallon	41
8997PG2	PG II Sharps Disposal Container, 8 Gallon, Hinged Lid	44
8997S	Large Volume Container with Sliding Lid, 8 Gallon	41
8997SPG2	PG II Sharps Disposal Container, 8 Gallon, Sliding Lid	44
8998	Large Volume Container with Hinged Lid, 18 Gallon	41
8998PG2	PG II Sharps Disposal Container 18 Gallon, Hinged Lid	44
8998S	Large Volume Container with Sliding Lid, 18 Gallon	41
8998SPG2	PG II Sharps Disposal Container 18 Gallon, Sliding Lid	44

Notes

Notes



COVIDIEN, COVIDIEN with logo and the Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG.
TM* Trademarks of their respective owners.
Other brands are trademarks of a Covidien company.
© 2010 Covidien AG.
H6200NSHP 7.5M 09/10