

FLEXON™

Temporary Cardiac Pacing Lead

FLEXON™ Temporary Cardiac Pacing Lead

Description

FLEXON™ temporary cardiac pacing lead is stainless steel twisted, multistrand wire coated with orange or white PTFE poly (terafluoroproethylene) or clear FEP poly (tetrafluoroethylene-co-hexafluoropropylene).

Characteristics

- Exceptional strength with great flexibility and handling characteristics.
- Non-magnetic and electro-passive in tissue fluids.

Indications

FLEXON™ temporary cardiac pacing lead is indicated for use in temporary atrial and ventricular pacing and sensing during and after cardiac surgery.

Contraindications

When permanent pacing or monitoring is required, the use of **FLEXON™** temporary cardiac pacing lead is contraindicated. The use of this product is contraindicated in patients with known sensitivities or allergies to its components. The particular medical condition and anatomy of the patient may dictate the lead systems and implantation procedure to be used.

Actions

FLEXON™ temporary cardiac pacing lead elicits minimal acute inflammatory reaction in tissues and is not absorbed.

Tensile Strength

Permanent

Sterilization

Ethylene oxide




Packaging

FLEXON™ temporary cardiac pacing lead is available sterile and in various sizes, lengths and quantities using ATRAUMATIC™ or BreakAway™ needles.



FLEXON™ temporary cardiac pacing lead are available coated with orange or white PTFE or clear FEP.





TAPER POINT 		Length	Color	Units Per Box	USP	4-0	3-0	2-0	0
					M	1.5	2	3	3.5
Straight 44mm  SC	1/2 Circle 17mm CV-23	24	Clear	12			2591-43		
Straight 60mm  SC-2	1/2 Circle 17mm CV-23	24	Clear	12			2592-43		
Straight 88mm  SC-6	1/2 Circle 17mm CV-23	24	Clear	12			2586-43 ➡		
Straight 60mm  SC-2	1/2 Circle 20mm CV-24	24	Clear	12					2593-63



TAPER POINT				Length	Color	Units Per Box	USP M	4-0 1.5	3-0 2	2-0 3	0 3.5
Straight 60mm	1/2 Circle 26mm	24	Clear	12							2597-63
											
SC-2	V-20										
Straight 88mm	1/2 Circle 26mm	24	Clear	12							2589-63 ➡
		24	Clear	36							2589-61 ➡
		24	Orange	12					2617-53 ➡		
		2x24	White Orange	12					2623-53 ➡		
		24	White	12					2615-53 ➡		
											
SC-6	V-20										