



Obstetric & Neonatal Products

1-800-962-9888



COVIDIEN

The products shown in this catalog may not be licensed or registered in all parts of the world.
Please contact your sales representative for availability.

Table of Contents

Obstetric & Neonatal Products

Kendall™ Obstetric Products

INTRAUTERINE PRESSURE CATHETER (IUPC)

Kendall™ IUPC & Cables	4
------------------------------	---

FETAL SPIRAL ELECTRODE (FSE)

Kendall™ FSE System & Cables	5
FSE Attachment Pads	5

ABDOMINAL BELTS, BINDERS AND STRAPS

Belts	6
Binders	7
Straps	7

ULTRASOUND GEL

Ultrasound Gel	7
----------------------	---

DELIVERY KITS

Vaginal Delivery Kits	8
C-Section Delivery Kits	9-10

FETAL MONITORING PAPER

Monitoring Paper	11
------------------------	----

CIRCUMCISION PRODUCTS

Devon™ Circumcision Trays	12
Single Sterile Clamps	12

Argyle™ Neonatal Products

NEONATAL CATHETERS

Argyle™ Umbilical Vessel Catheters (UVCs)	13
Umbilical Catheter Anchor	13
Argyle™ Umbilical Vessel Catheter Trays	14
Argyle™ Peripherally Inserted Central Catheters PICC(s) & Trays	15
NICU/PICU Procedure Tray	16
Argyle™ Replogle Catheters	16

NEONATAL ACCESSORIES

Cord Clamp Clippers	17
Argyle™ CPAP Cannulas	17
Baby Board I.V. Limb Support Boards	17
Infant Heel Warmers	18
Baby Beanies	18
Temperature Probe Covers	18
Baby Shades Phototherapy Eye Protectors	19
Baby Adhesive Tape	19


NEONATAL PRE-WIRED ELECTRODES

Pre-Wired Clear Tape Series	20
Pre-Wired Foam 4100 Series	21
Pre-Wired 1050 Cloth Series	21
Limb Band Specialty Pre-Wired	21

Intrauterine Pressure Catheter



Kendall™ IUPC



Transducer Tip



Catheter Body



Zero/Rezero Button



Cable Check



Fluid Resistant Connector

Kendall™ IUPC

The patented Kendall™ IUPC was developed to provide reliable tracings and convenience helping the clinician deliver the best patient care

- LATEX FREE
- Kendall™ IUPC's small, rounded transducer tip allows for easy insertion
- The catheter body is flat to gently deflect from tissue, giving added confidence during insertion and reduce trauma
- Clear, flexible amnioinfusion line delivers 55 ml/mm at gravity feed for fast infusion of liquids
- Fluid resistant connector protects the connection from fluid splash minimizing erroneous readings
- Convenient Zero/Rezero Button on cable has an Indicator Light that provides additional confirmation during the zero
- Unique Cable Check feature allows clinicians to easily identify a cable issue

Code	Description	Ship Case
56300	Kendall™ IUPC	
10		
Reusable Cables		
56311	For most Philips Medical Monitors	1
56312	For Philips Medical Avalon™ Monitors	1
56321	For GE Medical Corometrics™ Monitors	1
56330	For Spacelabs Monitors	1

Fetal Spiral Electrode



Kendall™ FSE

Leg Plate



Safety handle



Catheter Body



FSE Attachment Pads

Kendall™ FSE System

- The Kendall™ FSE sets the standard for compliance, convenience and performance
- LATEX FREE

Compliant

- Fully protected lead wires with a safety cap that plugs securely into our new cabling system
- Conforms to FDA safety guidelines for all Fetal Spiral Electrode systems
- Spiral is recessed to reduce risk of needle sticks when removing from package

Convenient and Easy to Use

- One step removal of the drive and guide tube
- Flexible guide tube can be gently shaped for easier insertion
- Simple connection to leg plate, no orientation required
- Unique plate for secure connections
- Attaches to the patient via a self-adhesive grounding electrode

Code	Description	Fetal Monitoring System	Ship Case
FSE			
31479549	Kendall™ FSE 2100		50
Attachment Pads			
31481396	FCB 2000 (small 1 7/8")		50
50000095	FCB 2001 (large 3 3/4" x 1 7/8")		50
Reusable Cables			
Philips Medical System			
31479754	FCB 100	8040A	1
31479846	FCB 300	For most Philips monitors—2 ft	1
50000331	FCB 308	For most Philips monitors—8 ft	1
GE-Medical			
31479812	FCB 400	Monitors with round connector	1
50000107	FCB 700	Monitors with rectangular connector	1
Spacelabs			
50000106	FCB 600	MOM 94000 (with 700-0024-00 pigtail supplied by Spacelabs, IM 77, AMS-02-7700)	1

Abdominal Belts



Buttonhole Belt



Buttonhole Belt



Tab Belt



Mesh Belt

Abdominal Belts and Binders

- Soft materials minimize skin irritation
- Large variety of sizes and packaging configurations available
- Now with more latex-free elastic to retain shape and hold transducers snugly
- Available in pre-cut or roll form

CODE	Description	Ship Case
Buttonhole Style		
Knit Elastic		
31410270	LATEX FREE - 2-3/8" x 48", buttonholes every 1 1/4", finished ends, 1 pink & 1 blue belt per set	50 sets
40000007	LATEX FREE - 2-3/8" x 48", buttonholes every 1 1/4", finished ends 1 pink & 1 blue belt per set	100 sets
30966348	2-3/8" x 150' roll, buttonholes every 1 1/4", white, must be cut to length	1 roll
30950300	2-3/8" x 48", buttonholes every 1 1/4", finished ends, 1 pink striped blue belt per set	50 sets
30950292	2-3/8" x 48", buttonholes every 1 1/4", finished ends 1 pink & blue striped belt per set	100 sets
56101	LATEX FREE - 2-3/8" x 48", buttonholes every 1 1/4", finished ends, 1 pink & blue striped belt per set	50 sets
Woven Elastic		
30817608	2-3/8" x 52", buttonholes every 3/4", navy blue	50 sets
Accessories		
30952090	100 bobbins per package	1 pack
Tab Style		
Knit Elastic- Latex Free		
40000008	1-1/2" x 36", 1 pink & blue belt per set	50 sets
40000009	1-1/2" x 36", 1 pink & blue belt per set	100 sets
56102	1-1/2" x 42", 1 pink & blue striped belt per set	50 sets
Woven Elastic		
530975232	1-1/2" x 40", 1 pink belt and 1 blue belt per set	50 sets
30950169	1-1/2" x 40", 1 pink belt and 1 blue belt per set	100 sets
Open Weave Mesh- Latex Free		
30955689	2" x 52", SPANDEX mesh belt , finished ends, 2 white belts per set	50 sets
Foam- Latex Free		
31325965	1-1/2" x 48", foam belt with pink or blue closure tabs	150 ea
31325981	1-1/2" x 60", foam belt with pink or blue closure tabs	150 ea

Binders & C-Section Straps



Binders



Straps

CODE	Width	Description	Ship Case
Binders- Latex Free			
31364394	10"	Small, beige, Pre-Cut, 16" in length	50
31364410	12"	Medium, beige, Pre-Cut, 16" in length	50
31364428	14"	Large, beige, Pre-Cut 16" in length	50
31364436	16"	X-Large, beige, Pre-Cut 16" in length	50
31363362	10"	Small, beige, , Roll, 33' in length, must be cut to length	1 roll
31364360	12"	Medium, beige, Roll, 33' in length, must be cut to length	1 roll
31364378	14"	Large, beige, Roll, 33' in length, must be cut to length	1 roll
C-Section Straps- Latex Free			
30951696		3" x 84"	25

Ultrasound Gels



CODE	Description	Ship Case
4060	Ultrasound Gel (250 mL bottle)	12
4061	Ultrasound Gel (5 liter "cube-tainer" with dispenser bottle)	1

Vaginal & C-Section Kits



Vaginal Deluxe Delivery Kit

Vaginal & C-Section Kits - Latex Free

- Standard kits reduce cost and waste while increasing efficiency
- Commonly used items in a single, convenient package
- C-Section kits feature the well-known, quality Devon™ product components such as basins, needle counters, lite gloves and Valleylab ES pencils
- Plus Kit - additional OB and major surgical components
- Deluxe Kit - Plus Kit components with drape and gown

Code	Description	Ship Case
52000013	Vaginal Plus Delivery Kit	8
	Plastics/Surgical	
	1-Placenta Basin/Lid	1-Red Blood Tube
	1-Bulb Syringe	1-Purple Blood Tube
	1-Umbilical Cord Clamp	
	Soft Goods/Drapes/Gowns	
	1-Maternity Pad	1-Baby Blanket
	4-Cotton Towel	10-4" x 4" X-Ray Dect.
	1-Vaginal Packing Sponge	1-Under Buttock Drape
	Wraps	
	1-50" x 50" CSR Fold	
52000014	Vaginal Deluxe Delivery Kit	8
	Plastics/Surgical	
	1-Placenta Basin/Lid	1-Red Blood Tube
	1-Bulb Syringe	1-Purple Blood Tube
	1-Umbilical Cord Clamp	
	Soft Goods/Drapes/Gowns	
	1-Maternity Pad	1-Baby Blanket
	4-Cotton Towel	10-4" x 4" X-Ray Dect.
	1-Vaginal Packing Sponge	1-Under Buttock Drape
	1-Leggings	1-Abdominal Drape
	1-X-Large Gown	
	Wraps	
	1-50" x 50" CSR Fold	1-50" x 50" Wrap

Vaginal & C-Section Kits



C-Section Deluxe Kit

Code	Description	Ship Case
52000015	C-Section Plus Kit Plastics/Surgical 2-7 Liter Basins 1-Heart Bowl 1-Bulb Syringe 1-Red Blood Tube 1-Needle Counter 1-Yankauer Tip Soft Goods/Drapes/Gowns 6-Cotton Towels 20-Lap Sponges Wraps 2-50" x 70" BTC Fold	4 1-Placenta Basin/Lid 1-Medium Cup 1-Umbilical Cord Clamp 1-Purple Blood Tube 1-Suction Tubing 2-#10 Surgical Blades 1-Baby Blanket 1- Suture Bag
52000016	C-Section Deluxe Kit Plastics/Surgical 2-7 Liter Basins 1-Heart Bowl 1-Bulb Syringe 1-Red Blood Tube 1-Needle Counter 1-Yankauer Tip Soft Goods/Drapes/Gowns 6-Cotton Towels 20-Lap Sponges 1-C-Section Drape 3-X-Large Gowns Wraps 1-18"x18" Wrap	4 1-Placenta Basin/Lid 1-Medium Cup 1-Umbilical Cord Clamp 1-Purple Blood Tube 1-Suction Tubing 2-#10 Surgical Blades 1-Baby Blanket 1-Non-Woven Towel 1-Mayo Stand Cover 1-Suture Bag 1-54" x 90" BTC Fold
31354429	C-Section Plus Kit Plastics/Surgical 2-7 Liter Basin 1-Heart Bowls 1-Bulb Syringe 1-Red Blood Tube 1-Suction Tubing 2- EZ Handle Light Handles 2-#10 Surgical Blades Soft Goods/Drapes/Gowns 1-Maternity Pad 1-Baby Blanket 20-Lap Sponges Wraps 1-50" x 70" Valentine Fold	6 1-Placenta Basin/Lid 1-Medium Cup 1-Umbilical Cord Clamp 1-Needle Counter 1-Yankauer Tip 1- ES Pencil 1-#20 Surgical Blade 6-Cotton Towel 1-Baby Beanie 1-Suture Bag

Vaginal & C-Section Kits



C-Section Delivery Kit with ES Pencil

Code	Description	Ship Case
520000076	C-Section Deluxe Kit, with ES pencil	4
	Plastics/Surgical	
	2-7 Liter Basin,	1-Placenta Basin/Lid
	1-Heart Bowls	1-Medium Cup
	2-Bulb Syringe	1-Umbilical Cord Clamp
	1-Needle Counter	1-Suction Tubing
	1-Yankauer Tip	2-LiteGlove™ light handle covers
	1-ES Pencil	2-#10 Surgical Blades
	Soft Goods/Drapes/Gowns	
	6-Cotton Towel	1-Non-Woven Towel
	1-Baby Blanket	10-Lap Sponges
	1-Suture Bag	1-C-Section Drape
	1-Mayo Stand Cover	3-X-Large Gown
	1-Large Gown	
	Wraps	
	1-18" x 18" Wrap	
	1-54" x 90" BTC Fold	

Fetal Monitoring Paper

Covidien offers a full line of recording charts for antepartum and intrapartum monitoring. We manufacture charts to fit monitors from major manufacturers including GE Medical, Philips (Agilent), Spacelabs, Oxford and others.



Fetal Monitoring Paper

- Easy-to-read grids with extended timelines at one-minute intervals
- Pages numbered consecutively for easy referencing and documentation
- End of pack warning alerts clinicians when 30 minutes of recording time remains

Code	Compatible to	Grid	Fetal Monitor	Ship Case
To Fit GE-Medical				
30597226	4305AAO	Red	Corometrics 115/116/118	40 packs
	4305CAO		120/126/128/129/150/151	
30767589	4305BAO	Red	Corometrics 115/116/118	40 packs
	4305DAO		20/126/128/129/150/151	
30589132	4483AAO	Red	Corometrics 145	6 rolls
To Fit Philips				
10005156	9270-0484	Green	Series 8025,8030, 8040	40 packs
10643709	9270-0485	Red	Series 8025, 8030, 8040	40 packs
30748696	M1910A	Orange	Series 50, M1350A/B	40 packs
30791761	M1913A	Green	Series 50, M1350A/B	40 packs
30930088	M1914A	Orange	Series 50, M1350A/B	60 packs
	(M1910A + labels)			
To Fit Other Monitors				
01082320		Red	AMS AM66, IM76 Air-Shields Fetascan 1200,1400 Litton/PPG FM 660, FM 670, FM 671	40 packs
31329702		Orange	AMS IM77	40 packs

Label Sets

Code	Description	Ship Case
30926524	Chart Wrapper-100 per pack	1 pack
30926516	Mother/Baby labels-200 of each per pack	1 pack
30926490	Chart Wrapper and Mother/Baby Labels (60 of each per pack)	1 pack

Circumcision Products



Tray with clamp



Single Sterile Clamp



Tray with straps, no clamp



Straps

Devon™ Circumcision Trays

- LATEX FREE
- Standardized procedure set-up
- Save time and labor
- Disposable; sterile
- Trays contain items most commonly used in one convenient, sterile tray
- Available with Povidone-Iodine (PVP) Swabsticks


Most trays include the following basic components:

Gauze sponges, 4" x 4"	Clamps, Mosquito, curved
Applicators 3"	Clamps, Mosquito, straight
Disposable scalpel #10	Absorbent drape 19" x 23"
Safety pins	Fenestrated drape
Eye probe 6"	Cotton balls
Scissors, sharp/blunt	3-compartment tray
Scissors, Iris straight	Wrap - 24" x 24"
Adson forceps with teeth	Retractable disposable scalpel
Circumcision straps	Povidone - Iodine (PVP) Swabsticks

Code	Description	Ship Case
Basic Circumcision Tray		
52000005	7380 Circumcision Tray - Basic (without straps)	20/case
52000006	7384 Circumcision Tray - Plus (with straps)	20/case
56455	Circumcision Tray with PVP- Basic (without straps)	16/case
56456	Circumcision Tray with PVP- Plus (with straps)	16/case
Deluxe Circumcision Tray with Clamps		
52000031	7381 clamps, 1.1 (without straps)	20/case
52000033	7383 clamps, 1.3 (without straps)	20/case
52000081	7391 clamps, 1.1 (with straps and retractable scalpel)	20/case
52000083	7393 clamps, 1.3 (with straps and retractable scalpel)	20/case
56408	7390 clamps, 1.45 (with straps)	20/case
Single Sterile Clamps		
56421	1.1 single sterile clamps	10/carton
56423	1.3 single sterile clamps	10/carton
56242	1.45 single sterile clamps	10/carton
Baby Board Straps - sold separately		
31309480	Four 10" foam straps per set	25 sets
31368858	Four 10" straps per set (made of ester/tricot)	25 sets

Made-to order Circumcision Trays are also available to meet your specific needs.

Argyle™ Umbilical Vessel Catheters



Single-Lumen



Double-Lumen



Triple-Lumen

Argyle™ Umbilical Vessel Catheters (UVCs)

- LATEX FREE
- DEHP FREE
- Single, dual and triple lumen catheter configurations
- Polyurethane formulation
- Clear, lipid-resistant Luer Lock hubs
- Completely encapsulated X-ray opaque lines
- Numerical depth markings every 1cm
- Smooth, rounded atraumatic tip
- Color-coded I.V. extensions on triple and dual lumen catheters
- Sterile, multiple sizes and lengths

Reorder Number	French Size	Primary/Secondary/Tertiary Lumen Size	Length	Ship Case
Polyurethane Single-Lumen UVC with Luer Lock Hubs				
8888160325	2.5		12"	10
8888160333	3.5		15"	10
8888160341	5.0		15"	10
Polyurethane Dual-Lumen UVC with Luer Lock Hubs				
8888160531	3.5	20ga / 23ga	15"	5
8888160556	5	18ga / 21ga	15"	5
Polyurethane Triple-Lumen UVC with Luer Lock Hubs				
8888160648	5	18ga / 21ga / 21ga	15"	5
8888160663	8	15ga / 18ga / 18ga	15"	5

Umbilical Catheter Anchor



Umbilical Catheter Anchor

- Provides a fast, convenient, secure way to anchor UVCs
- Forms a strain-relief bridge that holds catheter in place
- Replaces time-consuming, messy and complicated tape bridges
- With its two adhesive arms, the umbilical catheter anchor can hold two separate lines securely in place

Code	Description	Ship Case
MI30000	Umbilical catheter anchor	50

Argyle™ Umbilical Vessel Catheter Insertion Trays



Argyle UVC Insertion Tray

Contents-LATEX FREE

PREP TRAY (top section)

- 1 - Umbilical vessel catheter (optional) in pre-packaged sterile catheter system
- 1 - Pack povidone iodine swab sticks
- 1 - Fenestrated drape
- 2 - Towels
- 6 - Gauze pads 2" x 2"

PROCEDURE TRAY (bottom section)

- 1 - CSR wrap
- 2 - Filter aspiration needles
- 1 - 3-way stopcock
- 1 - Stainless steel scissors
- 1 - Stainless steel needle holder, 5"
- 1 - Stainless steel eye dressing forceps, half-curved
- 1 - Stainless steel eye dressing forceps, fully curved
- 2 - Stainless steel hemostats, curved
- 1 - Stainless steel hemostat, straight
- 1 - Blade scalpel #11
- 1 - Stainless steel Adson tissue forceps
- 1 - Stainless steel vessel probe dilator
- 1 - Monoject™ 1cc syringe
- 2 - Monoject™ 6cc syringes
- 1 - Monoject™ 3cc syringe
- 1 - 3.0 Silk suture with curved needle
- 1 - Umbilical tape 15"
- 1 - Syringe tip cap
- 6 - Gauze pads 4" x 4"
- 1 - Monoject™ 20ga blunt needle**
- 1 - Monoject™ 18ga blunt needle**
- 1 - Safety scalpel #11***
- 2 - Monoject™ 20ga x 1.5 safety needle***

**Only included in trays without catheter - 8888160408 and 160409

***Only included in tray - 43201

Code	Description	Ship Case
8888160119	3.5 and 5 fr. polyurethane UVC with luer lock hub	5
8888160424	3.5 fr. polyurethane UVC	5
8888160432	5 fr. polyurethane UVC	5
8888160838	3.5 fr. dual lumen polyurethane UVC	5
8888160853	5 fr. dual lumen polyurethane UVC	5
8888160408	UVC tray without catheter and with fine tip forceps	5
160409	UVC tray without catheter and with smooth tip forceps	5
43201	UVC tray without catheter and with safety components	5

Argyle™ Neonatal-Pediatric PICC



1.9 Fr Double Lumen PICC

Argyle™ Neonatal-Pediatric Peripherally Inserted Central Catheter (PICC)

- LATEX FREE
- DEHP FREE
- Dual and Single Lumen Products
- Polyurethane material for maximum strength and flow rates
- Available with procedure trays
- Radiopaque material for verification of proper placement
- Easy to read graduation marks
- Soft and smooth butterfly stabilizing wing
- Clear pigtails for visualization of fluid
- 30cm catheter length
- 0.6mm outer diameter



Procedure Tray

The Dual Lumen PICC has all of the above features and benefits, PLUS

- Independent lumens throughout the entire catheter to prevent mixing of incompatible infusates
- Color coded hubs to differentiate primary and secondary lumens
- Staggered pigtails for ease of securement and maintenance



1.9 Fr Single Lumen PICC

Code	Description	Ship Case
Single Lumen Products		
43303	1.9 Fr Single Lumen PICC	10
43309	1.9 Fr Single Lumen PICC Procedure Tray	5
Dual Lumen Products		
43304	1.9 Fr Dual Lumen PICC	10
43311	1.9 Fr Dual Lumen PICC Procedure Tray	5
Procedure Tray Only		
43308	PICC Procedure Tray, No Catheter	5
Introducer		
43320	Peel Away Cannula Introducer	5



Peel Away Cannula Introducer

Product sold in case quantities only.

NICU/PICU Procedure Tray



NICU/PICU Procedure Tray

- LATEX FREE
- All of the necessary components for common procedures such as: chest tube insertion, UVC insertion, cut down and many more

Tray components:

Forceps, Adson w/teeth	Scissors, iris straight
Needle holder, Webster	Clamp, Mosquito curved
Mosquito straight	Syringe, 1ml luer-lock tip
Probe with eye, 6"	Eye dressing forceps, semi-curved
Eye dressing forceps curved 4"	3-0 silk suture needle
Stopcock	60 cm ruler
Half curve smooth tip forceps	Full curve smooth tip forceps
Scalpel, handle w/#11 blade	Fenestrated drape 18"x 26"
Syringe 3cc luer lock	Syringe, 5cc luer lock
Gauze sponge, 2'x 2" 12-ply	Towel
Sponge, 4" x 4"	Tray, rectangular
Medium cup , 60cc	Wrap, 24" x 24"

Code	Description	Ship Case
50007250	7250 NICU/PICC Procedure Tray	18

Argyle™ Replogle Suction Catheters



Argyle™ Replogle Suction Catheters

- LATEX FREE
- DEHP FREE
- Clear dual-lumen tube
- X-ray opaque Sentinel Line™ tubing with Sentinel Eye™ catheter
- Integral tapered connector can be cut to form funnel
- Sterile, three smooth eyes

Code	French Size	Length	Ship Case
8888256545	6	24"	10
8888256529	8	24"	10
8888256503	10	24"	10

Neonatal Accessories



Cord Clamp Clipper

Cord Clamp Clipper

- LATEX FREE
- Cuts hinge to remove umbilical cord clamp from infant
- Compatible with most clamps
- Disposable, non-sterile
- Individually packaged

Code	Description	Ship Case
52000008	21333 Cord Clamp Clipper	60/carton



CPAP
Nasal Cannula

Argyle™ CPAP Nasal Cannula

CPAP

- LATEX FREE
- All silicone material
- Carefully angled nares
- Color coded, 15mm connectors

CPAP Kit (with cap)

- Includes CPAP nasal cannula
- Cushioned fabric
- Hook & loop straps

Code	Description	Size	Ship Case
8888162024 (CPAP)	For babies under 1000 grams	X-Small	10/case
8888162016 (CPAP)	For babies 1000-1500 grams	Small	10/case
8888162008 (CPAP)	For babies over 1500 grams	Large	10/case
8888163055 (Kit)	For babies under 1000 grams	X-Small	10/case
8888163030 (Kit)	For babies 1000-1500 grams	Small	10/case
8888163014 (Kit)	For babies over 1500 grams	Large	10/case



CPAP Kit

Baby Board I.V. Limb Support Boards

- LATEX FREE
- Pre-shaped, bendable but firm board supports and cushions elbow, leg or heel
- Provides maximum immobilization of infant's limb when securing I.V. lines
- Includes two optional soft hook & loop straps; 3 sizes



Code	Description	Ship Case
MI00880	I.V. Limb Support Board 1" x 3" (preemie)	50/carton
MI00775	I.V. Limb Support Board 2" x 4" (small)	50/carton
MI00850	I.V. Limb Support Board 2" x 5" (large)	50/carton

Neonatal Accessories



Heel Warmers

Infant Heel Warmer

- LATEX FREE
- Provides a fast, safe and economical means to apply just the right amount of heat to infant's heel for blood sampling
- Simply squeeze the pouch for instant activation
- Heel warmer quickly reaches the proper temperature and maintains this temperature for several minutes
- MH00002T has an adhesive tab for convenient attachment to the infant

Code	Description	Ship Case
MH00002T	Heel Warmer 3-1/2 x 5" (attachment pad)	100/case
MH00002N	Heel Warmer 3-1/2 x 5" (no attachment pad)	100/case



Baby Beanies

Baby Beanies

- LATEX FREE
- Baby beanies keep the head of the preemie or newborn warm
- Soft material; pink and blue striped
- Select from single ply baby beanie or 2-ply for added thickness and warmth

Code	Description	Ship Case
30953130	BBN4813 Baby Beanies (single-ply)	50/carton
52000010	BBN4814 Baby Beanies (2 ply)	50/carton



Temperature Probe Covers

Temperature Probe Covers

- LATEX FREE
- Gentle, non-irritating clear adhesive or hydrogel securely anchors temperature probe to infant's skin
- May be removed and repositioned
- Reflective mylar acts as a barrier against ambient and radiant heat
- Individually packaged in resealable packs

Code	Description	Ship Case
Hydrogel		
31424784	Arbo™ 2 (small)	200/case
31424792	Arbo™ 3 (large)	200/case
MI30002	CB-810 (3/4" x 11/4") Cover Babies	800/case
Adhesive		
43101	Arbo™ 2A (small)	200/case
43102	Arbo™ 3A (large)	200/case

Neonatal Accessories



Baby Shades

Baby Shades Phototherapy Eye Protectors

- LATEX FREE
- Used in the treatment of infantile jaundice
- Soft fabric shades gently adhere to infant's temples with our biocompatible hydrogel adhesive or optional strap

Code	Description	Ship Case
MI00575	6002 1-1/8" x 4-1/2" (2.8cm x 10.8cm) (small) for infants less than 1500 grams	50/carton
MI00625	6003 1-1/2" x 5-1/2" (3.4cm x 13.9cm) (large) for infants more than 1500 grams	50/carton



Baby Tape

Baby Adhesive Tape

- LATEX FREE
- Clear adhesive hydrogel formula gently yet firmly secures catheters, I.V. lines and lead wires to infant's skin
- Tapes may be repositioned and reapplied without irritating or damaging skin
- Available in individually packaged tape strips

Code	Description	Ship Case
MI00677	6100 Baby 3/4" x 2" Tape strips (50 strips/pack)	10 packs/case

Neonatal Pre-Wired Electrodes

Neonatal Pre-Wired Electrodes

- LATEX FREE
- DEHP FREE
- MERCURY FREE
- Patented conductive adhesive hydrogel provides firm adhesion while minimizing irritation to delicate newborn skin
- Color coded leadwires with industry standard .060" (1.5mm) safety socket connector
- Radiolucent leadwires available to minimize the need to remove electrodes during X-Ray
- Variety of shapes, sizes and substrates to address most neonatal patients



1041PTS
Radiolucent Wires

Pre-Wired Clear Tape Series

- Clear tape substrate is comfortable, allows visibility to skin, does not absorb fluids, and slows hydrogel moisture loss under warming lights
- 24" lead wire length
- 5 lead wire length
- Two shapes: Round 1" diameter, Square 1" x 1-1/4"



1042PTS
Radiolucent Wires

Code	Description	Quantity/Pack	Ship Case
31424780	1041PTS	3	300
31424785	1042PTS	3	300
31424786	1042PTS54	3	300

Neonatal Pre-Wired Electrodes



4103
Non-radiolucent



4105
Non-radiolucent

Pre-Wired Foam 4100 Series

- Foam substrate repels fluid
- 18" leadwire length
- Two Shapes: Round 7/8" diameter, Oval 7/8" x 1-3/16"

Code	Description	Quantity/Carton	Ship Case
Foam (3 per pack)			
30710001	4103	75	300
31086575	4105	75	300



1050NPSM
Non-radiolucent
wire



1051NPSM
Radiolucent wire

Pre-Wired 1050 Cloth Series

- Soft Spunlace cloth substrate
- 24" leadwire length
- Available with or without radiolucent wires
- Two Shapes: 1" diameter, 1-3/16" diameter

Code	Description	Ship Case
Cloth (3 per pack)		
31424768	1050NPSM	300
31424743	1051NPSM radiolucent wire	300
31424735	1052NPSM	300
31424719	1053NPSM radiolucent wire	300



1052NPSM
Non-radiolucent
wire



1051NPSM
Radiolucent wire

Limb Band Specialty Pre-Wired Electrodes

- Designed for wrapping around neonate limbs
- Reduces pressure on and increases access to the infant's torso
- Constructed of soft white spunlace cloth
- 18" lead wires with 4-1/2" limb bands

Code	Description	Quantity/Pack	Ship Case
EP30016	L36-LB3	3	300



EP30018 Limb Band

INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
01082320	Fetal Monitoring Paper	11	31479754	Reusable Cables, Philips Medical System FCB 100	5
10005156	Fetal Monitoring Paper To Fit Philips	11	31479812	Reusable Cables, GE-Medical FCB 400	5
10643709	Fetal Monitoring Paper To Fit Philips	11	31479846	Reusable Cables, Philips Medical System FCB 300	5
160409	Argyle™ UVC Insertion Trays w/o catheter, with smooth tip forceps	14	31481396	FSE Attachment Pads FCB 2000 (small 1-7/8")	5
30589132	Fetal Monitoring Paper To Fit GE-Medical	11	40000007	Abdominal Belts Buttonhole Style - Knit Elastic 2-3/8" x 48"	6
30597226	Fetal Monitoring Paper To Fit GE-Medical	11	40000008	Abdominal Belt Tab Style Knit Elastic 1-1/2" x 36"	6
30710001	Electrodes FOAM - 4103	21	40000009	Abdominal Belt Tab Style Knit Elastic 1-1/2" x 36"	6
30748696	Fetal Monitoring Paper To Fit Philips	11	43101	Temperature Probe Covers Arbo™ 2A	18
30767589	Fetal Monitoring Paper To Fit GE-Medical	11	43102	Temperature Probe Covers Arbo™ 2A	18
30791761	Fetal Monitoring Paper To Fit Philips	11	43201	Argyle™ UVC Insertion Trays w/o catheter, with safety components	14
30817608	Abdominal Belts Buttonhole Woven Elastic 2-3/8" x 52"	6	43303	1.9 Fr Single Lumen PICC	15
30926490	Chart Wrapper	11	43304	1.9 Fr Dual Lumen PICC	15
30926516	Mother/Baby labels	11	43308	PICC Procedure Tray, No Catheter	15
30926524	Chart Wrapper	11	43309	1.9 Fr Single Lumen PICC Procedure Tray	15
30930088	Fetal Monitoring Paper To Fit Philips	11	43311	1.9 Fr Dual Lumen PICC Procedure Tray	15
30950169	Abdominal Belt Tab Style Woven Elastic 1-1/2" x 40"	6	43320	Peel Away Cannula Introducer	15
30950292	Abdominal Belts Buttonhole Style - Knit Elastic 2-3/8" x 48"	6	50000095	Attachment Pads FCB 2001 (large 3-3/4" x 1 7/8")	5
30950300	Abdominal Belts Buttonhole Style - Knit Elastic 2-3/8" x 48"	6	50000106	Attachment Pads Spacelabs, FCB 600	5
30951696	C-Section Straps 3" x 84	7	50000107	Reusable Cables, GE-Medical, FCB 700	5
30952090	Abdominal Belt Bobbins	6	50000331	Reusable Cables, Philips Medical System™, FCB 308	5
30953130	BBN4813 Baby Beanies (single-ply)	18	50007250	7250 NICU/PICU Procedure Tray	16
30955689	Abdominal Belt Tab Style Open Weave Mesh 2" x 52"	6	52000005	7380 Circumcision Tray	12
30966348	Abdominal Belts Buttonhole Style - Knit Elastic 2-3/8" x 150' roll	6	52000006	7384 Circumcision Tray	12
31086575	Electrodes FOAM - 4105	21	52000008	21333 Cord Clamp Clipper	17
31309480	Baby Board Straps Four 10" foam straps per set	12	52000010	BBN4814 Baby Beanies (2 ply)	18
31325965	Abdominal Belt Tab Style Foam 1-1/2" x 48"	6	52000013	Vaginal Plus Delivery Kit	8
31325981	Abdominal Belt Tab Style Foam 1-1/2" x 60"	6	52000014	Vaginal Deluxe Delivery Kit, with drape & gown	8
31329702	Fetal Monitoring Paper	11	52000015	C-Section Plus Kit	9
31354429	C-Section Plus Kit	9	52000016	C-Section Deluxe Kit, with drape & gowns	9
31363362	Binders, Small, beige, Roll, 33' in length, must be cut to length	7	52000031	Deluxe Circumcision Tray with 7381 clamps	12
31364360	Binders, Medium, beige, Roll, 33' in length, must be cut to length	7	52000033	Deluxe Circumcision Tray with 7383 clamps	12
31364378	Binders, Large, beige, Roll, 33' in length, must be cut to length	7	52000076	C-Section Deluxe Kit, with ES pencil	10
31364394	Binders, Small, beige, Pre-Cut, 16" in length	7	52000081	Deluxe Circumcision Tray with 7391 clamps	12
31364410	Binders, Medium, beige, Pre-Cut, 16" in length	7	52000083	Deluxe Circumcision Tray with 7393 clamps	12
31364428	Binders Large, beige, Pre-Cut 16" in length	7	530975232	Abdominal Belt Tab Style Woven Elastic 1-1/2" x 40"	6
31364436	Binders X-Large, beige, Pre-Cut 16" in length	7	56101	Abdominal Belts Buttonhole Style - Knit Elastic -2-3/8" x 48"	6
31368858	Baby Board Straps Four 10" straps per set	12	56102	Abdominal Belt Tab Style Knit Elastic 1-1/2" x 42"	6
31410270	Abdominal Belts Buttonhole Style - Knit Elastic 2-3/8" x 48"	6	56242	1.45 single sterile clamps	12
31424719	Electrodes Cloth - 1053NPSM-Radiolucent	21	56300	Kendall™ IUPC	4
31424735	Electrodes Cloth - 1052NPSM	21	56311	REUSABLE CABLES For most Philips Medical™ Monitors	4
31424743	Electrodes Cloth - 1051NPSM-Radiolucent	21	56312	REUSABLE CABLES For Philips Medical Avalon™ Monitors	4
31424768	Electrodes Cloth - 1050NPSM	21	56321	REUSABLE CABLES For GE Medical Corometrics™ Monitors	4
31424780	Electrodes Clear Tape - 1041PTS Radiolucent	20	56330	REUSABLE CABLES For Spacelabs Monitors™	4
31424784	Temperature Probe Covers Arbo™ 2 (small)	18	56408	Deluxe Circumcision Tray with 7390 clamps	12
31424785	Electrodes Clear Tape - 1041PTS Radiolucent	20	56421	1.1 single sterile clamps	12
31424792	Temperature Probe Covers Arbo™ 3 (large)	18	56423	1.3 single sterile clamps	12
31479549	Kendall™ FSE 2100	5	56455	Circumcision Tray with PVP	12

INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
56456	Circumcision Tray with PVP	12	8888163014	Argyle™ Kit Nasal Cannula for babies over 1500 grams	17
8888160119	Argyle™ UVC Insertion Trays 3.5 and 5 fr. polyurethane	14	8888163030	Argyle™ Kit Nasal Cannula for babies 1000-1500 grams	17
8888160325	Polyurethane Single-Lumen UVC with Luer Lock Hubs	13	8888163055	Argyle™ Kit Nasal Cannula for babies under 1000 grams	17
8888160333	Polyurethane Single-Lumen UVC with Luer Lock Hubs	13	8888256503	Argyle™ Replogle Suction Catheters 10 FR	16
8888160341	Polyurethane Single-Lumen UVC with Luer Lock Hubs	13	8888256529	Argyle™ Replogle Suction Catheters 8 FR	16
8888160408	Argyle™ UVC Insertion Trays w/o catheter, with fine tip forceps	14	8888256545	Argyle™ Replogle Suction Catheters 6 FR	16
8888160424	Argyle™ UVC Insertion Trays 3.5 fr. polyurethane UVC	14	EP30016	Limb Band L36-LB	21
8888160432	Argyle™ UVC Insertion Trays 5 fr. polyurethane UVC	14	MH00002N	Heel Warmer 3-1/2 x 5" (no attachment pad)	18
8888160531	Polyurethane Dual-Lumen UVC with Luer Lock Hubs	13	MH00002T	Heel Warmer 3-1/2 x 5" (attachment pad)	18
8888160556	Polyurethane Dual-Lumen UVC with Luer Lock Hubs	13	MI00575	Baby Shades Phototherapy Eye Protectors - small	19
8888160648	Umbilical Vessel Catheters ,Polyurethane Triple-Lumen	13	MI00677	Baby Tape Adhesive Tape 6100 Baby 3/4" x 2" Tape strips	19
8888160663	Umbilical Vessel Catheters Polyurethane Triple-Lumen	13	MI00625	Baby Shades Phototherapy Eye Protectors - large	19
8888160838	Argyle™ UVC Insertion Trays 3.5 fr. dual lumen polyurethane	14	MI00775	I.V. Limb Support Board 2" x 4" (small)	17
8888160853	Argyle™ UVC Insertion Trays 5 fr. dual lumen polyurethane UVC	14	MI00850	I.V. Limb Support Board 2" x 5" (large)	17
8888162008	Argyle™ CPAP Nasal Cannula for babies over 1500 grams	17	MI00880	I.V. Limb Support Board 1" x 3" (preemie)	17
8888162016	Argyle™ CPAP Nasal Cannula for babies 1000-1500 grams	17	MI30000	Umbilical catheter anchor	13
8888162024	Argyle™ CPAP Nasal Cannula for babies under 1000 grams	17	MI30002	Temperature Probe Covers CB-810 (3/4" x 1-1/4") Cover Babies	18

TM * Avalon is a trademark of Koninklijke Philips Electronics, NV;
Conmetrics is a trademark of GE Medical Systems Info Tech., Inc.
COVIDIEN, COVIDIEN with logo and TM marked brands are trademarks of Covidien AG or its affiliate.
© 2009 Covidien AG. All rights reserved.

H-62000NEOB XM 1/09

