


# Sterile Trocar Catheter



## Trocar Catheter with SENTINEL LINE (X-Ray Opaque Line)

- Thermosensitive, non-pyrogenic, and inert in prolonged contact with body tissue.
- Sterile.
- Individually packaged.
- Medical grade catheter over aluminum trocar.

Code	French Size	O.D.	Length	Qty/Case
8888560805	8 Fr	2.7mm	9" (22.9cm)	1/Rigidpak 10/case
8888561019	10 Fr.	3.3mm	9" (22.9cm)	1/Rigidpak 10/case
8888561027	12 Fr.	4.0mm	9" (22.9cm)	1/Rigidpak 10/case
8888561035	16 Fr.	5.3mm	10" (25.4cm)	1/Rigidpak 10/case
8888561043	20 Fr.	6.7mm	16" (41cm)	1/Rigidpak 10/case
8888561050	24 Fr.	8.0mm	16" (41cm)	1/Rigidpak 10/case
8888561068	28 Fr.	9.3mm	16" (41cm)	1/Rigidpak 10/case
8888561076	32 Fr.	10.7mm	16" (41cm)	1/Rigidpak 10/case



## ARGYLE TURKEL Safety Thoracentesis System

- Spring-loaded blunt safety cannula extends beyond needle tip.
- One-way valve in needle housing.
- Color change indicator to identify blunt cannula position.
- Ball valve in catheter housing.
- Polyurethane catheter-over-needle design.
- 1cm depth marks on catheter shaft.
- Multiple side holes on catheter shaft.
- Sterile.

Code	Catheter Size	Catheter O.D.	Catheter Length	Qty/Case
566034	8 Fr.	2.7mm	3½" (9cm)	20