

**Cardio Thoracic Products**

Chest Drainage Products and Accessories

Thoracic Catheters

Specialty Products



## Table of Contents

### Cardio Thoracic Products

#### Chest Drainage Units

Argyle™ Altitude™ Dry Suction Chest Drainage Unit 1	
Argyle™ Aqua-Seal™ Chest Drainage Unit.....	1-3
Argyle™ Sentinel Seal™ Chest Drainage Unit .....	3-4
Argyle™ Thora-Seal™ Chest Drainage Unit.....	5-7
Argyle™ Double Seal™ Chest Drainage Unit.....	7
Clinical Products .....	7

#### Chest Drainage Accessories

Argyle™ Autotransfusion Accessory .....	8
Argyle™ Emerson™ Pump Adapter .....	8
Argyle™ Thora-Seal™ Emerson Pump/ Thoracic Regulator Adapter .....	8
Argyle™ Tubing Extension Sets .....	9

#### Thoracic Catheters

Argyle™ Thoracic Catheter .....	10
Argyle™ Ferguson™ Thoracic Catheter.....	10
Argyle™ Silicone Thoracic Catheter .....	11
Argyle™ Soft PVC Catheter .....	11
Argyle™ Silicone with Pre-packaged Catheter Connector System .....	12
Argyle™ Silicone Mediastinal Catheters .....	12
Argyle™ Trocar Catheters .....	13
Argyle™ Trocar Catheter Kits and Thoracostomy Trays .....	13

#### Specialty Products

Argyle™ Left Ventricular Sump Vent Catheters ....	14
Argyle™ Thi Aortic Perfusion Cannulas .....	14
Argyle™ Parallel “Y” Connector .....	15
Argyle™ Carotid Artery Shunt Kits .....	15
Argyle™ Vascular Tourniquet Kit .....	16
Argyle™ Sentinel Loops™ Silicone Vascular Tapes..	16
Argyle™ Penrose™ Tubing .....	17
Argyle™ Saratoga™ Sump Drains .....	17

Index .....	18-19
-------------	-------

## Chest Drainage Units

\*Autotransfusion Accessories are featured on page 8—9



### Argyle™ Altitude™ Dry Suction Autotransfusion Chest Drainage Unit (Autotransfusion-Ready)\*

- Precision suction regulator
- Source suction indicator
- Patient vacuum indicator
- Needle-free ports for collection chamber and water seal chamber
- Latex-free patient tube and components
- Automatic positive pressure relief valve
- Manual negative pressure relief valve
- Large fill spout
- Blue colored water seal
- Pediatric and adult collection chamber graduations identify drainage volume
- Easy-to-read graphics
- Set-up instructions printed on unit
- Wide, rotational floorstand and hangers
- Color-coded locking connectors
- Patient-tube slide clamp
- Clear visibility for vacuum set point
- Simple Autotransfusion hook and Velcro™\* attachment
- Sterile

CODE	DESCRIPTION	SHIP CASE
8888571370	Argyle Altitude Dry Suction Chest Drainage Unit	5



### Argyle™ Aqua-Seal™ All Purpose Chest Drainage Unit (Autotransfusion-Ready or Post-Evac)\*

- Three convenient post-operative Autotransfusion options including continuous, post-collection batch evacuation and bypass collection
- Out-the-top Autotransfusion Infusion Line for clean, convenient blood retrieval
- Exclusive Redi-Prime Valve for needle-free access to Autotransfusion Infusion Line
- Needle-free ports for collection chamber and water seal chamber
- Latex-free patient tube and components
- Clot collection chamber
- Manual and automatic negative pressure relief valves
- Blue-colored water seal
- Flow control valve for convenient suction flow control
- Wide, rotational floorstand and hangers
- Color-coded locking connectors
- Patient-tube slide clamp
- Simple Autotransfusion hook and Velcro™\* attachment
- Sterile

CODE	DESCRIPTION	SHIP CASE
8888571364	Aqua-Seal All Purpose Autotransfusion Chest Drainage Unit	5

## Chest Drainage Units

\*Autotransfusion Accessories are featured on page 8—9



### Argyle™ Aqua-Seal™ Basic Chest Drainage Unit

- Compact, clear, 3-chamber chest drainage unit
- Needle-free ports for collection chamber and water seal chamber
- Latex-free patient tube and components
- Automatic positive pressure relief valve
- Manual and automatic negative pressure relief valves
- Blue colored water seal
- Pediatric and adult collection chamber graduations identify drainage volume
- Easy-to-read graphics
- Set-up instructions printed on unit
- Wide, rotational floorstand and hangers
- Sterile

CODE	DESCRIPTION	SHIP CASE
8888571299	Aqua-Seal Basic Chest Drainage Unit	5



### Argyle™ Aqua-Seal™ Autotransfusion Chest Drainage Unit (Autotransfusion-Ready)\*

- Compact, clear, 3-chamber Autotransfusion-Ready Chest Drainage Unit
- Compatible with Argyle Autotransfusion Accessory Unit
- Needle-free ports for collection chamber and water seal chamber
- Latex-free patient tube and components
- Automatic positive pressure relief valve
- Manual and automatic negative pressure relief valves
- Blue colored water seal
- Pediatric and adult collection chamber graduations identify drainage volume
- Easy-to-read graphics
- Set-up instructions printed on unit
- Wide, rotational floorstand and hangers
- Color-coded locking connectors
- Patient-tube slide clamp
- Simple Autotransfusion hook and Velcro™ attachment
- Sterile

CODE	DESCRIPTION	SHIP CASE
8888571315	Aqua-Seal Autotransfusion Chest Drainage Unit	5

## Chest Drainage Units

\*Autotransfusion Accessories are featured on page 8—9



### Argyle™ Aqua-Seal™ Dual Drain Chest Drainage Unit

- Compact, clear, 3-chamber chest drainage unit
- Two independent collection chambers and two collection tubes
- Needle-free ports for collection chamber and water seal chamber
- Latex-free patient tube and components
- Automatic positive pressure relief valve
- Manual and automatic negative pressure relief valves
- Blue colored water seal
- Pediatric and adult collection chamber graduations identify drainage volume
- Easy-to-read graphics
- Set-up instructions printed on unit
- Wide, rotational floorstand and hangers
- Sterile

CODE	DESCRIPTION	SHIP CASE
8888571406	Aqua-Seal Dual Drain Chest Drainage Unit	5



### Argyle™ Sentinel Seal™ Dual Drain Chest Drainage Unit

- Two independent collection chambers and two collection tubes
- Compact, clear, 4-chamber chest drainage unit
- Needle-free ports for collection chamber and water seal chamber
- Latex-free patient tube and components
- Precision suction regulator
- Blue-colored Patient Assessment Manometer monitors actual negative intrathoracic pressure
- Automatic positive pressure relief valve
- Manual negative pressure relief valve
- Pediatric and adult collection chamber graduations identify drainage volume
- Large water seal fill spout
- Easy-to-read graphics
- Set-up instructions printed on unit
- Wide, rotational floorstand and hangers
- Sterile

CODE	DESCRIPTION	SHIP CASE
8888571513	Sentinel Seal Dual Drain Chest Drainage Unit	5

## Chest Drainage Units

\*Autotransfusion Accessories are featured on page 8—9



### Argyle™ Sentinel Seal™ Chest Drainage Unit

- Compact, clear, 4-chamber chest drainage unit
- Needle-free ports for collection chamber and water seal chamber
- Latex-free patient tube and components
- Precision suction regulator
- Blue-colored Patient Assessment Manometer monitors actual negative intrathoracic pressure
- Automatic positive pressure relief valve
- Manual negative pressure relief valve
- Pediatric and adult collection chamber graduations identify drainage volume
- Large water seal fill spout
- Easy-to-read graphics
- Set-up instructions printed on unit
- Wide, rotational floorstand and hangers
- Sterile

CODE	DESCRIPTION	SHIP CASE
8888571562	Sentinel Seal Chest Drainage Unit	5



### Argyle™ Sentinel Seal™ Autotransfusion Chest Drainage Unit (Autotransfusion-Ready)\*

- Compact, clear, 4-chamber Autotransfusion-Ready Chest Drainage Unit
- Compatible with Argyle™ Autotransfusion Accessory Unit
- Needle-free ports for collection chamber and water seal chamber
- Latex-free patient tube and components
- Precision suction regulator
- Blue-colored Patient Assessment Manometer monitors actual negative intrathoracic pressure
- Automatic positive pressure relief valve
- Manual negative pressure relief valve
- Pediatric and adult collection chamber graduations identify drainage volume
- Large water seal fill spout
- Easy-to-read graphics
- Set-up instructions printed on unit
- Wide, rotational floorstand and hangers
- Color-coded locking connectors
- Patient-tube slide clamp
- Simple Autotransfusion hook and Velcro™ attachment
- Sterile

CODE	DESCRIPTION	SHIP CASE
8888571489	Sentinel Seal Autotransfusion Chest Drainage Unit	5

## Chest Drainage Units

\*Autotransfusion Accessories are featured on page 8—9



### Argyle™ Thora-Seal™ III Chest Drainage Unit

- Compact, clear 3-chamber chest drainage unit
- Removable collection chamber
- Latex-free patient tube and components
- Automatic positive pressure and negative pressure relief valve
- Large fill openings for quick set-up
- Unique baffle system to prevent mixing of fluids
- Easy-to-read graphics
- Set-up instructions printed on unit
- Wide, rotational floorstand and hangers
- Sterile

CODE	DESCRIPTION	SHIP CASE
8884713308	Thora-Seal III Chest Drainage Unit	4



### Argyle™ Thora-Seal™ III Autotransfusion Chest Drainage Unit (Autotransfusion-Ready)\*

- Compact, clear 3-chamber Autotransfusion-Ready Chest Drainage Unit
- Compatible with Argyle™ Autotransfusion Accessory Unit
- Removable collection chamber
- Latex-Free patient tube and components
- Automatic positive pressure and negative pressure relief valve
- Large fill openings for quick set-up
- Easy-to-read graphics
- Set-up instructions printed on unit
- Wide, rotational floorstand and hangers
- Color-coded locking connectors
- Simple Autotransfusion hook and Velcro™\* attachment
- Sterile

CODE	DESCRIPTION	SHIP CASE
8884713176	Thora-Seal III Autotransfusion Chest Drainage Unit	4

## Chest Drainage Units

\*Autotransfusion Accessories are featured on page 8—9



### Argyle™ Thora-Seal™ III Autotransfusion Chest Drainage Unit (Autotransfusion-Ready)\*

- Complete Pre-Attached Autotransfusion System
- Compact, clear 3-chamber Autotransfusion Chest Drainage Unit
- Compatible with Argyle™ Autotransfusion Accessory Unit
- Removable collection chamber
- Latex-Free patient tube and components
- Automatic positive pressure and negative pressure relief valve
- Large fill openings for quick set-up
- Easy-to-read graphics
- Set-up instructions printed on unit
- Wide, rotational floorstand and hangers
- Color-coded locking connectors
- Simple Autotransfusion hook and Velcro™\* attachment
- Sterile

CODE	DESCRIPTION	SHIP CASE
8884713192	Thora-Seal III Autotransfusion Chest Drainage System (Pre-Attached Autotransfusion Accessory Unit)	3



### Argyle™ Thora-Seal™ III Replacement Collection Chamber

- For use with all Thora-Seal III Chest Drainage and Autotransfusion Chest Drainage Units
- Twist-on cap included
- Sterile

CODE	DESCRIPTION	SHIP CASE
8884713900	Replacement Chamber	6



## Chest Drainage Units

\*Autotransfusion Accessories are featured on page 8—9



### Argyle™ Thora-Seal™ I Chest Drainage Unit

- Disposable, single bottle chest drainage system
- Compact, clear, break-resistant plastic
- Can be used with gravity or regulated suction
- 5' clear PVC tubing connects to thoracic catheter
- Self-contained floor stand
- Sterile

CODE	DESCRIPTION	GRADUATIONS	TUBING L D.	SHIP CASE
8884713100	2000cc Plastic Bottle	50cc	9/32"	6



### Argyle™ Double Seal™ Chest Drainage Unit

- Transparent, 4-chamber Chest Drainage Unit
- Exclusive Patient Assessment Chamber monitors actual negative intrathoracic pressure
- Automatic positive pressure relief valve
- Automatic negative pressure relief system
- Blue-colored water seal
- Includes floorstand and hangers
- Sterile

CODE	DESCRIPTION	SHIP CASE
8888571547	Argyle Double Seal Chest Drainage Unit	5



### Clinical Products One and Two Bottle Chest Drainage Systems

- Crystal clear glass bottles
- Accurate graduations
- 5' clear PVC tubing connects to thoracic catheter
- Sterile

CODE	DESCRIPTION	GRADUATIONS	TUBING L D.	SHIP CASE
One Bottle Systems				
8884726100	Standard 2000cc Bottle	50cc	9/32"	6
8884726500	Dual Drainage 2000cc Bottle	50cc	9/32"	6
Two Bottle Systems				
8884726200	Standard 2000cc Bottles	50cc	9/32"	4


## Chest Drainage Accessories



### Argyle™ Autotransfusion Accessory Unit

- Compatible with all Altitude™, Aqua-Seal™, Sentinel Seal,™ and Thora-Seal™ III Autotransfusion-Ready Chest Drainage Units
- Fits easily into standard pressure infusion cuff
- Autotransfusion receptacle for collection and reinfusion of blood
- Automatic blood overflow to chest drainage unit
- Graduations printed on the autotransfusion bag
- Can be used with gravity or continuous reinfusion
- Needle-free injection port for addition of anti-coagulants
- Latex-free tubing and components
- Sterile

CODE	DESCRIPTION	SHIP CASE
8884713184	Argyle Autotransfusion Accessory Unit	5



### Argyle™ Emerson Pump Adapter

- For use with glass bottles and Thora-Seal™ I Chest Drainage Unit
- Clean, ready to use
- Sterile

CODE	DESCRIPTION	SHIP CASE
8884723000	Adapter for Emerson Pump	24



### Argyle™ Thora-Seal™ III Emerson Pump/Thoracic Regulator Adapter

- For use with all Thora-Seal™ III Chest Drainage Units
- Clean, ready to use
- Sterile

CODE	DESCRIPTION	SHIP CASE
8884723500	Argyle Thora-Seal III Emerson Pump/Thoracic Regulator Adapter	24

## Chest Drainage Accessories




### Argyle™ Tubing Extension Sets

- 18" tubing
- Croppable end connectors
- Latex-free patient tubing
- Sterile

CODE	DESCRIPTION	SHIP CASE
8884714200	For use with dual thoracic catheter	24

## Thoracic Catheters



### Argyle™ Thoracic Catheters

- Clear thermosensitive PVC catheter
- Radiopaque Sentinel Line™ tubing with Sentinel Eye™ catheter
- Smooth side eyes and open distal tip
- Numeric depth marks at 2cm increments
- Beveled connector end can be cropped to form integral funnel connector
- Rigid, Escon tube package
- Sterile

CODE	FRENCH SIZE	O.D. MM	SIDE EYES	LENGTH	SHIP CASE
<b>Straight</b>					
8888570507	12	4.0	4	20" (51cm)	10
8888570515	16	5.3	4	20" (51cm)	10
8888570523	20	6.7	5	20" (51cm)	10
8888570531	24	8.0	6	20" (51cm)	10
8888570549	28	9.3	6	20" (51cm)	10
8888570556	32	10.7	6	20" (51cm)	10
8888570564	36	12.0	6	20" (51cm)	10
8888570572	40	13.3	6	20" (51cm)	10
<b>Right Angle</b>					
8888571018	16	5.3	4	20" (51cm)	10
8888571026	20	6.7	5	20" (51cm)	10
8888571034	24	8.0	6	20" (51cm)	10
8888571042	28	9.3	6	20" (51cm)	10
8888571059	32	10.7	6	20" (51cm)	10
8888571067	36	12.0	6	20" (51cm)	10
8888571075	40	13.3	6	20" (51cm)	10



### Argyle™ Ferguson Thoracic Catheters

- Malleable wire insert
- Clear thermosensitive PVC catheter
- Radiopaque Sentinel Line™ tubing with Sentinel Eye™ catheter
- Smooth side eyes and open distal tip
- Numeric depth marks at 2cm increments
- Beveled connector end can be cropped to form integral funnel connector
- Rigid, Escon tube package
- Sterile

CODE	FRENCH SIZE	O.D. MM	SIDE EYES	LENGTH	SHIP CASE
8888572040	28	9.3	5	20" (51cm)	10
8888572057	32	10.7	5	20" (51cm)	10
8888572065	36	12.0	5	20" (51cm)	10

## Thoracic Catheters



### Argyle™ Silicone Thoracic Catheters

- Biocompatible, thrombo-resistant silicone catheter
- Soft, kink-resistant
- Radiopaque Sentinel Line™ tubing with Sentinel Eye™ catheter
- Smooth side eyes and open distal tip
- Beveled end
- Rigid, Escon tube package
- Sterile

CODE	FRENCH SIZE	O.D. MM	SIDE EYES	LENGTH	SHIP CASE
Straight					
8888572537	20	6.7	5	20" (51cm)	10
8888572545	24	8.0	6	20" (51cm)	10
8888572552	28	9.3	6	20" (51cm)	10
8888572560	32	10.7	6	20" (51cm)	10
8888572578	36	12.0	6	20" (51cm)	10
8888572586	40	13.3	6	20" (51cm)	10
Right Angle					
8888573048	28	9.3	6	20" (51cm)	10
8888573055	32	10.7	6	20" (51cm)	10
8888573063	36	12.0	6	20" (51cm)	10



### Argyle™ Soft PVC Catheters

- Clear thermosensitive PVC catheter
- Better ability to manipulate for insertion
- Radiopaque Sentinel Line™ tubing with Sentinel Eye™ catheter
- Numeric depth marks every 2cm
- Beveled connector end can be cropped to form integral funnel connector
- Flexible package
- Sterile

CODE	FRENCH SIZE	O.D. MM	SIDE EYES	LENGTH	SHIP CASE
Straight					
8888570341	28	9.3	6	20" (51cm)	10
8888570358	32	10.7	6	20" (51cm)	10
8888570366	36	12.0	6	20" (51cm)	10
Right Angle					
8888570408	16	5.3	6	20" (51cm)	10
8888570416	20	6.7	6	20" (51cm)	10
8888570432	28	9.3	6	20" (51cm)	10
8888570440	32	10.7	6	20" (51cm)	10
8888570457	36	12.0	6	20" (51cm)	10

## Thoracic Catheters



### Argyle™ Silicone Thoracic Catheters with Pre-packaged Catheter Connector System

- Pre-packaged catheter connector system assures secure connection to any size drainage tube
- Biocompatible, thrombo-resistant silicone catheter
- Soft, kink-resistant
- Radiopaque Sentinel Line™ tubing with Sentinel Eye™ catheter
- Smooth side eyes and open distal tip
- Beveled end
- Sterile

CODE	FRENCH SIZE	O.D. MM	SIDE EYES	LENGTH	SHIP CASE
Right Angle					
8888573022	32	10.7	6	20" (51cm)	10



### Argyle™ Silicone Mediastinal Catheters

- Biocompatible, thrombo-resistant silicone catheter
- Flat drainage section with six eyes on one side
- Round section exits wound
- Beveled end
- Sterile

CODE	FRENCH SIZE	O.D. MM	SIDE EYES	LENGTH	SHIP CASE
8888573550	28	9.3	20" (51cm)	11.0mm	10
8888573568	34	11.3	20" (51cm)	13.4mm	10

## Thoracic Catheters



### Argyle™ Trocar Catheters

- Clear thermosensitive PVC catheter over precision aluminum trocar
- Surgically sharp trocar tip
- Radiopaque Sentinel Line™ tubing with Sentinel Eye™ catheter
- Numeric depth marks at 2cm increments
- Integral funnel end
- Sterile

CODE	FRENCH SIZE	O.D. MM	LENGTH	SHIP CASE
8888560805	8	2.7	9" (23cm)	10
8888561019	10	3.3	9" (23cm)	10
8888561027	12	4.0	9" (23cm)	10
8888561035	16	5.3	10" (25cm)	10
8888561043	20	6.7	16" (41cm)	10
8888561050	24	8.0	16" (41cm)	10
8888561068	28	9.3	16" (41cm)	10
8888561076	32	10.7	16" (41cm)	10



### Argyle™ Trocar Catheter Kits and Thoracostomy Procedure Trays

#### Sterile Contents:

#### Prep Side

- (1) Pair Latex-free Gloves
- (2) Prep Swabs
- (2) Prep Cups
- (1) Fenestrated Drape with Self-Adhesive Attachment Tabs (22" x 26")
- (1) Roll Tape (3" x 36")
- (1) Medicine Cup (50ml)

#### Procedure Side

- (2) Monoject™ 12cc Syringes with Caps
- (1) Monoject™ 25ga. x 5/8" Needle
- (1) Monoject™ 20ga. x 3 1/2" Needle
- (2) Monoject™ 20ga. x 1 1/2" Needle
- (1) Scalpel with #11 Blade
- (1) 2 1/2" #96 Keith Cutting Needle with Attached 0 Silk Suture
- (4) Gauze Sponges (4" x 4")
- (2) Telfa™ Island Gauze Dressings, precut (4" x 4")
- (1) Argyle™ Trocar Catheter (except in Thoracostomy Tray)

CODE	FRENCH SIZE	O.D. MM	LENGTH	SHIP CASE
8888565028	12	4.0	9" (23cm)	10
8888565036	20	6.7	10" (25cm)	10
8888565044	28	9.3	10" (25cm)	10
8888565507	Thoracostomy Procedure Tray*			10

\*Tray does not include catheter

## Specialty Products



### Argyle™ Left Ventricular Sump Vent Catheters

- Double lumen sump tube
- Clear PVC catheter
- Bullet tip
- With Cath-Lok retention rings
- X-ray opaque Sentinel Line
- Four side eyes
- Integral tapered connector end
- Blue pigtail funnel connector on secondary lumen
- Sterile

CODE	FRENCH SIZE	O.D. MM	LENGTH	SHIP CASE
8888590208	14	4.7	20" (51cm)	10
8888590216	18	6.0	20" (51cm)	10
8888590224	20	6.7	20" (51cm)	10



### Argyle™ Thi Aortic Perfusion Cannulas

- Clear PVC catheter
- Blue orientation line runs the length of the cannula
- Smooth polished beveled tip
- Integral funnel connector
- 20" (51cm) cannula length
- With Cath-Lok retention ring
- Sterile

CODE	FRENCH SIZE	O.D. MM	LENGTH OF TIP BEYOND RING	PROXIMAL LUMEN	SHIP CASE
Angled Tip					
8888591008	10	3.3	3/4"(19mm)	1/4"(6.4mm)	10
8888591123	12	4.0	13/16"(20mm)	1/4"(6.4mm)	10
8888591024	14	4.7	7/8"(22mm)	1/4"(6.4mm)	10
8888591040	18	6.0	1"(25mm)	1/4"(6.4mm)	10
8888591065	21	7.0	1 1/4"(31mm)	3/8"(9.5mm)	10
8888591081	24	8.0	1 1/2"(38mm)	3/8"(9.5mm)	10
Straight Tip					
8888591057	18	6.0	1 1/16"(17mm)	1/4"(6.4mm)	10
8888591073	21	7.0	7/8"(22mm)	3/8"(9.5mm)	10
8888591099	24	8.0	1"(25mm)	3/8"(9.5mm)	10



## Specialty Products



### Argyle™ Parallel "Y" Connector

- For use with Vena Caval Catheters
- Durable, clear polycarbonate construction
- 22mm parallel leg separation
- Smooth, contoured interior walls
- Ribbed connector ends
- Sterile

CODE	SIZE	SHIP CASE
8888584912	3/8" x 3/8" x 1/2" (9.5 x 9.5 x 13mm)	5



### Argyle™ Carotid Artery Shunt Kits

- Four French sizes in each package
- Precision-extruded medical grade, clear
- PVC shunt
- X-ray opaque Sentinel Line
- Polished tapered ends
- Sterile

CODE	CONTENTS	LENGTH	SHIP CASE
8888577700	4 Tubes	11" (28cm)	10
8888577775	4 Tubes	6" (15cm)	10

\*Each Kit Contains 4 Tubes: 8 Fr. (2.6mm), 10 Fr. (3.3mm) 12 Fr. (4.0mm), 14 Fr. (4.6mm)

\* Tridodecylmethylammonium Chloride chemically bonded to Heparin. A process developed at Battelle Memorial Laboratories.

Specialty Products



**Argyle™ Vascular Tourniquet Kit**

- Includes red rubber sheath with retaining plug, an obturator and umbilical tape
- Used to secure umbilical tape, suture, or Sentinel Loops™ tapes
- Sterile

CODE	CONTENTS	SIZE	SHIP CASE
8888585000	2 vascular tourniquet sheaths with plugs and obturators and 2 cut lengths of umbilical tape (1/8" x 24")	5" (13cm)	10



**Argyle™ Sentinel Loops™ Silicone Vascular Tapes**

- Soft, silicone formulation
- Radiopaque
- Color-coded
- Elliptical cross section
- Smooth surface
- Sterile

CODE	COLOR	CROSS SECTION	LENGTH	PER BOX	SHIP CARTON
Regular Size					
8888585208	Blue	1.0 x 2.5mm	15" (38.0cm)	10	100
8888585216	Red	1.0 x 2.5mm	15" (38.0cm)	10	100
8888585224	Yellow	1.0 x 2.5mm	15" (38.0cm)	10	100
8888585232	White	1.0 x 2.5mm	15" (38.0cm)	10	100
Small Size					
8888585240	Blue	0.8 x 1.5mm	14" (35.5cm)	10	100
8888585257	Red	0.8 x 1.5mm	14" (35.5cm)	10	100
8888585265	Yellow	0.8 x 1.5mm	14" (35.5cm)	10	100
8888585273	White	0.8 x 1.5mm	14" (35.5cm)	10	100

## Specialty Products



### Argyle™ Penrose Tubing

- Latex construction
- Radiopaque
- Sterile

CODE	FLAT SIZE	DIAMETER	LENGTH	SHIP CASE
8888513002	7/16" (11mm)	1/4" (6mm)	12" (30cm)	50
8888513200	1/2" (13mm)	5/16" (8mm)	12" (30cm)	50
8888513408	5/8" (16mm)	3/8" (10mm)	12" (30cm)	50
8888513606	3/4" (19mm)	1/2" (13mm)	12" (30cm)	50
8888513804	1" (25mm)	5/8" (16mm)	12" (30cm)	50
8888514000	1 1/4" (32mm)	3/4" (19mm)	12" (30cm)	50
8888514208	1 3/8" (35mm)	7/8" (22mm)	12" (30cm)	50
8888514406	1 1/2" (38mm)	1" (25mm)	12" (30cm)	50
8888514604	7/16" (11mm)	1/4" (6mm)	18" (46cm)	50
8888514802	1/2" (13mm)	5/16" (8mm)	18" (46cm)	50
8888515007	5/8" (16mm)	3/8" (10mm)	18" (46cm)	50
8888515205	3/4" (19mm)	1/2" (13mm)	18" (46cm)	50
8888515403	1" (25mm)	5/8" (16mm)	18" (46cm)	50
8888515601	1 1/4" (32mm)	3/4" (19mm)	18" (46cm)	50
8888515809	1 1/2" (38mm)	1" (25mm)	18" (46cm)	50



### Argyle™ Saratoga™ Sump Drains

- Clear thermosensitive PVC catheter
- Two sump vents
- Radiopaque Sentinel Line™ tubing with Sentinel Eye™ catheter
- Beveled connector end can be cropped to form integral funnel connector
- Sterile

CODE	FRENCH SIZE	O.D. MM	EYE DISTRIBUTION	LENGTH	SHIP CASE
8888530519	16	5.3	From	12" (30cm)	10
8888530527	20	6.7	Distal Tip	12" (30cm)	10
8888530535	28	9.3	to 2" (5cm)	12" (30cm)	10
8888531020	20	6.7	From	20" (51cm)	10
8888531038	24	8.0	Distal Tip	20" (51cm)	10
8888531046	28	9.3	to 4" (10cm)	20" (51cm)	10

## INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
8884713100	Argyle Thora-Seal I Chest Drainage Unit	7	8888570366	Argyle Soft PVC Catheter	11
8884713176	Thora-Seal III Autotransfusion Chest Drainage Unit	5	8888570408	Argyle Soft PVC Catheter	11
8884713184	Argyle Autotransfusion Accessory Unit	8	8888570416	Argyle Soft PVC Catheter	11
8884713192	Thora-Seal III Autotransfusion Chest Drain System	6	8888570432	Argyle Soft PVC Catheter	11
8884713308	Thora-Seal III Chest Drainage Unit	5	8888570440	Argyle Soft PVC Catheter	11
8884713900	Replacement Chamber	6	8888570457	Argyle Soft PVC Catheter	11
8884714200	Argyle Tubing Extension Set	9	8888570507	Argyle Thoracic Catheter	10
8884723000	Adapter for Emerson Pump	8	8888570515	Argyle Thoracic Catheter	10
8884723500	Argyle Thora-Seal III Emerson Pump/Thor Reg Adapt	8	8888570523	Argyle Thoracic Catheter	10
8884726100	Standard 2000cc Bottle Chest Drainage System	7	8888570531	Argyle Thoracic Catheter	10
8884726200	Standard 2000cc Bottle Chest Drainage System	7	8888570549	Argyle Thoracic Catheter	10
8884726500	Dual Drainage 2000cc Bottle Chest Drain System	7	8888570556	Argyle Thoracic Catheter	10
8888577700	Argyle Carotid Artery Shunt Kit	15	8888570564	Argyle Thoracic Catheter	10
8888577775	Argyle Carotid Artery Shunt Kit	15	8888570572	Argyle Thoracic Catheter	10
8888584912	Argyle Parallel "Y" Connector	15	8888571018	Argyle Thoracic Catheter	10
8888513002	Argyle Penrose Tubing	17	8888571026	Argyle Thoracic Catheter	10
8888513200	Argyle Penrose Tubing	17	8888571034	Argyle Thoracic Catheter	10
8888513408	Argyle Penrose Tubing	17	8888571042	Argyle Thoracic Catheter	10
8888513606	Argyle Penrose Tubing	17	8888571059	Argyle Thoracic Catheter	10
8888513804	Argyle Penrose Tubing	17	8888571067	Argyle Thoracic Catheter	10
8888514000	Argyle Penrose Tubing	17	8888571075	Argyle Thoracic Catheter	10
8888514208	Argyle Penrose Tubing	17	8888571299	Aqua-Seal Basic Chest Drainage Unit	2
8888514406	Argyle Penrose Tubing	17	8888571364	Aqua-Seal All Purpose Autotrans Chest Drain Unit	1
8888514604	Argyle Penrose Tubing	17	8888571370	Argyle Altitude Dry Suction Chest Drainage Unit	1
8888514802	Argyle Penrose Tubing	17	8888571406	Aqua-Seal Dual Drain Chest Drainage Unit	3
8888515007	Argyle Penrose Tubing	17	8888571489	Sentinel Seal Autotransfusion Chest Drainage Unit	4
8888515205	Argyle Penrose Tubing	17	8888571513	Sentinel Seal Dual Drain Chest Drainage Unit	3
8888515403	Argyle Penrose Tubing	17	8888571547	Argyle Double Seal Chest Drainage Unit	7
8888515601	Argyle Penrose Tubing	17	8888571562	Sentinel Seal Chest Drainage Unit	4
8888515809	Argyle Penrose Tubing	17	8888572040	Argyle Ferguson Thoracic Catheter	10
8888530519	Argyle Saratoga Sump Drain	17	8888572057	Argyle Ferguson Thoracic Catheter	10
8888530527	Argyle Saratoga Sump Drain	17	8888572065	Argyle Ferguson Thoracic Catheter	10
8888530535	Argyle Saratoga Sump Drain	17	8888572537	Argyle Silicone Thoracic Catheter	11
8888531020	Argyle Saratoga Sump Drain	17	8888572545	Argyle Silicone Thoracic Catheter	11
8888531038	Argyle Saratoga Sump Drain	17	8888572552	Argyle Silicone Thoracic Catheter	11
8888531046	Argyle Saratoga Sump Drain	17	8888572560	Argyle Silicone Thoracic Catheter	11
8888560805	Argyle Trocar Catheter	13	8888572578	Argyle Silicone Thoracic Catheter	11
8888561019	Argyle Trocar Catheter	13	8888572586	Argyle Silicone Thoracic Catheter	11
8888561027	Argyle Trocar Catheter	13	8888573022	Argyle Silicone Thoracic Catheter	12
8888561035	Argyle Trocar Catheter	13	8888573048	Argyle Silicone Thoracic Catheter	11
8888561043	Argyle Trocar Catheter	13	8888573055	Argyle Silicone Thoracic Catheter	11
8888561050	Argyle Trocar Catheter	13	8888573063	Argyle Silicone Thoracic Catheter	11
8888561068	Argyle Trocar Catheter	13	8888573550	Argyle Silicone Mediastinal Catheter	12
8888561076	Argyle Trocar Catheter	13	8888573568	Argyle Silicone Mediastinal Catheter	12
8888565028	Argyle Trocar Catheter Kit, Thoracostomy Proc Tray	13	8888585000	Argyle Vascular Tourniquet Kit	16
8888565036	Argyle Trocar Catheter Kit, Thoracostomy Proc Tray	13	8888585208	Argyle Sentinel Loops Silicone Vascular Tape	16
8888565044	Argyle Trocar Catheter Kit, Thoracostomy Proc Tray	13	8888585216	Argyle Sentinel Loops Silicone Vascular Tape	16
8888565507	Thoracostomy Procedure Tray	13	8888585224	Argyle Sentinel Loops Silicone Vascular Tape	16
8888570341	Argyle Soft PVC Catheter	11	8888585232	Argyle Sentinel Loops Silicone Vascular Tape	16
8888570358	Argyle Soft PVC Catheter	11	8888585240	Argyle Sentinel Loops Silicone Vascular Tape	16

**INDEX**

<b>CODE</b>	<b>PRODUCT</b>	<b>PAGE</b>
8888585257	Argyle Sentinel Loops Silicone Vascular Tape	16
8888585265	Argyle Sentinel Loops Silicone Vascular Tape	16
8888585273	Argyle Sentinel Loops Silicone Vascular Tape	16
8888590208	Argyle Left Ventricular Sump Vent Catheter	14
8888590216	Argyle Left Ventricular Sump Vent Catheter	14
8888590224	Argyle Left Ventricular Sump Vent Catheter	14
8888591008	Argyle Thi Aortic Perfusion Cannula	14
8888591024	Argyle Thi Aortic Perfusion Cannula	14
8888591040	Argyle Thi Aortic Perfusion Cannula	14
8888591057	Argyle Thi Aortic Perfusion Cannula	14
8888591065	Argyle Thi Aortic Perfusion Cannula	14
8888591073	Argyle Thi Aortic Perfusion Cannula	14
8888591081	Argyle Thi Aortic Perfusion Cannula	14
8888591099	Argyle Thi Aortic Perfusion Cannula	14
8888591123	Argyle Thi Aortic Perfusion Cannula	14





The products shown in this catalog may not be licensed or registered in all parts of the world. Please contact your local sales representative for availability

™\*VELCRO is a trademark of Velcro Industries B.V.  
COVIDIEN, COVIDIEN with Logo and ™ marked brands are trademarks of Covidien AG or an affiliate.  
©2008 Covidien. All rights reserved.



H6200CT-1M-1208