



**Enteral Feeding, Access & GI Products**

**1-800-962-9888**

**[www.SafeEnteralConnections.com](http://www.SafeEnteralConnections.com)**



**COVIDIEN**

## Table of Contents

### Enteral Feeding, Access & GI Products

#### ENTERAL FEEDING

##### Enteral Feeding Pumps/Pump Sets

Kangaroo ePump™ .....	1
Kangaroo Joey™ .....	2
Kangaroo™ 924 .....	3

##### Nestlé®\* SpikeRight™\* System Compatible

##### Feeding Sets

Kangaroo Feeding Sets for use with SpikeRight System .....	4
--	---

##### Pump Sets

Kangaroo DEHP-Free .....	5
Kangaroo Anti-Free Flow .....	5-6
Kangaroo Pump Sets .....	6
Kangaroo Gravity.....	7

##### Enteral Irrigation Kits

Kangaroo I.V. Pole Pouch and Kits .....	8
---	---

#### ENTERAL ACCESS

##### Nasogastric Feeding Tubes

##### Pediatric

Argyle™ Indwell™ .....	9
Argyle PVC with Sentinel Line™ .....	10

##### Adult

Kangaroo (non-weighted) .....	11
Entriflex™ .....	11
Dobhoff™ .....	12

##### CO<sub>2</sub> Detector

Co <sub>2</sub> nfirm Now™ .....	13
----------------------------------	----

#### Specialty Tubes

Dobhoff Naso-Jejunal Feeding and Gastric Decompression Tubes .....	13
Endo-Tube™ Feeding Tubes .....	14
EntriStar™ Gastro-Jejunostomy Tubes .....	14
Dobhoff Jejunal Feeding/Gastric Decompression System .....	14

#### Safety P.E.G. Systems

Dobhoff Safety P.E.G. Kits .....	15
EntriStar Safety P.E.G. Kits .....	16
EntriStar/Dobhoff Y-Port-G-Tube/P.E.G. Adapter .....	16

#### Gastrostomy and Skin Level Tubes

Kangaroo Gastrostomy Tubes .....	17
NutriPort™ Skin Level Kits.....	18-19

#### GI PRODUCTS

##### Salem Sump Tubes

GiEntri™ Port.....	20
Silicone GiEntri Port .....	20
Argyle PVC Salem Sump .....	21
Argyle PVC Salem Sump w/Anti-Reflux Valve .....	21
Argyle Silicone Salem Sump .....	21
Argyle Silicone Salem Sump w/Anti-Reflux Valve .....	22
Argyle Salem Sump Anti-Reflux Valve .....	21

##### Stomach Tubes

Argyle Stomach Tube (Levin type).....	23
Stomach Tube (Levin type).....	22
Argyle Edlich Gastric Lavage Trays .....	23

## Enteral Feeding Pumps/Pump Sets



### Kangaroo ePump™ Enteral Feeding Pump

- Attitude independent
- "Smart" pump set technology
- "ATM-style" ease of operation
- Feeding history up to 72 hours
- DEHP-free pump sets
- Anti-free-flow pump sets
- Flushing capabilities
- Features the safe enteral connection system

#### SPECIFICATIONS

- 15 hour battery
- Weight – 2.34 lbs without pole clamp  
2.69 lbs with pole clamp
- Feed Rate – 1-400 in 1 ml increments
- Flush Volume – 10-500 in 1 ml increments
- Flush Interval – 1-24 in 1 hr increments
- Multilingual – 15 languages

CODE	DESCRIPTION	SHIP CASE
382400	Kangaroo ePump with Pole Clamp	1
382491	Power Adapter	1
382492	Pole Clamp	1
F010484	Kangaroo ePump Battery	1
770000	Kangaroo ePump Backpack	1
770001	Kangaroo ePump Shoulder Pack	1
<b>Pump sets</b>		
774655	Proximal Spike Set	30
772055	500 ml Pump Set	30
773656	1000 ml Pump Set	30
776255	600 ml Rigid Container Pump Set	30
776250	1200 ml Rigid Container Pump Set	30
774669	Proximal Spike Feed and 1000 ml Flush Set	30
773662	1000 ml Feed Set w/ 1000 ml Flush Set	30
672062	Spike Feed and 500 ml Flush Set	30
674668	500 ml Feed Set with 500 ml Flush Set	30
714650	Sterile Proximal Spike Set	30
712054	Sterile 500 ml Pump Set	30
713650	Sterile 1000 ml Pump Set	30
716154	Sterile 100 ml Burette Set	30
<b>Recertification set</b>		
776150	100 ml Burette Recertification Set	5
<b>Kangaroo Feeding Sets for use with Nestlé®* SpikeRight™* System</b>		
775655	Kangaroo ePump Safety Screw Spike	30
775662	Kangaroo ePump Safety Screw Spike w/1000 ml Flush Bag	30
775762	Safety Screw Spike Adapter Cap	30

## Enteral Feeding Pumps/Pump Sets



### Kangaroo Joey™ Enteral Feeding Pump

- Unique ambulatory feeding and flushing capabilities
- Attitude Independent
- “Smart” pump set technology
- “ATM-style” ease of operation
- Feeding history up to 72 hours
- DEHP-free pump sets
- Anti-free-flow pump sets
- Fully portable design with back-pack
- Compact, lightweight construction
- Easy, intuitive interface
- Features the safe enteral connection system
- New ambulatory-specific feed set design

#### SPECIFICATIONS

- 18 hour battery
- Weight – 1.7 lbs., with pole clamp – 2.3 lbs
- Feed Rate – 1-400 in 1 ml increments
- Flush Volume – 10-500 in 1 ml increments
- Flush Interval – 1-24 in 1 hr intervals
- Memory – 72 hours
- Multilingual – 15 languages
- Tubing length: 9" (from pump to feeding container).

CODE	DESCRIPTION	SHIP CASE
383400	Kangaroo Joey Enteral Feeding Pump w/ Pole Clamp	1
383493	Kangaroo Joey Pole Clamp	1
383494	Kangaroo Joey Flexible Pole Clamp	1
770024	Small Kangaroo Joey Backpack	1
383491	Kangaroo Joey Power Cord	1
F010506	Kangaroo Joey Battery	1
<b>Pump sets</b>		
764655	Kangaroo Joey Proximal Spike Set	30
762055	Kangaroo Joey 500 ml Pump Set	30
763656	Kangaroo Joey 1000 ml Pump Set	30
764669	Kangaroo Joey Spike Set w/ 1000 ml Flush Bag	30
763662	Kangaroo Joey 1000 ml Feed Set with 1000 ml Flush Set	30
<b>Recertification set</b>		
776150	100 ml Burette Recertification Set	5
<b>Kangaroo Feeding Sets for use with Nestlé®* SpikeRight™* System</b>		
765555	Kangaroo Safety Screw Spike for Kangaroo Joey	30
765662	Kangaroo Joey Safety Screw Spike w/1000 ml Flush Bag	30
775762	Safety Screw Spike Adapter Cap	30

## Enteral Feeding Pumps/Pump Sets



### Kangaroo™ 924 Enteral Feeding Pump

- Easy to use controls
- Preset volume and over-infusion safeguard
- Audio and visual alarms
- Highly visible LED display and built-in pole clamp
- Drip chamber technology
- DEHP-free pump sets
- Lightweight and compact
- Features the safe enteral connection system

#### SPECIFICATIONS

- 16 hour battery
- Weight – 3.89 lbs
- Operating Range – 1-300 ml/hr
- Volume to be Delivered – 1-2000 ml in 1ml/hr increments
- Memory – 24 hours
- Multilingual – 13 languages



CODE	DESCRIPTION	SHIP CASE
392455	Kangaroo 924 Enteral Feeding Pump with Pole Clamp	1
<b>Pump sets</b>		
773621	1000 ml Bag Set (Indiv. Pack)	30
703621	1000 ml Bag Set (Bulk Pack)	30
774625	Spike Set (Indiv. Pack)	30
772025	500 ml Bag Set	30
<b>Kangaroo Feeding Sets for use with Nestlé®** SpikeRight™** System</b>		
775756	Kangaroo Safety Screw Spike for Kangaroo 924 Pump	30
775762	Safety Screw Spike Adapter Cap	30

## Nestlé®\* SpikeRight™\* Compatible Feeding Sets



775655

### Kangaroo™ Feeding Sets for use with SpikeRight™\* System

- Covidien is the first provider of a safety spike compatible with the Nestlé SpikeRight System
- The new enteral-specific SpikeRight port is threaded and connects with screw-top Enteral spike sets



775756

CODE	DESCRIPTION	SHIP CASE
765555	Kangaroo Joey Safety Screw Spike	30
765662	Kangaroo Joey Safety Screw Spike w/1000 ml Flush Bag	30
775655	Kangaroo ePump Safety Screw Spike	30
775662	Kangaroo ePump Safety Screw Spike w/1000 ml Flush Bag	30
775756	Kangaroo Safety Screw Spike for Kangaroo 924 Pump	30



775756



775762

### Kangaroo™ Adapter Cap

- For use with a Kangaroo Safety Spike Set
- Provides compatibility with an Abbott™\* ready-to-hang feeding container

CODE	DESCRIPTION	SHIP CASE
775762	Safety Screw Spike Adapter Cap	30

## Pump Sets for use with all Kangaroo™ 224, Kangaroo 324, Kangaroo PET™ and Kangaroo Control™ Enteral Feeding Pumps



### Kangaroo™ DEHP-Free Pump Sets

- New PVC material eliminates the use of the compound/chemical DEHP

CODE	TYPE	DESCRIPTION	SHIP CASE
773630	Non-sterile	1000 ml DEHP-Free Pump Set	30
713006	Sterile	1000 ml DEHP-Free Pump Set	36
702035	Non-sterile	500 ml DEHP-Free Pump Set	30
704645	Non-sterile	DEHP-Free Proximal Spike Set	30
714004	Sterile	DEHP-Free Proximal Spike Set	36
713100	Sterile	DEHP-Free 100 ml Burette Pump Set	36
702520	Non-sterile	1000 ml DEHP-Free Gravity Set	30



### Kangaroo™ DEHP-Free Pump Sets with Anti-Free-Flow (AFF)

- New PVC material eliminates the use of the compound/chemical DEHP
- Incorporates the internal Anti-Free-Flow (AFF) Device for continuous protection for the patient from the possibility of a free flow incident

CODE	TYPE	DESCRIPTION	SHIP CASE
773620	Non-sterile	1000 ml DEHP-Free with AFF Pump Set	30
702025	Non-sterile	500 ml DEHP-Free with AFF Pump Set	30
704625	Non-sterile	DEHP-Free with AFF Proximal Spike Set	30



### Kangaroo™ Anti-Free-Flow Pump Sets

- Anti-Free-Flow (AFF) Device protects the patient from the possible free flow of formula that may result from the pump set becoming dislodged or loaded improperly
- Unique internal design ensures the safety mechanism may not be tampered with by either the patient or caregiver, providing continuous protection

CODE	TYPE	DESCRIPTION	SHIP CASE
773610	Non-sterile	1000 ml AFF Pump Set	30
702015	Non-sterile	500 ml AFF Pump Set	30
704615	Non-sterile	AFF Proximal Spike Set	30



## Pump Sets for use with all Kangaroo™ 224, Kangaroo 324, Kangaroo PET™ and Kangaroo Control™ Enteral Feeding Pumps



8884773600



8884704605



8884706204



8884716100



8884705008

### Kangaroo™ Pump Sets

- Overinfusion safeguard (when used with Kangaroo 224, Kangaroo 324, Kangaroo PET, or Kangaroo Control Pumps)
- Non-IV compatible connectors

CODE	TYPE	DESCRIPTION	SHIP CASE
8884712004	Sterile	500 ml Bag with Easy-Cap Closure and Ice Pouch	36
8884702005	Non-sterile	500 ml Bag with Easy-Cap Closure	30
8884702105	Non-sterile	Deluxe 500 ml Graduated Bag with Easy-Cap Closure	30
8884713600	Sterile	1000 ml Bag with Easy-Cap Closure and Ice Pouch	36
8884773600	Non-sterile	1000 ml Bag with Easy-Cap Closure and Ice Pouch	30
8884712100	Sterile	Deluxe 500 ml Graduated Bag with Easy-Cap Closure	36
8884714300	Sterile	Deluxe 1200 ml Graduated Bag with Easy-Cap Closure	36
8884704305	Non-sterile	Easy-Cap 1200 ml	30
8884714800	Sterile	Deluxe 1600 ml Graduated Bag with Easy-Cap Closure	36
8884704805	Non-sterile	1600 ml Bag with Easy-Cap Closure	30
8884714600	Sterile	Proximal Spike Pump Set	36
8884704605	Non-sterile	Proximal Spike Pump Set	30
8884704705	Non-sterile	Screw Cap Pump Set — Screw Cap Designed to Fit Glass Enteral Formula Bottles	30
8884706204	Non-sterile	1200 ml Rigid Container with Preattached Pump Set	30
8884706212	Non-sterile	1200 ml Rigid Container without Pump Set	30
8884706220	Non-sterile	600 ml Rigid Container with Preattached Pump Set	30
8884706238	Non-sterile	600 ml Rigid Container without Pump Set	30
8884716100	Sterile	100 ml Rigid Burette	36

### Kangaroo Pump Set Accessories

8884705008	Non-sterile	Y-site Extension Set	30
8884700207	Non-sterile	Pump Set 4 ft. Tubing Extension	30



8884700207



## Gravity Sets



### Kangaroo™ Gravity Sets

CODE	TYPE	DESCRIPTION	SHIP CASE
702505	Non-sterile	1000 ml Large Bore Gravity Set	30
702605	Non-sterile	500 ml Intermittent Feeding Set with Proximal Spike and Preattached Gravity Set	30
702520	Non-sterile	1000 ml DEHP-free, Large Easy-Cap Opening	30
8884702500	Non-sterile	1000 ml Gravity Bag	30
8884712300	Sterile	1000 ml Gravity Bag	24

## Kangaroo™ I.V. Pole Pouch and Kits




### Kangaroo™ I.V. Pole Pouch and Kits

- Convenient tear-to-open packaging, resealable flap
- Transparent, polyethylene bag hangs easily, clearly labeled
- Non-sterile



CODE	DESCRIPTION	SHIP CASE
84064	I.V. Pole Pouch, includes 60cc Flat-Top Syringe	30
84065	I.V. Pole Pouch, includes 60cc Piston Syringe	30
84066	I.V. Pole Kit, includes 60cc Piston Syringe w/500cc Container	30
8884700108	60cc Piston Syringe w/500cc Container and Basin	20
8884700116	60cc Piston Syringe w/500cc Container	30
8884700124	60cc Piston Syringe	30

## Nasogastric Feeding Tubes — Pediatric



### Purple Argyle™ Indwell™ Polyurethane Feeding Tubes with Safe Enteral Connections

- Enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert†
- Purple color to denote enteral feeding — oral tip syringe only
- Transparent tube with Purple radiopaque Sentinel Line™ & Sentinel Eye™ tubing
- Sterile
- DEHP-free

CODE	MARKINGS	FRENCH SIZE	LENGTH	SHIP CASE
461355	From 5cm to 25cm in 1cm increments	3.5 (1.2mm)	12" (31cm)	10
461412		5 (1.7mm)	20" (51cm)	10
461503		5 (1.7mm)	36" (91cm)	10
461420		6.5 (2.2mm)	20" (51cm)	10
461438		6.5 (2.2mm)	36" (91cm)	10
461701		8 (2.7mm)	20" (51cm)	10
461800		8 (2.7mm)	42" (107cm)	10
461042		10 (3.3mm)	42" (107cm)	10



### Argyle™ Indwell™ Polyurethane Feeding Tubes

- Polyurethane material remains flexible for prolonged use
- Transparent
- Colored radiopaque Sentinel Line™ & Sentinel Eye™ tubing
- Proximal end accepts male luer connector
- Attached cap to cover proximal end
- Sterile

CODE	MARKINGS	FRENCH SIZE	LENGTH	SHIP CASE
888261355	From 5cm to 25cm in 1cm increments	3.5 (1.2mm)	12" (31cm)	10
888261412		5 (1.7mm)	20" (51cm)	10
888261503		5 (1.7mm)	36" (91cm)	10
888261420		6.5 (2.2mm)	20" (51cm)	10
888261438		6.5 (2.2mm)	36" (91cm)	10
888261701		8 (2.7mm)	20" (51cm)	10
888261800		8 (2.7mm)	42" (107cm)	10
888261042		10 (3.3mm)	42" (107cm)	10

## Nasogastric Feeding Tubes — Pediatric



8888260802

### Argyle™ PVC Feeding Tubes

- Thermosensitive material with Slide-Tex finish
- Colored radiopaque Sentinel Line™ tubing
- Two eyes
- Proximal end accepts male Luer connector
- Funnel connector
- Attached cap to cover proximal end
- Sterile

CODE	MARKINGS	FRENCH SIZE	LENGTH	SHIP CASE
8888261206	From 5cm to 25cm in 1cm increments	3.5 (1.2mm)	12" (31cm)	10
8888260802		5 (1.7mm)	16" (41cm)	50
8888261008		5 (1.7mm)	36" (91cm)	50
8888260711		6.5 (2.2mm)	16" (41cm)	50
8888260604		8 (2.7mm)	16" (41cm)	50
8888260406		8 (2.7mm)	42" (107cm)	50
8888260208		10 (3.3mm)	42" (107cm)	50

### Curity™ PVC Feeding Tubes

155720*	8" (20cm)	5 (1.7mm)	15" (38cm)	50
155721*	8" (20cm)	5 (1.7mm)	36" (91cm)	50
155722*	8" (20cm)	8 (2.7mm)	15" (38cm)	50
155723*	8" (20cm)	8 (2.7mm)	42" (107cm)	50

\* Does not contain radiopaque Sentinel Line or cm markings.

## Nasogastric Feeding Tubes — Adult



### Kangaroo™ Non-Weighted Feeding Tubes

- Clear polyurethane material
- Proximal Y-site connector with catheter tip adapter
- Large exit holes and open distal end
- Radiopaque stripe
- DEHP-free

CODE	FRENCH SIZE	LENGTH	SHIP CASE
8884711568	8	36" (91cm)	10
8884711550	10	36" (91cm)	10
8884711501	12	36" (91cm)	10
8884711519	14	36" (91cm)	10



### Entriflex™ Nasogastric Feeding Tubes with Safe Enteral Connections

- Dual port enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert<sup>†</sup>
- Purple color to denote enteral feeding — oral/cath tip syringe only
- Weighted feeding tubes with and without Flow-Through stylets
- Radiopaque polyurethane material
- CM markings
- MRI conditional once stylet is removed
- DEHP-free

CODE	INCLUDES STYLET	RIGID OUTLET PORT	WEIGHT	FRENCH SIZE	LENGTH	SHIP CASE
8884720817	No	Yes	3g	8	36" (91cm)	10
8884720841	No	Yes	3g	8	43" (109cm)	10
8884720825	Yes	Yes	3g	8	36" (91cm)	10
8884720858	Yes	Yes	3g	8	43" (109cm)	10
8884720855	Yes	Yes	5g	8	55" (140cm)	10
8884721088	Yes	Yes	5g	10	43" (109cm)	10
8884721096	No	Yes	5g	10	43" (109cm)	10
8884721055	Yes	Yes	5g	10	55" (140cm)	10
8884721211	No	No	5g	12	36" (91cm)	10
8884721237	No	No	5g	12	43" (109cm)	10
8884721252	Yes	Yes	5g	12	43" (109cm)	10
8884721255	Yes	Yes	5g	12	55" (140cm)	10

## Nasogastric Feeding Tubes — Adult



### Dobbhoff™ Nasogastric Feeding Tubes with Safe Enteral Connections

- Dual port enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert†
- Purple color to denote enteral feeding
- Weighted feeding tubes with and without Flow-Through stylets
- Radiopaque polyurethane material
- CM markings
- MRI conditional once stylet is removed
- DEHP-free

CODE	INCLUDES STYLET	RIGID OUTLET PORT	WEIGHT	FRENCH SIZE	LENGTH	SHIP CASE
8884710867	Yes	Yes	5g	8	43" (109cm)	10
8884710842	No	Yes	7g	8	43" (109cm)	10
8884710859	Yes	Yes	7g	8	43" (109cm)	10
8884710826	Yes	Yes	7g	8	55" (140cm)	10
8884711006	Yes	Yes	7g	10	43" (109cm)	10
8884711246	No	No	7g	12	43" (109cm)	10
8884711253	Yes	Yes	7g	12	43" (109cm)	10

## CO<sub>2</sub>nfirm Now™ CO<sub>2</sub> Detector



777702

### CO<sub>2</sub> Detector For Placement of Gastric Tubes

- Proprietary colorimetric paper technology
- Detects presence of CO<sub>2</sub> within seconds
- Reduces unnecessary waiting time and minimizes the number of X-rays needed for procedure
- Ends the guesswork about proper gastric tube placement
- Convenient bellows device:
  - Clears obstructions from distal tip
  - Prevents CO<sub>2</sub> from entering tube prematurely
  - Facilitates airflow for quickest results
- Easily attaches to either a bellows or syringe
- Requires no change in insertion technique
- Universally adapts to all feeding and NG tubes

CODE	DESCRIPTION	SHIP CASE
777702	CO <sub>2</sub> nfirm Now CO <sub>2</sub> Detector, with Bellows	10/Box

## Specialty Tubes



8884752513

### Dobbhoff™ Naso-Jejunal Feeding and Gastric Decompression Tube

- For “push” insertion over guidewire (kit included)
- Radiopaque polyurethane material
- Unique adjustable length feeding tube design
- Spiral feeding outlet holes
- Jejunal tube can be retracted within gastric portion
- MRI conditional once stylet is removed

CODE	OVERALL DIAMETER	JEJUNAL TUBE DIAMETER	OVERALL LENGTH	SHIP CASE
8884752513	16Fr.	9Fr.	67" (170cm)	2



## Specialty Tubes



### Endo-Tube™ Feeding Tubes

- Designed for precision placement in the small bowel using either endoscopy and/or fluoroscopy (includes guidewire and transfer cannula)
- Unique hollow bolus weighted tip
- Radiopaque polyurethane material
- MRI conditional once stylet is removed

CODE	WEIGHT	FRENCH SIZE	LENGTH	SHIP CASE
8884752505	7g	12	60" (152cm)	2



### EntriStar™ Gastro-Jejunostomy Tubes

- For use with 20 Fr. EntriStar Single Pass P.E.G. and Dobbhoff™ Safety P.E.G. kits
- Guidewire insertion
- Simultaneous gastric decompression and jejunal feeding through a P.E.G.
- Sterile
- MRI conditional once stylet is removed

CODE	FRENCH SIZE OF JEJUNAL TUBE	FRENCH SIZE OF P.E.G.	LENGTH OF JEJUNAL TUBE	SHIP CASE
8884750095	9	20	35" (89cm)	2



### Dobbhoff™ Jejunal Feeding/Gastric Decompression System

- For G/J tube placement using surgery, endoscopy, and/or fluoroscopy
- Strategically placed aspiration holes and circular lumen on gastrostomy tube for decompression
- Rigid outlet port for feeding
- Unique weight system provides an alternative to the suture tip “pull” method by incorporating an over-the-wire “push” technique
- Extended intragastric tip can be trimmed
- Sterile
- MRI conditional once stylet is removed

CODE	FRENCH SIZE OF JEJUNAL TUBE	FRENCH SIZE GASTROSTOMY TUBE	LENGTH OF JEJUNAL TUBE	SHIP CASE
8884790224	9	22	35" (89cm)	2 kits
8884790240	9	24	35" (89cm)	2 kits

## Safety P.E.G. Systems



### Dobbhoff™ Safety P.E.G. Kits with Safe Enteral Connections

- Enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert†
- Purple color to denote enteral feeding
- Featuring the Grant 90° Bend which straightens easily for P.E.G. placement and insertion of accessory feeding tubes
- Endoscopic removal
- Radiopaque polyurethane material
- Accessory jejunal feeding/gastric decompression tube is available for 20 Fr. P.E.G.
- Sterile
- DEHP-free

CODE		SIZE	PROCEDURAL AIDS	SHIP CASE
<b>Safety Percutaneous Endoscopic Gastrostomy “Pull” Kits</b>				
761610	<u>Kit Contains:</u> 21ga x 1 1/2” Monoject Magellan™ Safety Needle, Dual Port Enteral Connector, 6cc Syringe, Pinch Clamp, #11 Safety Scalpel, 10 ml Plastic Ampoule Lidocaine, Catheter Assembly w/90 Grant Bend, 18ga Safety Introducer Needle, External Retention Band, Hemostat, PVP Iodine Swabstick, Obturator Removal Kit, Gauze Sponges, Looping Wire/Guidewire, Lubrication Jelly Packet, Fenestrated Drape, Disposable Snare	16 Fr.	Yes	2
762010	<u>Kit Contains:</u> 21ga x 1 1/2” Monoject Magellan™ Safety Needle, Dual Port Enteral Connector, 6cc Syringe, Pinch Clamp, #11 Safety Scalpel, 10 ml Plastic Ampoule Lidocaine, Catheter Assembly w/90 Grant Bend, 18ga Safety Introducer Needle, External Retention Band, Hemostat, PVP Iodine Swabstick, Obturator Removal Kit, Gauze Sponges, Looping Wire/Guidewire, Lubrication Jelly Packet, Fenestrated Drape, Disposable Snare	20 Fr.	Yes	2

## Safety P.E.G. Systems



### EntriStar™ Safety Single Pass P.E.G. Kits with Safe Enteral Connections

- Enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert†
- Purple color to denote enteral feeding
- Radiopaque internal retention bolster
- Non-endoscopic removal
- Easy to read placement marks
- Hydromer™\* hydrophilic grafting on the “Pull” dilator tip
- Accessory jejunal feeding/gastric decompression tube is available for 20 Fr. P.E.G.
- Available in both “Pull” and “Push” over the guidewire methods
- Sterile
- DEHP-free

CODE		SIZE	PROCEDURAL AIDS	SHIP CASE
<b>Safety Percutaneous Endoscopic Gastrostomy “Pull” Kits</b>				
761600	<u>Kit Contains:</u> 21ga x 1 1/2" Monoject Magellan™ Safety Needle, Dual Port Enteral Connector, 6cc Syringe, Pinch Clamp, #11 Safety Scalpel, 10 ml Plastic Ampoule Lidocaine, Catheter Assembly, 18ga Safety Introducer Needle, External Retention Band, Hemostat, PVP Iodine Swabstick, Obturator Removal Kit, Gauze Sponges, Looping Wire/Guidewire, Lubrication Jelly Packet, Fenestrated Drape, Disposable Snare	16 Fr.	Yes	2
762000	<u>Kit Contains:</u> 21ga x 1 1/2" Monoject Magellan™ Safety Needle, Dual Port Enteral Connector, 6cc Syringe, Pinch Clamp, #11 Safety Scalpel, 10 ml Plastic Ampoule Lidocaine, Catheter Assembly, 18ga Safety Introducer Needle, External Retention Band, Hemostat, PVP Iodine Swabstick, Obturator Removal Kit, Gauze Sponges, Looping Wire/Guidewire, Lubrication Jelly Packet, Fenestrated Drape, Disposable Snare	20 Fr.	Yes	2
<b>Safety Guidewire Percutaneous Endoscopic Gastrostomy “Push” Kits</b>				
762001	<u>Kit Contains:</u> 20 Fr. Catheter Assembly, External Retention Disk, Dual Port Enteral Connector, 5cc Ampoule Xylocaine, 21ga. 1 1/2" Needle, 25ga. 5/8" Needle, Package Swabsticks, 6cc Syringe, #11 Scalpel, 16ga. Needle Catheter, Looping Wire, Lubricating Jelly Packet, Fenestrated Drape, Retention Band, Obturator Removal Kit, Hemostat, Packet PVP-Iodine Swabstick, Scissor, and 4-Gauze Sponges	20 Fr.	Yes	2



### EntriStar™/Dobbhoff™ Y-Port-G-Tube/P.E.G. Adapter with Safe Enteral Connections

- Dual port enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert†
- Purple color to denote enteral feeding
- For use with the EntriStar or Dobbhoff P.E.G. Kits
- The 20 Fr. Y-port can be used as the replacement Y-port on all Kangaroo gastrostomy tubes

CODE	DESCRIPTION	SHIP CASE
8884751622	16 Fr. Y-port/P.E.G. Adapter, Non-sterile	20
8884752042	20 Fr. Y-port/P.E.G. Adapter, Non-sterile	20

## Gastrostomy and Skin Level Tubes



### Kangaroo™ Gastrostomy Feeding Tubes with Y-Ports with Safe Enteral Connections

- Dual port enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert<sup>†</sup>
- Purple color to denote enteral feeding — oral/cath tip syringe only
- Medical-grade silicone tubing
- Rounded skin disc for external retention with feet
- Graduated shaft
- Gastric balloon for internal retention
- Rounded flush tip with open distal end
- DEHP-free

CODE	FRENCH SIZE	BALLOON SIZE	SHIP CASE
8884715122	12 Fr.	5cc.	5
8884715148	14 Fr.	5cc.	5
8884720163	16 Fr.	20cc.	5
8884720189	18 Fr.	20cc.	5
8884720205	20 Fr.	20cc.	5
8884720221	22 Fr.	20cc.	5
8884720247	24 Fr.	20cc.	5
8884720265	26 Fr.	20cc.	5
8884720285	28 Fr.	20cc.	5

## Gastrostomy and Skin Level Tubes



### NutriPort™ Skin Level Gastrostomy Kits with Safe Enteral Connections

- Enteral connector engineered to reduce risk of tubing misconnections
- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert†
- Purple color to denote enteral feeding
- 100% silicone material
- Durable balloon construction
- Raised feet design
- ClipStar™ enteral security clip
- Large range of sizes available
- DEHP-free

#### Contents of NutriPort Kit:

Skin Level Gastrostomy Device, Step-Up Shim, Bolus Feeding Set, Continuous Feeding Set with ClipStar Enteral Security Clip, 60cc Irrigation Syringe, 12cc Syringe, Gauze Pads

CODE	DESCRIPTION	SHIP CARTON	CODE	DESCRIPTION	SHIP CARTON
712080	Balloon SLG 12 Fr. x 0.8cm	1	714450	Balloon SLG 14 Fr. x 4.5cm	1
712100	Balloon SLG 12 Fr. x 1.0cm	1	714500	Balloon SLG 14 Fr. x 5.0cm	1
712120	Balloon SLG 12 Fr. x 1.2cm	1	716080	Balloon SLG 16 Fr. x 0.8cm	1
712150	Balloon SLG 12 Fr. x 1.5cm	1	716100	Balloon SLG 16 Fr. x 1.0cm	1
712170	Balloon SLG 12 Fr. x 1.7cm	1	716120	Balloon SLG 16 Fr. x 1.2cm	1
712200	Balloon SLG 12 Fr. x 2.0cm	1	716150	Balloon SLG 16 Fr. x 1.5cm	1
712230	Balloon SLG 12 Fr. x 2.3cm	1	716170	Balloon SLG 16 Fr. x 1.7cm	1
712250	Balloon SLG 12 Fr. x 2.5cm	1	716200	Balloon SLG 16 Fr. x 2.0cm	1
712270	Balloon SLG 12 Fr. x 2.7cm	1	716230	Balloon SLG 16 Fr. x 2.3cm	1
712300	Balloon SLG 12 Fr. x 3.0cm	1	716250	Balloon SLG 16 Fr. x 2.5cm	1
712350	Balloon SLG 12 Fr. x 3.5cm	1	716270	Balloon SLG 16 Fr. x 2.7cm	1
712400	Balloon SLG 12 Fr. x 4.0cm	1	716300	Balloon SLG 16 Fr. x 3.0cm	1
712450	Balloon SLG 12 Fr. x 4.5cm	1	716350	Balloon SLG 16 Fr. x 3.5cm	1
712500	Balloon SLG 12 Fr. x 5.0cm	1	716400	Balloon SLG 16 Fr. x 4.0cm	1
714080	Balloon SLG 14 Fr. x 0.8cm	1	716450	Balloon SLG 16 Fr. x 4.5cm	1
714100	Balloon SLG 14 Fr. x 1.0cm	1	716500	Balloon SLG 16 Fr. x 5.0cm	1
714120	Balloon SLG 14 Fr. x 1.2cm	1	718080	Balloon SLG 18 Fr. x 0.8cm	1
714150	Balloon SLG 14 Fr. x 1.5cm	1	718100	Balloon SLG 18 Fr. x 1.0cm	1
714170	Balloon SLG 14 Fr. x 1.7cm	1	718120	Balloon SLG 18 Fr. x 1.2cm	1
714200	Balloon SLG 14 Fr. x 2.0cm	1	718150	Balloon SLG 18 Fr. x 1.5cm	1
714230	Balloon SLG 14 Fr. x 2.3cm	1	718170	Balloon SLG 18 Fr. x 1.7cm	1
714250	Balloon SLG 14 Fr. x 2.5cm	1	718200	Balloon SLG 18 Fr. x 2.0cm	1
714270	Balloon SLG 14 Fr. x 2.7cm	1	718230	Balloon SLG 18 Fr. x 2.3cm	1
714300	Balloon SLG 14 Fr. x 3.0cm	1	718250	Balloon SLG 18 Fr. x 2.5cm	1
714350	Balloon SLG 14 Fr. x 3.5cm	1	718270	Balloon SLG 18 Fr. x 2.7cm	1
714400	Balloon SLG 14 Fr. x 4.0cm	1	718300	Balloon SLG 18 Fr. x 3.0cm	1

## Gastrostomy and Skin Level Tubes

### NutriPort™ Skin Level Gastrostomy Kits (continued)

CODE	DESCRIPTION	SHIP CARTON	CODE	DESCRIPTION	SHIP CARTON
718350	Balloon SLG 18 Fr. x 3.5cm	1	724080	Balloon SLG 24 Fr. x 0.8cm	1
718400	Balloon SLG 18 Fr. x 4.0cm	1	724100	Balloon SLG 24 Fr. x 1.0cm	1
718450	Balloon SLG 18 Fr. x 4.5cm	1	724120	Balloon SLG 24 Fr. x 1.2cm	1
718500	Balloon SLG 18 Fr. x 5.0cm	1	724150	Balloon SLG 24 Fr. x 1.5cm	1
720080	Balloon SLG 20 Fr. x 0.8cm	1	724170	Balloon SLG 24 Fr. x 1.7cm	1
720100	Balloon SLG 20 Fr. x 1.0cm	1	724200	Balloon SLG 24 Fr. x 2.0cm	1
720120	Balloon SLG 20 Fr. x 1.2cm	1	724230	Balloon SLG 24 Fr. x 2.3cm	1
720150	Balloon SLG 20 Fr. x 1.5cm	1	724250	Balloon SLG 24 Fr. x 2.5cm	1
720170	Balloon SLG 20 Fr. x 1.7cm	1	724270	Balloon SLG 24 Fr. x 2.7cm	1
720200	Balloon SLG 20 Fr. x 2.0cm	1	724300	Balloon SLG 24 Fr. x 3.0cm	1
720230	Balloon SLG 20 Fr. x 2.3cm	1	724350	Balloon SLG 24 Fr. x 3.5cm	1
720250	Balloon SLG 20 Fr. x 2.5cm	1	724400	Balloon SLG 24 Fr. x 4.0cm	1
720270	Balloon SLG 20 Fr. x 2.7cm	1	724450	Balloon SLG 24 Fr. x 4.5cm	1
720300	Balloon SLG 20 Fr. x 3.0cm	1	724500	Balloon SLG 24 Fr. x 5.0cm	1
720350	Balloon SLG 20 Fr. x 3.5cm	1			
720400	Balloon SLG 20 Fr. x 4.0cm	1			
720450	Balloon SLG 20 Fr. x 4.5cm	1			
720500	Balloon SLG 20 Fr. x 5.0cm	1			
			<b>Accessories</b>		
			<b>CODE</b>	<b>DESCRIPTION</b>	<b>SHIP CASE</b>
			8884741805	Stoma Measuring Device, Sterile	5
			8884741813	Bolus Feeding Set, Non-sterile	10
			8884741821	Continuous Feeding Set, Non-sterile	10

## Salem Sump™ Tubes



### GiEntri™ Port with Safe Enteral Connections

- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert†
- Purple color to denote enteral feeding
- All in one closed port system
- Allows for Feeding, Suction, Irrigation and Medication
- Internal ARV design
- DEHP-free
- Reduce clinician exposure
- Improve clinician efficacy
- Radiopaque

CODE	DESCRIPTION	SHIP CASE
7771010	10 Fr. Salem Sump with GiEntri Port	10
7771210	12 Fr. Salem Sump with GiEntri Port	10
7771410	14 Fr. Salem Sump with GiEntri Port	10
7771610	16 Fr. Salem Sump with GiEntri Port	10
7771810	18 Fr. Salem Sump with GiEntri Port	10
771111	GiEntri Port with single lumen adaptor	10



### Silicone GiEntri™ Port with Safe Enteral Connections

- Addresses 2006 JCAHO Tubing Misconnections Sentinel Event alert†
- Purple color to denote enteral feeding
- All in one closed port system
- Allows for Feeding, Suction, Irrigation and Medication
- Internal ARV design
- DEHP-free
- Reduce clinician exposure
- Improve clinician efficacy
- Radiopaque

CODE	DESCRIPTION	SHIP CASE
7771011	10 Fr. Silicone Salem Sump w/Anti-Reflux Valve Irrigation, GiEntri Port	10
7771211	12 Fr. Silicone Salem Sump w/Anti-Reflux Valve Irrigation, GiEntri Port	10
7771411	14 Fr. Silicone Salem Sump w/Anti-Reflux Valve Irrigation, GiEntri Port	10
7771611	16 Fr. Silicone Salem Sump w/Anti-Reflux Valve Irrigation, GiEntri Port	10
7771811	18 Fr. Silicone Salem Sump w/Anti-Reflux Valve Irrigation, GiEntri Port	10
771111	GiEntri Port with single lumen adaptor	10



## Salem Sump™ Tubes



8888266106

### Argyle™ Salem Sump™ Tube with Anti-Reflux Valve

- Gastric reflux prevention valve
- Double lumen PVC tube – one for suction drainage; one for sump vent
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Integral funnel connector with removable 5-in-1 adapter
- Integral irrigation funnel for irrigating through sump vent lumen
- Funnel end of vent lumen can be used to cap 5-in-1 adapter
- Sterile

CODE	FRENCH SIZE	LENGTH	SHIP CASE
8888266106	10 (3.3mm)	36" (91cm)	50
8888266114	12 (4.0mm)	48" (122cm)	50
8888266122	14 (4.7mm)	48" (122cm)	50
8888266130	16 (5.3mm)	48" (122cm)	50
8888266148	18 (6.0mm)	48" (122cm)	50



8888266197

### Argyle™ Salem Sump™ Anti-Reflux Valve

- Gastric reflux prevention valve
- Allows for maximum air flow to assist with stomach ulcer prevention
- Sterile

CODE	DESCRIPTION	SHIP CASE
8888266197	Valve only	10



8888268060

### Argyle™ Salem Sump™ Tube

- Double lumen PVC tube – one for suction drainage; one for sump vent
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Integral funnel connector with removable 5-in-1 adapter
- Integral irrigation funnel for irrigating through sump vent lumen
- Funnel end of vent lumen can be used to cap 5-in-1 adapter
- Sterile

CODE	FRENCH SIZE	LENGTH	SHIP CASE
8888268060	6 (2.0mm)	24" (61cm)	10
8888268086	8 (2.7mm)	24" (61cm)	10
8888264911	10 (3.3mm)	36" (91cm)	50
8888264929	12 (4.0mm)	48" (122cm)	50
8888264945	14 (4.7mm)	48" (122cm)	50
8888264960	16 (5.3mm)	48" (122cm)	50
8888264986	18 (6.0mm)	48" (122cm)	50

## Salem Sump™ Tubes



8888266213

### Argyle™ Silicone Salem Sump™ Tube with Anti-Reflux Valve

- Gastric reflux prevention valve
- Double lumen silicone tube – one for suction drainage; one for sump vent
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Removable 5-in-1 connector
- Capping system with sump vent lumen cap, suction lumen adapter and suction lumen adapter cap
- Sump vent lumen can also be used for irrigation
- Sterile

CODE	FRENCH SIZE	LENGTH	SHIP CASE
8888266213	10 (3.3mm)	36" (91cm)	10
8888266221	12 (4.0mm)	48" (122cm)	10
8888266247	14 (4.7mm)	48" (122cm)	10
8888266262	16 (5.3mm)	48" (122cm)	10
8888266288	18 (6.0mm)	48" (122cm)	10



8888265108

### Argyle™ Silicone Salem Sump™ Tubes

- Double lumen silicone tube – one for suction drainage; one for venting action
- Radiopaque Sentinel Line™ with Sentinel Eye™ tubing
- Removable 5-in-1 connector
- Capping system with sump vent lumen cap, suction lumen adapter and suction lumen adapter cap
- Sump vent lumen can also be used for irrigation
- Sterile

CODE	FRENCH SIZE	LENGTH	SHIP CASE
8888265108	10 (3.3mm)	36" (91cm)	10
8888265116	12 (4.0mm)	48" (122cm)	10
8888265124	14 (4.7mm)	48" (122cm)	10
8888265132	16 (5.3mm)	48" (122cm)	10
8888265140	18 (6.0mm)	48" (122cm)	10

## Stomach Tubes



### Argyle™ Stomach Tubes (Levin Type)

- Clear, thermosensitive PVC tube with four eyes
- Radiopaque Sentinel Line™ with Sentinel Eye™
- Open distal end
- Integral funnel connector
- Markings 18" (46cm), 22" (56cm), 26" (66cm) and 30" (76cm) from distal end
- Non-sterile

CODE	FRENCH SIZE	LENGTH	SHIP CASE
8888265314	12 (4.0mm)	48" (122cm)	50
8888265322	16 (5.3mm)	48" (122cm)	50
8888265330	18 (6.0mm)	48" (122cm)	50

### Stomach Tubes (Levin Type)

- Tube has four staggered eyes and center opening to reduce trauma
- Depth marks begin 20" from distal tip with three additional marks at 4" intervals
- Smooth, rounded tip reduces insertion trauma
- A variety of French sizes, 48" long
- Non-sterile

CODE	FRENCH SIZE	SHIP CASE
155710	12 Fr.	50
155711	14 Fr.	50
155712	16 Fr.	50
155713	18 Fr.	50



### Argyle™ Edlich Gastric Lavage Trays

- Large 140cc (5 oz.) graduated Monoject™ syringe
- Transparent PVC gastric tube with closed end and 4 eyes
- Depth marks at 18" (46cm), 22" (56cm), 26" (66cm) and 30" (76cm)
- Clamp

CODE	FRENCH SIZE	LENGTH	SHIP CASE
8888750018	Large, 34Fr. (11.3mm)	36" (91cm)	20
8888750513	Medium/Pediatric, 24Fr. (8.0mm)	36" (91cm)	20

## INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
155710	Stomach Tubes, PVC (Levin Type)	23	713650	Kangaroo ePump Sterile 1000 ml Pump Set	1
155711	Stomach Tubes, PVC (Levin Type)	23	714004	Kangaroo Sterile DEHP-Free Proximal Spike Set	5
155712	Stomach Tubes, PVC (Levin Type)	23	714080	NutriPort Skin Level Gastrostomy Kit	18
155713	Stomach Tubes, PVC (Levin Type)	23	714100	NutriPort Skin Level Gastrostomy Kit	18
155720	Curity PVC Feeding Tubes	10	714120	NutriPort Skin Level Gastrostomy Kit	18
155721	Curity PVC Feeding Tubes	10	714150	NutriPort Skin Level Gastrostomy Kit	18
155722	Curity PVC Feeding Tubes	10	714170	NutriPort Skin Level Gastrostomy Kit	18
155723	Curity PVC Feeding Tubes	10	714200	NutriPort Skin Level Gastrostomy Kit	18
382400	Kangaroo ePump with Pole Clamp	1	714230	NutriPort Skin Level Gastrostomy Kit	18
382491	Kangaroo ePump Power Adapter	1	714250	NutriPort Skin Level Gastrostomy Kit	18
382492	Kangaroo ePump Pole Clamp	1	714270	NutriPort Skin Level Gastrostomy Kit	18
383400	Kangaroo Joey Enteral Feeding Pump	2	714300	NutriPort Skin Level Gastrostomy Kit	18
383491	Kangaroo Joey Power Cord	2	714350	NutriPort Skin Level Gastrostomy Kit	18
383493	Kangaroo Joey Pole Clamp	2	714400	NutriPort Skin Level Gastrostomy Kit	18
383494	Kangaroo Joey Flexible Pole Clamp	2	714450	NutriPort Skin Level Gastrostomy Kit	18
392455	Kangaroo 924 Enteral Feeding Pump, Pole Clamp	3	714500	NutriPort Skin Level Gastrostomy Kit	18
461042	Purple Argyle Indwell Feeding Tubes	9	714650	Kangaroo ePump Sterile Proximal Spike Set	1
461355	Purple Argyle Indwell Feeding Tubes	9	716080	NutriPort Skin Level Gastrostomy Kit	18
461412	Purple Argyle Indwell Feeding Tubes	9	716100	NutriPort Skin Level Gastrostomy Kit	18
461420	Purple Argyle Indwell Feeding Tubes	9	716120	NutriPort Skin Level Gastrostomy Kit	18
461438	Purple Argyle Indwell Feeding Tubes	9	716150	NutriPort Skin Level Gastrostomy Kit	18
461503	Purple Argyle Indwell Feeding Tubes	9	716154	Kangaroo ePump Sterile 100 ml Burette Set	1
461701	Purple Argyle Indwell Feeding Tubes	9	716170	NutriPort Skin Level Gastrostomy Kit	18
461800	Purple Argyle Indwell Feeding Tubes	9	716200	NutriPort Skin Level Gastrostomy Kit	18
672062	Kangaroo ePump Spike Feed and 500 ml Flush Set	1	716230	NutriPort Skin Level Gastrostomy Kit	18
674668	Kangaroo ePump 500 ml Feed Set, Flush Set	1	716250	NutriPort Skin Level Gastrostomy Kit	18
702015	Kangaroo 500 ml AFF Pump Set	5	716270	NutriPort Skin Level Gastrostomy Kit	18
702025	Kangaroo 500 ml DEHP-Free with AFF Pump Set	5	716300	NutriPort Skin Level Gastrostomy Kit	18
702035	Kangaroo 500 ml DEHP-Free Pump Set	5	716350	NutriPort Skin Level Gastrostomy Kit	18
702505	Kangaroo 1000 ml Large Bore Gravity Set	7	716400	NutriPort Skin Level Gastrostomy Kit	18
702520	Kangaroo 1000 ml DEHP-Free Gravity Set	5	716450	NutriPort Skin Level Gastrostomy Kit	18
702520	Kangaroo 1000 ml DEHP-Free, Large Opening	7	716500	NutriPort Skin Level Gastrostomy Kit	18
702605	Kangaroo Intermittent Feeding Set	7	718080	NutriPort Skin Level Gastrostomy Kit	18
703621	Kangaroo 924 1000 ml Bag Set (Bulk Pack)	3	718100	NutriPort Skin Level Gastrostomy Kit	18
704615	Kangaroo AFF Proximal Spike Set	5	718120	NutriPort Skin Level Gastrostomy Kit	18
704625	Kangaroo DEHP-Free with AFF Proximal Spike Set	5	718150	NutriPort Skin Level Gastrostomy Kit	18
704645	Kangaroo DEHP-Free Proximal Spike Set	5	718170	NutriPort Skin Level Gastrostomy Kit	18
712054	Kangaroo ePump Sterile 500 ml Pump Set	1	718200	NutriPort Skin Level Gastrostomy Kit	18
712080	NutriPort Skin Level Gastrostomy Kit	18	718230	NutriPort Skin Level Gastrostomy Kit	18
712100	NutriPort Skin Level Gastrostomy Kit	18	718250	NutriPort Skin Level Gastrostomy Kit	18
712120	NutriPort Skin Level Gastrostomy Kit	18	718270	NutriPort Skin Level Gastrostomy Kit	18
712150	NutriPort Skin Level Gastrostomy Kit	18	718300	NutriPort Skin Level Gastrostomy Kit	18
712170	NutriPort Skin Level Gastrostomy Kit	18	718350	NutriPort Skin Level Gastrostomy Kit	19
712200	NutriPort Skin Level Gastrostomy Kit	18	718400	NutriPort Skin Level Gastrostomy Kit	19
712230	NutriPort Skin Level Gastrostomy Kit	18	718450	NutriPort Skin Level Gastrostomy Kit	19
712250	NutriPort Skin Level Gastrostomy Kit	18	718500	NutriPort Skin Level Gastrostomy Kit	19
712270	NutriPort Skin Level Gastrostomy Kit	18	720080	NutriPort Skin Level Gastrostomy Kit	19
712300	NutriPort Skin Level Gastrostomy Kit	18	720100	NutriPort Skin Level Gastrostomy Kit	19
712350	NutriPort Skin Level Gastrostomy Kit	18	720120	NutriPort Skin Level Gastrostomy Kit	19
712400	NutriPort Skin Level Gastrostomy Kit	18	720150	NutriPort Skin Level Gastrostomy Kit	19
712450	NutriPort Skin Level Gastrostomy Kit	18	720170	NutriPort Skin Level Gastrostomy Kit	19
712500	NutriPort Skin Level Gastrostomy Kit	18	720200	NutriPort Skin Level Gastrostomy Kit	19
713006	Kangaroo Sterile 1000 ml DEHP-Free Pump Set	5	720230	NutriPort Skin Level Gastrostomy Kit	19
713100	Kangaroo Sterile DEHP-Free Burette Pump Set	5	720250	NutriPort Skin Level Gastrostomy Kit	19

## INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
720270	NutriPort Skin Level Gastrostomy Kit	19	7771010	10 Fr. Salem Sump with GiEntri Port	20
720300	NutriPort Skin Level Gastrostomy Kit	19	7771011	Silicone GiEntri Port, Safe Enteral Connections	20
720350	NutriPort Skin Level Gastrostomy Kit	19	7771210	12 Fr. Salem Sump with GiEntri Port	20
720400	NutriPort Skin Level Gastrostomy Kit	19	7771211	Silicone GiEntri Port, Safe Enteral Connections	20
720450	NutriPort Skin Level Gastrostomy Kit	19	7771410	14 Fr. Salem Sump with GiEntri Port	20
720500	NutriPort Skin Level Gastrostomy Kit	19	7771411	Silicone GiEntri Port, Safe Enteral Connections	20
724080	NutriPort Skin Level Gastrostomy Kit	19	7771610	16 Fr. Salem Sump, GiEntri Port	20
724100	NutriPort Skin Level Gastrostomy Kit	19	7771611	Silicone GiEntri Port, Safe Enteral Connections	20
724120	NutriPort Skin Level Gastrostomy Kit	19	7771810	18 Fr. Salem Sump with GiEntri Port	20
724150	NutriPort Skin Level Gastrostomy Kit	19	7771811	Silicone GiEntri Port, Safe Enteral Connections	20
724170	NutriPort Skin Level Gastrostomy Kit	19	777702	CO <sub>2</sub> firm Now CO <sub>2</sub> Detector	13
724200	NutriPort Skin Level Gastrostomy Kit	19	84064	Kangaroo I.V. Pole Pouch, 60cc Flat-Top Syringe	8
724230	NutriPort Skin Level Gastrostomy Kit	19	84065	Kangaroo I.V. Pole Pouch, 60cc Piston Syringe	8
724250	NutriPort Skin Level Gastrostomy Kit	19	84066	Kangaroo I.V. Pole Kit; 60cc Piston Syringe, Container	8
724270	NutriPort Skin Level Gastrostomy Kit	19	8884700108	60cc Piston Syringe w/500cc Container and Basin	8
724300	NutriPort Skin Level Gastrostomy Kit	19	8884700116	60cc Piston Syringe w/500cc Container	8
724350	NutriPort Skin Level Gastrostomy Kit	19	8884700124	60cc Piston Syringe	8
724400	NutriPort Skin Level Gastrostomy Kit	19	8884700207	Kangaroo Pump Set	6
724450	NutriPort Skin Level Gastrostomy Kit	19	8884702005	Kangaroo 500 ml Bag with Easy-Cap Closure	6
724500	NutriPort Skin Level Gastrostomy Kit	19	8884702105	Kangaroo Deluxe 500 ml Graduated Bag	6
761600	EntriStar Safety Single Pass P.E.G. "Pull" Kit	16	8884702500	Kangaroo 1000 ml Gravity Bag	7
761610	Dobhoff Safety P.E.G. Kit	15	8884704305	Kangaroo Easy-Cap 1200 ml	6
762000	EntriStar Safety Single Pass P.E.G. "Pull" Kit	16	8884704605	Kangaroo Proximal Spike Pump Set	6
762001	EntriStar Safety Single Pass P.E.G. "Push" Kit	16	8884704705	Kangaroo Screw Cap Pump Set for Glass Bottles	6
762010	Dobhoff Safety P.E.G. Kit	15	8884704805	Kangaroo 1600 ml Bag with Easy-Cap Closure	6
762055	Kangaroo Joey 500 ml Pump Set	2	8884705008	Kangaroo Y-site Extension Set	6
763656	Kangaroo Joey 1000 ml Pump Set	2	8884706204	Kangaroo 1200 ml Rigid Container, Pump Set	6
763662	Kangaroo Joey 1000 ml Feed Set	2	8884706212	Kangaroo 1200 ml Rigid Container	6
764655	Kangaroo Joey Proximal Spike Set	2	8884706220	Kangaroo 600 ml Rigid Container, Pump Set	6
764669	Kangaroo Joey Spike Set	2	8884706238	Kangaroo 600 ml Rigid Container	6
765555	Kangaroo Joey Safety Screw Spike	2, 4	8884710826	Dobhoff Nasogastric Feeding Tube	12
765662	Kangaroo Joey Safety Screw Spike, Flush Bag	2, 4	8884710842	Dobhoff Nasogastric Feeding Tube	12
770000	Kangaroo ePump Backpack	1	8884710859	Dobhoff Nasogastric Feeding Tube	12
770001	Kangaroo ePump Shoulder Pack	1	8884710867	Dobhoff Nasogastric Feeding Tube	12
770024	Small Kangaroo Joey BackPack	2	8884711006	Dobhoff Nasogastric Feeding Tube	12
771111	GiEntri Port with single lumen adaptor	20	8884711246	Dobhoff Nasogastric Feeding Tube	12
772025	Kangaroo 924 500 ml Bag Set	3	8884711253	Dobhoff Nasogastric Feeding Tube	12
772055	Kangaroo ePump/Joey 500 ml Pump Set	1	8884711501	Kangaroo Non-weighted Feeding Tubes	11
773610	Kangaroo 1000 ml AFF Pump Set	5	8884711519	Kangaroo Non-weighted Feeding Tubes	11
773620	Kangaroo 1000 ml DEHP-Free, AFF Pump Set	5	8884711550	Kangaroo Non-weighted Feeding Tubes	11
773621	Kangaroo 924 1000 ml Bag Set (Indiv. Pack)	3	8884711568	Kangaroo Non-weighted Feeding Tubes	11
773630	Kangaroo 1000 ml DEHP-Free Pump Set	5	8884712004	Kangaroo Sterile Bag, Easy-Cap, Ice Pouch	6
773656	Kangaroo ePump/Joey 1000 ml Pump Set	1	8884712100	Kangaroo Sterile Deluxe Bag, Easy-Cap	6
773662	Kangaroo ePump Feed set with Flush Set	1	8884712300	Kangaroo Sterile Gravity Bag	7
774625	Kangaroo 924 Spike Set (Indiv. Pack)	3	8884713600	Kangaroo Sterile Bag, Easy-Cap, Ice Pouch	6
774655	Kangaroo ePump Proximal Spike Set	1	8884714300	Kangaroo Sterile Deluxe Bag, Easy-Cap	6
774669	Kangaroo ePump Proximal Spike Feed and Flush Set	1	8884714600	Kangaroo Sterile Proximal Spike Pump Set	6
775655	Kangaroo ePump Safety Screw Spike	1	8884714800	Kangaroo Sterile Deluxe Bag, Easy-Cap	6
775662	Kangaroo ePump Safety Screw Spike w/ Flush Bag	1, 4	8884715122	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	17
775756	Kangaroo 924 Pump Safety Screw Spike	3, 4	8884715148	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	17
775762	Safety Screw Spike Adapter Cap	1, 2, 3, 4	8884716100	Kangaroo Sterile 100 ml Rigid Burette	6
776150	100 ml Burette Recertification Set	1, 2	8884720163	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	17
776250	Kangaroo ePump 1200 ml Rigid Container Pump Set	1	8884720189	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	17
776255	Kangaroo ePump 600 ml Rigid Container Pump Set	1	8884720205	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	17

INDEX

CODE	PRODUCT	PAGE	CODE	PRODUCT	PAGE
8884720221	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	17	8888261412	Argyle Indwell Feeding Tubes	9
8884720247	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	17	8888261420	Argyle Indwell Feeding Tubes	9
8884720265	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	17	8888261438	Argyle Indwell Feeding Tubes	9
8884720285	Kangaroo Gastrostomy Feeding Tubes, Y-Ports	17	8888261503	Argyle Indwell Feeding Tubes	9
8884720817	Entriflex Nasogastric Feeding Tubes	11	8888261701	Argyle Indwell Feeding Tubes	9
8884720825	Entriflex Nasogastric Feeding Tubes	11	8888261800	Argyle Indwell Feeding Tubes	9
8884720841	Entriflex Nasogastric Feeding Tubes	11	8888264911	Argyle Salem Sump Tube	21
8884720855	Entriflex Nasogastric Feeding Tubes	11	8888264929	Argyle Salem Sump Tube	21
8884720858	Entriflex Nasogastric Feeding Tubes	11	8888264945	Argyle Salem Sump Tube	21
8884721055	Entriflex Nasogastric Feeding Tubes	11	8888264960	Argyle Salem Sump Tube	21
8884721088	Entriflex Nasogastric Feeding Tubes	11	8888264986	Argyle Salem Sump Tube	21
8884721096	Entriflex Nasogastric Feeding Tubes	11	8888265108	Argyle Silicone Salem Sump Tube	22
8884721211	Entriflex Nasogastric Feeding Tubes	11	8888265116	Argyle Silicone Salem Sump Tube	22
8884721237	Entriflex Nasogastric Feeding Tubes	11	8888265124	Argyle Silicone Salem Sump Tube	22
8884721252	Entriflex Nasogastric Feeding Tubes	11	8888265132	Argyle Silicone Salem Sump Tube	22
8884721255	Entriflex Nasogastric Feeding Tubes	11	8888265140	Argyle Silicone Salem Sump Tube	22
8884741805	Stoma Measuring Device, Sterile	19	8888265314	Argyle Stomach Tube	23
8884741813	Bolus Feeding Set	19	8888265322	Argyle Stomach Tube	23
8884741821	Continuous Feeding Set	19	8888265330	Argyle Stomach Tube	23
8884750095	EntriStar Gastro-Jejunostomy Tubes	14	8888266106	Argyle Salem Sump Tube w/ARV	21
8884751622	16 Fr. Y-port/P.E.G. Adapter	16	8888266114	Argyle Salem Sump Tube w/ARV	21
8884752042	20 Fr. Y-port/P.E.G. Adapter	16	8888266122	Argyle Salem Sump Tube w/ARV	21
8884752505	Endo-Tube Feeding Tube	14	8888266130	Argyle Salem Sump Tube w/ARV	21
8884752513	Dobhoff Naso-Jejunal Tube	13	8888266148	Argyle Salem Sump Tube w/ARV	21
8884773600	Kangaroo 1000 ml Bag, Easy-Cap, Ice Pouch	6	8888266197	Argyle Salem Sump Anti-Reflux Valve	21
8884790224	Dobhoff Jejunal System	14	8888266213	Argyle Silicone Salem Sump Tube w/ARV	22
8884790240	Dobhoff Jejunal System	14	8888266221	Argyle Silicone Salem Sump Tube w/ARV	22
8888260208	Argyle PVC Feeding Tubes	10	8888266247	Argyle Silicone Salem Sump Tube w/ARV	22
8888260406	Argyle PVC Feeding Tubes	10	8888266262	Argyle Silicone Salem Sump Tube w/ARV	22
8888260604	Argyle PVC Feeding Tubes	10	8888266288	Argyle Silicone Salem Sump Tube w/ARV	22
8888260711	Argyle PVC Feeding Tubes	10	8888268060	Argyle Salem Sump Tube	21
8888260802	Argyle PVC Feeding Tubes	10	8888268086	Argyle Salem Sump Tube	21
8888261008	Argyle PVC Feeding Tubes	10	8888750018	Argyle Edlich Gastric Lavage Trays	23
8888261042	Argyle Indwell Feeding Tubes	9	8888750513	Argyle Edlich Gastric Lavage Trays	23
8888261206	Argyle PVC Feeding Tubes	10	F010484	Kangaroo ePump Battery	1
8888261355	Argyle Indwell Feeding Tubes	9	FO10506	Kangaroo Joey Battery	2

The products shown in this catalog may not be licensed or registered in all parts of the world. Please contact your local sales representative for availability

† Sentinel Event Alert, Joint Commission on Accreditation of Healthcare Organizations, Issue 36, April 3, 2006.

®\*/™\* Nestlé and SpikeRight are trademarks of the Société des Produits Nestlé S.A., Vevey, Switzerland.

COVIDIEN, COVIDIEN with Logo and ™ marked brands are trademarks of Covidien AG or an affiliate.

©2008 Covidien AG or an affiliate. All rights reserved.



H-6200ENAF-10M-0908