Non-Reinforced AAMI Level 3 Surgical Gowns

CardinalHealth

Fssential to care™

Protection you need. Value you want.

* Scrub Nurse Gowns are bulk non-sterile (BNS) gowns intended for procedure kits.

The gown and hand towel are wrapped in SMS material for placement in the kit.



Material attributes:

SMS gown base fabric – Blue spunbond-meltblown-spunbond (SMS) nonwoven fabric. Spunbond and meltblown composition is polypropylene. White Cuff – Polyester

Drawing is for illustrative purposes only.

Biocompatibility:

Cardinal Health™ Surgical Gowns are tested for and meet the toxicity requirements contained in the ISO 10993 standards for biological evaluation and biocompatibility testing of medical devices.

Flammability standard:

Cardinal Health™ Surgical Gown materials meet the "Standard for the Flammability of Clothing Textiles" for disposable surgical drapes and gowns required by the Consumer Product Safety Commission (16 CFR Part 1610, formerly CS 191-53). The base fabric comprising Cardinal Health™ Surgical Gowns meets Class 1: Normal Flammability, which is the rating required for consumer use. Cardinal Health promotes the safe use of our products and does not recommend placing any ignition or heat source in contact with our surgical gowns.

