



COVIDIEN

positive results for life™

Valleylab™ Bipolar Forceps Designed for surgical procedures

These foot-switching bipolar forceps are precision crafted of the highest grade stainless steel and titanium. They are designed to provide superior control, safety and surgical efficacy at minimum output power settings.



Providing precision to bipolar electrosurgery



Benefits


- Provide the degree of precision and safety vital to bipolar electrosurgery
- Lightweight, less fatiguing for long procedures
- Facilitates tactile sensitivity
- Reduced glare
- Delivers consistent pressure over a wide range of tissue thickness

Key Features

- Foot switching
- Insulated
- Titanium option , 40 percent lighter than stainless steel

ORDER INFORMATION

CATALOG NUMBER	DESCRIPTION	ORDER QUANTITY
E4051CT – E4062CT	Reusable, insulated stainless steel, foot switching, bipolar coagulation forceps	1 each
E4073CT	Reusable, insulated titanium, foot switching, bipolar coagulation forceps	1 each
E0509	Sterile, single use 3.6 m (12 ft) bipolar coagulation forceps cord; not for use with Force Triad™ energy platform	50/case
E0512	Sterile, single use 3.6 (12 ft) bipolar coagulation forceps cord with molded plug	50/case
E0020V	Reusable 4.6 m (14.76 ft) silicone, bipolar cord for two-pin U.S.-style bipolar coagulation forceps for use with any Covidien electrosurgical generator	1 each
E0021S	Reusable 4.6 m (14.76 ft) silicone, bipolar cord for select Karl Storz bipolar devices for use with any Covidien electrosurgical generator	1 each
E0022W	Reusable 4.6 m (14.76 ft) silicone, bipolar cord for select Richard Wolf devices for use with any Covidien electrosurgical generator	1 each



TECHNICAL SPECIFICATIONS

E4051CT	Scoville-Greenwood bayonet forceps	19.7 cm (7-3/4 in.) 1.5 mm
E4052CT	Cushing forceps	15.2 cm (6 in.) 0.7 mm
E4053CT	Cushing forceps	17.8 cm (7 in.) 1.5 mm
E4054CT	Cushing bayonet forceps	19.1 cm (7-1/2 in.) 2.0 mm
E4055CT	Jewelers forceps, straight tip	10.2 cm (4 in.) 0.4 mm
E4057CT	Gerald bayonet forceps	19.1 cm (7-1/2 in.) 0.7 mm
E4058CT	Cushing bayonet forceps	19.1 cm (7-1/2 in.) 0.7 mm
E4059CT	Hardy bayonet forceps with stops	20.9 cm (8-1/4 in.) 0.5 mm
E4060CT	Semkin forceps with stops	14.0 cm (5-3/4 in.) 0.5 mm
E4062CT	Curved iris forceps	8.9 cm (3-1/2 in.) 0.5 mm
E4073CT	Titanium bayonet forceps, straight medium tip	22.2 cm (8-3/4 in.) 1.0 mm

COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.

©2013 Covidien



US130025 07/13

5920 LONGBOW DRIVE
BOULDER, CO
80301

303-530-2300 [T]
800-722-8772 [T]

WWW.COVIDIEN.COM