



COVIDIEN

*positive results for life™*



## LigaSure Atlas™ Hand Switching Instruments

For open and laparoscopic procedures



Energized by the ForceTriad™ energy platform, LigaSure™ instruments provide the following to potentially reduce procedure time:

- Hand activation for sterile-field control
- 2-4 second fusion cycle
- Fuses vessels and tissue bundles up to and including 7 mm

## FORCETRIAD™ ENERGY PLATFORM

### ORDER INFORMATION

Catalog Number	Order Quantity
FORCETRIAD	1 each

Contact your local energy expert for a quote



## LIGASURE ATLAS™ HAND SWITCHING INSTRUMENTS

### ORDER INFORMATION

Catalog Number	Order Quantity
LS1037	6 each
LS1020	6 each

### SPECIFICATIONS

Shaft Length	LS1037	37 cm
	LS1020	20 cm
Shaft Diameter	10 mm	
Shaft Rotation	359°	
Jaw Angle	Straight	
Jaw Texture	Smooth surface with ceramic stops	
Seal Length	22 mm	
Seal Width	6 mm (2 x 3 mm)	
Cut Length	20 mm	



COVIDIEN, COVIDIEN with logo and "positive results for life" are trademarks of Covidien AG.  
© 2008 Covidien AG. All rights reserved.

R0005298 Rev. 2008/12