



# Diagnostic procedure trays

# More of what you need. Less of what you don't.

**When it comes to diagnostic procedure trays,  
you want only the items you need—nothing more, nothing less.**

## Introducing the new Cardinal Health diagnostic procedure trays

### **The right products for the right procedure**

With our customers' feedback in mind, we designed a line of standard procedure trays filled with the essentials that you need. We've reduced the clutter and added more of what you want -- such as 4 x 4 in. gauze and ChloroPrep®.

### **Safety matters**

To you AND your patients. We offer the brand names you know and trust when it comes to safety. Like the Bard-Parker® safety scalpel and BD SafetyGlide® hypodermic needles.

### **Save time and money**

When you have all the items you need in a prepackaged kit, say for a myelogram, you'll save time by using the items directly from the kit. No more wasting time putting together a procedure tray or accidentally grabbing the wrong item. And, because you only have the items you need, there is nothing that gets thrown away unused.

### **Custom kits for custom procedures**

With 29 standard kits to choose from, you are likely to find the perfect kit to fit your needs. If not, no worries, we can custom design a kit for any procedure. In fact, we've been custom designing kits for years, so you can trust you'll get exactly what you need.

### **Superior service**

Our experienced, dedicated representatives are committed to meeting — and exceeding — your needs. With a single call or email, we'll make sure you have what you need when you need it.

### **Learn more**

For more information contact your Presource® products and services sales representative.

# Table of contents

## Diagnostic procedures

|                             |     |
|-----------------------------|-----|
| Lumbar puncture trays ..... | 2-3 |
| Myelogram trays .....       | 4   |
| Arthrogram trays .....      | 5   |
| Amniocentesis trays .....   | 6   |

## Acute drainage

|  |     |
|--|-----|
| Paracentesis/Thoracentesis trays ..... | 7-8 |
|--|-----|

## Bone marrow biopsy


|   |      |
|---|------|
| Bone marrow aspiration and biopsy trays ..... | 9-10 |
|---|------|

## Soft tissue biopsy

|                                |    |
|--------------------------------|----|
| Soft tissue biopsy trays ..... | 11 |
|--------------------------------|----|

## Basic

|                              |    |
|------------------------------|----|
| Basic diagnostic trays ..... | 12 |
|------------------------------|----|



## Lumbar puncture trays

Cardinal Health lumbar puncture trays have been designed to assist with any lumbar puncture procedure. Choose from conventional and safety versions and options for adults or pediatrics. Additionally, we offer a tray that includes a BD™ Whitacre Spinal Needle, providing an option designed to help reduce the occurrence of post-dural puncture headache<sup>1</sup>, a potential side effect of lumbar punctures.

Looking for a lumbar puncture tray without a manometer? You have come to the right place, we have that option too!

All of our trays feature needles — standard, safety and spinal — that are manufactured by BD Medical, a name you know and trust. Our offering provides the components you want such as impregnated swabsticks to facilitate patient prep and 4 x 4 in. gauze. That's right – we have larger gauze, and more of it!



### Common components .....Qty

|                                       |    |
|---------------------------------------|----|
| Adhesive bandage, 1 x 1 ½ in. ....    | 1  |
| Drape, fenestrated with tape .....    | 1  |
| Extension tubing, 5 in. ....          | 1  |
| Gauze, 8-ply, 4 x 4 in. ....          | 10 |
| Lidocaine 1% ampoule, 5 mL .....      | 1  |
| Medication label sheet .....          | 1  |
| PVP-I swabstick, 3/pk.....            | 1  |
| Foam needle block .....               | 1  |
| Stopcock, 3-way .....                 | 1  |
| Syringe, luer-lock, 5 mL .....        | 1  |
| Towel, 17 x 19 in. ....               | 1  |
| Utility marker .....                  | 1  |
| Vial with cap, clear, #1, 10 mL ..... | 1  |
| Vial with cap, clear, #2, 10 mL ..... | 1  |
| Vial with cap, clear, #3, 10 mL ..... | 1  |
| Vial with cap, clear, #4, 10 mL ..... | 1  |
| Wrap, 30 x 30 in. ....                | 1  |

<sup>1</sup>[www.bd.com/anesthesia/products/spinal.asp](http://www.bd.com/anesthesia/products/spinal.asp)

### Accessory set

#### Vial set, 10 mL

### Cat. no. SPP99DPTVL ..... Qty

|                                       |   |
|---------------------------------------|---|
| Vial with cap, clear, #1, 10 mL ..... | 1 |
| Vial with cap, clear, #2, 10 mL ..... | 1 |
| Vial with cap, clear, #3, 10 mL ..... | 1 |
| Vial with cap, clear, #4, 10 mL ..... | 1 |



All products listed on this page are not made with natural rubber latex.

### Conventional offering

#### Lumbar puncture with spinal needle, 20g x 3 1/2 in.

| Cat. no. 26-LP1CO                                       | Qty |
|---|-----|
| BD Precision Glide® Needle, 22g x 1 1/2 in. ....        | 1   |
| BD Precision Glide® Needle, 25g x 5/8 in. ....          | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 1/2 in. .... | 1   |
| BD™ Quincke Spinal Needle, 20g x 3 1/2 in. ....         | 1   |
| Manometer, two-piece.....                               | 1   |

#### Lumbar puncture with spinal needle, pediatric, 22g x 1 1/2 in.

| Cat. no. 26-LP2CO                                       | Qty |
|---|-----|
| BD Precision Glide® Needle, 22g x 1 1/2 in. ....        | 1   |
| BD Precision Glide® Needle, 25g x 5/8 in. ....          | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 1/2 in. .... | 1   |
| BD™ Quincke Spinal Needle, 22g x 1 1/2 in. ....         | 1   |
| Manometer, one-piece.....                               | 1   |

#### Lumbar puncture with spinal needle, 22g x 3 1/2 in.

| Cat. no. 26-LP3CO                                      | Qty |
|--|-----|
| BD Precision Glide® Needle, 22g x 1 1/2 in. ....       | 1   |
| BD Precision Glide® Needle, 25g x 5/8 in. ....         | 1   |
| BD Precision Glide® Filter Needle 19g x 1 1/2 in. .... | 1   |
| BD™ Quincke Spinal Needle, 22g x 3 1/2 in. ....        | 1   |
| Manometer, two-piece.....                              | 1   |

#### Lumbar puncture, atraumatic with spinal needle, 22 g x 3 1/2 in.

| Cat. no. 26-LP4CO  | Qty |
|--|-----|
| BD Precision Glide® Needle, 22 g x 1 1/2 in. ....        | 1   |
| BD Precision Glide® Needle, 25 g x 5/8 in. ....          | 1   |
| BD Precision Glide® Filter Needle, 19 g x 1 1/2 in. .... | 1   |
| BD™ Whitacre Spinal Needle, 22 g x 3 1/2 in. ....        | 1   |
| Manometer, two-piece.....                                | 1   |

#### Lumbar puncture, no manometer

| Cat. no. 26-LP5CO                                       | Qty |
|---|-----|
| BD Precision Glide® Needle, 22g x 1 1/2 in. ....        | 1   |
| BD Precision Glide® Needle, 25g x 5/8 in. ....          | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 1/2 in. .... | 1   |
| BD™ Quincke Spinal Needle, 20g x 3 1/2 in. ....         | 1   |

#### Lumbar puncture with spinal needle, 18g x 3 1/2 in.

| Cat. no. 26-LP6CO                                       | Qty |
|---|-----|
| BD Precision Glide® Needle, 22g x 1 1/2 in. ....        | 1   |
| BD Precision Glide® Needle, 25g x 5/8 in. ....          | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 1/2 in. .... | 1   |
| BD™ Quincke Spinal Needle, 18g x 3 1/2 in. ....         | 1   |
| Manometer, two-piece.....                               | 1   |

### Safety offering

#### Safety lumbar puncture with spinal needle, 20g x 3 1/2 in.

| Cat. no. 26-LP1SF                              | Qty |
|--|-----|
| BD SafetyGlide® Needle, 22g x 1 1/2 in. ....   | 1   |
| BD SafetyGlide® Needle, 25g x 5/8 in. ....     | 1   |
| BD Quincke Spinal Needle, 20g x 3 1/2 in. .... | 1   |
| FILTER STRAW®, 5 Micron, 1 3/4 in. ....        | 1   |
| Manometer, two-piece.....                      | 1   |

#### Safety lumbar puncture with spinal needle, pediatric, 22g x 1 1/2 in.

| Cat. no. 26-LP2SF                              | Qty |
|--|-----|
| BD SafetyGlide® Needle, 22g x 1 1/2 in.....    | 1   |
| BD SafetyGlide® Needle, 25g x 5/8 in. ....     | 1   |
| BD Quincke Spinal Needle, 22g x 1 1/2 in. .... | 1   |
| FILTER STRAW®, 5 Micron, 1 3/4 in. ....        | 1   |
| Manometer, one-piece.....                      | 1   |

#### Safety lumbar puncture with spinal needle, 22g x 3 1/2 in.

| Cat. no. 26-LP3SF                              | Qty |
|--|-----|
| BD SafetyGlide® Needle, 22g x 1 1/2 in. ....   | 1   |
| BD SafetyGlide® Needle, 25g x 5/8 in. ....     | 1   |
| BD Quincke Spinal Needle, 22g x 3 1/2 in. .... | 1   |
| FILTER STRAW®, 5 Micron, 1 3/4 in. ....        | 1   |
| Manometer, two-piece.....                      | 1   |

#### Safety lumbar puncture, atraumatic with spinal needle, 22g x 3 1/2 in.

| Cat. no. 26-LP4SF                                | Qty |
|--|-----|
| BD SafetyGlide® Needle, 22g x 1 1/2 in. ....     | 1   |
| BD SafetyGlide® Needle, 25g x 5/8 in. ....       | 1   |
| BD™ Whitacre Spinal Needle, 22g x 3 1/2 in. .... | 1   |
| FILTER STRAW®, 5 Micron, 1 3/4 in. ....          | 1   |
| Manometer, two-piece.....                        | 1   |

#### Safety lumbar puncture with spinal needle, 18g x 3 1/2 in.

| Cat. no. 26-LP6SF                              | Qty |
|--|-----|
| BD SafetyGlide® Needle, 22g x 1 1/2 in. ....   | 1   |
| BD SafetyGlide® Needle, 25g x 5/8 in. ....     | 1   |
| BD Quincke Spinal Needle, 18g x 3 1/2 in. .... | 1   |
| FILTER STRAW®, 5 Micron, 1 3/4 in. ....        | 1   |
| Manometer, two-piece.....                      | 1   |



## Myelogram trays

Aiding in the imaging of the spinal column, Cardinal Health myelogram trays include the essentials you need like impregnated PVP swabsticks for prep solution application. Choose from a BD™ Quincke Spinal Needle or the BD™ Whitacre Spinal Needle featuring a pencil point designed to provide precise access and control for greater patient comfort and safety.<sup>1</sup> Conventional and safety trays are available.



| Common components                         | Qty |
|---|-----|
| Adhesive bandage, 1 x 1 ½ in.....         | 1   |
| BD™ Blunt Fill Needle, 18g x 1 ½ in. .... | 1   |
| Drape, fenestrated with tape .....        | 1   |
| Extension tubing, 20 in. ....             | 1   |
| Extension tubing, 5 in. ....              | 1   |
| Gauze, 8-ply, 4 x 4 in. ....              | 10  |
| Lidocaine 1% ampoule, 5 mL .....          | 1   |
| Medication label sheet .....              | 1   |
| PVP-I swabstick, 3/pk.....                | 1   |
| Foam needle block.....                    | 1   |
| Syringe, luer-lock, 5 mL .....            | 1   |
| Syringe, luer-lock, 20 mL .....           | 1   |
| Towel, 17 x 19 in.....                    | 1   |
| Utility marker .....                      | 1   |
| Vial with cap, clear, #1, 10 mL .....     | 1   |
| Vial with cap, clear, #2, 10 mL .....     | 1   |
| Wrap, 30 x 30 in. ....                    | 1   |

### Conventional offering

#### Myelogram with spinal needle, 22g x 3 ½ in.

| Cat. no. 27-MY1CO                                     | Qty |
|---|-----|
| BD Precision Glide® Needle, 22g x 1 ½ in. ....        | 1   |
| BD Precision Glide® Needle, 25g x ¾ in. ....          | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 ½ in. .... | 1   |
| BD™ Quincke Spinal Needle, 22g x 3 ½ in. ....         | 1   |

#### Myelogram, atraumatic with spinal needle, 22g x 3 ½ in.

| Cat. no. 27-MY2CO                                     | Qty |
|---|-----|
| BD Precision Glide® Needle, 22g x 1 ½ in. ....        | 1   |
| BD Precision Glide® Needle, 25g x ¾ in. ....          | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 ½ in. .... | 1   |
| BD™ Whitacre Spinal Needle, 22g x 3 ½ in. ....        | 1   |

### Safety offering

#### Safety myelogram with spinal needle, 22g x 3 ½ in.

| Cat. no. 27-MY1SF                             | Qty |
|---|-----|
| BD SafetyGlide® Needle, 22g x 1 ½ in. ....    | 1   |
| BD SafetyGlide® Needle, 25g x ¾ in. ....      | 1   |
| FILTER STRAW®, 5 Micron, 1 ¾ in. ....         | 1   |
| BD™ Quincke Spinal Needle, 22g x 3 ½ in. .... | 1   |

<sup>1</sup>www.bd.com/anesthesia/products/spinal.asp



All products listed on this page are not made with natural rubber latex.

# Arthrogram trays

Cardinal Health arthrogram trays include the key components to successfully facilitate joint imaging. With a spinal needle, Chloraprep® one-step preoperative skin prep and 4 x 4 in. gauze, our arthrogram trays will meet your joint imaging needs.



| Common components   | Qty |
|---|-----|
| Adhesive bandage, 1 x 1 ½ in. ....                        | 1   |
| BD™ Blunt Fill Needle, 18g x 1 ½ in. ....                 | 1   |
| BD™ Quincke spinal needle, 22g x 3 ½ in. ....             | 1   |
| Chloraprep® applicator, Hi-Lite Orange® Tint, 3.0 mL .... | 1   |
| Drape, fenestrated with tape ....                         | 1   |
| Extension tubing, 20 in. ....                             | 1   |
| Gauze, 8-ply, 4 x 4 in. ....                              | 10  |
| Lidocaine 1% ampoule, 5 mL ....                           | 1   |
| Medication label sheet ....                               | 1   |
| Foam needle block ....                                    | 1   |
| Syringe, luer-lock, 5 mL ....                             | 1   |
| Syringe, luer-lock, 20 mL ....                            | 1   |
| Towel, 17 x 19 in. ....                                   | 1   |
| Utility marker ....                                       | 1   |
| Wrap, 30 x 30 in. ....                                    | 1   |

## Conventional offering

### Arthrogram with spinal needle, 22g x 3 ½ in.

| Cat. no. 28-AR1CO                                      | Qty |
|--|-----|
| BD Precision Glide® Needle, 22 g x 1 ½ in. ....        | 1   |
| BD Precision Glide® Needle, 25 g x 5/8 in. ....        | 1   |
| BD Precision Glide® Filter Needle, 19 g x 1 ½ in. .... | 1   |

## Safety offering

### Safety arthrogram with spinal needle, 22g x 3 ½ in.

| Cat. no. 28-AR1SF                           | Qty |
|---|-----|
| BD SafetyGlide® Needle, 22 g x 1 ½ in. .... | 1   |
| BD SafetyGlide® Needle, 25 g x 5/8 in. .... | 1   |
| FILTER STRAW®, 5 Micron, 1 ¾ in. ....       | 1   |

 All products listed on this page are not made with natural rubber latex.

## Amniocentesis trays



Designed to assist with sampling amniotic fluid, Cardinal Health amniocentesis trays allow you to complete your procedure efficiently and effectively. Components include a BD™ Ultra-vue echogenic procedure needle to allow visualization of the needle under ultrasound, three amber vials for fluid collection, ChloraPrep® one-step preoperative skin prep, and 4 x 4 in. gauze. Conventional and safety trays available.

### Common components

### Qty

|   |    |
|---|----|
| Adhesive bandage, 1 x 1 ½ in. ....                            | 1  |
| Cap, male/female, luer-lock adapter .....                     | 1  |
| ChloroPrep® applicator, Hi-Lite Orange® Tint, 3.0 mL .....    | 1  |
| Drape, fenestrated with tape .....                            | 1  |
| Extension tubing, 5 in. ....                                  | 1  |
| Gauze, 8-ply, 4 x 4 in. ....                                  | 10 |
| Lidocaine 1% ampoule, 5 mL .....                              | 1  |
| Medication label sheet .....                                  | 1  |
| BD™ Ultra-vue echogenic procedure needle, 20 g x 3 ½ in. .... | 1  |
| Ruler, 6 in. ....   | 1  |
| Foam needle block.....  | 1  |
| Syringe, luer-lock, 5 mL .....                                | 1  |
| Syringe, luer-lock, 20 mL .....                               | 1  |
| Towel, 17 x 19 in. ....                                       | 1  |
| Utility marker .....  | 1  |
| Vial with cap, amber, #1, 10 mL .....                         | 1  |
| Vial with cap, amber, #2, 10 mL .....                         | 1  |
| Vial with cap, amber, #3, 10 mL .....                         | 1  |
| Wrap, 30 x 30 in. ....  | 1  |

### Conventional offering

#### Amniocentesis with Ultra-vue echogenic needle, 20g x 3 ½ in.

### Cat. no. 29-AM1CO

### Qty

|   |   |
|---|---|
| BD Precision Glide® Needle, 22g x 1 ½ in. ....        | 1 |
| BD Precision Glide® Needle, 25g x 5/8 in. ....        | 1 |
| BD Precision Glide® Filter Needle, 19g x 1 ½ in. .... | 1 |

### Safety offering

#### Safety amniocentesis with Ultra-vue echogenic needle, 20g x 3 ½ in.

### Cat. no. 29-AM1SF

### Qty

|  |   |
|--|---|
| BD SafetyGlide® Needle, 22g x 1 ½ in. .... | 1 |
| BD SafetyGlide® Needle, 25g x 5/8 in. .... | 1 |
| FILTER STRAW®, 5 Micron, 1 ¾ in. ....      | 1 |



All products listed on this page are not made with natural rubber latex.



## Paracentesis/Thoracentesis trays

Cardinal Health paracentesis/thoracentesis trays have been designed to provide components facilitating in acute drainage of the pleural space and peritoneal cavity. One tray, two procedures. Our eight french trays feature the Turkel® safety fluid drainage system that includes a blunt safety cannula housed within a sharp beveled hollow needle and a safety color change indicator. Other tray features include: Bard-Parker® safety scalpel, ChloroPrep® one-step preoperative skin prep and 4 x 4 in. gauze. Conventional and safety trays available.



| Common components  | Qty |
|--|-----|
| Adhesive bandage, 1 x 1 ½ in. ....   | 1   |
| Bard-Parker® safety scalpel, #11 .....   | 1   |
| Cap, male/female, luer-lock adapter .....  | 1   |
| ChloroPrep® applicator, Hi-Lite Orange® Tint, 3.0 mL .....                           | 1   |
| Drape, fenestrated with tape .....   | 1   |
| Gauze, 8-ply, 4 x 4 in. ....   | 10  |
| Lidocaine 1% ampoule, 5 mL .....   | 2   |
| Medication label sheet .....   | 1   |
| Foam needle block .....  | 1   |
| Syringe, luer-lock, 5 mL .....   | 2   |
| Syringe, luer-lock, 30 mL .....  | 1   |
| Towel, 17 x 19 in. ....  | 1   |
| Transparent dressing, 4 x 4 ¾ in. ....   | 1   |
| Turkel® extension tubing with roller<br>clamp and 16 g x 1 ½ in. needle, 36 in. .... | 1   |
| Utility marker .....   | 1   |
| Wrap, 30 x 30 in. ....   | 1   |

### Accessory set

#### Centesis drainage bag set

| Cat. no. <b>SPP99CEBAA</b>    | Qty |
|-------------------------------|-----|
| Drainage Bag .....            | 1   |
| Drainage Bag Tubing Set ..... | 1   |



### Conventional offering

#### Centesis, 8 french x 5 in.

| Cat. no. <b>30-CE1CO</b>                                      | Qty |
|---|-----|
| BD Precision Glide® Needle, 22g x 1 ½ in. ....                | 1   |
| BD Precision Glide® Needle, 25g x 5/8 in. ....                | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 ½ in. ....         | 1   |
| Turkel® safety fluid drainage system<br>8 french x 5 in. .... | 1   |

#### Centesis with aspirating needle, 17g x 6 in.

| Cat. no. <b>30-CE2CO</b>   | Qty |
|--|-----|
| BD Precision Glide® Needle, 22g x 1 ½ in. ....                               | 1   |
| BD Precision Glide® Needle, 25g x 5/8 in. ....                               | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 ½ in. ....                        | 1   |
| Curity® Thoracentesis needle,<br>14g Catheter over 17 g x 6 in. needle ..... | 1   |
| Turkel® extension tubing with stopcock, 3 in. ....                           | 1   |

#### Centesis, no catheter / needle

| Cat. no. <b>30-CE3CO</b>                              | Qty |
|---|-----|
| BD Precision Glide® Needle, 22g x 1 ½ in. ....        | 1   |
| BD Precision Glide® Needle, 25g x 5/8 in. ....        | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 ½ in. .... | 1   |
| Adapter .....   | 1   |
| Drainage Bag .....                                    | 1   |
| Drainage Bag Tubing Set .....                         | 1   |

### Safety offering

#### Safety centesis, 8 french x 5 in.

| Cat. no. <b>30-CE1SF</b>                                      | Qty |
|---|-----|
| BD SafetyGlide® Needle, 22g x 1 ½ in. ....                    | 1   |
| BD SafetyGlide® Needle, 25g x 5/8 in. ....                    | 1   |
| FILTER STRAW®, 5 Micron, 1 ¾ in. ....                         | 1   |
| Turkel® safety fluid drainage system<br>8 french x 5 in. .... | 1   |

#### Safety centesis with aspirating needle, 17g x 6 in.

| Cat. no. <b>30-CE2SF</b>  | Qty |
|---|-----|
| BD SafetyGlide® Needle, 22g x 1 ½ in. ....                                    | 1   |
| BD SafetyGlide® Needle, 25g x 5/8 in. ....                                    | 1   |
| Curity® Thoracentesis Needle,<br>14 g Catheter over 17 g x 6 in. needle ..... | 1   |
| FILTER STRAW®, 5 Micron, 1 ¾ in. ....   | 1   |

#### Safety centesis, no catheter / needle

| Cat. no. <b>30-CE3SF</b>                   | Qty |
|--|-----|
| BD SafetyGlide® Needle, 22g x 1 ½ in. .... | 1   |
| BD SafetyGlide® Needle, 25g x 5/8 in. .... | 1   |
| FILTER STRAW®, 5 Micron, 1 ¾ in. ....      | 1   |
| Adapter .....                              | 1   |
| Drainage Bag .....                         | 1   |
| Drainage Bag Tubing Set .....              | 1   |




All products listed on this page are not made with natural rubber latex.

## Bone marrow aspiration and biopsy trays

Cardinal Health bone marrow aspiration and biopsy trays provide many options to help you aspirate and capture bone marrow samples. Choose the J-type needle, I-type needle or both. Don't want a bone marrow needle inside your kit? We have that option too!

The J-type bone marrow biopsy needle is designed for both bone marrow aspirations and biopsies. It also comes with an extraction cannula to help sever and capture high-quality samples. The I-type bone marrow aspiration needle is designed for sternal or iliac crest aspiration with an adjustable depth guard helping with proper needle placement.

Other features of these trays include: Bard-Parker® safety scalpel, ChloroPrep® one-step preoperative skin prep and 4 x 4 in. gauze. Conventional and safety trays available.



| Common components  | Qty |
|--|-----|
| Adhesive bandage, 1 x 1 ½ in. ....                         | 1   |
| Bard-Parker® safety scalpel, #11 .....                     | 1   |
| ChloroPrep® applicator, Hi-Lite Orange® Tint, 3.0 mL ..... | 1   |
| Drape, fenestrated with tape .....                         | 1   |
| Gauze, 8-ply, 4 x 4 in. ....                               | 10  |
| Lidocaine 1% ampoule, 5 mL .....                           | 2   |
| Medication label sheet .....                               | 1   |
| Foam needle block .....                                    | 1   |
| Syringe, luer-lock, 10 mL .....                            | 1   |
| Syringe, luer-lock, 20 mL .....                            | 1   |
| Towel, 17 x 19 in. ....                                    | 1   |
| Utility marker .....                                       | 1   |
| Wrap, 30 x 30 in. ....                                     | 1   |

| Single needle offering   |             |               |
|--|-------------|---------------|
| Cat. no.   | Needle type | Needle size   |
| SPP99BM8A  | J-type      | 8g x 4 in.    |
| SPP99BM11A   | J-type      | 11g x 4 in.   |
| SPP99BM16A   | J-type      | 11g x 6 in.   |
| SPP99BM13A   | J-type      | 13g x 2 ½ in. |
| <i>J-type needles are packaged sterile with extraction cannula, probe and probe guide.</i> |             |               |
| SPP99BM15A   | I-type      | 15g           |
| <i>I-type needles are packaged sterile with adjustable depth guard.</i>                    |             |               |

 All products listed on this page are not made with natural rubber latex.

## Conventional offering

### Bone marrow biopsy with J-type needle, 11g x 4 in.

| Cat. no. 31-BM1CO                                     | Qty |
|---|-----|
| BD Precision Glide® Needle, 22g x 1 ½ in. ....        | 1   |
| BD Precision Glide® Needle, 25g x ⅝ in. ....          | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 ½ in. .... | 1   |
| Bone marrow biopsy extraction cannula .....           | 1   |
| Bone marrow biopsy probe .....                        | 1   |
| Bone marrow biopsy probe guide .....                  | 1   |
| Needle, bone marrow biopsy, J-Type, 11g x 4 in. ....  | 1   |

### Bone marrow biopsy with J-type and I-type needles

| Cat. no. 31-BM2CO   | Qty |
|---|-----|
| BD Precision Glide® Needle, 22g x 1 ½ in. ....                          | 1   |
| BD Precision Glide® Needle, 25g x ⅝ in. ....                            | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 ½ in. ....                   | 1   |
| Bone marrow biopsy extraction cannula .....                             | 1   |
| Bone marrow biopsy probe .....  | 1   |
| Bone marrow biopsy probe guide .....                                    | 1   |
| Needle, bone marrow biopsy,<br>I-type, with adjustable guard, 15g ..... | 1   |
| Needle, bone marrow biopsy, J-Type, 11g x 4 in. ....                    | 1   |

### Bone marrow biopsy with I-type needle, 15g

| Cat. no. 31-BM3CO  | Qty |
|--|-----|
| BD Precision Glide® Needle, 22g x 1 ½ in. ....                               | 1   |
| BD Precision Glide® Needle, 25g x ⅝ in. ....                                 | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 ½ in. ....                        | 1   |
| Needle, bone marrow biopsy,<br>I-type with adjustable depth guard, 15g ..... | 1   |

### Bone marrow biopsy, no needle

| Cat. no. 31-BM4CO                                     | Qty |
|---|-----|
| BD Precision Glide® Needle, 22g x 1 ½ in. ....        | 1   |
| BD Precision Glide® Needle, 25g x ⅝ in. ....          | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 ½ in. .... | 1   |

### Bone marrow biopsy with J-type and I-type needles

| Cat. no. 31-BM6CO   | Qty |
|---|-----|
| BD Precision Glide® Needle, 22g x 1 1/2 in. ....                        | 1   |
| BD Precision Glide® Needle, 25g x 5/8 in. ....                          | 1   |
| BD Precision Glide® Filter Needle, 19g x 1 1/2 in. ....                 | 1   |
| Bone marrow biopsy extraction cannula .....                             | 1   |
| Bone marrow biopsy probe .....  | 1   |
| Bone marrow biopsy probe guide .....                                    | 1   |
| Needle, bone marrow biopsy,<br>I-type, with adjustable guard, 15g ..... | 1   |
| Needle, bone marrow biopsy, J-Type, 8g x 4 in. ....                     | 1   |

## Safety offering

### Safety bone marrow biopsy with J-type needle, 11g x 4 in.

| Cat. no. 31-BM1SF                                    | Qty |
|--|-----|
| BD SafetyGlide® Needle, 22g x 1 ½ in. ....           | 1   |
| BD SafetyGlide® Needle, 25g x ⅝ in. ....             | 1   |
| FILTER STRAW®, 5 Micron, 1 ¾ in. ....                | 1   |
| Bone marrow biopsy extraction cannula .....          | 1   |
| Bone marrow biopsy probe .....                       | 1   |
| Bone marrow biopsy probe guide .....                 | 1   |
| Needle, bone marrow biopsy, J-type, 11g x 4 in. .... | 1   |

### Safety bone marrow biopsy, J-type and I-type needles

| Cat. no. 31-BM2SF  | Qty |
|--|-----|
| BD SafetyGlide® Needle, 22g x 1 ½ in. ....                                   | 1   |
| BD SafetyGlide® Needle, 25g x ⅝ in. ....                                     | 1   |
| FILTER STRAW®, 5 Micron, 1 ¾ in. ....  | 1   |
| Bone marrow biopsy extraction cannula .....                                  | 1   |
| Bone marrow biopsy probe .....   | 1   |
| Bone marrow biopsy probe guide .....   | 1   |
| Needle, bone marrow biopsy,<br>I-type with adjustable depth guard, 15g ..... | 1   |
| Needle, bone marrow biopsy, J-type, 11g x 4 in. ....                         | 1   |

### Safety bone marrow biopsy, no needle

| Cat. no. 31-BM4SF                          | Qty |
|--|-----|
| BD SafetyGlide® Needle, 22g x 1 ½ in. .... | 1   |
| BD SafetyGlide® Needle, 25g x ⅝ in. ....   | 1   |
| FILTER STRAW®, 5 Micron, 1 ¾ in. ....      | 1   |

### Safety bone marrow biopsy with J-type and I-type needles

| Cat. no. 31-BM6SF  | Qty |
|--|-----|
| BD SafetyGlide® Needle, 22g x 1 ½ in. ....                                   | 1   |
| BD SafetyGlide® Needle, 25g x ⅝ in. ....                                     | 1   |
| FILTER STRAW®, 5 Micron, 1 ¾ in. ....  | 1   |
| Bone marrow biopsy extraction cannula .....                                  | 1   |
| Bone marrow biopsy probe .....   | 1   |
| Bone marrow biopsy probe guide .....   | 1   |
| Needle, bone marrow biopsy,<br>I-type with adjustable depth guard, 15g ..... | 1   |
| Needle, bone marrow biopsy, J-type, 8g x 4 in. ....                          | 1   |



All products listed on this page are not made with natural rubber latex.

## Soft tissue biopsy trays

Cardinal Health soft tissue biopsy trays come equipped with the necessary components to facilitate your next soft tissue biopsy procedure. Conventional and safety trays are available with options that include: Bard-Parker® safety scalpel, ChloroPrep® one-step preoperative skin prep and 4 x 4 in. gauze.



| Common components                                    | Qty |
|--|-----|
| Adhesive bandage, 1 in. x 1 ½ in.                    | 1   |
| Bard-Parker® safety scalpel, #11                     | 1   |
| BD™ Blunt Fill Needle, 18 g x 1 ½ in.                | 1   |
| Cap, male/female, luer-lock adapter                  | 1   |
| ChloroPrep® applicator, Hi-Lite Orange® Tint, 3.0 mL | 1   |
| Drape, fenestrated with tape                         | 1   |
| Gauze, 8-ply, 4 x 4 in.                              | 10  |
| Lidocaine 1% ampoule, 5 mL                           | 1   |
| Medication label sheet                               | 1   |
| Ruler, 6 in.   | 1   |
| Foam needle block                                    | 1   |
| Syringe, luer-lock, 3 mL                             | 1   |
| Syringe, luer-lock, 5 mL                             | 1   |
| Syringe, luer-lock, 20 mL                            | 1   |
| Towel, 17 x 19 in.                                   | 1   |
| Utility marker                                       | 1   |
| Vial with cap, clear, #1, 10 mL                      | 1   |
| Vial with cap, clear, #2, 10 mL                      | 1   |
| Wrap, 30 x 30 in.                                    | 1   |

### Conventional offering

#### Soft tissue biopsy

| Cat. no. 32-ST1CO                                 | Qty |
|---|-----|
| BD Precision Glide® Needle, 22 g x 1 ½ in.        | 1   |
| BD Precision Glide® Needle, 25 g x 5/8 in.        | 1   |
| BD Precision Glide® Filter Needle, 19 g x 1 ½ in. | 1   |

### Safety offering

#### Safety soft tissue biopsy

| Cat. no. 32-ST1SF                      | Qty |
|--|-----|
| BD SafetyGlide® Needle, 22 g x 1 ½ in. | 1   |
| BD SafetyGlide® Needle, 25 g x 5/8 in. | 1   |
| FILTER STRAW®, 5 Micron, 1 ¾ in.       | 1   |



All products listed on this page are not made with natural rubber latex.

## Basic diagnostic trays

Cardinal Health basic diagnostic trays include supplies that can be used in almost every diagnostic procedure. The tray's unique design allows you to use your preferred specialty components while standardizing the basics. With conventional and safety trays available, components include the Bard-Parker® safety scalpel, ChloraPrep® one-step preoperative skin prep and 4 x 4 in. gauze.



| Common components   | Qty |
|---|-----|
| Adhesive bandage, 1 x 1 ½ in.....                         | 1   |
| Bard-Parker® safety scalpel, #11 .....                    | 1   |
| ChloraPrep® applicator, Hi-Lite Orange® Tint, 3.0 mL..... | 1   |
| Drape, fenestrated with tape .....                        | 1   |
| Gauze, 8-ply, 4 x 4 in. ....                              | 10  |
| Lidocaine 1% ampoule, 5 mL .....                          | 1   |
| Medication label sheet .....                              | 1   |
| Foam needle block.....                                    | 1   |
| Syringe, luer-lock, 5 mL .....                            | 1   |
| Towel, 17 x 19 in. ....                                   | 1   |
| Utility marker .....                                      | 1   |
| Wrap, 30 x 30 in. ....                                    | 1   |

### Conventional offering

| Basic diagnostic tray                                  | Qty |
|--|-----|
| Cat. no. <b>33-BA1CO</b>                               |     |
| BD Precision Glide® Needle, 22 g x 1 ½ in. ....        | 1   |
| BD Precision Glide® Needle, 25 g x 5/8 in. ....        | 1   |
| BD Precision Glide® Filter Needle, 19 g x 1 ½ in. .... | 1   |

### Safety offering

| Safety diagnostic tray                      | Qty |
|---|-----|
| Cat. no. <b>33-BA1SF</b>                    |     |
| BD SafetyGlide® Needle, 22 g x 1 ½ in. .... | 1   |
| BD SafetyGlide® Needle, 25 g x 5/8 in. .... | 1   |
| FILTER STRAW®, 5 Micron, 1 ¾ in. ....       | 1   |







BD, PRECISION GLIDE and SAFETYGLIDE are trademarks or registered trademarks of Becton Dickinson & Co.  
CURITY and TURKEL are registered trademarks of Covidien AG Corp.  
CHLORAPREP and HI-LITE ORANGE are registered trademarks of CareFusion, 2200 Inc.  
BARD-PARKER is a registered trademark of Aspen Surgical Products.  
FILTER STRAW is a registered trademark of B. Braun Medical Inc.

© 2013 Cardinal Health. All rights reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE and PRESOURCE are trademarks or registered trademarks of Cardinal Health. All other marks are the property of their respective owners. Lit. No. 2PRE12177 (08/2013)

[cardinalhealth.com/diagnostictrays](http://cardinalhealth.com/diagnostictrays)

Cardinal Health  
Presource® Products & Services  
7000 Cardinal Place  
Dublin, Ohio 43018