



CART-MATE
CABINETS™

You need your space, so make the most of it with Cart-Mate Cabinets by Clinton.

Because today's exam rooms are small, clinicians need to maximize the available space. Cart-Mate makes that possible by adding a mobile cart to the stationary cabinetry, making it easy to have needed instruments and accessories close at hand without cluttering the exam room space with extra furnishings. Just position the mobile cart for convenience during treatment and tuck it nicely away within the confines of the stationary cabinet when not in use.

Get affordable Cart-Mate cabinets today and save space—and money.

All Cart-Mate Cabinets Feature:

- Easy-clean surfaces
- Unit is pre-assembled to save installation time
- Fully adjustable, soft-close hinges
- 3.25" backsplash included with top
- Standard cart with 4-drawers, 3" swivel casters, plastic bumper and one-piece, vacuum-formed, gray-haircell, textured, top rimmed on all sides
- Smooth glide drawers with one-piece metal sides hold up to 75 lbs.
- Satin finish pulls
- 5-year limited warranty
- Optional -022 stainless steel sink with chrome faucet set available
- Optional door and drawer locks available
- Optional cart configurations
- Choice of Classic Clinton laminates or Fashion Finish surfaces



48" Cart-Mate Cabinet

DIMENSIONS: 48"W X 18"D X 35"H

- 2 doors and 2 adjustable shelves



48048L Shown with Optional -022 Sink



48048R Shown with Optional -022 Sink

60" Cart-Mate Cabinet

DIMENSIONS: 60"W X 18"D X 35"H

- 2 doors, 2 adjustable shelves, 2 drawers



48060L Shown with Optional -022 Sink



48060R Shown with Optional -022 Sink



48060ML Shown with Optional -022 Sink



48060MR Shown with Optional -022 Sink

! WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer and reproductive harm.

For more information go to www.P65Warnings.ca.gov.

! WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer and reproductive harm.

For more information go to www.P65Warnings.ca.gov.

66" Cart-Mate Cabinet

- 3 doors, 2 adjustable shelves, 2 drawers

DIMENSIONS: 66"W X 18"D X 35"H



48066L Shown with Optional -022 Sink



48066R Shown with Optional -022 Sink



48066MR Shown with Optional -022 Sink



48066ML Shown with Optional -022 Sink

72" Cart-Mate Cabinet

- 4 doors, 3 adjustable shelves, 2 drawers

DIMENSIONS: 72"W X 18"D X 35"H



48072L Shown with Optional -022 Sink



48072R Shown with Optional -022 Sink



48072MR Shown with Optional -022 Sink



48072ML Shown with Optional -022 Sink

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer and reproductive harm.

For more information go to www.P65Warnings.ca.gov.

72" Cart-Mate Split Cabinet

- 4 doors, 3 adjustable shelves, 2 drawers

DIMENSIONS: 72"W X 18"D X 35"H



48072SL Shown with Optional -022 Sink



48072SR Shown with Optional -022 Sink

72" Dual, Cart-Mate Cabinet

- 4 doors, 3 adjustable shelves, 2 drawers

DIMENSIONS: 72"W X 18"D X 35"H



48072D Shown with Optional -022 Sink



LAMINATE COLORS



1GR—GRAY



1SG—SLATE GRAY



1MP—MAPLE



1DC—DARK CHERRY

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer and reproductive harm.

For more information go to www.P65Warnings.ca.gov.

48" Fashion Finish, Cart-Mate Cabinet

DIMENSIONS: 48"W X 18"D X 35"H

- 2 doors and 2 adjustable shelves
- Laminate countertop



58048L Shown with Optional -022 Sink



58048R Shown with Optional -022 Sink

60" Fashion Finish, Cart-Mate Cabinet

DIMENSIONS: 60"W X 18"D X 35"H

- 2 doors, 2 adjustable shelves, 2 drawers
- Laminate countertop



58060L Shown with Optional -022 Sink



58060R Shown with Optional -022 Sink



58060ML Shown with Optional -022 Sink



58060MR Shown with Optional -022 Sink

! WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer and reproductive harm.

For more information go to www.P65Warnings.ca.gov.

66" Fashion Finish, Cart-Mate Cabinet

DIMENSIONS: 66"W X 18"D X 35"H

- 3 doors, 2 adjustable shelves, 2 drawers
- Laminate countertop



58066L Shown with Optional -022 Sink



58066R Shown with Optional -022 Sink



58066ML Shown with Optional -022 Sink



58066MR Shown with Optional -022 Sink

72" Fashion Finish, Cart-Mate Cabinet

DIMENSIONS: 72"W X 18"D X 35"H

- 4 doors, 3 adjustable shelves, 2 drawers
- Laminate countertop



58072L Shown with Optional -022 Sink



58072R Shown with Optional -022 Sink



58072ML Shown with Optional -022 Sink



58072MR Shown with Optional -022 Sink

! WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer and reproductive harm.

For more information go to www.P65Warnings.ca.gov.

72" Fashion Finish, Cart-Mate, Split Cabinet

DIMENSIONS: 72"W X 18"D X 35"H

- 4 doors and 3 adjustable shelves, 2 drawers
- Laminate countertop



58072SL Shown with Optional -022 Sink



58072SR Shown with Optional -022 Sink

72" Fashion Finish, Dual, Cart-Mate Cabinet

DIMENSIONS: 72"W X 18"D X 35"H

- 2 doors and 2 adjustable shelves
- Laminate countertop



58072D Shown with Optional -022 Sink

FASHION FINISH COLORS

NOTE: When selecting Fashion Finish, choose a countertop color to complement your Fashion Finish color.



1—ARCTIC WHITE



2—TWILIGHT



4—METROPOLIS GRAY

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer and reproductive harm.

For more information go to www.P65Warnings.ca.gov.

CART

DIMENSIONS: 20"W X 18.5"D X 29"H

Cart

- Cabinet mounted on four, 3" locking, swivel casters
- One-piece, vacuum-formed, gray-haircell, textured, top rimmed on all sides
- ClintonClean™ solid plastic bumper base
- Satin finish pulls
- Smooth glide drawers with one-piece metal sides hold up to 75 lbs.
- Finish matches cabinet selection
- Choose the configuration to fit your needs
 - Standard Configuration—4 drawers
 - Option 042—1 drawer, 1 door and 1 shelf
 - Option 043—2 drawers and 1 door



CART-MATE OPTIONS




042—1 drawer & 1 door



043—2 drawers & 1 door



022 Stainless Steel Sink and Gooseneck Chrome Faucet set




055 Door & Drawer Locks



066 Door Locks & Inside Latch

FASHION FINISH, LAMINATE COUNTERTOP COLORS



P1—White Carrera



P2—Meteorite



P3 Black Alicante

⚠ WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer and reproductive harm.

For more information go to www.P65Warnings.ca.gov.



See and download all of Clinton's product catalogs online.

- Medical Tables
- Medical Casework
- Power Tables
- Seating & Medical Accessories
- Phlebotomy Equipment
- Pediatric Equipment
- Physical Therapy Equipment

525 East Market St., York, PA 17403

800-441-9131 • 717-848-3519

Toll-Free Fax 866-522-1233

Fax 717-843-5871

www.clinton-ind.com